
**Page Denied**


Next 2 Page(s) In Document Denied

STAT

Vakuum pumpen  
LP

STAT

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


## **Aqua destillata**

**HIGHER QUALITY  
AND LESS EXPENSIVE**

In the past distilled water was used exclusively in scientific laboratories for the preparation of reagent solutions, for the production of medicaments, and similar purposes.

Today, wherever water free of

salts is adequate for the required purpose water is used purified by the ionexchange method, i. e. by contact with active substances only.

The most effective exchangers, however, can neither remove organic substances nor neutralize bacteria contained in the crude water when applying this method. Therefore, wherever water free of all inorganic and organic substances, or, in addition, also free of bacteria is required it is necessary to apply the classic method of distillation which does not mean, however, also the use of classic apparatus with metal vessels heated by a direct flame.

### **SIMAX MORE SUITABLE THAN METAL**

Today this method is replaced by more economic utilization of heat in electrically heated apparatus of completely new design made of high-quality glass and it is certainly not by chance that these highly efficient and most suitable apparatus are produced in the

land of the world is Bohemian glass - Czechoslovak SIMAX bo-  
rosilicate glass is the best material for the construction of distillation apparatus as it  
excels by its high chemical resistance (1st hydrolytic class) and also  
by its high thermal resistance.

### CAPACITIES ACCORDING TO REQUIREMENTS

The new distillation apparatus are compact types of original construction and  
are applied in the following capacities:  
of 10 l, 20 l, 50 l, 100 l of distilled water per hour.

### ABOVE ALL ECONOMY OF ENERGY

The double boiling vessel made of SIMAX glass is directly connected to the condenser.  
In the inner boiling vessel the electric heating coil is entirely submerged in water  
to ensure complete utilization of heat.

Through the outer boiling vessel vapour stream to the condenser (also made of SIMAX  
glass) and condensed there. The drops separator is of original design - a dephlegmator -  
and a special sludge trap systematically drains part of the water from the boiling vessel  
so that the concentration of salts impurities in the capacity of the apparatus remains within  
the limits ensuring economical operation of the apparatus.

A too low salt content would, however, prevent quick attainment of the full capacity of the  
apparatus and therefore when working with especially soft water it is necessary to add  
a small quantity of sodium chloride ( $\text{NaCl}_2$ ) to the boiling vessel before starting distillation.  
Any other substance increasing conductivity may also be added.

GLASS EXPOR  
STAT

## **AUTOMATIC OPERATION**

Compact design makes possible a very short connection of the condenser and feeding

preheated water from the condenser and the maintaining of a constant water level in the boiling vessel being safeguarded by regulating the overflow tube.


## **COMPLETE SAFETY**

The electrode heating elements are made of non-corrosive metals, apparatus of a capacity of 4 litres per hour being produced for single-phase current (they can be connected to the electric mains) and apparatus of greater capacities for three-phase current. To ensure safety of operation a special safety device is supplied as complementary equipment. This safety device switches off the current automatically as soon as all three phases are not under current. Thus the formation of foding is prevented which could reach high intensity during operation of the apparatus on two or one phase, leak into the water mains and seriously endanger the safety of the operators of the apparatus. The operating switch and glow discharge tube signalling the operation of the apparatus are fitted in a special connection. The entire device is housed in a metal case and forms the control box of the apparatus:


## **EASY ASSEMBLY**

The described distillation apparatus are of compact design and require a relatively small space. They are supplied with metal brackets for fixing to the wall at a height enabling easy access and easy manipulation with receiver vessels underneath. On special request the apparatus are also supplied with supporting structures; apparatus of 4 litres per hour capacity can thus be placed directly on the laboratory table. Assembly is very easy:

CLASSIFIED REPORT

STAT

the glass parts are joined by flanges. The feed and the boiler tubes of the assembled apparatus is attached by flanges.

### **PRACTICALLY NO ATTENTION REQUIRED**

Servicing is restricted merely to putting the apparatus into operation and to filling: the vessel is filled with water and the electric current is switched on (never switch the apparatus unless it is filled with water) and the water level in the boiling vessel does not reach the upper part of the electric heating elements. Then the inflow of condenser water is adjusted to ensure that the distillate flows out adequately cool (20 to 25° C). The temperature of the cooling water flowing out of the condenser should be about 70° C. After this adjustment the apparatus operates quite automatically, servicing being restricted only to periodic checks for correct operation and to rinsing the vessels after 8 to 10 hours of operation. Once a week the boiling vessel should be rinsed with a solution of hydrochloric or acetic acid to remove sediments. SIMAX glass is quite resistant to these

### **HIGH QUALITY OF DISTILLED WATER**

Apart from economic operation and the possibility of automatic control of distillation, the main advantage of these distillation apparatus is the high quality of the distillate with a very small entrainment of 2.5 mg up to 4.5 per 1 liter.

This high-quality distillate is obtained due to almost complete prevention of any overflow by a special dephlegmator and mainly due to the material used for the construction of the apparatus, viz. SIMAX glass, with its high chemical resistance. A further important element is the fact that as a result of the suitable design of the apparatus vapours and distillate practically do not come during distillation into contact with any foreign substances

which could contaminate the distilled water (and as far as apparatus of a 4 litres per hour capacity are concerned no such substances whatsoever are met).

### CONCLUSION

Finally it can be said that electrically heated SIMAX borosilicate distillation apparatus made in Czechoslovakia are highly efficient and very economical in operation, work automatically and, require practically no servicing, the distillate attaining almost the quality limits set for re-distilled water.

### TECHNICAL DATA


Type	4	10	30	50
Capacity (quantity of distillate in litres per hour)	4	10	30	50
Consumption of condenser water in litres per hour	50-70	80-100	300-350	450-450
Electric energy consumption in kW	2.4-2.8	6.8	13.8	19.5
Operating voltage in V	220	380	380	380
Fuse connection in Amp	15	25	35	50
Water content of boiling vessel in litres	1.2	4	4	8.5
Height of apparatus in cm	80	120	120	150
Width of apparatus in cm	25	40	40	50
Connecting cable (cross section) in sq. mm	3x2.5	3x4	3x10	3x16

GLASS REPORT

STAT

STAT

STAT

**Page Denied**

Next 14 Page(s) In Document Denied