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
MOSCOW

CONTENTS

Biological Microscopes	
Biological Microscope, Model M-10	Page 3
School Microscope, Model MY	7
Meat Inspection Microscope, Model MHC-7	9
Textile Microscope, Model MHC-9	11
Biological Microscope, Model MBH-1	13
Biological Research Microscope, Model MBH-2	18
Travelling Biological Microscope, Model MBH-4	23
Accessories for Biological Microscopes	
Binocular Attachment, Model AY-12	26
Universal Photomicrographic Attachment, Model MΦH-4	28
Micrometer Eyepiece, Model AM-9-2	30
Drawing Apparatus, Model PA-4	31
Phase-Contrast Equipment, Model KΦ-1	32
Mechanical Stage, Model CT-12	34
Normal Incident Illuminator, Model OH-1	35
Dark Ground Condenser, Model OH-13	37
Microscopy Lamp, Model OH-7	39
Polarizing Microscopes	
Polarizing Microscope, Model MHH-4	40
Petrological Microscope, Model MHH-5	44
Ore Microscope, Model MHH-6	46
Metallographic Microscopes	
Horizontal Metallographic Microscope, Model MHM-8	49
Vertical Metallographic Microscope, Model MHM-6	56
Electron Microscope, Model 3M-3	59

BIOLOGICAL MICROSCOPES

As a rule, microscopes designed for the examination of translucent objects, illuminated by transmitted light, are called biological microscopes.

Translucent specimens may be examined in either a bright or dark field.

The biological microscope may also be employed for the examination of opaque objects in low power work using an epicondenser.

The fields of application of biological microscopes are exceptionally wide. Botany, zoology, biochemistry, medicine, agriculture, industrial laboratories are a far from complete list of those branches of science and economy in which biological microscopes have found wide-spread use.

In accordance with the field of application and the type of research, the biological microscopes may be furnished with special accessories and devices enlarging their operating capacity.

BIOLOGICAL MICROSCOPE, MODEL M-10

The M-10 Microscope (Fig. 1) is a medium type biological microscope, designed for use in biological and other laboratories, clinics, universities and scientific research institutions.

This microscope satisfies the chief requirements of the biologist, physician or agronomist. Magnifications achieved by the various combinations of objectives and eyepieces range from 56x to 600x.

Fig. 2 shows a diagrammatical sectional view of the microscope.

The base 1 of the stand is of horse-shoe form and has three supporting pads which provide a stable position of the microscope on the table.

The weight of the base is such that the microscope is kept from upsetting even when the limb 2 is in a horizontal position.

To protect the stand from falling over from accidental side blows, two additional lugs are provided, underneath, between the main supporting pads.


Fig. 1
M-10 Biological Microscope

The vertical extension of the base is fitted into the slot of the limb 2. The latter has the form of a segment.

The limb is pivoted on the base with the pin 3 and can be inclined to any angle convenient for observation.

By the aid of a special wrench, furnished with the microscope, the effort required for tilting the limb can be regulated. In this way the desirable position of the body 4 in relation to the observer can be ensured.


Fig. 2
Diagrammatical sectional view of M-10 Microscope

The total travel of the body due to the fine focusing mechanism equals 2.4—2.5 mm (24—25 revolutions of the heads). This travel is limited by a special restricting device. The extreme positions of the body are determined by two scratches on the limb. The reference mark on the movable part of the body indicates its position in relation to the limb under the action of the fine focusing mechanism. The fine motion mechanism displaces the body together with the coarse focusing mechanism.

The coarse focusing mechanism comprises a rack 8 fastened to the body and a pinion 9 engaging the rack.

The stop screws 5 provide for setting the body in a precise vertical or horizontal position.

The cut-out portion in the center of the limb allows the latter to be used as a convenient handle for carrying the microscope from place to place.

The space, provided by this same cut-out portion, allows objects having a large lateral size (for instance, flat plates or pans) to be set up on the stage 6.

The upper part of the limb carries the fine focusing mechanism 7.

The fine focusing mechanism comprises a system of gears and a lever. It is actuated by two heads arranged on the right- and left-hand side of the limb. One revolution of the heads, on which shaft the driving gear of the mechanism is mounted, advances the body 0.1 mm.

The left-hand head is furnished with a graduated drum having 50 divisions on its circumference. One drum scale division corresponds to a body movement of 0.002 mm.

Heads 10 are mounted on the pinion shaft 9. Rotating the heads raises or lowers the body. One full revolution of the heads 10 provides a body travel of 20 mm. Rotating the heads 10 away from the observer lowers the body and vice versa.

The tube is of the telescopic type. By extending the draw tube 11, the mechanical tube length can be varied in a range from 150 to 200 mm. A scale inscribed on the draw tube indicates the regulated tube length.

Interchangeable eyepieces 12 are inserted into the upper part of the tube. The standard microscope set contains Huyghenian eyepieces, designed in such a manner, that a change of eyepiece power does not affect the focusing adjustment. The initial magnification is engraved on each eyepiece.

The lower part of the body carries the revolving nosepiece 13 for holding and rapidly changing the objectives 14. Rotating the body of the nosepiece brings one of the objectives in line with the tube opening. The correct position of the objective is fixed by a spring click.

The revolving nosepiece and its objective holes are aligned with the microscope axis with such accuracy that on changing the higher powers, the specimen set in the center of the field of view with the low power objective always remains in the field of view of the consequently replaced high power objective. Besides this, when changing from one objective to another, the specimen always remains visible and only a small adjustment of the fine focusing heads is necessary to obtain a sharp image.

The standard microscope set includes achromatic objectives computed for a mechanical tube length of 160 mm and a cover-glass thickness of 0.17 mm. Each objective is stored in a plastics case on which the initial magnification and numerical aperture of the objective are engraved.

An angle shaped bracket 15, carrying the rack 16, is fastened by screws to the tail of the limb. This rack engages a pinion 17 mounted in the movable part 18 of the bracket.

The movable part of the bracket or substage carries the spring fitting 19 into which the mount 20 of the two-lens condenser 21 and the swing-out holder 22 of the filter 23 are inserted.

The pinion 17 and the rack 16 provide for vertical adjustment of the substage together with the condenser.

The design of the condenser adjustment provides a smooth movement with sufficient friction to prevent the substage 18 from running down under its own weight.

The upward movement of the condenser is limited by a stop screw in such manner, that in the extreme upper position of the substage a space remains, between the front lens of the condenser and the object stage, to introduce an oil-immersion fluid used in some cases.

The two-lens condenser 21 has a maximum numerical aperture of 1.2. The condenser is furnished with an iris diaphragm and operates in conjunction with the plain and concave mirror 24.

The fork-shaped mirror holder 25 allows either mirror surface to be directed toward the source of light.

The object stage 6 of the microscope is fastened to the stationary part of bracket 15. The stage is circular and its upper part has rotary and centring movements actuated by special screws.

The stage clips 26, inserted into holes provided in the upper part of the stage, serve to clamp the specimen.

If required, the condenser in its mount 20 may be removed from the spring fitting 19 and replaced by some other condenser or by cylindrical diaphragms.

The microscope is stored in a wooden case having the form of a cabinet with locking door and a handle on top for carrying.

There are sliding holders in the case for storing objectives and eyepieces, as well as a box for clips, filters, spare wrench and other accessories.

SPECIFICATIONS

Range of total magnification	from 56 \times to 600 \times
Achromatic objectives:	8 \times 40 \times
initial magnification	8 \times 40 \times
numerical aperture	0.20 0.65
focal length, mm	18.2 4.35
working distance, mm	8.91 0.6
field of view with 10 \times eyepiece, mm	1.75 0.35
Huyghenian eyepieces:	7 \times 10 \times 15 \times
magnification	7 \times 10 \times 15 \times
focal length, mm	36 25 17
linear field of view, mm	18 14 8
Objective changer	double revolving nosepiece
Body tube variable length, mm	from 150 to 200
Focusing adjustments:	
coarse	by rack and pinion
fine	micrometric mechanism, reading to 0.002 mm
Illuminating system	reversible mirror plane and concave two-lens condenser 1.2 N.A. with iris diaphragm and interchangeable filter
Object stage	circular, with centring and rotating adjustments
Overall dimensions of microscope (height \times length \times width), mm	315 \times 190 \times 170
Overall dimensions of case, mm	350 \times 250 \times 195
Weight of microscope, kg	7.9
Weight of microscope in case, kg	11.8

MICROSCOPE SET

- Microscope stand comprising horse-shoe base, limb, body with coarse and fine focusing adjustments, revolving nosepiece, circular rotary object stage and illuminating system
- Achromatic objective, 8 \times \times 0.20, in case
- Achromatic objective, 40 \times \times 0.65, in case
- Huyghenian eyepiece, 7 \times
- Huyghenian eyepiece, 10 \times
- Huyghenian eyepiece, 15 \times
- Opal glass
- Blue filter
- Clips for holding specimens (2 pieces)
- Wrenches for adjustment (2 pieces)
- Flannel napkin
- Squirrel-hair brush
- Case for microscope
- Box for microscope accessories
- Description and instruction manual
- Certificate

SCHOOL MICROSCOPE, MODEL MV

The MV School Microscope (Fig. 3) is a medium type biological microscope constructed on the basis of the M-10 microscope.

The School Microscope is designed for the use of students of biological, medical and agricultural institutions as well as for studies in public and technical schools. Model MV is considerably simplified in comparison with model M-10.

The revolving nosepiece for changing objectives is absent in the simplified model. A bushing, of the same height as the revolving nosepiece, is screwed into the body instead of it. Either of the two objectives, furnished with the microscope, can be screwed into the bushing.

A circular stationary object stage with two clips for holding specimens is substituted on the school microscope for the rotary stage. The sub-stage vertical adjustment is absent on this microscope and the condenser has been replaced by interchangeable cylindrical diaphragms of various diameter.

Two achromatic objectives and two Huyghenian eyepieces are furnished with the instrument.

The microscope is stored in a wooden case having the form of a cabinet with locking door and a handle on top for carrying.

If desired, various accessories may be furnished on special order with the school microscope to provide more universal application.

Experience has shown that the school microscope meets the requirements of schools and colleges and provides the possibility of familiarizing students and pupils with the principles of microscopy and the chief methods of its procedure.


Fig. 3
MV School microscope

SPECIFICATIONS

Range of total magnification	from 80 \times to 600 \times
Achromatic objectives	8 \times and 40 \times (optical data same as for Model M-10 objectives)
Huyghenian eyepieces	10 \times and 15 \times (optical data same as for Model M-10 eyepieces)
Object stage	circular, stationary, with clips
Focusing adjustments:	
coarse	by rack and pinion
fine	micrometric mechanism without graduated drum
Illuminating system	reversible mirror plane and concave and three interchangeable diaphragms of 1, 3 and 6 mm diameter
Overall dimensions of microscope (height \times length \times width), mm	315 \times 190 \times 170

Overall dimensions of case, mm 350 x 250 x 195
 Weight of microscope, kg 7.75
 Weight of microscope in case, kg 11.65

MICROSCOPE SET

- Microscope stand comprising horse-shoe base, limb, body with coarse and fine focusing adjustments, stationary object stage with clips and mirror
- Achromatic objective, 8x x 0.20, in case
- Achromatic objective, 40x x 0.65, in case
- Huyghenian eyepiece, 10x
- Huyghenian eyepiece, 15x
- Interchangeable diaphragms (3 pieces)
- Clips for holding specimens (2 pieces)
- Wrench for adjustment
- Flannel napkin
- Squirrel-hair brush
- Case for microscope
- Box for microscope accessories
- Description and instruction manual
- Certificate

MEAT INSPECTION MICROSCOPE, MODEL MHC-7

The MHC-7 Meat Inspection Microscope (Fig. 4) is a medium type biological microscope constructed on the basis of the M-10 microscope.

The Meat Inspection Microscope is designed for the inspection of meat to determine the presence of trichina and finds wide application in the food-stuffs industry and in sanitary inspection.

In comparison with the M-10 microscope the Meat Inspection Microscope has a number of simplifications in design concerning the object stage, limb, body, revolving nosepiece and illuminating system.

The instrument has a straight body of invariable tube length. The body is supplied with a coarse focusing mechanism. The revolving nosepiece is replaced by a bushing into which the objectives are screwed.

The illuminating system — without condenser and adjustable substage — comprises a reversible plane and concave mirror fastened on a pivoted fork-shaped holder and three interchangeable cylindrical diaphragms of various diameter.

The stationary object stage is rectangular in form and the compressing slide (Fig. 5) is moved about on its surface.

The compressing slide of the Meat Inspection Microscope comprises two plates of thick glass fastened together by two screws. It serves to flatten out the specimens of meat and sinew which are being examined for the presence of trichina.

The compressing slide can be moved by hand on the stage along a straight-edge. This allows all of the specimens, between the glass plates, to be examined. Rectangles with numerals from 1 to 28 are engraved on both the lower and upper plates. The inspected specimens are arranged inside these rectangles.

The microscope is furnished with two achromatic objectives and one Huyghenian eyepiece.

The microscope is stored in a wooden case having the form of a cabinet with locking door and a handle on top for carrying.


Fig. 4
 MHC-7 Meat Inspection
 Microscope


Fig. 5
 Compressing Slide

SPECIFICATIONS

Range of total magnification	from 25.9 \times to 56 \times
Achromatic objectives:	3.7 \times 8 \times
initial magnification	3.7 \times 8 \times
numerical aperture	0.11 0.20
focal length, mm	33.1 18.2
working distance, mm	27.7 6.91
Huyghenian eyepiece	7 \times (optical data same as for Model M-10 eyepiece)
Object stage	stationary, rectangular, size: 100 \times 170 mm
Focusing adjustment	by rack and pinion
Illuminating system	reversible plane and concave mirror and three interchangeable diaphragms of 1, 3 and 6 mm diameter
Size of compressing slide, mm	50 \times 220
Overall dimensions of microscope (height \times length \times width), mm	315 \times 190 \times 170
Overall dimensions of case, mm	350 \times 250 \times 195
Weight of microscope, kg	7.55
Weight of microscope in case, kg	11.65

MICROSCOPE SET

- Microscope stand comprising horse-shoe base, limb, body with coarse focusing adjustment, rectangular object stage with straight-edge and mirror
- Achromatic objective, 3.7 \times \times 0.11, in case
- Achromatic objective, 8 \times \times 0.20, in case
- Huyghenian eyepiece, 7 \times
- Interchangeable diaphragms (3 pieces)
- Clips for holding specimens (2 pieces)
- Compressing slide in case
- Wrench for adjustment
- Flannel napkin
- Squirrel-hair brush
- Case for microscope
- Box for microscope accessories
- Description and instruction manual
- Certificate

TEXTILE MICROSCOPE, MODEL MHC-9

The MHC-9 Textile Microscope (Fig. 6) is a medium type biological microscope constructed on the basis of the M-10 microscope.

The Textile Microscope has found wide application in agricultural laboratories and in inspection stations of cotton manufacturing and textile factories.

The Textile Microscope differs from the M-10 model, in that, it has a simplified body, object stage and illuminating system.

The instrument has a straight body of invariable tube length. The microscope is furnished with a coarse focusing mechanism. The circular object stage is of the stationary type.

The revolving nosepiece is replaced by a bushing into which objectives are screwed.

Instead of a condenser, this microscope uses interchangeable cylindrical diaphragms. Two achromatic objectives and two eyepieces are furnished with the microscope.

The microscope is provided with a mechanical stage for shifting the inspected cotton fibre in a longitudinal direction.

The microscope is stored in a wooden case having the form of a cabinet with locking door and a handle on the top for carrying.


Fig. 6
MHC-9 Textile Microscope

SPECIFICATIONS

Range of total magnification	from 56 \times to 300 \times
Achromatic objectives:	8 \times 20 \times
initial magnification	8 \times 20 \times
numerical aperture	0.20 0.40
Huyghenian eyepieces	7 \times and 15 \times (optical data same as for Model M-10 eyepieces)
Object stage	circular, stationary type
Illuminating system	reversible plane and concave mirror and three interchangeable diaphragms, of 1, 3 and 6 mm diameter
Overall dimensions of microscope (height \times length \times width), mm	315 \times 190 \times 170
Overall dimensions of case, mm	350 \times 250 \times 195
Weight of microscope, kg	7.8
Weight of microscope in case, kg	11.5

Transverse movement of mechanical stage, mm	75
Reading accuracy of mechanical stage scale, mm	0.1
Overall dimensions of mechanical stage case, mm	170 × 121 × 50
Weight of mechanical stage, kg	0.15

MICROSCOPE SET

Microscope stand comprising horse-shoe base, limb, body with coarse focusing adjustment, circular stationary object stage and mirror
 Achromatic objective, 8× × 0.20, in case
 Achromatic objective, 20× × 0.40, in case
 Huyghenian eyepiece, 7×
 Huyghenian eyepiece, 15×
 Interchangeable diaphragms (3 pieces)
 Mechanical stage in case
 Wrenches for adjustment (2 pieces)
 Clips for holding specimens (2 pieces)
 Flannel napkin
 Squirrel-hair brush
 Case for microscope
 Box for microscope accessories
 Description and instruction manual
 Certificate

BIOLOGICAL MICROSCOPE, MODEL MBH-1

The MBH-1 Microscope (Fig. 7) is a medium type biological microscope and is a further development and an improvement of the M-10 microscope.

This instrument provides for research microscopy with a magnification from 56× to 1350×; the latter value being the maximum useful power for optical microscopes.

A further increase in magnification is of no avail and even detrimental as it leads to a decrease in illumination without gain in raising the limit of the resolving power of the microscope, that is, the possibility of examining more minute details of the specimen.

The MBH-1 Microscope differs from the M-10 microscope chiefly in its form which possesses a number of distinct features in comparison with the earlier arrangement. It is lower in height than the M-10 model and has an inclined eyepiece tube which allows the observer to be comfortably seated at a table during operation.

The object stage is always horizontal. This is necessary when conducting research on liquids.

The coarse and fine focusing mechanisms are arranged at the lower end of the limb. This allows both arms to rest on the table during microscopy and decreases the fatigue of the observer.

The optical system (Fig. 8) of the MBH-1 Microscope differs from the M-10 model in that a prism 1, changing the path of the rays to an angle of 45° to the horizontal plane, is arranged between the objective and the eyepiece.

The main parts of the instrument are: base 2 of the stand, housing with fine focusing mechanism 3, object stage 4, limb 5, revolving nosepiece 6 on guides, inclined monocular body 7, condenser substage bracket 8, condenser 9, objectives 10 and eyepiece 11.

The base of the stand has a horse-shoe form and imports high stability to the microscope.

The housing 3 is fastened by screws to the base. On one side it has guides for the condenser substage bracket 8 and on the other side, guides for the limb 5. The fine focusing mechanism is mounted inside the housing.

One revolution of the fine focusing head on whose shaft the driving gear of the mechanism is mounted, advances the body 0.1 mm. The total adjustment of the body, by means of the fine focusing mechanism, equals 2.2—2.4 mm.

The extreme positions of the body are determined by scratches on the base. A reference mark on the movable part and two scratches on the stationary part


Fig. 7
MBH-1 Biological Microscope

correspond to the extreme positions of the fine focusing adjustment. The fine focusing mechanism displaces the limb together with the coarse focusing mechanism.

The limb has the form of an arc. This allows large-sized objects to be placed on the microscope stage and also facilitates carrying the microscope from place to place.

The lower end of the limb has guides and a fitting with two heads 12, which serve for coarse focusing adjustment.


Fig. 8
Diagrammatical sectional view of the MBH-1 Microscope

The dimensions of the guides have been selected to provide a body travel of 50mm. One full revolution of the heads 12 corresponds to a movement of 20 mm.

The heads 12 have been designed so that, rotating one head in relation to the other regulates the smoothness of the movement to suit the observer.

The mounting 13 is fastened to the upper end of the limb. It has dovetail guides for the revolving nosepiece 6 and a socket for fastening the inclined monocular tube 7. A straight tube or a binocular attachment, the latter not included in the standard microscope set, but furnished on special order, can also be fastened in this mounting.

The inclined monocular tube can be rotated about a vertical axis and fixed in any position desired by the observer.

14

The lower end of the body comprises a spherical housing in which the prism 1 is arranged. The tube, in which the eyepiece 11 is inserted, is screwed into this housing in inclined position.

All the components of the inclined (or straight) tube are designed to provide a mechanical tube length of 160 mm.

The standard set contains Huyghenian eyepieces designed in a series so that a change in eyepiece power does not affect the focusing adjustments. The initial magnification is engraved on each eyepiece.

The straight tube, necessary for microphotography and other work, is included in the microscope set.

The revolving nosepiece, for holding and rapidly changing objectives, has four threaded holes for inserting objectives arranged in a spherical dish. The upper part of the revolving nosepiece is provided with dovetail guides for attachment to the mounting 13 on the limb.

The correct position of the revolving nosepiece in reference to the body axis is fixed by the screw 14. The spherical dish of this nosepiece can be rotated about its axis so that any of the four holes with its objective can be aligned with the body axis. A spring click device inside the revolving nosepiece centers each objective in reference to the optical axis of the microscope body. The alignment accuracy is such that on changing to higher powers the specimen set in the center of the field of view with the low power-objective always remains in the field of view of the consequently replaced high power objective. When changing from one objective to another, the specimen always remains visible and only a small adjustment of the fine focusing head is necessary to obtain a sharp image.

The only exception, to the above, is the oil-immersion objective with a magnification of 90x. It has a working distance somewhat higher than for dry systems.

The standard set of the instrument includes three achromatic objectives computed for a mechanical tube length of 160 mm and a cover glass thickness of 0.17 mm. Each objective is stored in a plastic case to protect it from dust. The initial magnification and the numerical aperture are engraved on the mount and on the bottom of the case of each objective.

The condenser substage bracket 8 can be adjusted up to 20 mm by a rack and pinion. The substage carries a cylindrical spring fitting into which the condenser 9, in its mount, is fastened by means of screw 15.

A washer with two holes is arranged on the right-hand side of the condenser adjustment mechanism. Rotating this washer with a special wrench, provided with the microscope, adjusts the friction of the substage movement to prevent it from running down under its own weight. This arrangement is especially important when using a heavy condenser with a phase contrast device.

The two-lens condenser, with a numerical aperture of 1.2 is provided with iris diaphragm and a swing-out holder 16 for the light filter. The condenser operates on conjunction with the plane and concave mirror 17.

The fork-shaped mirror holder is reversible and allows either mirror surface to be directed toward the source of light and be properly positioned.

The upper front lens of the condenser may be removed. This decreases the numerical aperture to 0.5 and is necessary when working with low power objectives as, for instance, with an objective 8x.

15

The condenser aperture is equal to 1.2 only when oil-immersion is used between the front lens and the object glass. Without oil-immersion, the condenser aperture is approximately unity.

The upward movement of the condenser substage is limited by a stop so that in the extreme upper position, a clearance of about 0.1 mm remains between the condenser front lens and the cover glass.

By special order a dark ground condenser OII-13 may be supplied with the microscope. It can be inserted in the fitting instead of the standard type.

The object stage is fastened on a special bracket. The latter, in turn, is fastened to the housing 3. The upper part of the stage may be rotated by means of the knurled ferrule. Besides this, the stage may be adjusted 8 mm in the longitudinal and cross directions by means of two screws and a spring to bring any point of the specimen to the center of the field of view.

There are seven threaded holes on the top of the object stage. The four central holes serve for fastening the spring clips 18 for holding the specimen while the three holes at the sides are for fastening a superimposed mechanical stage, not included in the standard microscope set, but furnished on special order.

The microscope is stored in a wooden case having the form of a cabinet with locking door and a handle on top for carrying.

There are sliding holders in the case for storing objectives and eyepieces, as well as a box for microscope accessories.

SPECIFICATIONS

Range of total magnification	from 56 \times to 1350 \times		
Achromatic objectives:	8 \times	40 \times	90 \times
initial magnification	8 \times	40 \times	90 \times
numerical aperture	0.20	0.65	1.25
focal length, mm	18.2	4.35	1.96
working distance, mm	8.91	0.61	0.19
field of view with 10 \times eyepiece, mm	1.75	0.35	0.15
Huyghenian eyepieces:	7 \times	10 \times	15 \times
magnification	7 \times	10 \times	15 \times
focal length, mm	36	25	17
linear field of view, mm	18	14	8
Objective changer	quadruple revolving nosepiece		
Number of body tubes	two (straight and inclined monocular tubes)		
Focusing adjustments:			
coarse	by rack and pinion		
fine	micrometric mechanism, reading to 0.002 mm		
Illuminating system	reversible plane and concave mirror, two-lens condenser 1.2 N.A. with iris diaphragm and interchangeable filter		
Object stage	circular with centring and rotating adjustments		
Overall dimensions of microscope (height \times length \times width), mm	285 \times 210 \times 190		
Overall dimensions of case, mm	365 \times 200 \times 245		
Weight of microscope, kg	4.6		
Weight of microscope in case, kg	9.5		

16

MICROSCOPE SET

Microscope stand comprising horse-shoe base, limb, inclined monocular body tube, coarse and fine focusing mechanisms, revolving nosepiece, circular rotary object stage and illuminating system

Achromatic objective, 8 \times \times 0.20, in case

Achromatic objective, 40 \times \times 0.65, in case

Achromatic objective, 90 \times \times 1.25 (oil-immersion), in case

Huyghenian eyepiece, 7 \times

Huyghenian eyepiece, 10 \times

Huyghenian eyepiece, 15 \times

1.2 N.A. Aplanatic condenser

Straight tube

Opal glass

Blue filter

Clips for holding specimens (2 pieces)

Flannel napkin

Squirrel-hair brush

Wrenches (2 pieces)

Case for microscope

Vial with immersion, oil in case

Box for microscope accessories

Sliding holders for storing objectives and eyepieces

Description and instruction manual

Certificate

17

BIOLOGICAL RESEARCH MICROSCOPE, MODEL MBH-2

The MBH-2 Biological Microscope (Fig. 9) is an ideal instrument for detailed and comprehensive research work.

This microscope has found wide application in biological, bacteriological, biochemical, medical and other scientific research institutions.

The MBH-2 Microscope has an inclined eyepiece tube which allows the observer to work in a comfortable seated position.

The low position of the coarse and fine focusing mechanism is especially convenient during prolonged research work when it is necessary to focus on some layer of the specimen for a long time.

The corresponding position of the coarse and fine focusing heads allows the arms to rest on the table during microscopy and decreases the fatigue of the observer.

The main binocular tube, furnished with the microscope, provides for observation with both eyes. This decreases eye fatigue during prolonged observation by creating more natural working conditions. Besides this, the binocular tube provides for some stereoscopic relief which increases the research potentialities of the microscope.

The MBH-2 Microscope with a binocular tube produces an upright image of the object. An inclined monocular tube or a straight tube intended for photographic purposes may be substituted for the inclined binocular tube. All these tubes are interchangeable and may be attached without preliminary fitting operations.

The microscope has a highly perfected illuminating system mounted in the base of the instrument.

The condensers are arranged in a triple revolving substage device operating in conjunction with the pancreatic system. The pancreatic system allows the numerical aperture of the condenser to be varied from 0.16 to 1.4.

Changing condensers by means of the revolving substage device has considerable advantages over the usual method of changing condensers by inserting them into the spring fitting.

Microscopy may be carried out with transmitted light on either a bright or dark ground. The illuminating system provides sufficient light for dark ground operation or for microphotography.

The object stage of the microscope is always horizontal which is convenient for examining liquids.

This microscope is furnished with a set of high-quality apochromatic objectives and compensating eyepieces.


Fig. 9
MBH-2 Biological Research
Microscope

The optical system of the microscope (Fig. 10) comprises the light source — an electric lamp 1, lens 2, aperture unit consisting of prism 3, lens 4 and aperture iris diaphragm 5, pancreatic system 6, interchangeable condensers 7, objective 8, prism 9 and eyepiece 10.

The specimen 11 is placed between the object and cover glasses arranged between the condenser 7 and the objective 8.

The base 1 of the microscope (Fig. 11) has the form of a circle with cut-off segments.

The bracket-housing 2 is fastened to the base. On one side it has guides for the condenser substage bracket 3 and on the other side, guides for the limb 4.

The fine focusing mechanism is mounted inside the housing. It is actuated by rotating the heads 5 arranged on the right- and left-hand sides. On the right-hand


Fig. 10
Optical System of the MBH-2 Microscope

Fig. 11
Main parts of the MBH-2 Microscope

side a graduated drum with 50 divisions is mounted on the shaft of the heads. One revolution of the drum corresponds to a body adjustment of 0.1 mm. The drum is graduated 1 div. = 0.002 mm.

The total movement of the body from the fine focusing mechanism equals 2.3—2.5 mm and is determined by two scratches on the right-hand side of the housing. A reference mark on the movable part of the housing and two scratches on the stationary part correspond to the extreme positions of the body during adjustment. The fine focusing mechanism displaces the limb together with the coarse focusing mechanism.

The limb has the form of an arc. This allows large sized objects to be placed on the microscope stage and also facilitates carrying the microscope from place to place.

The coarse focusing mechanism, actuated by revolving the heads 6, is mounted in the lower part of the limb. The design of the coarse focusing mechanism allows the smoothness of the movement to be regulated to suit the observer. This mechanism provides a vertical body movement of 50 mm (two and one half turns of the heads).

The mounting 7, with dovetail guides, is fastened to the upper end of the limb.

The revolving nosepiece 8 is inserted into these guides and properly positioned by screw 9. The revolving nosepiece serves for holding and rapidly changing objectives. It has four threaded holes for inserting objectives arranged in a spherical dish. The rotation of the spherical dish is indexed in the four positions by a spring device arranged inside the nosepiece. The alignment accuracy is such that on changing to higher powers the specimen set in the center of the field of view, with the low power objective always remains in the field of view of the consequently replaced high power objective. When changing from one objective to another, the specimen always remains visible and only a small adjustment of the fine focusing head is necessary to obtain a sharp image.

The only exception, to the above, is the oil-immersion objective with a magnification of 90x. It has a working distance somewhat higher than for dry systems.

The standard set of the instrument includes four apochromatic objectives computed for a mechanical tube length of 160 mm and a cover glass thickness of 0.17 mm.

The oil-immersion objective, 60x x 1.0, is furnished with an iris diaphragm for varying the objective aperture. This objective is used, chiefly, in conjunction with cardioid-condenser. The aperture of the objective is changed by rotating a knurled ring on the mount of the objective.

The initial magnification and the numerical aperture are engraved on the mount and on the bottom of the storage cases of each objective.

The mounting 7 has a socket on top for fastening the interchangeable body tubes. An annular tapered neck at the base of each tube enters the socket of the mounting.

The inclined monocular and binocular tubes 10 can be rotated about a vertical axis and fixed in any position by means of screw 11 on the right-hand side of the limb mounting.

The mechanical tube length of any of the MBH-2 Microscope body tubes equals 160 mm.

The length of the straight tube is variable and is adjusted by extending a draw-tube to a reading taken on the scale on the stationary section of the straight tube (Fig. 12).

Five types of compensating eyepieces are included in the standard microscope set. The eyepieces are designed so that after interchanging them, only a small adjustment of the fine focusing heads is required to obtain a sharp image. The type and magnification of the eyepieces are engraved on the mount of the eye lens.

By using suitable objectives and eyepieces furnished with the microscope total magnifications from 75x to 1350x can be obtained with the binocular body (including the initial magnification of the tube equal to 1.5x) and a magnification from 50x to 1350x with the monocular body tube.

The microscope is furnished with a rotary centring object stage having a mechanism for crosswise movement of the specimen. The stage is moved in a longitudinal

direction by rotating the screw 12 (Fig. 11). Cross movement is accomplished by aid of the mechanical stage 13, fastened to the movable part of the object stage. It is actuated by rotating the screw 14.

The amount of movement in either direction can be read on the scales and verniers to an accuracy of 0.1 mm. The upper part of the object stage can be rotated by releasing screw 15 on the left-hand side of the stage.

The illuminating system is mounted in a cylindrical recess inside the base 1. The lampholder and bulb are inserted and aligned in an eccentric bushing. The aperture unit is aligned by means of screws 16.

To obtain oblique illumination in any direction, the aperture diaphragm can be eccentrically positioned by means of screw 17.

The opening of the diaphragm is adjusted by means of the ring 18 which, at the same time, serves as a mount for the interchangeable light filters.

The panoramic system and the revolving condenser substage are arranged on the bracket 3 and fastened by the screw 19. This bracket is raised and lowered by rotating the head 20.

The condensers are mounted in a detachable fixture fastened to the bracket.

The following condensers are held in the revolving substage device: aplanatic for apertures from 0.16 to 1.4, condenser for low power objectives and a cardioid-condenser for dark ground operation.

The value of the aperture set on the aplanatic condenser can be read on a scale engraved on ring 21 whose rotation changes the aperture.

The microscope is stored in a wooden case having the form of a cabinet with locking door and a handle on top for carrying.


Fig. 12
The straight body tube and its installation in the limb mounting of the MBH-2 Microscope

SPECIFICATIONS

Range of total magnification:					
with binocular body	from 75x to 1350x				
with monocular body tube	from 50x to 1350x				
Apochromatic objectives:	10x	20x	60x	90x	
initial magnification	10x	20x	60x	90x	
numerical aperture	0.30	0.65	1.0	1.30	
focal length, mm	16.1	8.43	3.0	2.0	
working distance, mm	5.17	0.94	0.44	0.32	
field of view with 10x eyepiece, mm	1.3	0.65	0.22	0.15	
Compensating eyepieces:	5x	7x	10x	15x	20x
magnification	5x	7x	10x	15x	20x
focal length, mm	49.8	36	25.2	16.8	12.6
linear field of view, mm	23	18	13	12	9

Illuminating system	electric bulb 8V, 20W; aperture unit with iris diaphragm; paracentric system and three interchangeable condensers: aplanatic with variable aperture from 0.10 to 1.4, low power objective and cardioid-condenser
Objective changer	quadruple revolving nosepiece
Number of body tubes	two (straight monocular and inclined binocular)
Focusing adjustments:	
coarse	by rack and pinion
fine	micrometric mechanism, reading to 0.002 mm
Range of condenser substage adjustment, mm	20
Object stage	circular with centring and rotating adjustments, with crosswise movement, reading to 0.1 mm
Overall dimensions of microscope (height x length x width), mm	370 x 260 x 165
Overall dimensions of case, mm	400 x 235 x 280
Weight of microscope, kg	7.3
Weight of microscope in case, kg	13.7

MICROSCOPE SET

Microscope stand comprising a box-shaped foot containing the illuminating equipment, limb, coarse and fine focusing mechanisms, body tube mounting and crosswise type object stage, quadruple revolving nosepiece for objectives

Binocular attachment, type AY-12

Revolving substage with three condensers

Straight monocular body tube

Apochromatic objective, 10x x 0.30, in case

Apochromatic objective, 20x x 0.65, in case

Apochromatic objective 60x (1.0-0.7 oil-immersion) with iris diaphragm, in case

Apochromatic objective, 90x x 1.3 (oil-immersion), in case

Compensating eyepieces, 5x (2 pieces)

Compensating eyepieces, 7x (2 pieces)

Compensating eyepieces, 10x (2 pieces)

Compensating eyepiece, 15x

Compensating eyepiece, 20x

Opal glass

Green filter

Blue filter

Smoke-coloured filter

Lampholder with bulb and plug

Spare bulbs (8V, 20W - 2 pieces)

Transformer 127/220/8V with rheostat

Set of wrenches and screwdrivers (4 pieces)

Protecting caps (2 pieces)

Vial with immersion oil, in case

Flannel napkin

Squirrel-hair brush

Case for microscope

Description and instruction manual

Certificate

TRAVELLING BIOLOGICAL MICROSCOPE, MODEL MBH-4

The MBH-4 Microscope (Fig. 13) is an up-to-date in design and is furnished with a similar set of objectives and eyepieces as the MBM-1 Microscope. It is, however, somewhat different in form and in the size of the stand.

The Travelling Microscope is designed for use in field and travelling conditions. It can be used for the examination of translucent objects with transmitted light

on a bright or dark ground. By means of a normal incident illuminator, the microscope can be used for examining opaque objects with reflected light.

The Travelling Microscope is unsurpassed for research work in expeditions for botanists, zoologists, bacteriologists and other scientific workers.

This microscope has a small height and the coarse and fine focusing heads are arranged on the lower part of the microscope, so that the arms may be rested on the table during operation.

In the optical system (Fig. 14) of the microscope, a prism 3, changing the path of the rays to an angle of 45° to the horizontal plane is arranged between the objective 1 and the eyepiece 2. The prism 3 is designed so that the mechanical tube length remains 160 mm.

The base 4 of the instrument is a rectangular plate with three supporting pads and the screw 5 which serves to keep the microscope stable on an uneven surface.

The intermediate housing 6 is a rectangular parallelepiped screwed to the base.

Two sides of the housing 6 carry guides; on one side for the limb 7 and on the other side for the condenser substage bracket 8.

The fine focusing mechanism is mounted inside the housing. It is actuated by rotating the heads 9 arranged on the right- and left-hand side of the housing. A drum with a scale on the left-hand head allows the body to be adjusted with an accuracy of 0.002 mm.

Coarse focusing is provided by a rack and pinion actuated by rotating the heads 10.

The limb has an arched form and, on its upper end, carries the mounting 11 with dovetail guides for fastening the revolving nosepiece 12.

The revolving nosepiece for holding and rapidly changing objectives has three threaded holes for inserting objectives and a spring click device which assures proper centring of each of the objectives with the optical axis of the body.

The mounting 11 has a socket on top for inserting the inclined monocular tube 13.

The inclined tube may be rotated to any position about a vertical axis and can be fixed in the required position by the screw 14.

The condenser substage has a spring fitting for holding the condenser 15. This

Fig. 13
MBH-4 Travelling Biological
Microscope

substage, together with the inserted condenser, may be adjusted in height by means of the rack and pinion 16.

Three achromatic objectives and three Huyghenian eyepieces are furnished with the microscope. Magnifications achieved by the various combinations of objectives and eyepieces range from 56 \times to 1350 \times .


Fig. 14
Diagrammatical sectional view of the MBH-4 Microscope

A two-lens aplanatic condenser, with an aperture of 1.2, having an iris diaphragm and a swing-out mount for the light filter 17, is furnished with the microscope. The condenser operates in conjunction with the mirror 12.

The microscope with all the necessary accessories is stored in a convenient travelling case with a handle.

SPECIFICATIONS

Range of total magnification	from 56 \times to 1350 \times
Achromatic objectives	8 \times \times 0.20, 40 \times \times 0.65 and 90 \times \times 1.25 (optical data same as for Microscope MBH-1)

Huyghenian eyepieces	7 \times , 10 \times and 15 \times (optical data same as for Microscope MBH-1)
Objective changer	triple revolving nosepiece
Focusing adjustments:	
course	by rack and pinion
fine	micrometric mechanism, reading to 0.002 mm
Illuminating system	reversible mirror plane and concave; 1.2 N.A. two-lens condenser with iris diaphragm and interchangeable filter
Object stage	rectangular
Overall dimensions of microscope (height \times length \times width), mm	280 \times 213 \times 92
Overall dimensions of case, mm	268 \times 225 \times 115
Weight of microscope, kg	3.6
Weight of microscope in case, kg	5.7

MICROSCOPE SET

- Microscope stand comprising rectangular base, limb, revolving nosepiece, object stage and illuminating system
- Inclined monocular tube
- Achromatic objective, 8 \times \times 0.20, in case
- Achromatic objective, 40 \times \times 0.65, in case
- Achromatic objective, 90 \times \times 1.25 (oil-immersion), in case
- Huyghenian eyepiece, 7 \times
- Huyghenian eyepiece, 10 \times
- Huyghenian eyepiece, 15 \times
- 1.2 N.A. aplanatic condenser
- Opal glass
- Blue filter
- Clips for holding specimens (2 pieces)
- Tanned napkin
- Squirrel-hair brush
- Wrenches (2 pieces)
- Vial with immersion oil, in case
- Travelling case
- Description and instruction manual
- Certificate

ACCESSORIES FOR BIOLOGICAL MICROSCOPES

BINOCULAR ATTACHMENT, MODEL AV-12

When working with a microscope having a monocular tube, only one eye is used in observation. This differs from the natural conditions of using both eyes and sustained observation leads to eye fatigue.

The AV-12 Binocular Attachment (Fig. 15) allows the observer to use both eyes, providing more natural conditions and decreasing eye fatigue. Besides this, the resolving power of the microscope is increased by use of the Binocular Attachment.

The Binocular Attachment operates in the following manner:

After passing through the microscope objective, rays of light fall on the lens 1 (Fig. 16), which carries the image to the focal plane of the eyepieces.

The prism 2 diverts the passing rays and directs them, at an angle of 45° to the vertical, to the system of cemented prisms 3 comprising a rhombic and a triangular prism. The contact face of the rhombic prism is semi-chronium plated.


Fig. 15
The AV-12 Binocular Attachment


Fig. 16
Optical System of the AV-12
Binocular Attachment

26

This semi-chronium plated cementing surface passes 50% of the light rays and directs them through prism 4 to the eyepiece 5. The reflected 50% of the rays pass through the rhombical part of the prisms 3 and are deviated through an angle of 90° through the compensator 6 and further to the second eyepiece 7 of the attachment.

The optical system of the Binocular Attachment is designed to provide an upright image of the object examined.

The interchangeable eyepieces are inserted into the eyepiece tubes 1 (Fig. 17) screwed into the right- and left-hand bodies 2.

The dioptic mechanism 3 on the left eyepiece tube is adjusted by rotating the knurled ring 4 having a scale.

The attachment has an interocular adjustment which allows the distance between the eyepieces to be regulated to suit the distance between the pupils of the observer. The required interocular distance is set to a scale on the flange of the eyepiece tube.

Due to the spherical fitting 5, the Binocular Attachment may be rotated to any position and inclined to any angle convenient for the observer.

The Binocular Attachment is stored in a special case containing holders for the eyepieces.


Fig. 17
Main parts of the AV-12
Binocular Attachment

SPECIFICATIONS

Initial magnification of attachment.....	1.5X
Magnification of Huyghenian eyepieces	7X and 10X
Dioptic adjustment of eyepiece, diopters	± 5
Range of interocular adjustment, mm	from 55 to 75
Overall dimensions of attachment, mm	160 × 120 × 75
Overall dimensions of case, mm	230 × 140 × 90
Weight of attachment, kg	0.76
Weight of attachment in case, kg	1.15

BINOCULAR ATTACHMENT SET

Huyghenian eyepieces 7X (2 pieces)
Huyghenian eyepieces 10X (2 pieces)
Caps (2 pieces)
Case
Description and instruction manual
Certificate

27

UNIVERSAL PHOTOMICROGRAPHIC ATTACHMENT, MODEL MΦH-1

The MΦH-1 Photomicrographic Attachment (Fig. 18) serves in an universal manner for photographing objects through a microscope.

The attachment is designed to be supported by the microscope and may be attached to any instrument having an eyepiece tube of standard external diameter of 25 mm.

Either of two sizes of cameras may be furnished with the attachment in accordance with the Buyer's requirements. Type MΦK-1 Camera has a plate size of 6.5 x 9 cm while type MΦK-2 has a plate size of 9 x 12 cm.


Fig. 18
MΦH-1 Photomicrographic Attachment


Fig. 19
Diagrammatical sectional view of the MΦH-1 Photomicrographic Attachment

The Photomicrographic Attachment operates in the following manner:
Rays of light, passing through the microscope eyepiece 1 (Fig. 19) either directly to the photographic plate 2, through the shutter 3, or if the prism 4 is inserted, are reflected at an angle of 70° to the microscope axis and are directed to the visual observation tube.

Objective 5, of the visual observation tube, focus the image of the object in the plane of the reticule 6, arranged in the focal plane of the eyepiece 7.

The eyepiece reticule 6 is so positioned that a sharp image of the object appears simultaneously on both the reticule and the photographic plate.

The eyepiece 7 has a dioptric adjustment of ±5 diopters which allows the eyepiece to be focused to obtain a sharp image on the reticule to suit the observer.

The body 8, of the attachment, has an adaptor 9 for fastening it to the eyepiece tube 1. It is fixed by the screw 10.

The swing-out prism 4 is mounted in the body 8 and is interlocked by a lever with the control mechanism of the shutter 3. When the cable release button is depressed, the prism 4 is automatically swung aside out of the path of the rays. At the end of the exposure, the prism returns to its initial position.

The shutter provides for two exposures: "bulb"—"K" and "snap-shot"—"M". The "snap-shot" exposure is obtained by pressing the cable button and does not depend on the length of time the button is held depressed. When arranged for "bulb" exposure, the shutter remains open as long as the button is depressed.

After the shutter has operated and the button is released, the prism returns to its initial position and directs the rays of light to the visual observation tube.

The camera 11 is connected to the body 8 and may be easily removed after unscrewing the threaded ring 12.

The plate holder 13 with the photographic plate is inserted into the camera. To make cameras, types MΦK-1 and MΦK-2, fully interchangeable—and to eliminate adjustment when changing from one camera to the other, a lens is arranged in the lower flange of the MΦK-2 camera to compensate for the difference in length of the cameras.

The objective mounting 15, reticule mounting 16 and eyepiece mounting 17 are arranged inside the tubular body 14 of the visual observation tube.

The dioptric graduations are engraved on the external tube of the eyepiece mounting. A smoke-colored light filter 18 is placed on the eyepiece to protect the eye from the intense light necessary during photography.

SPECIFICATIONS

Plate size of camera, cm:		
MΦK-1	6.5 x 9	
MΦK-2	9 x 12	
Coefficient of magnification of photograph:		
MΦK-1	0.5	
MΦK-2	1	
Overall dimensions of attachment, mm:	with camera MΦK-1	153 x 155 x 126
	with camera MΦK-2	224 x 155 x 148
Overall dimensions of case, mm:	with camera MΦK-1	234 x 173 x 108
	with camera MΦK-2	268 x 185 x 126
Weight of attachment, kg:	with camera MΦK-1	0.55
	with camera MΦK-2	0.70
Weight of attachment in case, kg:	with camera MΦK-1	1.3
	with camera MΦK-2	1.7

PHOTOMICROGRAPHIC ATTACHMENT SET

- Camera MΦK-1 or MΦK-2
- Plateholders 6.5 x 9 or 9 x 12 (3 pieces either)
- Light filter in mounting
- Cable release
- Case for attachment
- Description and instruction manual
- Certificate

MICROMETER EYEPIECE, MODEL AM-9-2

The AM-9-2 Micrometer Eyepiece (Fig. 20) is an eyepiece with a fixed and a movable reticule in its field of view.

The Micrometer Eyepiece is designed for accurate measurement of objects being studied under the microscope.

The instrument can be used on all models of biological microscopes to replace the standard eyepiece.

The body 1 (Fig. 21) of the Micrometer Eyepiece has a fitting 2 by means of which it is fastened on the microscope tube with the aid of screw 3.


Fig. 20
AM-9-2 Micrometer Eyepiece


Fig. 21
Main parts of the AM-9-2 Micrometer Eyepiece

The compensating eyepiece 4, containing the stationary reticule, is inserted into the upper part of the body. The reticule has a millimeter scale 8 mm in length. The eyepiece has a dioptric adjustment to allow the scale images to be properly focused.

The graduated drum 5, connected with the tube 6, has a circular uniform scale graduated to 100 divisions. A line engraved on the tube serves as a reference for drum readings.

One revolution of the drum displaces the movable reticule 1 mm. A rotation through one division of the drum corresponds to a displacement of 0.01 mm.

The Micrometer Eyepiece is stored in a polished wooden case.

SPECIFICATIONS

Eyepiece magnification	15×
Total movement of reticule, mm	8
Value of drum scale divisions, mm	0.01
Overall dimensions of instrument	
(length × width × height), mm	84 × 42 × 58
Overall dimensions of case, mm	120 × 90 × 48
Weight of instrument, kg	0.21
Weight of instrument in case, kg	0.38

MICROMETER EYEPIECE SET

- Case
- Description and instruction manual
- Certificate

DRAWING APPARATUS, MODEL PA-4

The PA-4 Drawing Apparatus (Fig. 22) is used for tracing the outlines of objects examined under a microscope, as well as for checking the magnification provided by the optical system of a microscope.

This apparatus is applicable to all biological microscopes constructed on the basis of the M-10 microscope, as well as to the MBH-1 and MBH-2 microscopes when using a straight monocular body.

The main part of the instrument is the prism-cube 1 (Fig. 23). It comprises two rectangular prisms cemented together on their hypotenuse faces. The hypotenuse face of one prism is semi-silvered.

This face of the prism-cube directs rays of light from the drawing paper 2, reflected by the mirror 3, to the observer's eye. It also passes approximately 50% of the light from the object.

In this way the observer simultaneously sees both the object and the paper and pencil located 250 mm from the eye.

To even the backgrounds, and consequently, to provide better visibility of both the specimen and the paper and pencil, two systems with neutral and blue light filter 4 are furnished in the instrument.

The prism-cube is mounted in a swing-back fitting which carries the sector and the drum with light filters. Both the sector and the drum have four openings. One opening is free, two have neutral filters of different density while the fourth opening carries a blue filter.

The fitting has an opening on the side facing the mirror through which rays pass from the paper to the prism-cube.

The mirror mount is fastened to a bar whose other end is fixed in a hole in an upright member. The latter is rigidly fastened to the fitting by means of which the apparatus is attached to the microscope eyepiece tube.

The Drawing Apparatus is stored in a wooden case.

SPECIFICATIONS

Overall dimensions of the instrument	
(length × width × height), mm	215 × 90 × 40
Overall dimensions of case, mm	240 × 120 × 55
Weight of instrument, kg	0.23
Weight of instrument in case, kg	0.64


Fig. 22
PA-4 Drawing Apparatus


Fig. 23
Optical System of the PA-4 Drawing Apparatus

PHASE-CONTRAST EQUIPMENT, MODEL KΦ-1

As a rule, microscopy is concerned with low contrast specimens which must be stained before observation under the microscope. Living specimens usually perish after staining.

The simplest method of increasing the contrast by diaphragming the aperture diaphragm of the condenser, decreases the revolving power and the illumination without increasing the contrast in any considerable degree.

Dark ground observation provides good contrast but it is reversed. Light parts of the specimen are shown dark and vice versa. This method, besides, only allows the outline of the specimen to be determined without revealing its internal structure.


Fig. 24
KΦ-1 Phase-Contrast Equipment

The most perfected method of research is observation by means of the phase-contrast method. This method provides for observation of unstained noncontrasting specimens. A contrasting image is obtained in which the lighter or darker portions correspond to various thickness or optical densities of the specimen.

In this way, the phase-contrast method opens wide possibilities for research on living unstained specimens and finds extensive application in bacteriology, biology, medicine and other branches of science.

The equipment for the phase-contrast method of observation can be used on the biological microscopes, models M-10, MB11-1 and MB14-4.

The KΦ-1 Equipment (Fig. 24) comprises the following main parts:

1. Objectives for phase-contrast observation.
2. Phase-contrast condenser.
3. Auxiliary microscope.

The objectives in the equipment differ from usual achromatic objectives only in that a phase plate has been placed in the plane of the exit lens. This phase plate changes the phase of the zero maximum by 90° and decreases its intensity.

Besides the usually engraved figures indicating the initial magnification and the aperture, these objectives have the letter "Φ" engraved on the mount and case indicating that they are of the phase type.

The condenser included in the equipment does not differ from the usual 1.2 N. A. condenser except that it has a revolving disc with annular diaphragms arranged in its focal plane.

Each objective requires its annular diaphragm. Numbers appear in an opening of the revolving disc which indicate the objective that should be used in conjunction with the indexed annular diaphragm.

For observation in the usual manner, the condenser is provided with an iris diaphragm and the revolving disc with a free opening for passing light rays. The condenser is inserted into the microscope condenser holder and is fastened by a screw in the usual manner.

The auxiliary microscope included in the equipment serves for checking the centring of the annular diaphragm in reference to the phase plate of the objective. The auxiliary microscope is inserted into the microscope tube instead of the eyepiece and after centring is completed, it is replaced by the usual eyepiece. The microscope comprises a fitting in which the objective is mounted and into which the eyepiece is inserted. The latter can be adjusted in the fitting and fixed in the required position by a screw.

The equipment is stored in a wooden case.

SPECIFICATIONS

Objectives	special with phase plate — 10× × 0.30; 20× × 0.40; 40× × 0.65 and 90× × 1.25
Condenser	aperture 1.2; with iris diaphragm and revolving disc holding four annular diaphragms
Magnification of auxiliary microscope	20×
Overall dimensions of equipment, mm	120 × 130 × 60
Overall dimensions of case, mm	220 × 150 × 80
Weight of equipment, kg	0.51
Weight of equipment in case, kg	1.68

PHASE-CONTRAST EQUIPMENT SET

- Objective 10× × 0.3
- Objective 20× × 0.4
- Objective 40× × 0.65
- Objective 90× × 1.25
- 1.2 N.A. condenser, with revolving disc
- Auxiliary microscope
- Case
- Description and instruction manual
- Certificate

MECHANICAL STAGE, MODEL CT-12

The CT-12 Mechanical Stage (Fig. 25) has been designed for smooth movement of the specimen on the object stage of the microscope in two directions at right angles to each other.

The CT-12 Mechanical Stage is an attachment which can be fitted to any biological microscope.

The frame 1 of the stage (Fig. 26) has guides 2 along which slides, with clips 3 and 4, can be moved. The slides with clips can be fixed by the clamping screws 5 to suit the length of the object glass.


Fig. 25
CT-12 Mechanical Stage

Fig. 26
Main parts of the CT-12 Mechanical Stage

The specimen is shifted in a transverse direction by the head 6 and in the longitudinal direction by the head 7.

The scales 8 and 9 for the transverse and longitudinal movements are read by means of the verniers 10 and 11 fastened to the movable parts of the stage.

The Mechanical Stage is stored in a wooden case.

SPECIFICATIONS

Transverse movement scale range, mm	from 50 to 135
Longitudinal movement scale range, mm	from 0 to 35
Reading accuracy, mm	0.1
Overall dimensions of instrument (length × width × height), mm	150 × 135 × 30
Overall dimensions of case, mm	172 × 130 × 40
Weight of instrument, kg	0.278
Weight of instrument in case, kg	0.550

MECHANICAL STAGE SET

- Centring plate
- Stud
- Case
- Description and instruction manual
- Certificate

NORMAL INCIDENT ILLUMINATOR, MODEL OI-1

The OI-1 Normal Incident Illuminator is designed for the intensive illumination of opaque or semi-translucent objects seen under biological microscopes.

With the Normal Incident Illuminator (Fig. 27) the object is illuminated through the microscope objective. Consequently the specimen is seen on a bright background, that is, the portions of the specimen which reflect less light are seen as dark spots or bands in accordance to their form.

When using the illuminator it is necessary to use objectives in short mountings computed for a mechanical tube length of 190 mm and corrected for a specimen without cover-glass.


Fig. 27
OI-1 Normal Incident Illuminator


Fig. 28
Principle of operation of the OI-1 Illuminator

The Normal Incident illuminator is applicable on the M-10 Microscopes and all other models constructed on its basis.

When using the illuminator the draw tube of the M-10 Microscope should be set at the 160 mm division as the 30 mm length of the illuminator body supplements to the computed mechanical tube length of 190 mm.

Fig. 28 indicates the principle of operation of the Normal Incident Illuminator. The source of light is an electric bulb 1, supplied a. c. mains through a step-down transformer.

The rays of light pass through lens 2 and diaphragm 3 and fall on the thin glass plate 4, arranged between the objective 5 and the microscope eyepiece not shown on the diagram.

The glass plate is positioned at an angle of 45° to both the direction of the rays and to the optical axis of the microscope.

Part of the light passes through the plate while the other part is reflected through the objective 5 and concentrated by the latter on the observed surface of the specimen 6.

On being reflected from the surface, the rays again pass through the objective and falling on the plate 4 are again diverged. A part is reflected back to the illuminator while the rest, after passing through the plate, is directed through the eyepiece to the eye of the observer.

The cylindrical body 1 of the illuminator (Fig. 29) connects the microscope objective, which is screwed into the socket 2, with the microscope body which is screwed on the threaded ring 3.

The socket 5 with bulb is inserted into one end of the tube 4 while an iris diaphragm is mounted inside the other end. Rotating the handle 6 changes the opening from 0.5 mm to 7.5 mm and in this way varies the brightness of illumination of the object.

Rotating the head 7 revolves the glass plate through an angle of 45° to the axis of the body 1. This allows the plate to be adjusted to the most advantageous position for illuminating the object.

The necessary friction to prevent accidental rotation of the plate is provided by two spring washers.


Fig. 29
Main parts of the OI-1 Normal Incident Illuminator

The illuminator body has a port 8 on the side opposite to the illuminating tube. It serves to pass the light from the lamp, that passes through the glass plate. This is necessary to eliminate a glare that interferes with observation.

The Normal Incident Illuminator is stored in a special wooden case.

SPECIFICATIONS

Source of illumination	electric bulb 3.5 V, 0.16 A
Range of diameter variation of iris diaphragm, mm	from 0.5 to 7.5
Overall dimensions of instrument, mm	90 × 35 × 33
Overall dimensions of case, mm	120 × 90 × 48
Weight of instrument, kg	0.12
Weight of instrument in case, kg	0.285

ILLUMINATOR SET

- Pin (wrench)
- Case
- Description and instruction manual
- Certificate

DARK GROUND CONDENSER, MODEL OI-13

When examining weakly contrasting objects under the microscope, conditions are often met with in which observation on a bright ground, with the usual type of condenser, is impossible as the image of the object is invisible.

To increase the contrast of the image in reference to the general background of the microscope field of view, so-called dark ground illumination is used. This provides the possibility of revealing details of the object that were invisible before.

This type of illumination can be provided by means of the Dark Ground Condenser, Model OI-13 (Fig. 30).

The Dark Ground Condenser is applicable to all models of biological microscopes except MV and MHC-7.

Increasing the resolving power of the microscope by the aid of the Dark Ground Condenser is based on the diffraction of light and on obtaining a light contrast of an illuminated object on a dark background.

This effect is achieved by means of the illumination created by the Dark Ground Condenser much in the same way as particles of dust, invisible in the air under usual conditions, become distinctly visible in a dark room in rays of sunlight.

The optical system of the condenser is shown in Fig. 31.

A parallel pencil of light, from the microscope mirror, passes through the plane-concave lens 1, the plane-parallel plate 2 having a spherical recess internally silvered, the cardioid lens 3 and exits in an oblique direction.

In this way the object receives oblique illumination in an absolutely dark background as the direct rays, due to the spherical recess in the plate 2, do not reach the objective of the microscope.

The Dark Ground Condenser comprises the external ring 1 (Fig. 32) into which are mounted the internal ring with a cylinder 2 and the sleeve 3 containing the optical equipment 4.

Rotating the two screws 5 displaces the cylinder with the optical system in a plane perpendicular to the optical axis of the microscope. The screws are rotated with the aid of a special wrench 6.


Fig. 30
OI-13 Dark Ground Condenser


Fig. 31
Optical System of the OI-13 Dark Ground Condenser


Fig. 32
Main parts of the OI-13 Dark Ground Condenser

The Dark Ground Condenser can be used with either oil-immersion objectives or with dry system objectives. The dark ground is less satisfactory in the second case.

When working with the Dark Ground Condenser it is advisable to use a cover-glass from 0.8 to 1.2 in thickness.

The condenser is furnished in a special wooden case.

SPECIFICATIONS

Numerical aperture	1.2
Overall dimensions of condenser (diameter x height), mm	58 x 36
Overall dimensions of case, mm	70 x 70 x 40
Weight of condenser, kg	0.19
Weight of condenser in case, kg	0.48

CONDENSER SET

- Diaphragm for objective 90x x 1.25
- Wrenches for adjustment (2 pieces)
- Case
- Description and instruction manual
- Certificate

MICROSCOPY LAMP, MODEL OII-7

The OII-7 Microscopy Lamp (Fig. 33) is designed for illuminating the specimen. The lamp comprises a two-lens condenser with an iris diaphragm and a clamping device by means of which the lamp body is fastened to the vertical stand at the required height.

The source of light is a special electric bulb, 8V, 20W, supplied from an a. c. mains through step-down transformer 127/220/8V.


Fig. 33
OII-7 Microscopy Lamp

SPECIFICATIONS

Overall dimensions of instrument (length x width x height), mm	250 x 240 x 110
Overall dimensions of case, mm	310 x 200 x 150
Weight of lamp, kg	1.37
Weight of lamp in case, kg	4.6

MICROSCOPY LAMP SET

- Lamp on stand
- Connecting strip
- Transformer
- Electric bulbs, 8V, 20W (2 pieces)
- Case for lamp
- Description and instruction manual
- Certificate

POLARIZING MICROSCOPES

POLARIZING MICROSCOPE, MODEL MHH-4

The MHH-4 Polarizing Microscope (Fig. 34) has wide applications in mineralogy, petrography and mineralography. It can also be used in chemistry and biology.

This instrument is a large-size up-to-date polarizing microscope and its design features fully satisfy all requirements made to this class of microscopes.


Fig. 34
MHH-4 Polarizing Microscope

This microscope is designed for the inspection of translucent specimens in transmitted (ordinary and polarized) light with either conoscopic or orthoscopic optical systems.

Besides its usual applications, the design of the microscope provides for the examination of minerals by Prof. Fyodorov's method with reflected ordinary and polarizing light, as well as for the photography of specimens by means of the MΦH-1 photomicrographic attachment.

The Polarizing Microscope, in addition to the usual optical system of objective, eyepiece and illuminating system, has prisms for polarizing light and for analysing the polarized light. The microscope has an additional lens for conoscopic operation.

Fig. 35 illustrates the optical system of the microscope.

From the source of light 1, the rays pass through the condenser 2 to the mirror 3 and are reflected to the polarizer 4.

After passing through the diaphragm 5, the polarizing light enters one of three interchangeable condensers 6. Further on the rays pass through the lens 7 and the specimen 8 into the objective 9 and analyzer 10. They then pass either directly to the eyepiece 11 (for ortho-

scopic operation) or first through a Bertrand type lens, not shown on the diagram (for conoscopic operation).

To prevent the analyzer 10 from distorting the image, it is located between the lenses 12 creating a parallel pencil of rays.

A quartz plate 13 may be inserted between the objective and the analyzer. It compensates for the variation in the travel of the rays caused when the light passes through the examined specimen.

The microscope has the following constructive features:

The base 1 of the stand (Fig. 36) has a horse-shoe form with a boss having two lugs. The limb 2 is pivoted between the lugs. The limb is designed so that the object stage can accommodate the Fyodorov stage or other attachments that can be rotated through 360°. At the same time, the limb serves for carrying the instrument from place to place.

Coarse focusing adjustment is effected by a rack and pinion arranged in the upper part of the limb. They are operated by rotating the head 4.


Fig. 35
Optical System of the MHH-4 Microscope

Fig. 36
Main parts of the MHH-4 Microscope

The head 5 for fine focusing adjustment is located somewhat lower. The head has a graduated drum on which readings can be made to an accuracy of 0.002 mm.

The body 6 has a slot in the upper part for the Bertrand type lens mechanism. It has a handle 7 for controlling the variable diameter diaphragm and a head 8 for displacing the lens along the optical axis of the microscope.

The lower end of the tube has dovetail guides for the mechanism with the analyser prism mount. This mechanism provides for rotating the analyser to the required angle read on a special scale, for inserting the quartz compensating plate and for including or excluding the analyser in the optical system of the microscope.

The clamping device 9 serves to hold the objectives.

A drawtube 10 is screwed into the upper end of the body to accommodate interchangeable eyepieces or for fastening the photomicrographic camera.

A large set of special objectives, computed for operation in polarized light, and eyepieces are included in the instrument set. Part of the eyepieces are furnished with reticule or cross hairs for reading purposes.

The microscope object stage is of the circular rotary type. The working part of the stage is furnished with a circular scale graduated to 360° in 1° divisions. The angle of rotation of the object stage is read on a vernier with a reading accuracy of 6 minutes. The object stage has two holes for mounting standard spring clips, three holes for attaching the Fyodorov stage and three holes for fastening the mechanical stage 11 included in the microscope set.

The illuminating system of the microscope comprises the mirror 12 fastened in a fork-shaped bracket, the interchangeable condenser 13 and the polarizer.

The condenser together with the polarizer, can be adjusted vertically by means of the head 14.

The Polarizing Microscope is stored in a wooden case.

SPECIFICATIONS

Range of total magnification	from 1.85 ^x to 1350 ^x
Achromatic objectives:	3.7 ^x 8 ^x 20 ^x 60 ^x 90 ^x
initial magnification	3.7 ^x 8 ^x 20 ^x 60 ^x 90 ^x
numerical aperture	0.11 0.20 0.40 0.85 1.25
focal length, mm	33.1 18.2 8.4 2.99 1.96
working distance, mm	27.7 9.19 1.8 0.41 0.32
field of view with 5 ^x eyepiece, mm	6.24 2.85 1.15 0.38 0.25
Wide-angle Huyghenian eyepieces:	With reticule With cross and scale Γ-5 ^x hair Γ-8 ^x
magnification	5 ^x 8 ^x
focal length, mm	50.6 31.4
linear field of view, mm	23 21
Photographic eyepiece, 10 ^x :	
magnification	10 ^x
focal length, mm	25
linear field of view, mm	13.4
Symmetrical eyepiece, 15 ^x :	
magnification	15 ^x
focal length, mm	17
linear field of view, mm	12
Objective changing mechanism	clamping device
Illuminating system	reversible mirror, plane and concave; interchangeable condensers with 1.27; 0.94 and 0.22 N.A.; polarizer; interchangeable opal glass
Focusing adjustments:	
coarse	by rack and pinion
fine	micrometric mechanism, reading to 0.002 mm
Object stage	circular, with rotary movement and with attached mechanical stage
Range of condenser substage adjustment, mm	20
Angle of polarizer rotation	360°
Rotation angle reading accuracy	2.5°
Angle of analyser rotation	90°
Rotation angle reading accuracy	2.5°
Eyepiece scale reading, mm	0.1

42

Value of reticule squares, mm	0.5 × 0.5
Angle of rotation of object stage	360°
Object stage scale reading	1°
Overall dimensions of microscope (height × length × width), mm	350 × 250 × 155
Overall dimensions of case, mm	400 × 259 × 280
Weight of microscope, kg	5.9
Weight of microscope in case, kg	12.8

MICROSCOPE SET

Microscope stand comprising horse-shoe base, limb with coarse and fine focusing mechanisms, illuminating system and rotary object stage

Achromatic objective, 3.7^x × 0.11

Achromatic objective, 8^x × 0.20

Achromatic objective, 20^x × 0.40

Achromatic objective, 60^x × 0.85

Achromatic objective, 90^x × 1.25

Wide angle Huyghenian eyepiece, 5^x

Wide angle Huyghenian eyepiece, 8^x

Photographic eyepiece, 6.5^x

Photographic eyepiece, 10^x

Symmetrical eyepiece, 15^x

Stage micrometer, for transmitted light, in case

Quartz compensating wedge in mount

Quartz compensating plate in mount

1.27 N.A. interchangeable condenser

0.94 N.A. interchangeable condenser

0.22 N.A. interchangeable condenser

Opal glass in mount

Mechanical stage in case

Fitting for M4H-1 photomicrographic attachment

Adapter for standard eyepieces

Superimposed diaphragm for oscoscopic operation

Wrench for clamping device

Combination wrench

Watchmaker's screwdriver

Squirrel-hair brush

Flannel napkin

Cover

Case for microscope

Case for objectives

Case for eyepieces and accessories

Description and instruction manual

Certificate

43

PETROLOGICAL MICROSCOPE, MODEL MHH-5

The MHH-5 Petrological Microscope (Fig. 37) is designed for all types of petrological work, as well as for educational purposes.

This microscope has the same specifications as the MHH-4 Microscope, except that the polarizer and analyzer are manufactured of high-quality polyvinyl poluroid plastics, and the included set of objectives and eyepieces provide for magnification from 18.5x to 600x.

The microscope and its accessories are stored in a wooden case.


Fig. 37
MHH-5 Petrological Microscope

SPECIFICATIONS

Range of total magnification	from 18.5x to 600x			
Achromatic objectives:	3.7x	8x	20x	60x
initial magnification	3.7x	8x	20x	60x
numerical aperture	0.11	0.20	0.40	0.85
focal length, mm	33.1	18.2	8.4	2.99
working distance, mm	27.7	0.19	1.8	0.41
field of view with 5x eyepiece, mm	6.24	2.85	1.15	0.38
Wide-angle Huyghenian eyepieces:	With reticule and scale P-5x		With cross hair P-8x	
magnification	5x		8x	
focal length, mm	50.0		31.4	
linear field of view, mm	23		21	

Photographic eyepiece, 10x:	
magnification	10x
focal length, mm	25
linear field of view, mm	13.4
Objective changing mechanism	clamping device
Illuminating system	reversible mirror plane and concave; interchangeable condensers with 1.27, 0.94 and 0.22 N.A.; polarizer; interchangeable opal glass
Focusing adjustments:	
coarse	by rack and pinion
fine	micrometric mechanism, reading to 0.002 mm
Object stage	circular with rotary movement
Overall dimensions of microscope (height x length x width), mm	350 x 250 x 155
Overall dimensions of case, mm	400 x 259 x 280
Weight of microscope, kg	6.0
Weight of microscope in case, kg	12.6

MICROSCOPE SET

- Microscope stand comprising horse-shoe base, limb with coarse and fine focusing mechanisms, illuminating system and rotary object stage
- Achromatic objective, 3.7x x 0.11
- Achromatic objective, 8x x 0.20
- Achromatic objective, 20x x 0.40
- Achromatic objective, 60x x 0.85
- Wide-angle Huyghenian eyepiece, 5x
- Wide-angle Huyghenian eyepiece, 8x
- Photographic eyepiece, 10x
- Stage micrometer for transmitted light, in case
- Quartz compensating wedge in mount
- Interchangeable condenser, 1.27 N.A.
- Interchangeable condenser, 0.94 N.A.
- Interchangeable condenser, 0.22 N.A.
- Opal glass in mount
- Fitting for M411-1 photomicrographic attachment
- Adapter for standard eyepieces
- Superimposed diaphragm for conoscopic operation
- Wrench for clamping device
- Combination wrench
- Watchmaker's screwdriver
- Squirrel-hair brush
- Flannel napkin
- Cover
- Case for microscope
- Case for objectives
- Case for eyepieces and accessories
- Description and instruction manual
- Certificate

ORE MICROSCOPE, MODEL MIII-6

The MIII-6 Ore Microscope (Fig. 38) is designed for exact examination of opaque objects in ordinary and polarized light at magnifications from 42x to 750x. Besides its usual applications, the design of the microscope provides for the examination of translucent specimens in transmitted (ordinary or polarized) light, as well as for photography by means of the MΦH-1 photomicrographic attachment. This microscope is used for all types of mineralographic research. Fig. 39 illustrates the main parts of the Ore Microscope.


Fig. 38
MIII-6 Ore Microscope

Fig. 39
Main Parts of the MIII-6 Ore Microscope

The microscope stand comprises the base 1, limb 2 with head 3 for coarse focusing and drum 4 for fine micrometric focusing adjustment.

The microscope body 5 is of tubular form with a dovetail recess at its lower end. Into this recess is inserted the slide 6 with the polyvinyl polaroid plastics analyser. Transverse movement of the slide brings the analyser in or out of the microscope optical system.

The eyepiece tube 7 is screwed into the upper end of the body. It is designed to accommodate standard eyepieces or a photomicrographic attachment, model MΦH-1. The illuminating system of the microscope comprises the polyvinyl polaroid plastics polarizer in a mount 8, which can rotate through an angle of 360°, and the reversible plane and concave mirror 9.

The polarizer mount is fastened to the swing-out bracket 10.

The instrument includes the vertical polarized light illuminator 11, designed for the intensive illumination of the field of view of the examined surface of opaque and semi-translucent objects.

When using the vertical illuminator, the object is illuminated through the objective 12 which is inserted into the body by means of the clamping device 13.

The microscope is furnished with a rotary object stage 14 having a braking device 15.

The microscope is stored in a wooden case.

SPECIFICATIONS

Range of total magnification	from 42x to 750x				
Achromatic objectives:	4.2x	6x	10x	23.5x	35.3x 50.4x
initial magnification with additional lens, focal length — 141 mm	4.2x	6x	10x	23.5x	35.3x 50.4x
numerical aperture	0.12	0.17	0.30	0.65	0.75 1.25
focal length, mm	33.1	23.17	13.89	5.99	4.0 2.77
working distance, mm	18.9	6.2	5.71	0.96	0.32 0.44
field of view with 7x eyepiece, mm	4.2	3.0	1.8	0.75	0.36 0.52
Eyepieces:	Huyghenian 7x	Huyghenian 10x	Symmetrical 15x	Compensating 15x	
focal length, mm	36	25	17	16.8	
linear field of view, mm	18	14	12	12	
Objective changing mechanism	clamping device				
Illuminating system	reversible mirror, plane and concave; polarizer				
Focusing adjustments:	by rack and pinion				
coarse	micrometric mechanism, reading to 0.002 mm				
fine	circular with rotary movement				
Object stage	360°				
Range of object stage rotation	1°				
Object stage scale reading	350 x 250 x 155				
Overall dimensions of microscope (height x length x width), mm	400 x 259 x 280				
Overall dimensions of case, mm	5.0				
Weight of microscope, kg	14.2				
Weight of microscope in case, kg					

MICROSCOPE SET

Microscope stand comprising horse-shoe base, limb with coarse and fine focusing mechanisms, illuminating system, rotary, circular object stage and vertical illuminator

- Achromatic objective 6x x 0.17
- Achromatic objective 10x x 0.30
- Achromatic objective 4.2x x 0.12
- Achromatic objective 50.4x x 1.25
- Achromatic objective 35.3x x 0.75
- Achromatic objective 23.5x x 0.65
- Huyghenian eyepiece 7x with reticule and scale
- Huyghenian eyepiece 10x
- Symmetrical eyepiece 15x
- Compensating eyepiece 15x

Stage micrometer for transmitted light
Blue filter
Transformer 127/220/8V
Electric bulbs 8V, 8.5W (3 pieces)
Hand press
Wrench for removing clamping device
Wrenches for adjusting objective holder (2 pieces)
Vial with immersion oil
Combination wrench
Watchmaker's screwdriver
Flannel napkin
Squirrel-hair brush
Case for microscope
Case for objectives
Case for eyepieces
Description and instruction manual
Certificate

METALLOGRAPHIC MICROSCOPES

HORIZONTAL METALLOGRAPHIC MICROSCOPE, MODEL MHM-8

The MHM-8 Horizontal Metallographic Microscope (Fig. 40) provides for comprehensive research in the microstructure of metals and has a total effective magnification up to 1300 \times (for visual observation) and up to 3000 \times (for photography).


Fig. 40
MHM-8 Horizontal Metallographic Microscope

In its design features and optical characteristics, this microscope fully satisfies up-to-date requirements of metallography and provides for observing and photographing objects under the following conditions:

- a) bright field with normal and oblique illumination;
- b) dark field;
- c) polarized light.

Opaque objects are examined under the microscope in reflected light.

The optical system of the microscope (Fig. 41) comprises the following elements: light source — arc lamp 1 with clockwork mechanism, collector 2, heat filter 3 for

protecting the polarizer prisms from overheating, swing-out lens 4 used only for dark field observation, incandescent lamp 5 used for visual observation of objects with bright field (and sometimes for polarized light), light filters 6, polarizer 7 inserted into the system for polarized light observations, iris diaphragm 8 for adjusting the illumination of the object, illuminator lenses 9, 10 and 11, field and annular diaphragms 12 mounted in one strip, oblique illumination prism 13 for oblique illumination of the object, reflecting plate 14 for illuminating the object through the objective in bright field observation, annular mirror 15 and metallic condenser 16 for dark field illumination, objective 17, analyser 18 for observation of the object in


Fig. 41
Optical System of the MIM-8 Microscope

polarized light, achromatic lens 19 of the visual tube, prism 20 of the visual tube, prism 21 of the photographic tube, achromatic lens 22 of the photographic tube and Homal eyepiece 23.

With bright field observations the collector 2 forms an image of the source of light 1 on the aperture diaphragm 8.

An image of the aperture diaphragm is formed on the exit lens of the objective 17 by means of the system of three lenses (9, 10 and 11) and the plate 15. In this way, a pencil of light from the arc lamp 1, after passing through the aperture and field diaphragms, is thrown on plate 15. The latter reflects a portion of the light through the objective 17 to the object. The remainder of the light passes through the plate and does not take part in illuminating the object.

If the oblique illumination prism 13 is introduced in the optical system, all the light is reflected to the object, through the objective. For this reason, illumination by means of the prism illuminator is considerably brighter than when using the reflecting plate.

The metallic annular condenser 16 is used for dark ground observations. It has a reflecting parabolic mirror surface arranged around the objective 17. The swing-out lens 4 and the annular diaphragm are introduced when making dark field observations.

The hollow cylinder of light passes through the annular diaphragm to the annular mirror 15 in whose center the plate 14 is arranged. This mirror reflects rays

of light to the surface of the dark ground condenser 16 which converges them to the plane of the object.

Diffused rays of light, reflected from the object, as with light field observations, pass through the objective, reflecting plate and achromatic lenses to the visual tube or photographic camera and produce an image of the object.

The arc lamp is too bright for visual observation and consequently, for this purpose, an incandescent bulb 5 is used and located directly in front of the aperture diaphragm 8.

The microscope comprises four main parts: a) illuminating system; b) the microscope proper; c) the photographic camera arranged on a massive base and d) the special table.

The optical bench is a bar of prismatic cross-section with slots for erecting the separate parts of the instrument.

The arc lamp, mounted in a housing, is fastened to the bench with a special screw. The carbons are held in special clamps which are fed by the clockwork mechanism. The mechanism and the carbons are enclosed in a hood which can be lifted for changing carbons. The arc lamp is supplied from a 127 or 220 V. mains in series with a rheostat.

The housing of the incandescent lamp is fastened on a holder together with the opaque screen. The mechanism with the lamp can be brought out of the optical system simply by being swung downward. The screen serves to protect the observer from stray light.

On a holder, arranged directly after the arc lamp housing, the parabolic collector and the swing-out lens for dark field observation are fastened. The heat-absorbing device with distilled water is located after the collector.

The central part of the instrument comprises the microscope proper whose parts and mechanisms are mounted inside and outside of the casing fastened to the bench.

This part of the instrument comprises: the illuminating tube, central prism system, object stage, objective holder, coarse focusing mechanism, fine focusing mechanism, monocular or binocular attachments and the tube connecting the microscope with the photographic camera.

The coarse focusing mechanism is connected with the object stage. Besides, the latter is provided with crosswise movement and rotates about a vertical axis. In this way, the object stage is universal in type.

The fine focusing mechanism is connected with the objective holder.

Unlike the majority of other microscopes, in the MIM-8 Microscope the objectives are not screwed into the tube but are inserted into a special bearing ring of the objective holder. This considerably simplifies changing objectives.

Another distinctive feature of the objectives of this microscope is that they are computed and corrected to an infinitely distant image.

The photographic camera, fastened to the bench, comprises the extensible bellows, front and rear boards, shutter, scale and mirror.

The photographic camera is designed for 13×18 cm plates but may also be used for 9×12 cm plates as well. The camera is furnished with a wooden plate holder.

The microscope table comprises: the table top and two cabinets with drawers for holding the accessories.

SPECIFICATIONS

Range of total magnification	from 45× to 1350×
Image scale when photographing	from 100: 1 to 3000: 1
Achromatic objectives for light and dark field:	11× 18× 40×
initial magnification	11× 18× 40×
numerical aperture	0.17 0.30 0.65
focal length, mm	23.17 13.89 6.16
working distance, mm	6.2 5.71 0.87
Apochromatic objectives for bright and dark field:	15× 30×
initial magnification	15× 30×
numerical aperture	0.30 0.65
focal length, mm	15.70 8.37
working distance, mm	4.9 0.87
Objectives only for bright field:	Achromat 90× Apochromats 60× 90×
initial magnification	90× 60× 90×
numerical aperture	1.25 0.05 1.30
focal length, mm	2.77 4.30 2.79
working distance, mm	0.44 0.22 0.20
Apochromatic objective only for dark field:	90×
initial magnification	90×
numerical aperture	1.00
focal length, mm	2.77
working distance, mm	0.68
Compensating eyepieces:	Focal length, Linear field of view, mm
3×	83 20
5×	50 20
7×	36 18
10×	25.2 13
15×	16.8 12
20×	12.6 9
Huyghenian eyepieces:	Focal length, Linear field of view, mm
4×	62.8 20
7×	35.0 18
10×	25 14
15×	17 8
Homal eyepieces:	
Type II	70.36 15
Type IV	20.28 8
Type VI	37.61 13
Illuminating system	arc lamp or incandescent bulb with illuminator
Object stage	circular, with rotary and cross-wise movements
Focusing	object stage adjustment: coarse — by rack and pinion; fine — by micrometric mechanism, reading 0.002 mm
Plate size, cm	13 × 18 or 9 × 12
Set of light filters	yellow, green, orange, yellow-green, sky-blue and blue
Max. allowable load on object stage, kg	10
Overall dimensions of microscope (length × width × height) mm	1860 × 630 × 490
Overall dimensions of case, mm	2100 × 550 × 900
Weight of microscope, kg	103.6
Weight of microscope in case, kg	195.6

MICROSCOPE SET

Table with two cabinets and optical bench.
 Microscope with inclined tube for visual observation, object stage, illuminator with swing-out polarizer, coarse and fine focusing mechanisms, objective holder and adapter tube for photographic camera

Apochromatic objectives in cases: 15× × 0.30
 30× × 0.65
 60× × 0.95
 90× × 1.30
 90× × 1.00

Achromatic objectives in cases: 11× × 0.17
 18× × 0.30
 40× × 0.65
 90× × 1.25

Compensating eyepieces: 3×, 5× (2 pieces), 7× (3 pieces), 10×, 15×, 20×
 Huyghenian eyepieces: 4×, 7×, 10×, 15×
 Homal eyepieces: II, IV, VI
 Magnifying glass 8×
 Binocular attachment
 Stage micrometer, opaque, in case
 Lamp I with ring
 Lamp II with ring
 Lamp III with ring
 Additional ring for lamp
 Inserts for object stage (glass) (2 pieces)
 Inserts for object stage (metal) (3 pieces)
 Crosswire for centring stage
 Eyepiece attachment
 Holder for unstable objects
 Object stage clips (2 pieces)
 Vial with immersion oil
 Protecting cap for the objective socket
 Caps for photographic tube (2 pieces)
 Caps for the visual tube (3 pieces)
 Marking apparatus
 Arc lamp with clockwork mechanism
 Collector
 Rheostat, 3.5 ohm, 10A
 Opal electric bulb 40W, 120V (3 pieces)
 Swing-out lens in mount
 Hollow heat absorbing device
 Plate holder 13 × 18 cm (3 pieces)
 Holder with opaque screen, revolving discs with light filters and lamp housing
 Opal glass in frame
 Adapter
 Metal plate holders, 9 × 12 cm (6 pieces)
 Metallic diaphragms (3 pieces)
 Photographic camera
 Arc lamp carbons, 5 mm (100 pieces)
 Arc lamp carbons, 9 mm (100 pieces)
 Cord for connecting microscope illumination
 Squirrel-hair brush
 Flannel napkin
 Cambric or madapolam napkin
 Watchmaker's screwdriver (2 pieces)
 Can with lubricant
 Cover
 Description and instruction manual
 Certificate

VERTICAL METALLOGRAPHIC MICROSCOPE, MODEL MHM-6

The MHM-6 Microscope (Fig. 42) is a vertical metallographic microscope, designed for the examination of the microstructure of metals.

The Microscope provides for the examination of opaque objects with bright field or polarized light as well as for the photography of these objects.

The set of objectives and eyepieces provides for a total magnification up to 600 \times for visual observation and photography.

The design features of the microscope and its optical characteristics fully satisfy up-to-date requirements of metallographic research.

The Vertical Metallographic Microscope finds wide application in industrial, educational and scientific research laboratories.

As the Vertical Metallographic Microscope is designed for the examination of opaque objects, observation is carried out with reflected light.

Fig. 43 shows the optical system of the instrument.

Light from the lamp 1 passes through the condenser 2, aperture diaphragm 3, lens 4, field diaphragm 5, lenses 6 and 7 and falls on the plane glass plate 8 of the central illuminator. The plate is positioned at an angle of 45 $^{\circ}$ to the axis of the tube.

The object is illuminated, by reflected light, by means of this plate. A part of the light is reflected by the plate through the micro-objective 9 to the object 10. In

this way, the objective 9 is used not only for providing an image of the object but also as a part of the illuminating system.

Oblique illumination of the object can be effected by displacing and subsequently rotating the aperture diaphragm.

The condenser 2 forms an image of the light source 1 on the aperture diaphragm 3.

An image of the aperture diaphragm is formed on the end lens of the objective 9 by means of the lens system 6 and 7 and the plate 8.

Rays, reflected from the object, pass again through the objective 9 and plate 8 to the prism 11 from where they are directed to the eyepiece 12. The objective 9 forms an image in the focal plane of the eyepiece 12.

During photography the prism 11 is swung out of the microscope axis and the beam of rays passes through the photographic eyepiece 13 and is reflected by the


Fig. 42

MHM-6 Vertical Metallographic Microscope


Fig. 43

Optical System of the MHM-6 Microscope

mirror 14 on the opal glass 15 of the photographic camera. An image is formed on this glass.

For work with polarized light, the superimposed polarising filters are brought in the optical system. Filter 16, a polarizer, is inserted into the illuminating system before the illuminator lens 6 while the second filter 17, used as an analyzer, is superimposed on the eyepiece 12.

To create more uniform illumination of the object, the semi-opal plate 18 is arranged before the aperture diaphragm 3.

The interchangeable light filters 19 can be brought in the optical system to provide monochromatic illumination.

Fig. 44 shows the main parts of the Vertical Metallographic Microscope.

The instrument comprises three main assemblies: a) illuminating device 1 fastened on bracket 2; b) upper housing 3 with the illuminator tube 4, visual tube 5; object


Fig. 44

Main parts of the MHM-6 Microscope

stage 6 and the coarse and fine focusing mechanism 7 and 8; c) lower housing 9 with the photographic camera and the base 10.

The illuminating lamp is arranged inside a spherical housing 11 on which are mounted the condenser body 12 and the camera 13 with a set of light filters.

The lamp holder 14 is inserted into a sleeve of the centring ring 15. The screws 16 serve to center the lamp.

Button 17 adjusts the condenser along the optical axis to regulate the illumination of the object.

The handles 18 and 19 serve for adjusting the opening in the aperture and field diaphragms. The latter are iris diaphragms with openings that can be varied from 0.8 to 7.5 mm.

Head 20 swings the prism in the system for visual observation or swings it out for photography.

The instrument has provisions for separate adjustment of the object stage 6 by rotating the coarse focusing head 7 and of the illuminator and objective by rotating the fine focusing head 8.

The object stage is arranged so that the object is placed over the objective with the polished surface directed downward.

The stage is square and does not have a rotary movement but has cross-wise adjustment actuated by rotating the heads 21 and 22.

The object is held by the clips 23 and a special fixture for clamping unstable objects (not illustrated).

A clamping device, on the shaft of the coarse focusing head 7, is provided to prevent the object stage from running down of its own weight. This device is actuated by the handle 24.

The visual tube 5 is arranged at an angle of 75° to the upper housing.

The lower housing 9 contains the photographic camera for 9 × 12 cm plateholders and the opal glass covered by the lid 25.

The microscope is stored in a wooden case containing a full set of accessories.

SPECIFICATIONS

Range of total magnification:			
visual observation	up to 600×		
photography	up to 600×		
Achromatic objectives:			
initial magnification	9×	21×	40×
focal length, mm	18.2	8.4	4.6
numerical aperture	0.20	0.40	0.65
working distance, mm	8.7	1.9	0.66
Photographic eyepieces:			
magnification	6.5×	10×	
focal length, mm	38.28	25.0	
Symmetrical eyepiece for photography:			
magnification	15×		
focal length, mm	17.0		
Huyghenian eyepieces:			
magnification	7×	10×	15×
focal length, mm	36.17	24.76	16.9
Light source	electric bulb, 8V, 20W with vertical illuminator		

56

Focusing adjustments:	
coarse	by rack and pinion for object stage
fine	by micrometric mechanism moving the illuminator and objective, reading to 0.002 mm
Object stage	square with lateral and transverse movement and special circular attachment for rotating specimen
Plate size, cm	9 × 12
Overall dimensions of microscope (height × length × width), mm	480 × 410 × 300
Overall dimensions of case, mm	540 × 515 × 355
Weight of microscope, kg	17
Weight of microscope in case, kg	34

MICROSCOPE SET

Microscope stand with coarse and fine focusing mechanisms, object stage, illuminating system with a set of light filters, inclined tube for visual observation and photographic camera

Achromatic objective, 9× × 0.20, in case
 Achromatic objective, 21× × 0.40, in case
 Achromatic objective, 40× × 0.65, in case
 Huyghenian eyepiece for visual observation, 7×
 Huyghenian eyepiece for visual observation, 10×
 Huyghenian eyepiece for visual observation, 15×
 Eyepiece for photography, 6.5×
 Eyepiece for photography, 10×
 Eyepiece for photography, 15× (symmetric)
 Eyepiece attachment with light filter
 Holder for unstable objects
 Object stage clips (2 pieces)
 Glass stage micrometer in case
 Object stage inserts (3 pieces)
 Magnifier, 8×
 Bulb 8V, 20W (2 pieces)
 Can with lubricant
 Metal plateholder, 9 × 12 cm (3 pieces)
 Squirrel-hair brush
 Flannel napkin
 Transformer 110/127/220/8V
 Frame with opal glass
 Caps (2 pieces)
 Cover
 Case for accessories
 Description and instruction manual
 Certificate

57

ELECTRON MICROSCOPE, MODEL 9M-3

The best modern optical microscopes provide for a magnification of the object not exceeding 1400 \times . Magnification beyond this limit with an optical system does not reveal new details but merely make the pictures larger.

The resolving power of an optical microscope is limited to about one half of the wave length of the light employed.

With visible light, objects, having a size not less than 0.2 microns, can be separated.

Further penetration into the world of micro objects, inaccessible to the most highly perfected optical microscopes, became possible only with the advent of electron radiation in microscopy. Electron beams have wavelength many times shorter than light beams.

Instruments, by whose aid electron beams provide images of objects, have been called electron microscopes.

Modern electron microscopes provide an effective magnification up to 100000 \times . This allows objects, considerably more minute than in optical microscopes, to be observed and studied at ense.

The 9M-3 Electron Microscope (Fig 45) is constructed with full application of electromagnetic optics.

It provides for a magnification from 250 \times to 25000 \times for normal exposures of the object and when better exposures are obtained, up to 100000 \times .

By using one or another of the methods known for preparing the specimen, the electron microscope allows the observation of almost all of the objects studied under optic microscopes.

The Electron Microscope, Model 9M-3 can be successfully used in research on biological specimens, bacteria, viruses, colloidal solutions, dyes, silicates, thin metal films obtained by vaporization in a vacuum, rubbers, caoutchouc, plastics, structures of metals and their alloys, etc.

All of the component elements of the electron microscope are analogous to the corresponding elements in the optical microscope, but are replaced by electrical devices.

The source of light is replaced by a source of electrons and the glass lenses are replaced by magnetic ones.

The image, produced by the electrons, may be either viewed on a fluorescent screen, which is luminous under the action of the electrons, or it may be revealed by the photographic layers which are darkened by electrons in the same way as ordinary photographic plates are darkened by light rays.

Fig. 46 illustrates the electron optical ray path of the 9M-3 Electron Microscope.

The electron beam is produced by an "electron gun" comprising a tungsten filament 1 0.1 mm in diameter, heated by an electric current and the anode 2.

The electrons, leaving the filament, are accelerated by the electric field applied between the filament and the anode.

The electrons, which are further to be used in forming the image, pass through an aperture in the center of the anode.

Electrons, from the "electron gun" pass further on by inertia until they reach the field of the first magnetic condenser lens 3. The latter changes their direction so that they are focused on the specimen 4.

The origin of the contrasts in the image, as obtained in the Electron Microscope, is caused by the fact that various portions of the specimen diffuse the electrons in different manner.

Only those electrons, which are deviated through comparatively small angles when passing through the specimen, are focused on the first screen 5 located before the magnetic projection lens 6. The corresponding portions of the screen will be lighted.

If the electrons, when passing through the specimen, are deviated through large angles, they are, for the most part, cut off by the aperture diaphragm 7 inside the objective 8. The corresponding elements of the image on the screen will be dark.

The screen for the intermediate image has an opening in the center through which the electrons, corresponding to the part of the image in the zone of the opening, pass to the second magnifying electromagnetic lens.

The electrons are again focused in the magnetic field of the projector lens and a magnified image of the object is produced on the final image screen 9 located in the lower part of the microscope.

As the electrons are stopped by air molecules, all air must be evacuated from their path in the instrument. In this way, the Electron Microscope is a vacuum instrument with continuous air evacuation.

The Electron Microscope has been designed as a column erected on a special desk fastened to a cabinet behind the column.

The following main parts of the microscope are arranged in the column:

- illuminating system of the microscope comprising the electron source ("electron gun") and the magnetic condenser lens focusing the electrons on the object which is to be examined;

- specimen chamber with object stage; the design of the stage allows the object to be moved in two perpendicular directions, as well as to be tilted through an angle of 4° to either side for producing pairs of stereomicrographs;

- electromagnetic objective lens which provides the first magnification of 130 \times of the observed object on the first screen;

- electromagnetic projector lens which allows the image to be remagnified in a range from 2 \times to 192 \times . In this way, the total magnification of the microscope can be varied, by means of these lens from 250 \times to 25000 \times .

- photographic camera of the microscope with a fluorescent screen for obtaining the final image and the plate holder device located under the screen.

Fig. 45
9M-3 Electron
Microscope


Fig. 46
Electron Optical Ray Path

The photographic camera has three viewing ports for observation of the final image. They are arranged so as to provide convenient observation through any port.

During photography exposition is accomplished by means of the screen which, at the same time, fulfills the functions of a camera shutter.

Small viewing ports are arranged on the front of the column for watching the electron beam path in the microscope and for viewing the image on the first screen of the projector lens.

The objective lens and the photographic camera are furnished with special doors through which the specimen may be mounted on the object stage or the plate holders may be changed and removed from the microscope.

The vacuum equipment of the microscope consists of the microscope column, the preliminary vacuum pump, the oil-diffusion vacuum pump, the vacuum pipelines and a special distributing mechanism which provides the required vacuum arrangements during microscope operation.

The distributing mechanism is connected by vacuum pipelines to the microscope column and to both vacuum pumps. By means of the distributor, air entering the microscope first passes through a chemical dryer located under the desk.

The controls of the distributing mechanism is furnished with a special indicator which shows the vacuum conditions in the microscope.

All the elements of the electrical supply system of the microscope, except the resonant voltage stabilizer, are mounted in the microscope cabinet.

The resonant voltage stabilizer, whose operation is accompanied by strong stray magnetic fields, should be sufficiently distant from the microscope (it is usually erected in a neighbouring room).

The design and wiring of the electric supply system is arranged in five independent units: distribution switchboard, control desk, low voltage unit, high voltage unit and the high frequency transformer. The power supplies are divided into two main lines: low voltage and high voltage supplies.

The main part of the supply circuit is connected to a 220V three-phase line through the resonant voltage stabilizer.

The low voltage supplies feed the "electron gun" filament, the magnetic lens coils as well as the vacuum pump electric motor, the oil-diffusion pump heater, the thermo-couple vacuum gauge and all the electric bulbs for illuminating the scales of the instrument.

The high voltage supplies provide the Electron Microscope with a voltage of 50 kV in steps of 30, 40 and 50 kV.

The microscope cabinet is furnished with a special discharging device brought into action when any of the cabinet doors or the "electron gun" cowlings is opened. This device serves to protect the observer in accidental cases or erroneous operation.

SPECIFICATIONS

Optical characteristics of the instrument

Electronic-optical magnification:	
without projector lens pole pieces	from 250× to 1100×
with low magnification pole pieces	from 2000× to 10000×
with high magnification pole pieces	from 5000× to 25000×
Resolving power of the microscope	up to 20 Angstrom units
Useful photographic enlargement	up to 100000×

60

Power supply characteristics

Line voltage supplied to microscope	3-phase, 220V, 50 c. p. s.
Accelerating voltage of Electron Microscope, kV	first step — 30 second step — 40 third step — 50
High tension current, mA	up to 170
Total power consumption, kW	up to 2

Research Facilities

Visual observation and photography of the image of the object, observed in transmitted electron beam.

Stereoscopic photography at a stereoscopic angle of 8°. Electron diffraction patterns of the specimen area can be made with an electron diffraction unit.

Photographic Arrangement

Four exposures, mm	45 × 45
--------------------	---------

Vacuum Characteristics

Operating vacuum	1.10 ⁻¹ to 5.10 ⁻¹ mm Hg
Time for creating operating vacuum with normal cycle (air — preliminary vacuum — operating vacuum)	3—4 minutes after replacing specimen; 10—12 min. after replacing photographic plates

Overall Dimensions and Weight of Instrument

Height of microscope, mm	2080
Width of microscope (without preliminary vacuum pump), mm	630
Width of microscope with pump, mm	870
Length of microscope (cabinet and desk), mm	870
Length of microscope with pump, mm	1300
Weight of complete microscope set, kg	980

ELECTRON MICROSCOPE SET

Microscope column
Microscope desk with vacuum distributing mechanism
Microscope cabinet
Control panel
Low-voltage unit
High-voltage unit
High-frequency transformer
Distribution switchboard
Oil-diffusion vacuum pump
Electron diffraction unit
Preliminary vacuum pump
Resonant voltage stabilizer
Vibration insulators (7 pieces)
Plate holders (2 pieces)
Magnifiers in mounts (3 pieces)
Object stage with objective lens pole pieces
Low magnification pole pieces for projector lens
High magnification pole pieces for projector lens
Steel needle in mount
Brass needle in mount
Stage for drying specimens

61

Fixture for assembling "electron gun" housing
Box with caps for the specimens (9 pieces)
Ampoules with oil for oil-diffusion pump (3 pieces)
Fixture for welding filaments
Tank with oil for preliminary vacuum pump
Set of tools (wrenches) (15 pieces)
Diaphragms (10 pieces)
Holders for specimens with caps (5 pieces)
Filament assemblies (25 pieces)
Bushing with diaphragm
Large valve
Small valve
Resistors (2 pieces)
Spare bellows (8 pieces)
Sockets for electron diffraction unit (8 pieces)
Rubber gaskets and rings (30 pieces)
Fuses for 10A, 5A, 2A and 0.5A (10 pieces)
Screen axle collars (2 pieces)
Porcelain rods (2 pieces)
Straps for clamping rings (10 pieces)
Specimen grids (500 pieces)
Electric heater coils for oil-diffusion pump (2 pieces)
Diaphragms for low magnification pole shoe (2 pieces)
Diaphragms for high magnification pole shoe (2 pieces)
Screws and bolts (35 pieces)
Tungsten wire, 0.1 mm diameter (2 m)
Rubber hose (20 m)
Discharge tube rectifiers, type ИГ-129 (4 pieces)
Electronic tubes, type 6X17 (4 pieces)
Electronic tubes, type 6I13 (10 pieces)
Stabilivolts (2 pieces)
Kenotrons, type B 40/100 (3 pieces)
Neon tube, type MIL-7
Electric bulbs for illuminating scales (3 pieces)
Thermoelectric lamp, type ИТ-2
Automobile bulbs, type А-16 (3 pieces)
Boxes for accessories (3 pieces)
Storage cases (4 pieces)
Description and instruction manual
Certificate

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"STANKOIMPORT"

EXPORTS AND IMPORTS

Machine Tools
Woodworking Machinery
Metal Working Machinery (presses, hammers, shears, cold forming machines, punching machines)
Rolling Mills (imports)
Measuring Instruments and Apparatus (for metal industry)
Testing Machines and Instruments (for metals)
Optical Instruments and Equipment
Portable Electric and Pneumatic Tools (for metal and wood-working)
Metal and Wood Cutting Tools
Mechanic's Tools and Chucks
Sintered Carbide and Hard-Alloy Products
Abrasive Products
Ball and Roller Bearings
Microscopes of all types
Motion-Picture Equipment
Photographic Cameras
Binoculars
Magnifiers
Lenses
Crude Optical Glass Blocks and Blanks

All inquiries and correspondence to be forwarded to:
Vsesojusnoje Exportno-Importnoje Objedinenije
"Stankoimport"
32/34, Smolenskaja-Sennaja pl., Moscow, USSR.
For cables: Stankoimport Moscow

Design and specifications of microscopes illustrated herein are subject to change without notice

CABLE ADDRESS:
STANKOIMPORT
MOSCOW

UNCLASSIFIED

GEODETTIC INSTRUMENTS

UNCLASSIFIED

STAT

STANKOIMPORT

UNCLASSIFIED

1/3

STAT

UNCLASSIFIED

STAT

UNCLASSIFIED

STAT

Geodetic Instruments, produced in the U. S. S. R., are the result of many years of earnest study of the requirements of the engineering profession and the application of modern scientific methods in their manufacture.

The first class materials, perfect modern design and skillful workmanship provide accuracy and durability of the instruments.

Optical properties of Geodetic Instruments are famous and are maintained at the highest possible standard of excellence.

Their fine appearance is due to careful manufacture and decorative finish of various parts and surfaces of the instruments

UNIVERSAL ASTRONOMICAL THEODOLITE, MODEL AY 2/10

The Universal Astronomical Theodolite (Fig. 1) is a high-precision astronomical geodetic instrument designed for astronomical observations and for measuring horizontal and vertical angles at points of first order triangulation. In first order triangulation, not only the astronomical latitudes and longitudes


Fig. 1

are determined for intermediate astronomical points but also their astronomical azimuth in reference to one of the adjoining trigonometric points. The astronomical azimuth, in conjunction with the astronomical longitude of the point whose azimuth has been determined, is of vital importance for controlling

measurements in triangulation. Readings on the horizontal circle are made with the aid of two micrometer microscopes with drum scale divisions of 2". The vertical circle is read by means of two vernier microscopes having scale divisions of 10".

The main telescope is of the astronomical, broken-line, central type comprising a two-lens objective, a rectangular prism and a rotary eyepiece micrometer with two interchangeable eyepieces.

The Universal Astronomical Theodolite consists of the following main parts:

- a) Base of instrument (lower part) carrying the horizontal circle and the alidade with the two micrometer microscopes.
- b) Main telescope (upper part) with the vertical circle, vernier microscopes and the Talcott spirit level.
- c) Striding spirit level.
- d) Auxiliary telescope.

The high-quality materials used in manufacturing of the instrument (bronze, brass) ensure its anticorrosive properties. The vital parts, such as the limb, base, lower movement and horizontal axis undergo a special heat treatment to ensure stability of dimensions, one of the most important factors effecting the accuracy of the instrument.

The design of the instrument provides for reliable stability in operation and high accuracy in measurements.

The instrument is stored in two cases.


Fig. 2

The lower movement and base, striding spirit level, screw drivers, brush, studs, wrench and levelling shoes are stored in specially provided recesses of one case (Fig. 2).

The other case (Fig. 3) contains the main telescope (upper part of instrument), auxiliary telescope, blind, spare eyepiece, hand hammer and a special box with electrical accessories. All these items are fitted into specially provided recesses.


Fig. 3

SPECIFICATIONS

Main telescope

Aperture of objective	55 mm
Resolving power of objective	2.50"
Focal length of objective	450 mm ± 1%
Focal lengths of eyepieces	8 and 10 mm
Total magnification of telescope	56× and 45×
Diameters of exit apertures	1.0 and 1.2 mm
Field of view	0° 54'
Telescope focusing range	from 5 m to infinity
Distances to exit apertures	1.4 and 2.76 mm
Spider cross-hairs comprising 9 stationary hairs with an angular interval	90" ± 5"
and adjustable bisector and hairs; width of bisector	25" — 30"
Distance between bisector and parallel hair	115" ± 3"
The drum of the eyepiece micrometer has 100 divisions	
Value of micrometer drum scale divisions	1.0"
Angle of rotation of eyepiece micrometer together with cross-hairs	up to 90°
The angle of rotation is set to a positioning circle having a sector scale of 90° with 1° divisions	
Maximum elevation of telescope with striding level	73°

Auxiliary telescope

Auxiliary telescope — astronomical, straight type comprising a two-lens objective and an eyepiece furnished with a micrometer	
Aperture of objective	36 mm
Resolving power of objective	4"
Focal length of objective	360 mm ± 1%
Focal length of eyepiece	12 mm
Total magnification of telescope	30×
Diameter of exit aperture	1.2 mm
Field of view	1°
Telescope focusing range	from 7 m to infinity
Distance to exit aperture	5.0 mm

Spider cross-hairs consisting of three horizontal hairs at an angular interval..... $150^{\circ} \pm 10^{\circ}$ and adjustable vertical bisector with an angular distance..... $30^{\circ} - 35^{\circ}$
 The eyepiece micrometer drum has 100 divisions
 Value of drum divisions..... 1.4°

Reading micrometer microscopes for horizontal circle

Magnification of micrometer microscopes..... $40\times$
 Pitch of micrometer screws..... 0.25 mm
 Microscope field of view (visible portion of circle)..... 3.1 mm
 this corresponds to a limb reading of $1^{\circ}37'$
 Value of drum scale divisions..... 2°
 The drum has 60 divisions
 2.5 rotations of the screw correspond to a cross-hair bisector movement through 1 limb division
 The bisector consists of two pairs of spider cross-hairs
 Distance between axes of bisector..... $4' \pm 2^{\circ}$
 Width of bisectors..... $53^{\circ} \pm 2^{\circ}$
 The micrometer cases have engraved inscriptions:
 "A" and "B" for horizontal circle microscopes and "I" and "II" for vertical circle microscopes

Scale reading microscopes for vertical circle

Magnification of microscopes..... $28\times$
 Focal length of objective..... 29.96 mm
 Focal length of eyepiece..... 13.5 mm
 Microscope field of view (visible portion of circle)..... 4 mm
 this corresponds to a limb reading of $3^{\circ}10'$
 Distance to exit aperture..... 6.3 mm
 The vernier scale of the microscopes has 30 divisions corresponding to 29 circle divisions

Horizontal and vertical circles

Circle scale diameter (to outer ends of graduations):
 horizontal circle..... 220 mm
 vertical circle..... 135 mm
 Smallest division on circles..... $5'$

For approximate adjustment of the alidade section, one-degree divisions are engraved on the horizontal circle. Readings on the one-degree scale are the same as microscope readings with a tolerance of $\pm 5'$

Vertical axis is of the Repsold conical type
 The adjusting device, a knob, has a screw with a pitch of 0.25 mm and a head with 50 divisions

Spirit levels

Chamber-type vial
 Value of striding spirit level divisions..... $2'' - 2.5''$
 to 2 mm of arc
 Value of Talcott spirit level divisions..... $1.5''$ to 2 mm
 of arc

Overall dimensions and weight

Height of theodolite..... 500 mm
 Height of tripod..... 1100 mm
 Overall dimensions of storage case for lower part of instrument..... $428 \times 488 \times 428\text{ mm}$
 Overall dimensions of storage case for upper part of instrument..... $488 \times 478 \times 208\text{ mm}$
 Weight of tripod..... 13.2 kg
 Weight of theodolite..... 37.5 kg
 Weight of lower part of instrument in storage case..... 42.7 kg
 Weight of upper part of instrument in storage case..... 28.2 kg

Attachments

Interchangeable eyepiece for main telescope for a magnification of $56\times$
 Illuminators for field of view of main and auxiliary telescopes (2 pcs.)
 Caps for objectives of main and auxiliary telescopes (2 pcs.)
 Dark glasses for main telescope eyepiece (2 pcs.)
 Blind with 3 diaphragms of various size
 Shoes (3 pcs.)
 Cords with plugs and contacts (2 pcs.)
 Hand hammers (2 pcs.)
 Pencil illuminator
 Spare electric light bulbs (20 pcs.)
 Spare flash-light bulbs (4 pcs.)
 Large screw driver
 Small screw driver
 Watchmaker's screw driver
 Straight studs (2 pcs.)
 Bent studs (2 pcs.)
 Wrench for axial nut (adjustable)
 Wrench
 Soft brush
 Chamois leather, $200 \times 200\text{ mm}$
 Oil can with oil; in case
 Sectional-type tripod
 Circular spirit level for tripod
 Wrench for tripod nuts
 Storage case for lower part of instrument
 Storage case for upper part of instrument
 Storage case for tripod
 Box with tripod accessories (stored in tripod case)
 Packing cases for lower and upper parts of instrument with shock-absorber devices (2 pcs.)
 Spare keys for storage and packing cases of lower and upper parts of instrument (4 pcs. fitted in special recesses on the bottoms of the cases)
 Canvas covers for storage cases of lower and upper parts of instrument (2 pcs.)
 Cloth hood for the instrument
 Certificate and operating instructions for the instrument

TRIANGULATION THEODOLITE, MODEL TT 2/6


Fig. 4

The Triangulation Theodolite (Fig. 4) is a highly accurate geodetic instrument designed for measuring horizontal angles and zenith distances at first order triangulation points.

The survey of the stations of the geodetical base-line is carried out by triangulation.

Points are selected, in the region to be surveyed, that are located at a considerable distance from each other. When connected by straight lines (sighting lines), the system of points should form a number of triangles; each

triangle being made as nearly equilateral as the conditions will permit. With the aid of this Theodolite, which is an instrument unsurpassed in accuracy, all the angles of the triangles are measured. When laying out a system or chain of triangles for highly accurate triangulation, the geographic co-ordinates (the latitude and longitude) of the apexes of the triangles can be determined by the aid of this instrument, on the basis of astronomical-geodetic observations.

In design, the instrument is sufficiently stable. This is one of the factors ensuring its high accuracy.

The instrument is manufactured of high quality metals and non-ferrous alloys.

The lacquered finish of the instrument as well as the golden lacquer coating on the geodetic micrometers enhance its appearance.

Readings, on the horizontal circle, are made by the aid of two micrometer microscopes having drum scale division values of 2". On the vertical circle, the readings are taken by means of scale microscopes with division values of 6". The instrument is furnished with electrical illumination for use in night surveying.

The main telescope of the instrument provides for measuring angles between signals located at a distance up to 60 km from the observer.

The Triangulation Theodolite comprises the following main parts:

Base of instrument (lower part) carrying the horizontal circle and the alidade with the micrometer microscopes;

Main telescope (upper part) with vertical arc and scale microscopes;

Striding spirit level;

Auxiliary telescope.


Fig. 5

The instrument is stored in two cases. The lower movement and base, striding spirit level, screw drivers, brush, studs, wrench, shoes and oil can with oil are stored in specially provided recesses of one case (Fig. 5).

The other case (Fig. 6) contains the main telescope with the horizontal axis, auxiliary telescope, blind, detachable mirror, spare eyepiece, hand hammer, light filters and a special box with electrical accessories. All these items are fitted into special recesses.


Fig. 6

SPECIFICATIONS

Main telescope

Aperture of objective	65 mm
Resolving power of objective	2.15"
Focal length of objective	520 mm
Focal lengths of eyepieces	8 and 10 mm
Total magnification of telescope	52 \times and 65 \times
Diameters of exit apertures	1.0 and 1.25 mm
Fields of view	0°43' and 0°37'
Telescope focusing range	from 5 m to infinity
Distance to exit aperture	3.2 and 4.0 mm
Spider cross-hairs consisting of three horizontal hairs with an angular interval of 8'35" and adjustable vertical bisector with an angular hair interval of 25"—30"	
The drum of the eyepiece micrometer has 100 divisions	
Value of micrometer drum scale divisions	1.0"

Auxiliary telescope

Auxiliary telescope — astronomical, straight type comprising a two-lens objective and an eyepiece with a micrometer	
Aperture of objective	36 mm
Resolving power of objective	4"
Focal length of objective	360 mm
Focal length of eyepiece	12 mm
Total magnification of telescope	30 \times

10

Diameter of exit aperture	1.2 mm
Field of view	1°
Telescope focusing range	from 7 m to infinity
Distance to exit aperture	8.0 mm
Spider cross-hairs consisting of two horizontal hairs at an angular interval of	150"
and an adjustable vertical bisector with an angular interval of ..	30"—35"
The drum of the eyepiece micrometer has 100 divisions	
Value of micrometer drum scale divisions	1.4"

Reading micrometer microscopes for horizontal circle

Magnification of micrometer microscopes	40 \times
Focal length of objective	28.21 mm
Focal length of eyepiece	20.0 mm
Microscope field of view (visible portion of scale)	3.1 mm (1°37' on limb)
Distance to exit aperture	0.56 mm
Value of micrometer drum divisions	2"
The drum has 60 divisions	
2.5 rotations of the screw correspond to a cross-hair bisector adjustment of 1 limb division	
The bisector consists of two pairs of spider cross-hairs	
Distance between axes of bisector	4'
Width of bisectors	53"
Pitch of micrometer screws	0.25 mm

Scale reading microscopes for vertical arc

Magnification of microscopes	48 \times
Focal length of objective	15.6 mm
Focal length of eyepiece	10.0 mm
Microscope field of view (visible portion of scale)	3 mm (2°9' on vertical arc)
Distance to exit aperture	4.3 mm
The microscope scale has 10 divisions	
Ten scale divisions correspond to one arc scale division;	
reading accuracy	6"

Horizontal circle and vertical arc

Circle scale diameters (to outer ends of graduations):	
horizontal circle	220 mm
vertical arc	160 mm
Smallest division of horizontal circle	5'
Smallest division of vertical arc	10'
For approximate adjustment of the alidade section, one-degree divisions are engraved on the horizontal circle	
Readings on the one-degree scale are the same as microscope readings with a tolerance of $\pm 5'$	
Vertical axis is of the Repsold conical type	
The adjusting device, a knob, has a screw with a pitch of 0.25 mm and a head with 50 divisions	

11

Spirit levels

Value of striding spirit level divisions	2"—2.5" to 2 mm of arc
The striding spirit level vial is of the chamber type	
Value of vertical arc spirit level divisions	6"—10"

Overall dimensions and weight

Overall dimensions of storage case for lower part of instrument ..	420×425×400 mm
Overall dimensions of storage case for upper part of instrument ..	600×300×210 mm
Overall dimensions of packing case for lower part of instrument ..	660×660×605 mm
Overall dimensions of packing case for upper part of instrument ..	850×390×450 mm
Weight of lower part in storage case	27 kg
Weight of upper part in storage case	17 kg
Weight of lower part in storage and packing cases	59 kg
Weight of upper part in storage and packing cases	38 kg

Attachments

Interchangeable eyepiece for main telescope for a magnification of 65×	Adjustable wrench
Illuminators for field of view of main and auxiliary telescopes (2 pcs.)	Soft brush
Caps for objectives of main and auxiliary telescopes (2 pcs.)	Chamois leather, 200×200 mm
Mirror for vertical arc spirit level	Oil can with oil, in case
Blind with 3 diaphragms of various size	Storage cases for upper and lower parts of instrument (2 pcs.)
Shoes (3 pcs.)	Packing cases for upper and lower parts of instrument with shock-absorbing devices (2 pcs.)
Cords with plugs and contacts (2 pcs.)	Spare keys for storage and packing cases of upper and lower parts of instrument (4 keys fitted into special recesses on the outside of the bottoms of the cases)
Hand hammers (2 pcs.)	Canvas covers for storage cases of upper and lower parts of instrument (2 pcs.)
Pencil illuminator	Cloth hood for the instrument
Spare electric light bulbs (20 pcs.)	Certificate and operating instructions for the theodolite
Spare flash-light bulbs (4 pcs.)	
Large screw driver	
Small screw driver	
Watchmaker's screw driver	
Straight studs (2 pcs.)	
Bent studs (2 pcs.)	

OPTICAL THEODOLITE, MODEL OT-02

Fig. 7

The Optical Theodolite (Fig. 7) is designed for measuring horizontal and vertical angles at higher order triangulation and polygonometric points, as well as for astronomical observations.

The small size, hermetic design, and comparatively small weight of the instrument, in conjunction with the speed and convenience of its operation provide for the possibility of its wide application in geodetic surveys of difficult, mountainous or distant regions. Besides this, the instrument can be advantageous-

ly used in industry (machine-tool building, underground construction, etc.) where the accurate measurement of horizontal and vertical angles is required.

Readings on the horizontal and vertical limbs are combined by the aid of the optical systems of the horizontal and vertical circle microscope objectives into a single field of view of the reading microscope whose eyepiece is arranged side by side with the telescope eyepiece. Readings are taken with a single optical microscope having a seconds disc scale divisions value of 0.2".


The instrument can be used the whole year around, as it provides for normal operation at temperatures from -25°C to $+50^{\circ}\text{C}$ (from -13°F to $+122^{\circ}\text{F}$).

The instrument has electric illumination.

For making astronomical observations, the instrument is furnished with a removable prism attachment having a dark light filter.


The Optical Theodolite comprises three main parts: lower part of the theodolite, intermediate part and telescope with horizontal axis.

The telescope is the central, astronomical type with interior focusing. It consists of a telephotolens and three interchangeable eyepieces.

Limb reading — $14^{\circ}20'$
 1st second scale reading 20.07"
 2nd second scale reading 20.06"
 RESULT $14^{\circ}20' 40.13''$

Fig. 8

Limb reading — $156^{\circ}38'$
 1st second scale reading 40.27"
 2nd second scale reading 40.21"
 RESULT $156^{\circ}39' 20.48''$

Fig. 9

Figs. 8 and 9 illustrate the field of view of the reading microscope. Rotating the head actuates the optical micrometer for achieving an exact coincidence of the limb graduations. When the graduations coincide, the degrees and minutes are read on the upper scale. The lower (seconds) scale reading is taken and is multiplied by 2. An alternate method (to increase the accuracy of coincidence of the graduations) is to repeat the coinciding of the graduations and to add the two readings taken on the seconds scale.

14

SPECIFICATIONS

Main telescope

Visible magnification of telescope	24x, 30x and 40x
Aperture of telephotolens	60 mm
Focal length of telephotolens	350 mm
Fields of view	$1^{\circ}40'$; $1^{\circ}20'$ and $1^{\circ}0'$
Resolving power of telephotolens	2.4"
Focal lengths of eyepieces	8.6; 11.7 and 14.6 mm
Diameters of exit apertures	1.5; 2.0 and 2.5 mm
Telescope focusing range	from 5 m to infinity
Length of telescope	265 mm
Angular interval of bisector	35"
Thickness of cross-hairs	7—8 microns
Maximum inclination of telescope	65°

Reading devices:

for horizontal limb — one optical system of the horizontal circle microscope objective with magnification of 4x;

for vertical limb — one optical system of the vertical circle microscope objective with magnification of 3x;

for horizontal and vertical limbs — one reading microscope with a magnification of 9.3x; consequently the total magnification of the microscope is:

for the horizontal limb	37x
for the vertical limb	28x

The reading microscope has a field of view of 8×3.8 mm which corresponds to:

on horizontal limb	$1^{\circ}40'$
on vertical limb	$3^{\circ}20'$

Distance to exit apertures of reading microscope:

for horizontal limb	25.5 mm
for vertical limb	25.0 mm

Diameters of exit apertures of reading microscope:

for horizontal limb	1.7 mm
for vertical limb	1.6 mm

The horizontal and vertical limbs are made of optical glass

Scale circle diameters (to inner ends of graduations):

horizontal limb	135 mm
vertical limb	90 mm

Smallest division of limbs:

horizontal	4'
vertical	8'

Thickness of limb graduation lines:

horizontal limb	6—7 microns
vertical limb	7—8 microns

15

Axes

Vertical axis — cylindrical, self-adjusting type with a ball support on the tapered part of the box
Horizontal axis — cylindrical sectional type with internal support on cylindrical bearings.

Spirit levels

Value of horizontal limb alidade level scale divisions to 2 mm of arc... 6"—7"
Value of vertical limb alidade level scale divisions to 2 mm of arc... 10"—12"

Overall dimensions and weight

Height of theodolite 404 mm
Overall dimensions of storage case..... diam. 233×423 mm
Weight of theodolite 11 kg
Weight of theodolite in storage case 15 kg

Attachments

Illuminating mirrors (2 pcs.)	Wrench
Illuminating lamps (2 pcs.)	Stud for adjustment screw
Eyepiece attachment	Various screw drivers
Interchangeable eyepieces 24× and 40× (2 pcs.)	Studs (2 pcs.)
Socket plug	Adjustable wrench
Cap for objective	Brush
Centering plate	Flannel napkin
Tripod	Can of oil
Metal sheath	Electric light bulbs (10 pcs.)
Device for carrying the theodolite	Hoods (2 pcs.)
Hand hammer	Bag for plate and accessories
	Certificate and operating instructions for the theodolite

OPTICAL THEODOLITE, MODEL OT-10


Fig. 10

The Optical Theodolite (Fig. 10) is designed for various types of geodetic surveys and can be used for third order triangulation, second and third order polygonometry as well as for running theodolite traverses.

The horizontal and vertical limbs are made of optical glass.

Readings on the horizontal and vertical limbs are combined in the field of view of a micrometer microscope located on the standard opposite the vertical circle.

Readings on the limbs are taken by a single micrometer microscope whose smallest drum scale division equals 20".

The micrometer microscope tube and the telescope can be inverted through the zenith.

Figs. 11 and 12 illustrate the field of view of the reading microscope and the reading drum of the micrometer screw.

Horizontal circle readings are taken as follows (Fig. 11):


Fig. 11

1. On the upper image of the horizontal circle scale, in reference to the stationary index, read the degrees and full 20-minute intervals on the circle;

2. Rotate the micrometer drum until the vertical adjustable bisector coincides with the upper graduation of the circle to the left of the index and take the drum scale reading;

3. Rotate the micrometer drum until the vertical adjustable bisector coincides with the lower graduation of the circle to the left of the index and take the drum scale reading;

4. Add the circle and micrometer readings.

Limb reading	245° 40'
1st reading on drum scale	4' 25"
2nd reading on drum scale	4' 23"

RESULT 245° 48' 48"


Fig. 12

Vertical circle readings are taken as follows (Fig. 12):

1. Rotate the micrometer drum until the adjustable bisector coincides with the stationary index;

2. Using the horizontal bisector as an index, read the number of degrees and full 20-minute intervals (below the bisector) on the vertical circle;

3. Rotate the micrometer drum until the adjustable horizontal bisector coincides with the lower graduation on the vertical circle and take the micrometer drum scale reading;

4. Repeat by coinciding the bisector on the same circle graduation;

5. Add the readings.

Reading on vertical circle	2° 20'
1st reading on drum scale	6' 32"
2nd reading on drum scale	6' 30"

RESULT 2° 33' 02"

SPECIFICATIONS

Magnification of telescope	25.3x
Telescope field of view	1° 10'
Aperture of objective	40 mm
Diameter of exit aperture	1.5 mm
Focal length of objective	253 mm
Telescope focusing range	from 1.5 m to infinity
Value of micrometer drum scale divisions for vertical and horizontal circles	20"
Value of smallest division of vertical and horizontal circles	20"
Value of horizontal circle alidade level scale divisions to 2 mm of arc	40"—60"
Value of vertical circle alidade level scale divisions to 2 mm of arc	15"—20"

Overall dimensions and weight

Total height of theodolite	280 mm
Length of tripod with legs retracted	970 mm
Maximum length of tripod	1510 mm
Overall dimensions of theodolite case	340x280x340 mm
Overall dimensions of case with shock-absorber device for transporting the theodolite	410x380x520 mm
Weight of theodolite with levelling base	6.3 kg
Weight of theodolite in storage case with accessories	17 kg
Weight of tripod	5.8 kg
Weight of case with shock-absorber device for transporting the theodolite	13 kg
Weight of complete outfit	35.8 kg

Attachments

- | | |
|---|---|
| Eyepiece prism | Sun blind |
| Theodolite case with shoulder straps | Canvas case for theodolite |
| Extensible tripod with fastening screw | Oil can with oil |
| Case with shock-absorber device for transporting the theodolite | Adjustable wrench |
| Vertical circle spirit level | Tripod wrench |
| Horizontal circle spirit level | Screw driver with 4 blades |
| Plumb bob hook | Studs for screws (2 pcs.) |
| Plumb bob with counter-weight, hook and cord | Brush |
| | Napkin, 200×200 mm |
| | Certificate and operating instructions for the theodolite |

OPTICAL THEODOLITE, MODEL ТБ-1


Fig. 13

The Optical Theodolite (Fig. 13) is designed for measuring angles in second and third order triangulation, for astronomical observations as well as for measuring angles in first and second order polygonometry.

Readings on the horizontal and vertical limbs are combined, by the aid of the optical systems of the horizontal and vertical limb microscope objectives, into a single field of view of the reading microscope whose eyepiece is arranged side by side with the telescope eyepiece.

Readings are taken with a single optical micrometer to an accuracy of 1".

The theodolite has stadia hairs for measuring distances with the aid of horizontal or vertical stadia rods.

The instrument can be fitted for night operations.

The instrument is furnished with zenith attachments on the telescope and microscope for astronomical observations.

Due to its small size and comparatively small weight, as well as the ease in reading the horizontal and vertical limbs, this instrument provides for the possibility of achieving highly efficient performance.

Besides geodetic and astronomical operations, this instrument can be used for measuring horizontal and vertical angles in industrial enterprises, in the erection of various types of construction, etc.

SPECIFICATIONS

Telescope

Magnification of telescope	25.6x
Field of view	1° 30'
Diameter of exit aperture	1.5 mm
Distance of exit aperture from last surface of eyepiece	6.7
Minimum sighting distance	1.2 m
Resolving power	4"
Stadia constant coefficient	100
Magnification of eyepiece	25.5
Adjustment of eyepiece	± 5 diopters
Focal length of objective with focusing lens	240.7 mm
Focal length of eyepiece	9.8 mm

Microscopes

	Horizontal limb diam. 85 mm	Vertical limb diam. 75 mm
Magnification	47x	52.3x
Diameter of exit aperture	1.3 mm	1.7 mm
Distance of exit aperture	12 mm	12 mm
Adjustment of eyepiece	± 5 diopters	± 5 diopters

Optical plumbing device

Magnification	1.4x
Field of view	8° 17'
Diameter of exit aperture	4 mm
Distance of exit aperture	7 mm
Minimum sighting distance	0.7 m
Adjustment of eyepiece	± 5 diopters

Spirit levels

Value of horizontal limb alidade cylindrical level divisions	12"—20" to 2 mm of arc
Value of vertical limb cylindrical level divisions	17"—25" to 2 mm of arc
Value of spherical spirit level divisions	7"—19" to 2 mm of arc

Overall dimensions and weight

Height of instrument with tripod	2035 mm
Overall dimensions of case	285×210×400 mm
Overall dimensions of box with storage batteries	238×165×180 mm
Weight of theodolite without case	5.1 kg
Weight of theodolite in case and with accessories	9.1 kg
Weight of tripod with sighting rod	6.25 kg
Weight of box with storage batteries	5.9 kg

Attachments

Tubular surveying compass	Zenith attachment for telescope
Bridge for sighting rod	Zenith attachment for microscope
Brightening chamber	Storage battery with box and cable
Tripod with fastening screw	Set of accessories
Sighting rod	Certificate and operating instructions for using the theodolite
Centering device	

THEODOLITE TACHEOMETER, MODEL TT-50


Fig. 14

The Theodolite Tacheometer (Fig. 14) is an angle measuring instrument designed for measuring horizontal and vertical angles with a reading accuracy of 30" as well as for measuring distances by the stadia method using stadia rods.

This theodolite finds its widest application in geotopographic surveys.

In design the Theodolite Tacheometer is of the repeating type of theodolites. A surveying compass is provided with the theodolite for orientation in reference to the magnetic meridian.

In operation, the theodolite is set up on the tripod and fastened with the locking screw. The theodolite is centered with a plumb-line.

The optical system of the theodolite has coated ("blue") lenses. This increases the transmission of light and facilitates observation under conditions of poor illumination.

The Theodolite Tacheometer ensures normal operation at temperatures from -40°C to $+45^{\circ}\text{C}$ (from -40°F to $+113^{\circ}\text{F}$).

SPECIFICATIONS

Telescope magnification	25.3x
Telescope field of view	1° 10'
Stadia constant coefficient	100
Resolving power of objective, not over	4.5"
Telescope focusing range	from 1.5 m to infinity
Value of horizontal circle spirit level divisions	40"—60"
Value of vertical circle spirit level divisions	25"—40"
Reading accuracy on horizontal and vertical circles	30"

Overall dimensions and weight

Overall dimensions of theodolite case	340x260x340 mm
Total weight of the theodolite outfit with accessories, case and tripod ..	16.7 kg

Attachments

Surveying compass (striding type)	Metal screw driver with four blades
Tripod	Studs for screws (2 pcs.)
Theodolite case with shoulder straps	Tripod wrench
Plumb-bob with counter-weight, hook and cord	Conical center nut wrench
Light filter in mount	Oil can with oil
Blind	Brush
Canvas case for theodolite	Napkin
Adjustable wrench	Certificate and operating instructions for the theodolite

MINE THEODOLITE, MODEL TF-1


Fig. 15

The Mine Theodolite (Fig. 15) is designed for the measurement of horizontal and vertical angles. The arrangement of a spirit level on the telescope allows the instrument to be used for levelling operations.

Linear distances can be measured by using the stadia cross-hairs.

The theodolite is used for all underground and surface mine surveying, for various topographical work as well as for astronomical-geodetic observations of the corresponding accuracy. For the latter purpose, the theodolite is furnished with a prism and a dark glass (light filter).

The theodolite comprises the following main parts:

1. Lower part consisting of a limb with a silver ring, on which the divisions are engraved, and the levelling base which are connected together by a system of conical centers;
2. Upper part consisting of the telescope to which the vertical graduated limb is attached. They are connected by the horizontal axis of rotation of the telescope;
3. Horizontal axis of rotation of telescope which is supported in the bearings of the standard that connect the upper and lower parts of the theodolite;
4. Four spirit levels (striding, on the guard and standard, on the telescope, on the vertical circle alidade).

SPECIFICATIONS

Reading accuracy on horizontal circle	30"
Reading accuracy on vertical circle	30"
Focusing range	from 2 m to infinity
Telescope magnification	21 × ± 5 %
Telescope field of view	1.8° ± 5 %
Focal length of objective	210 mm
Aperture of objective	25 mm
Stadia constant coefficient	100
Value of spirit level divisions:	
striding	20" ± 2"
on guard and standard	50" ± 10"
on telescope	35" ± 5"
on vertical circle alidade	35" ± 5"
Diameter of horizontal circle	120 mm
Diameter of vertical circle	90 mm

Overall dimensions and weight

Height of theodolite without tripod	300 mm
Weight of theodolite:	
in packing	9.3 kg
without packing	4.3 kg
Weight of extensible tripod	6.1 kg

Attachments

Extensible tripod, with fastening screw	Shoes (3 pcs.)
Case for theodolite with accessories	Canvas case
Striding spirit level	Spare cross-hairs diaphragm in mount
Zenith prism	Studs (2 pcs.)
Dark glass in mount (light filter)	Can with oil
Blind	Wrench for axes
Reflector for illuminating cross-hairs	Napkin
Plumb-bob	Spare screws in cross-hair diaphragm mount (4 pcs.)
Brush	Certificate and operating instructions for the theodolite
Screw driver	

PILOT BALLOON THEODOLITE, MODEL III T


Fig. 16

The Pilot Balloon Theodolite (Fig. 16) is designed for determining the azimuth and height of pilot balloons during meteorologic observations.

One of the methods of investigating phenomena taking place in the upper layers of the atmosphere (as for instance: velocity, direction and constancy of

28

the wind, etc.) is the observation of the flight of pilot balloons (or sounding balloons). Special instruments are used to observe the flight of such balloons and to register data concerning their height, direction and velocity. They are also used for the solution of a number of problems concerning aircraft in general. The simplest of these special instruments is the pilot balloon theodolite.

For operation, these instruments are set up in two or even three locations in the region. The bearings of the locations and the distance between them is known beforehand. The position of an observed point in space will be determined if the angular values of the spherical co-ordinates or some other, for instance, graphical expression of these values is measured from the locations of the instruments. One such paired measurement is sufficient for stationary points. If the point is moving continuously, in each successive moment, the co-ordinates determining its position, will differ from the previous co-ordinates.

Consequently, in this case, it is necessary to make a number of such determinations and the paired observations should be made simultaneously.

The computed results of such data, registered by pilot balloon theodolites, in the case of a moving point will furnish a number of its successive positions in space. This allows the path of the point to be determined as well as its velocity if the time is known. The Pilot Balloon Theodolite ensures normal operation at temperatures from -40°C to $+45^{\circ}\text{C}$ (from -40°F to $+113^{\circ}\text{F}$).

SPECIFICATIONS

Telescope magnification	12x
Telescope field of view	3° 45'
Diameter of exit aperture	3.4 mm
Value of horizontal and vertical circle divisions	1°
Reading accuracy on circles	0.1°
Value of spirit level divisions	0' to 0.6 mm of arc

Overall dimensions and weight

Overall dimensions of theodolite case	200×250×340 mm
Overall dimensions of storage battery case	125×110×160 mm
Overall dimensions of levelling base case	120×130×110 mm
Length of folded tripod	950 mm
Weight of theodolite	3.8 kg
Weight of tripod	5.8 kg
Weight of complete outfit	15.8 kg

Attachments

Illuminating device	Case with storage batteries
Canvas case	Certificate and operating instructions for the theodolite
Storage case	
Tripod	

29

PRECISE LEVEL, MODEL HA-1


Fig. 17

The Precise Level (Fig. 17) is an instrument designed for determining the relative elevation of points in a locality. It is designed for carrying out first order levelling operations.

In comparison with other designs of levels for high-accuracy levelling, the Precise Level has a number of advantageous features:

1. The application in this instrument of an interior focusing telescope shortens the length of the telescope and eliminates the possibility of the penetration of dust and dirt into the optical parts.
2. The spirit level is of the contact type. The images of the ends of the spirit level bubble are seen on the telescope cross-hairs. This is of considerable convenience in operation and speeds up readings on the spirit level.
3. The plane-parallel plate of the reading mechanism is arranged before the telescope objective. This allows various methods of levelling to be used.
4. The vertical axis is of the cylindrical type. This ensures normal operation at a sufficiently large range of temperatures, requires no adjustments and lengthens the term of service of the instrument.

The difference in elevation of point in a locality is determined by the difference in readings on the foresight and backsight rods. During this operation the sighting axis of the telescope must be parallel to the axis of the cylindrical spirit level. Measurement as well as checking of the instrument should be commenced from 15 to 30 minutes after unpacking the level, when the instrument reaches the ambient temperature.

SPECIFICATIONS

Telescope magnification	44×
Vertical field of view of the telescope	60'
Horizontal field of view of the telescope	40'
Aperture of objective	55 mm
Diameter of exit aperture	1.25 mm
Distance to exit aperture	6 mm
Magnification of eyepiece	27×
Resolving power of telescope	3"
Total focal length of objective	411 mm
Sighting range	from 3.6 m to infinity
Stadia constant coefficient	100
Value of cylindrical spirit level divisions	10" to 2 mm of arc
Magnification of bubble image	2.5×
Value of cross-type spirit level divisions	2' to 2 mm of arc
Normal levelling distance	60—65 m
Mean systematic error in measurement for double-rod leveling per km	± 0.3 mm

Overall dimensions and weight

Length of telescope with attachment	400 mm
Height of instrument with levelling screws screwed-in	200 mm
Overall dimensions of case	400×240×170 mm
Length of tripod	1490 mm
Weight of instrument with attachment	5.8 kg
Weight of case	4.0 kg
Weight of tripod	6.5 kg

Attachments

Storage case	Stud
Tripod	Can with oil
Fastening screw	Napkin, 200×200 mm
Wrench for tripod	Sun blind
Screw driver	Certificate and operating instructions for the level

ENGINEER'S LEVEL, MODEL HT


Fig. 18

The Engineer's Level (Fig. 18) is of the type of levels in which the telescope can be lifted out of its supports, turned end for end and replaced. It is designed for determining the difference in elevation between points in a locality for technical levelling operations.

The Engineer's Level is used for engineering, as well as ordinary geodetic work requiring an accuracy of performance equalling a mean square error of ± 3 mm per station or ± 2 mm per rod-reading with average distances between the rods of about 100 m. The sighting axis can be levelled with an accuracy of 5"—10".

A direct reading is taken on the rod to 1 mm.

SPECIFICATIONS

Telescope magnification	31.4×
Telescope field of view	1°
Aperture of objective	34 mm
Diameter of exit aperture	1.1 mm
Distance to exit aperture	8.0 mm

Magnification of eyepiece	25×
Resolving power of telescope	4.5"
Sighting range	from 3 m to infinity
Stadia constant coefficient	100
Value of cylindrical spirit level divisions	17"—25" to 2 mm of arc
Value of circular spirit level divisions	7"—15" to 2 mm of arc

Overall dimensions and weight

Length of telescope	270 mm
Height of level	165 mm
Overall dimensions of case	200×170×300 mm
Length of tripod	1425 mm
Weight of level	2.45 kg
Weight of case with accessories	3.1 kg
Weight of tripod	4 kg
Weight of complete outfit	9.55 kg

Attachments

Fastening screw	Flannel napkin, 200×200 mm
Sun blind	Storage case with accessories
Tripod	Certificate and operating instructions for the level
Can of oil	

DUMPY LEVEL, MODEL HT


Fig. 19

The Dumpy Level (Fig. 19) is designed for determining the difference in elevation between points in a locality during levelling operations.

This Level finds wide applications in various engineering and geodetic work. This instrument differs from other designs in that the spirit level and telescope are rigidly fastened to the upper part of the instrument. The use of a system of prisms for observing the spirit level bubble doubles the accuracy of adjustment in levelling the telescope.

SPECIFICATIONS

Telescope magnification	31×
Telescope field of view	1°
Diameter of exit aperture	1.1 mm
Distance to exit aperture	7 mm
Resolving power of telescope	4.5"
Minimum limit of sighting	3 m
Stadia constant coefficient	100
Magnification of eyepiece	25.5×
Total focal length of objective system	314 mm

Value of cylindrical spirit level divisions	17"—25" to 2 mm of arc
Value of circular spirit level divisions	7"—15" to 2 mm of arc

Overall dimensions and weight

Length of telescope	270 mm
Height of instrument	150 mm
Overall dimensions of storage case	200×170×300 mm
Length of tripod	1500 mm
Weight of instrument	2.25 kg
Weight of case with accessories	3.1 kg
Weight of tripod	4 kg
Weight of complete outfit	9.35 kg

Attachments

Tripod	Storage case with necessities
Fastening screw	Certificate and operating instructions for the level
Sun blind	

MINE LEVEL, MODEL НП-1


Fig. 20

The Mine Level (Fig. 20) is an instrument used for determining the difference in elevation in a locality and is designed for surface and underground third and fourth order levelling operations.

Linear distances can be measured during levelling.

This model is an engineering level of the wye type with the spirit level fastened to the telescope. It is widely used for levelling operations in pits and mines.

SPECIFICATIONS

Visible telescope magnification	31x
Telescope field of view	1.5°
Focal length of the telescope optical system	316 mm
Aperture of objective	36 mm
Stadia constant coefficient	100
Value of spirit level divisions	15" to 2 mm of arc

Overall dimensions and weight

Length of telescope	242 mm
Height of level	155 mm
Weight of instrument without packing	2.7 kg
Weight of instrument in packing	5.2 kg
Weight of extensible tripod	5.2 kg

Attachments

Tripod with fastening screw	Can with oil
Case for level with accessories	Brush
Illuminator	Napkin
Screw driver	Certificate and operating instructions for the level
Studs (2 pcs.)	

OPTICAL ALIDADE, MODEL KB-1


Fig. 21

The Alidade (Fig. 21) complete with a metal plane table, tripod and drawing board are designed for topographical surveying of localities in scales of 1 to 100 and 1 to 10000, as well as for inscribing horizontal lines on photomaps.

The alidade has a glass vertical circle which, besides the usual circular scale, has special curves, engraved on its surface, to provide for reading directly elevations and horizontal distances without reading angles and carrying out tiresome computations. The working surface of the vertical circle is in the focal plane of the objective. Due to this, the observer can either read vertical angles or, using the curves, he can directly read elevations and horizontal distances, or, finally, both methods can be combined.

Direct reading on the curves speeds up surveying operations by 1.5—2 times.

The vertical circle of the alidade is hermetically enclosed and protected by a metal guard.

A box compass is arranged on the ruler. The ruler is also furnished with a parallelogram device which eliminates the necessity of simultaneously sighting the stadia rod and aligning the edge of the ruler with the point of the plane table.

Vertical angles can be read to an accuracy of 1 minute. The error in elevation readings on the curves does not exceed 50 mm at a distance of 100 m. The error in measuring horizontal distances on the curves does not exceed 0.5 %.

SPECIFICATIONS

Magnification of telescope	20×
Telescope field of view	1° 30'
Diameter of exit aperture	1.95 mm
Distance to exit aperture	8.2 mm
Resolving power	4.5"
Aperture of objective	40 mm
Value of telescope spirit level divisions	30"
Value of limb spirit level divisions	30"
Value of scale divisions	2 mm
Value of limb divisions	10'
Value of compass scale divisions	30'

Overall dimensions and weight

Length of alidade ruler	580 mm
Height of alidade	230 mm
Length of tripod	1570 mm
Size of drawing-board	600×600 mm
Weight of alidade and plane table in case	10.25 kg
Weight of drawing-board in canvas case	7 kg
Weight of tripod	5.4 kg
Weight of complete outfit	23 kg

Attachments

Plane table in case	Sun blind
Tripod	Accessories and spare parts
Drawing-board in canvas case	Certificate and operating instructions for the alidade
Box compass	

ALIDADE, MODEL KB


Fig. 22

The Alidade (Fig. 22) complete with a plane table, tripod and drawing-board is designed for topographical surveys of localities and for inscribing horizontal lines on photomaps.

The alidade has an open-type vertical circle and two reading glasses for opposite readings on the circle. Vertical angles can be read to an accuracy of one minute.

A box compass is furnished to orient the plane table in reference to the magnetic meridian.

SPECIFICATIONS

Magnification of telescope	25×
Aperture of objective	34 mm
Focal length of objective	380 mm
Focal length of eyepiece	15 mm

40

Telescope field of view	1°
Minimum sighting distance	5 m
Stadia constant coefficient	100
Value of smallest vertical circle divisions	30'
Value of vertical circle spirit level divisions	30"—50"
Value of ruler spirit level divisions	50"—80"
Smallest division of compass scale	30'

Overall dimensions and weight

Length of alidade rule	530 mm
Height of alidade (to horizontal axis)	205 mm
Length of tripod	1250 mm
Size of drawing-board	600×600 mm
Weight of alidade	5.0 kg
Weight of plane table	5.5 kg
Weight of drawing-board	4.5 kg
Weight of tripod	4.6 kg
Weight of complete outfit	38.2 kg

Attachments

Plane table in case	Box compass
Tripod	Sun blind
Drawing-board in canvas case	Certificate and operating instructions for the alidade

41

SURVEYING COMPASS, MODEL BC


Fig. 23

The Surveying Compass (Fig. 23) is an improved geodetic compass and is used as an independent field instrument for orientation in reference to points of the compass, as well as for measuring horizontal angles and azimuths.

This compass is used for finding the bearings of base points in reconnaissance operations; especially in heavily wooded localities and near rivers, as well as for sketching details located within the surveyed areas.

SPECIFICATIONS

Value of limb divisions	1°
Accuracy of reading angles and azimuths	5'
Distance between vertical sights	84 mm
Value of compass circle divisions	1°
Weight of instrument	0.8 kg

Attachments

- Plumb line
- Screw driver
- Storage case
- Certificate and operating instructions for the compass

GONIOMETER, MODEL GP


Fig. 24

The Goniometer (Fig. 24) is the simplest of angle measuring instrument for use in the field and it fully replaces a surveying compass and an optical square. It is used to measure horizontal angles in reference to compass points, azimuths and between sighting directions.

The angles between directions can be read, simultaneously, on the limb and, using the compass, by their bearings.

The goniometer is used for surveying and pegging out a locality.

SPECIFICATIONS

Value of limb divisions	1°
Reading accuracy of angles and azimuths	5'

Distance between sights 80 mm
 Value of compass circle divisions 1°
 Weight of instrument 1.2 kg

Attachments

Plumb line
 Screw driver
 Storage case
 Certificate and operating instructions for the goniometer

GEODETIC OPTICAL SQUARE, MODEL ӘГ-2


Fig. 25

The Optical Square (Fig. 25) is a geodetic instrument used for pegging out perpendicular directions to the line of sight. It is of the mirror type.

SPECIFICATIONS

Angle between mirrors 45° ± 2'
 Length of instrument 140 mm
 Width of instrument 50 mm
 Weight of optical square 0.11 kg
 Weight of canvas case 0.035 kg

Attachments

Canvas case
 Flannel napkin, 100 × 100 mm
 Certificate and operating instructions for the optical square

PRECISE CHECKERED LEVELLING RODS, MODEL 51-T-86

Precise Checkered Levelling Rods are used in first and second order levelling operations.

They are made of fine-grained pine-wood (aviation type).

The rods are graduated on both faces; one face has centimeter checker squares and half-centimeter divisions in black paint; while the opposite face has 11-millimeter squares in red paint and graduations analagical with the first face.

A circular spirit level is attached to one side of the rod while a plumb line may be fastened to the other side for checking the spirit level.

A steel plate is fastened to the lower end of the rod and is used for setting the rod up on the convex spherical projection of the levelling shoe.

The values of the divisions are inscribed as follows:

- on black face — at decimeter intervals (from 0 to 30),
- on red face — at 11-centimeter intervals (from 01 to 28).

Three control markings are provided on each face. They are arranged at a distance of one meter from each other.

Precise levelling rods are furnished in sets of two pieces.

SPECIFICATIONS

Overall dimensions and weight

Overall dimensions of rod	47×84×3000 mm
Overall dimensions of rod extension	28×72×1200 mm
Weight of rod.....	3.4 kg
Weight of rod extension.....	1.2 kg

Attachments

Circular spirit levels (2 on rods and 2 spares in mounts)	Screw drivers for spirit level adjusting screws (2 pcs.)
Plumb bobs with lines (2 pcs.)	Studs for attaching rod extension (4 pcs.)
Handles (rod supports) (4 pcs.)	Device for holding plumb lines and checking circular levels (on the rods)
Canvas cases for rods (2 pcs.)	Storage case (with lock) for precise rods and accessories for extension rod
Wooden extension rod	
Canvas case for extension rod	

C O N T E N T S

	Page
Universal Astronomical Theodolite, Model AY 2/10	3
Triangulation Theodolite, Model TT 2/0	8
Optical Theodolite, Model OT-02	13
Optical Theodolite, Model OT-10	17
Optical Theodolite, Model TB-1	21
Theodolite Tacheometer, Model TT-50	24
Mine Theodolite, Model TT-1	26
Pilot Balloon Theodolite, Model HIT	28
Precise Level, Model HA-1	30
Engineer's Level, Model HT	32
Dumpy Level, Model HP	34
Mine Level, Model HII-1	36
Optical Alidade, Model KB-1	38
Alidade, Model KB	40
Surveying Compass, Model BC	42
Goniometer, Model GP	43
Geodetic Optical Square, Model 3Г-2	45
Precise Checkered Levelling Rods, Model 51-T-86	46

VSESOJUZNOJE EXPORTNO-IMPORTNOJE OBJEDINENIJE

“STANKOIMPORT”

exports and imports:


Machine Tools
Woodworking Machinery
Metal Working Machinery (Presses, Hammers, Shears, Cold Forming
Machines, Punching Machines)
Rolling Mills (imports)
Measuring Instruments and Apparatus (for metal industry)
Testing Machines and Instruments (for metals)
Optical Instruments and Equipment
Portable Electric and Pneumatic Tools (for metal and woodworking)
Metal and Wood Cutting Tools
Mechanic's Tools and Chucks
Sintered Carbide and Hard-Alloy Products
Abrasive Products
Ball and Roller Bearings
Microscopes of all types
Motion-Picture Equipment and Accessories
Geodetic Instruments and Equipment
Photographic Cameras
Binoculars
Magnifiers
Lenses
Crude Optical Glass Blocks and Blanks

All inquiries and correspondence to be forwarded to:

Vsesojuznoje Exportno-Importnoje Objedinenije
“STANKOIMPORT”
32/34, Smolenskaja-Sennaja pl., Moscow, U. S. S. R.

For cables: STANKOIMPORT MOSCOW

Design and specifications of the instruments illustrated herein are subject
to change without notice.

CABLE ADDRESS:
STANKOIMPORT
MOSCOW

UNCLASSIFIED

MOTION
PICTURE
EQUIPMENT

UNCLASSIFIED

USSR

STANKOIMPORT

MOSCOW

STAT

MOTION
PICTURE
EQUIPMENT

VSESOJUZNOJE EXPORTNO-IMPORTNOJE OBJEDINENIJE

S T A N K O I M P O R T

U S S R

M O S C O W

otion picture equipment and accessories manufactured in the U.S.S.R. embody the latest achievements in this field, and are perfect in design and workmanship.

Superior materials, modern design and skilful workmanship ensure outstanding performance and durability of the motion picture equipment and accessories.

High operating merits are combined with a finely-styled outer appearance, excellent finish, light weight and convenience in handling.

CONTENTS

	Page
Motion Picture Cameras	5
"Rodina" (KCX) 35-mm Motion Picture Camera for Newsreel and Expeditionary Filming	7
"Moskva" (KC-32) 35-mm Motion Picture Camera for Synchronous Filming	10
KC-50E 35-mm Motion Picture Camera for Newsreel Filming	14
IIIKC-2 Camera Tripod	15
IIIC-3 Camera Tripod	17
Sound Film Recording Equipment	19
K3IIV Portable Sound Film Recorder	21
K3YC Stationary Sound Film Recorder	29
KII3-1 Film Re-recording Equipment	36
Laboratory Equipment	43
40 II-1 Automatic Developing Machine	45
KII3-2 Film Sample Printer	48
YKA Sound Film Printers	51
MM-11 Film Renovation Machine	55
PVII-1 Cutting Machine	57
KCH-3 Film Stapler	59
35-3MA-3 Sound Editor	60
35-MMIC-3 Subtitle Making Machine	62
35-CIA-2 Semi-Automatic Film Splicing Machine	65
35-CO-1 Synchronizer	67
16-IIICII 16-mm Film Splicer	68
35-IIICII-3 35-mm Film Splicer	69
ΦC-2 Film Inspection Bench	70
35-MOI-3 Horizontal Film Rewinder	71
Motion Picture Projectors and Equipment	73
"Ukraina" 16-mm Portable Sound-on-Film Projector	75
KIIIC-M 35-mm Portable Sound-on-Film Projector	79
KIIT-1 35-mm Stationary Sound-on-Film Projector	83
K3BT-3 Sound-reproducing Amplifying Device	88

	Page
КУСУ-52 Sound-reproducing Amplifying Device	91
РУ-65 Switching Rack	94
КАТ-14 Auto-Transformer	96
ТРД-50 Reactive Arc Transformer	98
ЭПД Portable Diffuse Reflection Screen	101
ТС-5 Light Dimmer	102
АЗС-9-10 Automatic Fire Shutters	104
ЛЗ-2 Automatic Screen Curtain Winch	106
ФС-5 Film Storage Container	106
"КИЕВ-1" Mobile Power Plant	107
КЭС-5 Mobile Power Plant	109
КЭС-4 Mobile Power Plant	111
Measuring and Inspection Instruments	113
16-MM-3 Film Length Counter	115
35-MM-3 Film Length Counter	116
СО-301-1 Film Inspection Magnifier	117
ЛШП-16 Perforation Pitch Inspection Rule	118
ЛШП-35 Perforation Pitch Inspection Rule	119
ПКП-2 Perforation Inspection Projector	120
РИФ-3 Sound Track Measuring Instrument	121
ИН-3 Inspection Kit	122
ДФЭ-2 Photoelectric Densitometer	124
ЭПП-4 Universal Exposure Meter	126
Lighting Apparatus	129
ПР-60 Apparatus for Diffused Light Illumination	131
РД-5 Apparatus for Diffused Light Illumination	133
КПД-15 Arc Lamp Spotlight	135
КПД-25 Arc Lamp Spotlight	137
КПД-М Arc Lamp Spotlight	139
КПЛ-25 Incandescent Lamp Spotlight	141
КПЛ-35 Incandescent Lamp Spotlight	143
КПЛ-50 Incandescent Lamp Spotlight	145

MOTION PICTURE CAMERAS

**"RODINA" (KCX) 35-mm MOTION PICTURE CAMERA
FOR NEWSREEL AND EXPEDITIONARY FILMING**

The "Rodina" Motion Picture Camera (Fig. 1) is designed to perform newsreel and expeditionary filming with black-and-white as well as colour 35-mm film.


Fig. 1. Motion Picture Camera "Rodina"

New design of the intermittent film mechanism ensures high-precision registration of separate picture frames within the film channel, and allows use of the camera for special and trick filming.

Special design of the friction clutch permits use of magazines of 120-metres and 300-metres film capacity (Fig. 2). Quick change lens mounts provide ease and facility in the interchange of lenses.

Picture composition and visual control during the filming process are achieved by means of a magnifying view-finder.

Control of filming speed is provided by a tachometer with dial divisions in picture-frames per second.

SPECIFICATIONS

Film size 35 mm
 Path of film travel within camera . . . in three planes
 Film gate dimensions 16 X 22 mm
 Shutter aperture angle 0° to 160°
 Magazines single compartment type of 120 and 300-meter capacity
 Film length and picture-frame counter drum type with zero setting


Fig. 2. Motion Picture Camera "Rodina" with 300-metre capacity magazines


Fig. 3. Motion Picture Camera "Rodina", right-side view

Camera-to-tripod attachment by means of 3/8" screw
 Interlocking devices to effect stoppage in case of film breakage, termination of film roll, or slackening of film tension
 View-finder parallaxless magnifier of 5.5-time magnifying power
 Set of coated photographic lenses in bayonet mounts 28, 35, 50, 75 and 100-mm focal length

Minimum focal length of photographic lens 28 mm
 Focusing of photographic lens by magnified image on view-finder screen; by focusing scales
 Tachometer readings 8, 16, 24, 32 and 48 picture frames per second
 Normal operation temperature conditions from + 40° to - 25° C
 Power drive 12 V, 60 W D. C. electric motor
 Camera run reversible
 Exposure speed 8 to 48 picture frames per sec
 Noise level 50 + 2 db
 Overall dimensions of camera in operating position with electric motor (length X height X width) . 650 X 280 X 260 mm
 Weight of camera (less tripod and film) 14 kg with hand drive 16.2 kg with motor drive
 Weight of camera set in carrying cases (less storage battery) 35 kg
 Weight of storage battery 11.5 kg

"MOSKVA" (KC-32) 35-mm MOTION PICTURE CAMERA FOR SYNCHRONOUS FILMING

The "Moskva" Camera (Fig. 4) is designed to perform synchronous filming of images and simultaneous sound-track recording-


Fig. 4. Motion Picture Camera "Moskva"

ing-by means of a recorder on 35-mm film and is adapted for studio and outdoor work.

Specially silenced to eliminate mechanical noises the camera permits simultaneous sound-recording, provided the microphone is stationed not less than 1 metre from the camera front.

High-precision performance of the intermittent film mechanism allows use of the camera for special (trick) filming.

Focusing of the photographic lens is achieved by ground glass, film, or focusing scales; simultaneous compensation of parallax and focusing of the view-finder lens are automatic.

A left-side view of the camera is given in Fig. 5, and a diagram of film threading in Fig. 6.

The camera is provided with an interlocking device to shut off the motor when normal operation of the camera is interfered with.

On-and-off switching of the camera is achieved by means of a switch mounted on the power-supply line.


Fig. 5. Motion Picture Camera "Moskva", left-side view


Fig. 6. Film Threading Diagram for "Moskva" Camera

Camera power supply is fed through a 220 V autotransformer from 220 or 380-volt 50-cycle A. C. mains. Voltage variations in the mains are compensated by a regulator switch in the autotransformer.

Lighter weight and smaller overall dimensions add greatly to the camera's operating merits over other known synchronous motion picture cameras.

Handy location of all operation and control elements on the back and right side-walls of the camera, and the design of the easily removable parts (intermittent film mechanism, sprocket assembly, etc.) as well as the interchangeability of the detachable parts in all cameras of this type ensure operation efficiency and convenience in inspection and cleaning.

SPECIFICATIONS

- Film size 35 mm
- Path of film within camera in one plane
- Exposure speed 24 picture frames per sec
- Noise level not over 29 db
- Accuracy of picture frame registration in the film channel 0.008 mm
- Film gate dimensions 16 X 22 mm
- Shutter aperture angle 0° to 170°
- View-finder detachable, with automatic parallax compensation and lens focusing
- Magazines single compartment type, in pairs of 300-metre capacity each
- Set of coated lenses 28, 35, 50, 75, 100 mm focal length
- Minimum focal length of photographic lens 24 mm
- Power drive three-phase 220 V, 150 W, 1,500 r.p.m. synchronous motor
- Speed of reduction gear shaft 1,440 r.p.m.
- Camera run reversible
- Interlocking devices to effect stoppage in case of film breakage, slackening of film tension or idling of intermittent film mechanism
- Focusing of photographic lens by ground glass, by film, by focusing scales
- Camera to tripod attachment by means of 1/8" screw
- Overall dimensions of camera in working position 1,120X585X470 mm
- Distance between optical axis and base of camera housing 169.7 mm
- Weight of camera (less tripod and film) 62 kg
- Weight of camera set in carrying cases 130 kg

The motion picture camera is supplied with the following accessories:

- Focusing magnifier; light-protective device with filter-holder;
- 6 magazines of 300-metre capacity; belt tensioning and switching mechanism; carrying cases for magazines and motor; switch.
- Photographic lenses (PO coated type) with 1:2 relative aperture and focal lengths of 28, 35, 50, 75 and 100 mm (5 lenses).
- IM-35 three-phase, 220 V, 150 W electric motor.

- 7C-II View-finder.
- KAT-24-I auto-transformer in jacket.
- Set of tools.
- Description of camera and instructions for maintenance.
- Carrying cases (Fig. 7).

Fig. 7. Motion Picture Camera "Moskva", packed in carrying cases

The motion picture camera "Moskva" is additionally supplied with a ИИСК-2 tripod on special order.

**KC-50 B 35-mm MOTION PICTURE CAMERA
FOR NEWSREEL FILMING**

The KC-50 B 35-mm Motion Picture Camera (Fig. 8) is designed for filming of newsreel and documentary films.

The camera has 3 interchangeable lenses mounted on a revolving turret, and a rotating view-finder; the focal lengths of the view-finder lenses are proportionate to their corresponding


Fig. 8. KC-50 B Motion Picture Camera

photographic lenses. Focusing of the lenses is achieved by means of a distance collar on each of the lens mounts.

A spring-drive mechanism enables to conduct continuous filming up to a 16.5-meter film run.

Filming can be performed with the camera being hand-held or with the camera set on a tripod.

The camera is equipped with a hand drive; one revolution of the hand crank corresponds to the exposure of eight picture frames.

The camera has a special carrying case for convenient transportation.

Accessories and necessary tools are provided with the camera.

SPECIFICATIONS

Number of photographic lenses	3
Lenses relative aperture	1:2
Focal length of lenses	35, 50 and 75 mm
Reel capacity	30-34 m
Operating speed	8, 12, 16, 24, 32 exp. per sec
Film length counter division	0.5 m
Shutter	of the constant 160° angular aperture allows the following exposures:
	at 8 frames per sec
	" 12 " " "
	" 16 " " "
	" 24 " " "
	" 32 " " "
	" " " "
Distance range	for each lens is from 1 m to infinity.
Overall dimensions of camera:	285×250×450 mm
Weight of camera	12 kg

IIIKC-2 CAMERA TRIPOD

The IIIKC-2 Tripod (Fig. 9) is designed for motion picture cameras of various types for studio and outdoor filming.


Fig. 9. IIIKC-2 Tripod

The tripod is adapted for use on different soils and rough surfaces, in premises with smooth and hard floors, as well as on special vehicles and trucks.

The tripod comes in a set, and the use of its components in different combinations gives the camera a height ranging from 300 mm to 1,800 mm.


Fig. 10. Tripod Head with Auxiliary Platform

The set includes:

- Tripod head (Fig. 10)
- Auxiliary platform (Fig. 10)
- Normal length tripod legs (Fig. 11)
- Short length tripod legs (Fig. 11)
- Low mount supports (Fig. 11)


Fig. 11. Tripod Legs

The tripod head allows vertical and horizontal panning. Special levers are provided to regulate the friction devices of the tripod head. Vertical panning is facilitated by spring compensators set in the tripod head.

Control of the tripod motions is achieved by means of a single handle (bar). For convenience of operation the handle may be attached to the right or left side of the tripod head, telescoped, and swung into any desired position.

Two lock levers, one for vertical tilt, the other for horizontal, are provided to secure the tripod head in any operating position.

The motion picture camera is secured to the tripod head with a 3/8" screw or by means of a special dovetailed auxiliary platform (Fig. 10).

The tips of the oaken telescopic tripod legs are fitted with double-spurred shoes. When telescoped, the legs are fitted in position by a double-action clamp.

For convenience in transportation the tripod set is provided with jackets.

SPECIFICATIONS

Height of tripod may be varied:

main tripod	1,000—1,800 mm
small tripod	600—1,050 mm
low mount support	300—360 mm

Motions of the tripod head:

horizontal motion	endless over 360°
downward tilting	to any angle up to 45°
upward tilting	to any angle up to 35°

Tripod is designed to support cameras weighing up to 100 kg

Weight of tripod:

head	11.5 kg
main tripod	10.5 kg
small tripod	9.3 kg
support	5.2 kg
auxiliary platform	2.3 kg
overall weight of tripod	38.8 kg

IIC-3 CAMERA TRIPOD

The IIC-3 Tripod (Fig. 12) is designed for both silent cameras and cameras provided with a sound recorder in filming indoors and outdoors.

The camera is secured to the tripod by means of an auxiliary platform and a special clamping device.

The tripod set includes:

- An inertia type head with an worm and gear device.
- Normal length legs.
- Short length legs.
- Leg-supporting shoes to ensure stability of the tripod on smooth and hard surfaces.
- Two tripod jackets.
- Two sets of tripod legs allow the wide height-range from 0.5 to 1.8 metres.
- The tripod head is equipped with a device ensuring inertia, worm and gear and free motions.

Besides free and variable inertia motion, the inertia-type head possesses two worm and gear drives for each.


Fig. 12. IHC-3 Tripod

A special brake regulates the movement and stoppage of the camera during inertia panoraming.
The tripod head design permits horizontal circular panoraming (360°) and vertical panoraming with an up-and-down range


Fig. 13. IHC-3 Tripod Head

of 40° each. The design ensures easy and quick adjustment of the tripod head to the legs.

The tripod is adapted to support cameras weighing up to 30 kg. Finely styled, light-inweight, packed in compact jackets, the tripod is convenient for transportation and storage.

Weight of the tripod — 18 kg.

**SOUND FILM RECORDING
EQUIPMENT**

K3IIV PORTABLE SOUND FILM RECORDER

The K3IIV Portable Sound Film Recorder is designed to perform synchronous sound recording on 35-mm film by the photographic method in expeditionary conditions.

Several models of portable sound recorders are available at present for recording on normal or double-area sound tracks, and enable connection of the recorder to different supply lines under varying local conditions.

The table below includes brief characteristics of the models available and their sets.

The K3IIV Sound Recorder contains:

33II Recording Device (Fig. 14) which is designed for recording on 35-mm film by means of the variable area method. Recording is of normal or push-pull noiseless sound track type of standard or double area.

High uniformity of film travelling speed is ensured by a rotary speed stabilizer.


Fig. 14. 33II Recording Device

The light modulation system of the recording device allows for recording with white and ultra-violet light.
High-efficiency recording is provided by the frequency range of the device.

The recording device is mounted and secured in a solid carrying case convenient for transportation.

Removal of the easily detachable case cover and installation of the motor and magazines are sufficient to bring the device into instant use.

Overall dimensions of recording device (height X length X width) 600 X 720 X 300 mm
Weight of recording device 54 kg

1 V 50-A Pre-amplifier (Fig. 15), which mixes and pre-amplifies incoming signals from two microphones.


Fig. 15. 1 V 50-A Pre-amplifier

The amplifying range of the device permits use of microphones of the lowest response.

The provided correction range of the frequency characteristics in the pre-amplifier and the main amplifier covers recording requirements in the studio and outdoors.

The electromechanical properties of the volume indicator mounted in the pre-amplifier enable to maintain complete visual control of the volume of signals recorded.

The electrical characteristics of the transformers in the pre-amplifier ensure a remarkably low noise level of the amplifying channel, and do not, therefore, limit the dynamic range of the recording.

The pre-amplifier is mounted in a metal portable case and is readied for use by a simple connection of its cables. Removal of the easily detachable front cover gives access to the control board.

Overall dimensions of pre-amplifier (height X length X width) 305 X 350 X 293 mm
Weight of pre-amplifier 15 kg

12 V-3 Main Amplifier (Fig. 16), which achieves the subsequent amplifying of signals coming in from the pre-amplifier and makes provisions for the possibility of compression of the output volume.

The main amplifier is mounted in a metal cabinet containing frames for recorder and intermediate amplifiers. Clamps on the cabinet side walls are provided to secure the cabinet in a truck.


Fig. 16. 12 V-3 Main Amplifier

The cabinet front cover is easily removable to give free access to the tubes. The frame is mounted on loops, and special hinges afford free access to the amplifier parts.

Overall dimensions of the amplifier (height X length X width) 712 X 535 X 365 mm
Weight of amplifier 41 kg

Power-supply Device, which supplies motion picture camera and sound-recorder electric motors with three-phase current and feeds D. C., high and low voltage circuits of the pre-amplifier, the main amplifier and the exposure lamp of the sound recorder. Four types of power supply devices — KIICV-1, KIICV-2, KIICV-3 and KIICV-4 are available. Regardless of the type of Power-supply device, stabilized supply is ensured for D. C. circuits of high- and low voltage.

a) The KIICV-2 Power-supply device (Fig. 17) is operated on high-capacity storage batteries and is used where no A. C. mains are available.

50-cycle frequency of the three-phase current is constantly maintained irrespective of battery discharge in the bounds of normal operating conditions. Battery capacity provides for 12 hour operation without reloading.

All parts of the equipment for converting and distributing the electric power supplied from the batteries are mounted in a durable aluminum frame.

To provide transportation convenience the frame is secured in a solid case of special design. Inspection and adjustment of the rotary converters is facilitated by the sliding guides on which they are mounted.


Fig. 17. KIICV-2 (KIICV-4) Power-supply Device

The side walls of the case may be swung aside to allow free access to the switch panels.

Overall dimensions of the power-supply device (height X length X width) 430 X 605 X 330 mm
Weight of device 55 kg

b) The KIICV-4 Power-supply device is intended for the same purposes as the KIICV-2 device with the difference that the alternating current it converts is of 60-cycle frequency.


Fig. 18. KIICV-1 (KIICV-3) Power-supply Device

c) The KIICV-1 Power-supply device (Fig. 18) is intended for use when the K3IIY sound recorder is supplied from 220 or 380 V three-phase 50-cycle mains.

To make it portable the KIICV-1 power-supply device is divided into two functionally related parts. The first contains the general stabilizer and exposure lamp rectifier, the second — the power-supply rectifier of the amplifying channel.

Both parts are of similar construction, styled of duraluminum frames, mounted in solid cases provided with hinged covers.

Overall dimensions of each of the parts (height X length X width) . . . 430 X 615 X 300 mm
Weight of device 94 kg

d) The KIICV-3 Power-supply device is intended for the same purposes as device KIICV-1, with the difference that it operates on 60-cycle frequency mains.

1 Y-101 Public Address Amplifier (Fig. 19), which intensifies orders relayed from the camera-crew during mass outdoor filming.


Fig. 19. 1 Y-101 Public Address Amplifier

The public address amplifier is styled as a light metal chassis covered with a casing. The hinged rear cover of the casing provides free access to the tubes.

Overall dimensions of the amplifier (height X length X width) 245 X 312 X 190 mm
Weight of amplifier 6 kg

25 A-1 Loudspeaker (Fig. 20), which reproduces instructions from the camera-crew during mass outdoor filming.

The loudspeaker consists of a head provided with a permanent magnet and mounted in a solid wooden case with a grip for transportation. A metal grid protects the head against mechanical damages and a hood of light fabric keeps out the dust.

Overall dimensions of loudspeaker (height X length X width) 430 X 432 X 201 mm
Weight of loudspeaker 12 kg

26

TABLE
of KЗПУ Sound Recorder Sets

Model		Main units of set							
Symbol	Characteristics of the model	Recording device	Preamplifier	Main amplifier	Power-supply device	Public address amplifier	Loud-speaker	Microphone	Communication phones
KЗПУ-2	Basic type. Designed for supply from storage batteries.	ЗЗП-1/Б	1 У 50-А	12 У-3	КПСУ-2	1 У-101	25 А-1	8 А-5	УНАФ-42 М
KЗПУ-3	Designed for supply from A. C. mains	ЗЗП-1/Б	1 У 50-А	12 У-3	КПСУ-1	1 У-101	25 А-1	8 А-5	УНАФ-42 М
KЗПУ-4	Recording device adapted to record sound tracks of double area. Power supply from storage batteries.	ЗЗП-2	1 У 50-А	12 У-3	КПСУ-2	1 У-101	25 А-1	8 А-5	УНАФ-42 М
KЗПУ-5	Recording device for double area sound tracks. Designed for supply from A. C. mains.	ЗЗП-2	1 У 50-А	12 У-3	КПСУ-1	1 У-101	25 А-1	8 А-5	УНАФ-42 М

KЗПУ-6	Special set of optics and spare parts allows change-over to double area sound tracks in studio. Power supply from storage batteries.	ЗЗП-1/Б	1 У 50-А	12 У-3	КПСУ-2	1 У-101	25 А-1	8 А-5	УНАФ-42 М
KЗПУ-7	Special set of optics and spare parts allows change-over to double area sound tracks in studio. Power supply from A. C. mains.	ЗЗП-1/Б	1 У 50-А	12 У-3	КПСУ-1	1 У-101	25 А-1	8 А-5	УНАФ-42 М
KЗПУ-8	Similar to KЗПУ-6 with power supply frequency being 60 cycles	ЗЗП-1/Б	1 У 50-А	12 У-3	КПСУ-4	1 У-101	25 А-1	8 А-5	УНАФ-42 М
KЗПУ-9	Similar to KЗПУ-7 with power supply from A. C. mains being 60 cycles	ЗЗП-1/Б	1 У 50-А	12 У-3	КПСУ-3	1 У-101	25 А-1	8 А-5	УНАФ-42 М

27

8 A-5 Dynamic Microphone, with an 11-A-8 type tripod, is used with the public address amplifier when instructions are being relayed. The solid design of the 8 A-5 dynamic microphone ensures reliable performance in expeditionary conditions.

Two УНАФ-42М Communication Telephones, which serve for communications between the sound-recording operator and the sound-recording room.

Any requirements liable to arise with the use of the K3ИY Sound Recorder in expeditions have been taken into consideration in the design of the units.


Fig. 20. 25A-1 Loudspeaker

Switching and inter-connecting of separate parts of the recorder are achieved by means of connector cables with connection plugs. Cable extension is provided to allow the preamplifier, the motion picture camera and parts of the public address channel to be moved to distances up to 100 metres from the remaining equipment (usually installed on a truck).

For railroad transportation and storage of spare parts and service tools the amplifying section of the recorder is provided with three sturdy cases.

Weight of set, complete with carrying cases — approx. 1,050 kg.

K3YC STATIONARY SOUND FILM RECORDER

The K3YC Stationary Sound Film Recorder performs synchronous sound recording on 35-mm film by the photographic method in studios.

In order to meet requirements of recording sound tracks of normal, as well as double area, and connection of the device to different supply lines under varying local studio conditions, several models of stationary sound film recorders are available.

The table below gives brief characteristics of available sound recorder models and their sets.

The K3YC Sound Recorder includes the following units:


Fig. 21. 3K Recording Device

3K Recording Device (Fig. 21), which records sound tracks on 35-mm film by the variable area method. Recording is of the normal or push-pull noiseless method with standard or double area sound track.

The speed stabilizer (magnetic drive) used with the device ensures unfailing uniformity of film speed past the exposure slit. The light modulation system of the recording device allows for recording with ultra-violet and white light.

Reliability of construction and the high-precision workmanship of its parts provide lasting service, the device being readied for operation by simple insertion of magazines and threading of film. Correct setting of exposure lamp and magnetic drive exciter coil is achieved by means of rheostats, which ensure smooth adjustment controlled by adequate instruments.

Overall dimensions of recording device (height X length X width) 640 X 560 X 550 mm
Weight of recording device 102.3 kg

12 Y-5, type 2 Preamplifier (Fig. 22), which mixes and pre-amplifies incoming signals from the microphones, and allows simultaneous connection of four microphones. The amplifying range permits use of microphones of any type, even under most unfavourable conditions.


Fig. 22. 12Y-5, type 2 Preamplifier

The adjustment range of the frequency characteristics in both preamplifier and main amplifier makes provisions for the requirements of recording in studios.

The electro-mechanical properties of the level indicator mounted in the preamplifier enable full visual control of the volume of signals recorded.

The extremely low level of interferences throughout the amplifying channel is achieved by the electric characteristics of the amplifier input transformers and the type of shields used with them, thus limiting the dynamic range of the recording only by the quality of the film used.

The preamplifier is styled as a panel supported by a movable table equipped with swivel rollers. The device is ready for operation after simple connection to the switchboard by means of connecting cables.

Overall dimensions of preamplifier (height X length X width) 1,010 X 890 X 600 mm
Weight of preamplifier 92.2 kg

12 Y-4, type 2 Main Amplifier (Fig. 23), which is intended for subsequent amplifying of signals from the preamplifier, and for compressing and limiting the output volume. The main amplifier includes a monitor amplifier channel and a loudspeaker, which enable perfect monitoring of the recording.


Fig. 23. 12Y-4, type 2 Main Amplifier

Control of the main amplifier circuit current is provided by a highly-sensitive instrument with push-button change-over control.

The main amplifier is styled as a metal cabinet accomodating the chassis of both recorder and intermediate amplifiers; the front panel of the cabinet is easily removable to give access to the tubes. Free access to the amplifier parts is provided by the hinged chassis and special rods, allowing to swing the chassis off.

Overall dimensions of the amplifier (height X length X width) 912 X 525 X 365 mm
Weight of amplifier 50 kg


Fig. 24. 12M-3 Power-supply Device

12M-3 Power-supply Device (Fig. 24), which supplies all high and low voltage D. C. circuits with stabilized rectified current. Power is fed to the device from standard 127/220 V, 50-cycle A.C. mains.

The power-supply device ensures the unflinching stabilization of rectified current in a wide range of line voltage variations.

The rectifiers and stabilizers of the power-supply device are mounted in a single cabinet.

Overall dimensions of the power-supply device (height \times length \times width) . . . 1,240 \times 640 \times 400 mm
Weight of device 151 kg

Overall dimensions of panel (height \times length \times width) . . . 490 \times 550 \times 550 mm
Weight of panel 22.6 kg

Switching Equipment, which consists of the 6K-50 main switching unit (Fig. 26), two 6K-18 studio panels (Fig. 27 and 28), a set of connecting cables and designed for centralized switching of preamplifier and recorder.


Fig. 25. 7K-7 Switchboard Panel

Overall dimensions of the 6K-50 main switching unit (height \times length \times width) 643 \times 563 \times 123 mm
Weight of main switching unit 20 kg
Dimensions of the 6K-18 panel 90 \times 482 \times 103 mm each
Weight of each panel 2.3 kg

Communication System, which is designed for telephone communications between studio and central recording-room.

25 A-7 Monitor Loudspeaker (Fig. 29), which serves for monitoring recording being made in central recording-room.

Overall dimensions of loudspeaker (height \times length \times width) 600 \times 700 \times 300 mm
Weight of loudspeaker 9.3 kg


Fig. 26. Switching Unit

25 A-6 Monitor Loudspeaker (Fig. 30), which enables high-standard audio checking in a special demonstration room thanks to the electroacoustical properties of the speaker.

Fig. 27. 6K-18 Switching Panel


Fig. 28. 6K-18 Switching Panel

Overall dimensions of loudspeaker (height \times length \times width) 860 \times 720 \times 430 mm
Weight of loudspeaker 41 kg

TABLE
of K3VC Sound Film Recorder Sets

Symbol	Characteristics of the model	Main Units of the Set						
		Recording device	Preamplifier	Main amplifier	Power supply device	Switching equipment	Communication system	Monitor loudspeaker
K3VC-5	Basic model	3K-1	12 Y-5, model 2	12 Y-4, model 2	12 M-3, model 2 with 7 K-7, model 2 switchboard	6 K-50 with two 6 K-18 panels	Available	25 A-7 and 25 A-6
K3VC-6	Recording device with optical system for double area sound track recording	3K-5	12 Y-5, model 2	12 Y-4, model 2	12 M-3, model 2 with 7 K-7, model 2 switchboard	6 K-50 with two 6 K-18 panels	Available	25 A-7 and 25 A-6
K3VC-7	Recording device has set of optics and spare parts, which allow change-over to double area sound tracks for studio work	3K-1	12 Y-5, model 2	12 Y-4, model 2	12 M-3, model 2 with 7 K-7, model 2 switchboard	6 K-50 with two 6 K-18 panels	Available	25 A-7 and 25 A-6
K3VC-8	Similar to K3VC-5, but supplied with power from 60-cycle frequency mains	3K-6	12 Y-5, model 2	12 Y-4, model 2	20 B-6 with 50 K-2 switchboard	6 K-50 with two 6 K-18 panels	Available	25 A-7 and 25 A-6
K3VC-9	Similar to K3VC-5, but supplied with power from 60-cycle frequency mains	3K-6	12 Y-5, model 2	12 Y-4, model 2	20 B-6 with 50 K-2 switchboard	6 K-50 with two 6 K-18 panels	Available	25 A-7 and 25 A-6

The main wiring of the recorder is of the stationary type. Detachable wiring consists of flexible cables provided with special connectors to ensure dependable connections.


Fig. 29. 25A-7 Monitor Loudspeaker


Fig. 30. 25A-6 Monitor Loudspeaker

Weight of the complete sound recorder set packed in carrying cases is about 1,100 kg.

KII3-1 FILM RE-RECORDING EQUIPMENT

The KII3-1 Re-recording Equipment is adapted for sound re-recording on 35-mm film from photographic sound tracks and magnetic tapes as well as for sound recording of announcer's speech.

The equipment provides for simultaneous high-standard sound re-recording from eight sound tracks recorded on 35-mm film.

The equipment provides a wide range of independent volume control and adjustment of signal frequency characteristics of the re-recorded sound from 50 to 8,000 cycles.


Fig. 31. General Scheme of KII3-1 Re-recording Equipment

The interlock drive system of the equipment ensures travelling of the film through the sound reproducing device and projector at a speed of 24 picture frames per sec with power frequency of 50 cycles.

The complete set includes the following units (Fig. 31):

- Eight sound reproducing devices.
- Re-recording control panel.
- Interlock drive system.

- Loudspeaker communication and command system.
- Recorder device.
- Sound reproducing unit and accessories.
- Projector.

4 P-1 Sound Reproducing Device (Fig. 32) of the KII3-1 re-recording equipment has a new design: besides the film-transporting mechanism with an optical scanning system it is provided with a photoelectronic multiplier and a stabilized power supply source for exposure lamp and amplifier.

The mechanism is power-driven by the 5 M-1 motor of the interlock system.

A smooth-surface drum of reduced diameter provides for a high degree of stabilization of film-travelling speed.

High-standard performance of the device and the possibility for sound reproducing from any kind of sound track are achieved by the system of optical scanning.

The device provides passage for film loops with a diameter from 2 to 15 metres.

Swift rewinding of 300-metre film roll in 1.5 minutes is an important service feature.


Fig. 32. 4P-1 Sound Reproducing Device

70 K-1 Re-recorder Control Panel (Fig. 33) centralizes the main units controlling re-recording operations, such as: adjustment of volume and frequency characteristics of 8 channels; switching of input circuits and their supplementary adjusters, amplifier adjusting devices and limiting filters, signal and loudspeaker communication system, etc.

Ease in transportation is a special feature of the panel design: it may be dismantled into three separate parts and easily assembled again on arrival.

Easy access is provided for inspection and repairs of panel parts. Highly identical channel characteristics and reduced response of adjuster link coils to exterior magnetic fields is achieved by the use of asyfferrous rings in the design of filters and coils.

The channel layout is of the two-way system using mixing transformers, and achieves high values of transitory channel fading.


Fig. 33. 70K-1 Re-recorder Control Panel

Interlock Electric Drive of sound-reproducing devices and the projector is effected by the 4M-4 generator-governor (Fig. 34) driven by a synchronous motor and provides a film-travelling speed of 24 pict. frame per sec, with a 50-cycle frequency of the power supply mains.


Fig. 34. 4M-4 Generator-Governor


Fig. 35. 7K-9 Main Terminal Board


Fig. 36. 80K-1 Control Desk

The drive is fed from the 7K-9 main terminal board (Fig. 35). Remote control of the drive is achieved by means of the 80K-1 control desk (Fig. 36) installed in the reproducer booth. Synchronous stoppage of the system may be effected from the re-recorder control panel and the signal panel in the re-recorder booth.

Loudspeaker Communication and Signal System consists of self-contained 11K-3 communication panels (Fig. 37) situated in adequate operation spots, and various units of communication and signals of the control panels 11K-2 (Fig. 38).


Fig. 37. 11K-3 Communication Panels


Fig. 38. 11K-2 Control Panels

The signal system is intended to relay general and local signals in accordance with studio routine.

The loudspeaker communication system maintains communications between the main channel control panel and corresponding operation locations, as well as direct two-way communications between any of these locations.

Recorder Device used with the КПЗ-1 re-recording equipment represents a modified set of the K3VC sound recorder serial. It consists of the following units:

- 3K-4 Sound recording apparatus (Fig. 39).
- 100Y-2 Main amplifier (Fig. 40).
- 12M-3 Power-supply device (Fig. 41).
- 6K-50 Switching unit (Fig. 42).
- 30A-3 Two-band speaker (Fig. 43).
- 7K-7 Control panel (Fig. 44).

The 7K-7 Control panel for 12M-3 power-supply unit is somewhat modified to meet requirements of supplying the equipment directly from the supply mains or through the main terminal board.

Sound Reproducing Unit and Accessories used with the КПЗ-1 re-recording equipment are of the KYCY-52 sound reproducing equipment type, detailed features of which are given in the catalogue.


Fig. 39. 3K-4 Sound Recording Apparatus


Fig. 40. 100Y-2 Main Amplifier


Fig. 41. 12M-3 Power-supply Device

The KIIT-1 Stationary Projector has been somewhat modified for use with the KIIT-3 re-recording equipment, the asynchronous electric projector motor being substituted by a 5M-2 self-synchronizing electric motor with a special reduction gear providing


Fig. 42. 6K-50 Switching Unit


Fig. 43. 30A-3 Two-band Speaker

precise conformity of the film-travelling speed within the projector and reproducing device. The projector is equipped with an automatic shutter which cuts the light beam when operation ceases.


Fig. 44. 7K-7 Control Panel

For the convenience of small studio outfits a somewhat modified set of re-recording equipment KИ3-2 is made available for simultaneous re-recording from four sound reproducing devices.



Fig. 45. General Scheme of KИ3-2 Re-recording Equipment

The design of the KИ3-2 equipment is similar to that of the KИ3-1, but it has only four 4P-1 sound reproducing devices and a 70K-2 re-recorder panel.

A general scheme of KИ3-2 equipment is shown in figure 45. Three-phase A.C. 220 V, 50-cycle mains provide the power supply.

Power consumption is approx. 10 kW.

LABORATORY EQUIPMENT

40П-1 AUTOMATIC DEVELOPING MACHINE

The 40П-1 Automatic Developing Machine (Fig. 46) performs photo-chemical processing of negative and positive 35-mm black-and-white film.

The film-threading mechanism of the machine is covered with a light-proof metal hood, allowing film processing in a lighted room.

The machine consists of two self-contained developing systems for both positive and negative films with a multi-loop arrangement of the films.


Fig. 46. 40П-1 Automatic Developing Machine, view from side of drier cabinet

The frame of the Developing Machine is conveniently furnished with four rollers, allowing easy moving of the equipment. Mechanical damage to film perforations during the processing is eliminated by the use of frictional methods of film transportation.

The Developing Machine is provided with an interlocking device for:

- a) automatically shutting off the machine drive when film breakage occurs in the wet part of the film path, and when the lower rollers of the threading magazine reach their extreme upper position;

b) sounding of warning signal when the feed reel is empty.
 A constant temperature of the developing solution at $20 \pm 0.3^\circ\text{C}$ is automatically maintained by an air-conditioning device (heater and cooling device).
 Uniform conditions of film drying are ensured by the smooth adjustment of heating temperatures and incoming air.
 The Automatic Developing Machine consists of the following parts:
 Framework on which all assemblies and units are mounted;
 Tanks;
 Wet part of the film-threading path, including film transportation mechanism;
 Drier cabinet with transportation and winder mechanisms;


Fig. 47 Developing Machine Control Panel

Air-conditioning system, consisting of: a) filter, b) blower fan with electric motor, c) baffle for controlling air blast entering the drier cabinet, d) air duct to squeegee, e) electric air-heating apparatus;

Plumbing system, consisting of: a) water conduit to washing tanks and thermostat, b) sewer line for drainage from all tanks, c) hypo line for pumping hypo from the second tank to the first one by means of electrically-driven pump, d) hypo recovery line from first tank of final washing, e) two independent developer lines (negative and positive) with pumps providing circulation of developing solution and thermostat for maintaining constant temperature, f) supply lines of developer systems with replenishment from metering tanks;

Driving mechanism which is an electric motor with variable speed reduction gear;

Control panel on which are mounted all instruments of control (Fig. 47);

Hoods for the replenisher metering tank, machine driving mechanism and solution conduits;
 Light-proof magazines.

SPECIFICATIONS

Output of the Automatic Developing Machine varies from 50 to 300 metres per hour, the developing time varying accordingly from 16 to 2.5 min.
 Power supply from three-phase 220 V, A. C. mains.
 Power consumption — 11 kW.

Time of film processing stages, in minutes:

Processing Stage	For Positive Film		For Negative Film
	Basic processing speed of 200 metres per hour	Increased processing speed: 300 metres per hour	
Developing	4	2.7	16
Intermediate rinsing	1.4	0.9	5.6
Fixing	4	2.7	16
Final washing	6.1	4	24.4
Drying	20	13.3	80
Total time	35.5	23.6	142

Overall dimensions of Automatic Developing Machine $2,900 \times 1,965 \times 1,100$ mm
 Weight of machine:
 without solutions 1,000 kg
 with solutions 2,060 kg

KII9-2 FILM SAMPLE PRINTER

The Sample Printer (Fig. 48) prints samples from 35-mm negatives on positive film, facilitating correct selection of the exposure and light number for release printing of motion picture films by continuous printing.


Fig. 48. KII9-2 Film Sample Printer

The Sample Printer provides for simultaneous printing of negative lengths of ten picture frames each (Fig. 49 and 50), showing uniform photographic density with ten different illumination intensities corresponding to those of the continuous printer; printing is achieved by contact method. This facilitates the selection of convenient exposure conditions for the printing of each negative film length by visual inspection of the positive sample and comparison of the ten images of different exposure.

Illumination of the printing gate is provided by a single light source. Interposition in the light path of the film-channel gate of

neutral-grey light filters allows for the establishment of different illumination stages.

The exposure time is adjusted by modifying the shutter slit dimensions.


Fig. 49. KII9-2 Printer Head

The head is lowered by hand.

Rewinding of positive film, shutter operations, and lifting of head are effected automatically by pressing a lever.

The Printer is equipped with adequate instruments providing for control of operating conditions and adjustment of printer lamp.


Fig. 50. KII9-2 Printer Head, rear view

SPECIFICATIONS

Capacity of film magazines:
for positive film 60 metres
for negative film 300 metres
Illuminating system fed with 110-120 V D.C.
Printer lamp 300 W, 110 V
Maximum exposure time 0.05 sec
Overall dimensions of printer (length
× width × height) 1,360×660×1,230 mm
Weight of printer 150 kg

YKA SOUND FILM PRINTERS

Available are three serial models of the YKA Sound Film Printer:

YKA-M — for mass release printing on 35-mm film from density-corrected duplicated negatives of black-and-white motion picture films.

YKA-II — for mass release printing on positive 35-mm multilayer colour film.

YKA-T — for routine printing on 35-mm film of colour and black-and-white motion pictures.

Model YKA-M

The YKA-M Printer (Fig. 51) is for mass release printing on 35-mm film of black-and-white motion pictures in motion-picture printing laboratories.

The machine prints from density-corrected duplicated negatives, each portion film roll being printed under constant illumination of the printer film gate. Changes of image printing illumination are achieved by diaphragms with different apertures, which are inserted into a special recess in the light beam path.

The YKA-M Printer is the basic model for all other printers.


Fig. 51. YKA-M Printer

Model YKA-II

The YKA-II Printer (Fig. 52) is for mass printing on 35-mm multilayer film of colour sound films in motion-picture printing laboratories.

The machine prints from multi-layer colour negatives or duplicated negatives. Change of image printing illumination is produced by the AP automatic device, which shifts a light control band with punched openings and light-filters. Fore filter holders mounted in special recesses in the path of the light beam are provided.


Fig. 52. YKA-II and YKA-T Printer

Model YKA-T

The YKA-T Printer (Fig. 52) is for regular printing on 35-mm film of black-and-white sound films in motion-picture printing laboratories.

Printing is done from negatives. Change of image printing illumination is achieved by the AP automatic device which shifts a light control band with punched openings. The machine may be adapted for regular printing of colour films if a control band similar to that of the YKA-II model and fore filters are fitted to it.

Reverse drive allows printing without rethreading the negatives.

All the models of the YKA Printer are designed on basically similar lines, which facilitates repairs and maintenance.

Spare lamps and necessary tools are supplied with the YKA Printers. On special order the YKA Printers may be provided with the following accessories:

- a) Jack-bogie for transportation (Fig. 53).


Fig. 53. Jack-bogie for transportation of YKA Printers

b) ПИМД-3 Light Control Band Puncher (Fig. 54) for punching circular apertures on a 35-mm perforated paper band when making light control bands for YKA-II and YKA-T Printers.


Fig. 54. ПИМД-3 Light Control Band Puncher

The Puncher has 20 circular punches of various diameters from 6.2 to 22 mm; the punching is done by hand.

c) The МБИ-3 Slot Cutting Machine (Fig. 55) for cutting lateral slots in negative films to achieve switchover impulses of the light control band in YKA-II and YKA-T Printers.


Fig. 55. МБИ-3 Slot Cutting Machine

d) The ПИМД-1 machine for stapling of light control bands used with YKA-II Printers, the stapling being done by means of metal staples.

SPECIFICATIONS

(common to all YKA models)

Path of film travel continuous in one plane
Direction of film travel reversible
Speed of film travel up to 1,640 metres per hour

Printing method by contact, over sprockets with bilateral compressed air contacting arrangement

Capacity of film magazine 300 metres

Printing lamp 300 W incandescent lamp

Power supply of printing lamps from 100 V D.C. line

Adjustment of power supply to printing lamp by means of rheostats

Control of power supply to printing lamp by means of measuring instruments on control board

Adjustment of sound track printing gate illumination by means of diaphragms

Electric drive 1,440 r.p.m., 220 V, 50 cycles asynchronous, three-phase electric motor

Control of electric drive by means of reversible electromagnetic starter

Local illumination by means of a 220/24 V transformer

Interlocking devices work when:

- a) printing lamps burn out,
- b) spool ends,
- c) pressing clamp of the printing assembly is not secured,
- d) negative film (in YKA-M model) tears.

Film cleaning by filtered air provided by the ventilation system

Overall dimensions of printer (length X width X height) 1,320X560X1,660 mm

Weight of Printer (in carrying cases) 1,000 kg

MM-11 FILM RENOVATION MACHINE

The MM-11 Film Renovation Machine (Fig. 56 and 57) eliminates mechanical damage from the film base, and plastifies it.

Elimination of damage is achieved by the method of superficial dissolving of the film base by highly-volatile solvents, followed by rolling of the base over a glass disk with a matte or glossy surface.


Fig. 56. MM-11 Film Renovation Machine

The essential parts of the machine are: frame, electric drive, winding device for 600 metres of film, cleansing device, post-cleansing drying chamber, pressure rollers, renovation assembly, take-off rollers, post-renovation drying chamber, winding device, mechanism for lifting and lowering of the solvent bath, solvent tank, solvent vapour-exhausting system, control panel.

The cleansing device effects damp cleansing of the film base and dry cleansing of the emulsion side.

Design of the pressure rollers ensures adjustment of the film pressure against the glass disk.
Solvents are fed automatically.


Fig. 57. MM-11 Film Renovation Machine, rear view

Quick lifting and lowering of the bath is achieved by a special mechanism, which also adjusts the humidifying rate of the glass disk and of the solvent level in the bath.

SPECIFICATIONS

- Speed of film travel 250 and 500 metres per hour
- Machine drive 0.52 kW explosion-proof asynchronous electric motor
- Film-threading chambers airtight
- Overall dimensions of machine 1,320×1,590×550 mm
- Weight of machine 400 kg

PVII-1 CUTTING MACHINE

The PVII-1 Cutting Machine (Fig. 58) cuts 32-mm film into two 16-mm film strips.

Machine drive is provided by an asynchronous electric motor with reduction gear.

The cutting mechanism consists of two knives fitted on parallel shafts and rotating in opposite directions.

The special design of the knives achieves high-precision film-cutting.

The upper and lower disk-shaped knives are protected by hoods for operation safety and prevention of winding of cut film around the knives.

A 165° portion of the lower knife circumference is encircled by the film held down by pressure roller.

Film transportation within machine is performed by the lower knife and two take-up devices.

The film-feed device is installed in the lower part of the machine, consisting of an open spool fastened on a freely-rotating spindle with an adjustable brake. The spool is designed for use with a standard 50-mm core or with a special 100 mm-diameter core.

Over the feeding device is fitted a guide roller which directs the film to a 20-tooth sprocket. The sprocket and the two upper take-up devices ensure the necessary film tension around the lower knife.


Fig. 58. PVII-1 Cutting Machine

The upper 20-tooth sprocket is intended to align tensioning on the take-up device of both film strips (after cutting) in order to prevent disalignment of the film on the knife.
The moving parts of the machine are mounted on ball bearings.

SPECIFICATIONS

Speed of film transportation 28.9 metres per min or 1,730 metres per hour
Maximum length of roll to be cut . 600 metres
Diameter of core:
with 120-metre reel 50 mm
with 600-metre reel 100 mm
Width of strip being cut 15.95 ± 0.05 mm
Power supply from three-phase, 127/220 V, 50-cycle mains
Electric motor drive 0.25 kW, 1,500 r.p.m.
Overall dimensions of the machine 1,635×570×655 mm
Weight of machine (in carrying cases) 320 kg

KCII-3 FILM STAPLER

The KCII-3 Stapler (Fig. 59) is designed for mechanical splicing of 35 mm motion-picture film ends.
Wire staples are used for the splicing. They are loaded into the stapler in 100-piece packages.


Fig. 59. KCII-3 Film Stapler

Pressing of the Stapler handles automatically releases one staple to fasten the film ends.
The Stapler is hand-operated and extremely easy to handle.

Overall dimensions of the Stapler . . 32×85×160 mm
Weight of Stapler 400 g
The staples are of 11 × 6 mm size

35-3MA-3 SOUND EDITOR

The 35-3MA-3 Sound Editor (Fig. 60) is designed for sound editing of 35-mm film and performs the following processes:
Synchronization of basic sound track, music and noises with image.

Synchronization of sound track and sound film during doubling and editing.

Synchronous visual and audible control of film from separate "picture record" and "sound record" films.

Visual and audible control of the edited film or parts of it.

The essential process of sound editing — synchronization of "picture record" and "sound record" films is achieved in two stages:

a) The "picture record" film is shifted by means of an electrical drive in forward or reverse direction at a speed of 5—30 picture frames per sec, while the "sound record" film is stationary.

b) Shifting of the "sound record" film relatively to the "picture record" film is accomplished by hand through a differential mechanism in the process of a simultaneous running of both films. This provides great convenience in editing sound films.

The Editor is adapted for simultaneous visual and audible control of image and sound track printed on one film. This makes the Editor desirable for foreign-language doubling.

The design of both Editor viewing and sound heads provides for easy threading of the film, and allows winding of film on spools or special hubs, as well as a free travel of the film from the device. Provision is made for editing of film lengths of 3—4 picture frames.


Fig. 60. 35-3MA-3 Sound Editor

Starting, stopping and reversal of sound-head and projector-head electric motors, adjustment of motor speed, switching of amplifier of sound and projector lamps are achieved by means of switches and rheostats conveniently located for operation.

The Editor's projector system does not require a darkened room.

The Editor's mechanism, its electric drive, amplifier and loudspeaker are mounted on a metal table equipped with four swivel wheels on rubber tires.

Light weight, small overall dimensions and ease in dismantling make the device convenient for use in studios and laboratories, and in outdoor conditions.

SPECIFICATIONS

Power supply from 110 V, 50 cycle A. C. mains
Power consumption 600 W

Projector (viewing) head. Film propelling is achieved by means of a maltese cross mechanism with a 16-tooth transporting sprocket. Precision workmanship of the mechanism parts ensures a noiseless performance.

The image is viewed through a picture-frame gate by means of a magnifier mounted in the gate of the film channel.

The image is magnified to 3 times its size.

Projection lamp — 15 W, 110 V with E-14 socket.

Power drive — 65 W, 110 V, 6,000 r.p.m. commutator motor.

The design of the picture-frame gate provides for easy setting of synchronous control marks on the film.

Sound head. Uniform propelling of the film is achieved by a 16-tooth sprocket.

Differential gearing allows quick and slow propelling of the "sound record" film. The head also permits use of 17.5-mm films with unilateral perforation (halved 35-mm film).

The design of the film channel gate provides for setting of synchronous control marks on the film.

A current of supersonic frequency is supplied to a 3 W, 4 V exciter lamp from a lamp generator.

The photoelectric cell — of type IIF-3.

Power drive — 50 W, 110 V, 50 cycle, 1,500 r.p.m. asynchronous single-phase condenser-type electric motor.

Film-winding device. Projector and sound heads are provided with reversible winding devices and dismountable reels of 300-metre film capacity. For small film rolls (up to 80 metres) the device has special hubs mounted on swinging brackets.

Amplifier and loudspeaker. The 1Y-7 Amplifier is installed in the metal box of the projector table.

Nominal output volume of the sound frequency is 2.5 W, the non-linear distortion factor not exceeding 3% at 1000-cycle frequency and 5% at 100-cycle frequency.

The band of reproduced frequency ranges from 50 to 8,000 cycles.

The amplifier power-supply voltage is 110 V, 50 cycles.

Electronic tubes used with the amplifier are: 6X7 (2 pcs.), 6I3 or 6J6, 6Y0186, 5I4ac.

The 4-A-18 loudspeaker is fitted on a special bracket above the sound head.

Various accessories and spare parts are supplied with the Editor.

Overall dimensions of the Editor 500×600×1,600 mm

Weight of Editor 100 kg

35-MMIC-3 SUBTITLE MAKING MACHINE

The Subtitle Making Machine prints subtitles on 35-mm film by the mechanical method with use of clichés.

Main advantages of the mechanical method of subtitle printing are: printing on ready-made films and simplicity of technical process allowing printing of subtitles at film-exchange offices.

The complete set includes:

35-MMIC-3 Subtitle Making Machine for Mechanical Printing (Fig. 61), which consists of the following units: assembly for


Fig. 61. 35-MMIC-3 Subtitle Making Machine for Mechanical Printing

dampening of the film-emulsion layer by special solvents; assembly for printing subtitles on softened film-emulsion layer by means of heated matrix-clichés; film drying assembly (after printing of subtitles) working on warm air from heating devices and fan; control panel for all elements of operation.

The subtitle printing assembly has a printing mechanism, a light-control band head and a projector device intended for checking quality of printed subtitles.

The design of the machine allows operation of the printing machine by hand, semi-automatically, or automatically.

Subtitle-printing speed ranges from 2 to 15 picture frames per sec.

Length of film rolls to be subtitled . . . 300 metres
Projection lamp 21 W, 12-16 V
Pilot lamp 8 W, 110 V

35-CIIC-1 Bench for Synchronized Printing (Fig. 62), which produces control bands for automatic subtitle printing on 35-mm film with the 35-MMIC-3 machine.


Fig. 62. 35-CIIC-1 Bench for Synchronized Printing

The bench consists of the following units: synchronizer, control band puncher, and rewinding device

Rewinder reel capacity 300 metres
Rotation of synchronizer achieved by film tenstioning
Film travelling in horizontal plane.

MPK-2 Cliché Cutting Machine (Fig. 63), which is a milling machine designed to cut ready 120x160-mm cliché plates into separate 8x22-mm matrixes.


Fig. 63. MPK-2 Cliché Cutting Machine

SPECIFICATIONS

Power supply	220/380 V, 50 cycle, A. C. mains
35-MMIC-3 machine total power consumption	2.7 kW
Overall dimensions:	
35-MMIC-3 machine	1,490×1,670×510 mm
35-CIIC-1 machine	1,300×1,250×600 mm
MPK-2 machine	970×670×630 mm
Total weight of 35-MMIC-3 sub-title making machine	approx. 650 kg

35-CIIA-2 SEMI-AUTOMATIC FILM SPLICING MACHINE

The 35-CIIA-2 Semi-Automatic Film Splicing Machine (Fig. 64) is designed for splicing 35-mm standard perforated or unperforated film. The Semi-Automatic Splicing Machine performs shearing, scraping, splicing, and pressing of spliced film ends.

The machine is for use in film-manufacturing factories, film-printing laboratories and film-exchange offices.


Fig. 64. 35-CIIA-2 Semi-Automatic Film Splicing Machine

The Semi-Automatic Splicing Machine offers speed, ease and high efficiency in the film-splicing process.

Speedy drying of the splices is ensured by a heating device. The machine has illuminators under the path of the film and above the splicing machine head.

Replaceable ground-glasses (white) and a set of nonactinic red light-filters are provided with the illuminating devices.

SPECIFICATIONS

Power supply from 110 V A. C. mains
 Straight splice, width of splice . . . 3 mm
 Film rewinding by hand or electric drive
 Geared rewind reel 1:3 ratio, acceleration type
 Commutator motor 110 V, 25 W, 3,600 r. p. m., sealed type
 Change of motor speed by pedal-operated rheostat
 Travel of motion picture in horizontal plane
 Power consumption approx. 150 W
 Height (from floor) of bench surface 780 mm
 Overall dimensions of Splicer 697 × 1,000 × 1,285 mm
 Weight of Splicer 95 kg

35-CO-1 SYNCHRONIZER

The Synchronizer (Fig. 65) is intended for synchronous alignment of two and three 35-mm motion picture films with sound and picture records during studio film-editing.

The synchronizer is a portable apparatus conveniently fitted on the Editor table without fastenings. If permanent installation is required special holes are provided to secure it to table with screws.


Fig. 65. 35-CO-1 Synchronizer

The Synchronizer's mechanism arranges the film rolls horizontally, and consists of a freely rotating 64-tooth sprocket and pressure rollers. The sprocket stops automatically when the rollers are swung aside.

The sprocket shaft of the Synchronizer is mounted on ball bearings and is easily rotated by the film edited. Easily rotating pressure rollers supported by special carriage effect

contact of film and sprocket. A spherical push-button on the top of the sprocket swings aside the pressure rollers to release the films.

Radial etches on the Synchronizer's upper crown provide perfect picture frame adjustment by alignment with interframe film spaces.

The Synchronizer has a cast support.

Performance reliability is a characteristic feature of the Synchronizer.

Overall dimensions of the Synchronizer 180 × 145 × 85 mm
 Weight of Synchronizer 3 kg

16-IIICII 16-mm FILM SPLICER

The portable, easy-to-handle 16-IIICII Film Splicer (Fig. 66) of the bench type is intended for splicing of 16-mm positive and negative film in film manufacturing factories, film printing laboratories, studios, film-exchange offices and motion picture theatres.


Fig. 66. 16-IIICII Film Splicer

A special device in the Splicer shears off damaged film ends when the film is fitted into the Splicer. Another special device scrapes emulsion off the film ends.

Quick replacement of the blades is facilitated by their construction. The design of the Splicer ensures standard and high-quality splicing.

Overall dimensions of Film Splicer 150 X 90 X 51 mm
Weight of Film Splicer 900 g

35-IIICII-3 35 mm FILM SPLICER

The 35-IIICII-3 Film Splicer (Fig. 67) is designed for splicing 35-mm positive and negative films in film manufacturing factories, film printing laboratories, studios, film-exchange offices and motion picture theatres.


Fig. 67. 35-IIICII-3 Film Splicer


The Splicer effects shearing and scraping of the emulsion off the film ends.

The design of the Splicer ensures standard and high-quality splicing.

A bottle of cement is fitted in the Splicer frame.

The blades are easily removed for grinding.

Overall dimensions of Splicer 215X130X100 mm
Weight of Splicer 4 kg

**MOTION PICTURE PROJECTORS
AND EQUIPMENT**

"UKRAINA" 16-mm PORTABLE SOUND-ON-FILM PROJECTOR

The Portable Sound-on-Film Projector "Ukraina" (Fig. 70) is adapted for showing 16-mm black-and-white and colour sound films in halls with a seating capacity of 200.


Fig. 70. Sound-on-Film Projector "Ukraina"

Its excellent performance and high technical features make it one of the best modern sound-on-film projectors in their class. Superb visual and sound projection, powerful light flood, noiseless operation, absolute fire safety, ease in film threading and

convenience in handling, low power consumption, facility in transportation, and a high degree of reliability — these and many more features characterize the unit.


Fig. 71. Set of Sound-on-Film Projector "Ukraina"

Bright illumination of a 5-sq. metre screen is ensured by an efficient light flood of 250 lumens achieved by a lamphouse of special design, lighting optics, high-quality coated lens, and a powerful projection lamp of 400 W.

Use of non-inflammable film allows installation of the Projector directly in the auditorium, eliminating the need for a special projection booth. This is an additional feature which makes the Projector desirable for schools, clubs, village motion picture theatres and home entertainment.

The Sound-on-Film Projector "Ukraina" includes the following units (Fig. 71):

III-16-1 Motion Picture Projector.

KIIV-50 Amplifier (Fig. 72).

KAT-14 Auto-transformer (Fig. 73).

ЭПН-2 Screen.

III-16-1 Motion Picture Projector. Power supply of Projector from single phase 110 V or 220 V, 50 cycle A. C. mains. Voltage variations in the mains are controlled by an auto-transformer.


Fig. 72. KIIV-50 Amplifier with Loudspeaker

Projector drive — asynchronous single-phase 35 W, 3,000 r.p.m. electric motor.

Film motion achieved by intermittent mechanism at speed of 24 picture frames per sec.

Light source — 400 W, 30 V projection lamp.

Coated anastigmatic projection lens of 1:1.2 relative aperture and 50-mm focal

length, providing brilliant illumination of 2.5-metre-wide screens. Lens with 35 or 65-mm focal length may also be installed in the Projector.

The efficient light flood with working shutter (without film) is approximately 250 lumens at 2,700-stilb brightness of projection lamp.


Fig. 73. KAT-14 Auto-transformer

The 3 W, 4 V exciter lamp is supplied from a selenium rectifier, mounted in the amplifier.

The optical scanning system is cylindrical; light slit dimensions are 1.9×0.018 mm.

Take-up and feed mechanisms ensure standard performance with 120 and 600-metre film reels.

Power consumption — 550 W.

Overall dimensions of projector (packed in carrying case) 255 × 450 × 515 mm

Weight of projector (packed in carrying case) 26 kg

KIIV-50 Amplifier with Loudspeaker. Power supply — single-phase 110 V, 50 cycle A. C. mains.

Power consumption — 100 W.

Rated output volume of sound frequency — 10 W, the nonlinear distortion factor not exceeding 3% on frequencies of 500-1,000 cycles, and not exceeding 6% on boundary frequencies.

The following electronic tubes and auxiliary lamps are used in the unit:

6ZK7, 6H9M, 6П3 (2 pcs.), 5И14С, MH3 (neon indicator), MH-15 (6.3 V, 0.28A incandescent lamp), ФЭУ-2 (photo-electronic multiplier).

The loudspeaker consists of two heads mounted in a grid-case. The sound coils of both heads are series-connected.

The frequency range of the loudspeaker extends from 100 to 6,000 cycles.

Overall dimensions of Amplifier (packed in carrying case) 295 X 410 X 560 mm
 Weight of Amplifier (packed in carrying case) 25 kg
 Overall dimensions of loudspeaker (packed in carrying case) 245 X 460 X 550 mm
 Weight of loudspeaker (packed in carrying case) 21 kg

KAT-14 Auto-transformer. The Auto-transformer permits adjustment of voltage variations of supply mains ranging from 65 to 130 V for 127 V mains, and from 165 to 230 V for 220 V mains.

The Auto-transformer is equipped with a voltmeter for control of voltage at the output terminals.

Power of Auto-transformer 750 W
 Overall dimensions of Auto-transformer 155 X 285 X 390 mm
 Weight of Auto-transformer 13 kg

ЭИИИ-2 Screen. Light and portable, the ЭИИИ-2 suspension 2,600 X 1,900 mm screen is most convenient for motion picture showings with the "Ukraina" Projector. The efficient screen surface is coated with a diffusing and reflecting barite layer which ensures a reflection factor of 0.7 to 0.75.


Fig. 74. Sound-on-Film Projector "Ukraina", packed in cases

Light in weight and simple in design the Screen can be quickly hung up for showings and as quickly taken down.

Overall dimensions of Screen (packed in jacket) 100 X 100 X 2,700 mm
 Weight of Screen 14 kg

A view of the complete "Ukraina" Sound-on-Film Projector set packed in carrying cases is shown in Fig. 74.

Total weight of the "Ukraina" Sound-on-Film Projector...99 kg.

KIIC-M 35-mm PORTABLE SOUND-ON-FILM PROJECTOR

The improved Portable KIIC-M Sound-on-Film Projector (Fig. 75) is designed for demonstrating 35-mm black-and-white and colour films.


Fig. 75. KIIC-M Sound-on-Film Projector

Reliability of the Projector design has been tested by operation in conditions of prolonged travel. This and splendid visual and sound projection plus ease in handling all mechanisms places the projector among the best in its class (Fig. 76 and 77).

A new design of the picture frame correction device has noticeably reduced wear of film and ensures steadiness of screen picture.

Use of ball bearings and high quality heat-treated material ensure lasting trouble-free operation of the projector.

The Projector is successfully employed in auditoriums with a seating capacity of 200 and is installed directly in the auditorium, eliminating the need of a special projection booth. This and other features make it particularly suitable for use in schools, clubs and village motion picture theatres.

SPECIFICATIONS

Power supply — from 110 or 220 V 50 cycle A. C. mains. Voltage variations in the mains are controlled by a special auto-transformer.

The Projector is supplied with two powerful specially-coated projection lenses with 52.5-mm diameter, relative apertures of 1:1.9 and 1:2.1, and focal lengths of 90 mm and 120 mm respectively.

Efficient illuminating power of Projector 250 to 300 lumens
 Light source 400 W, 30 V (type K-22) projection lamp


Fig. 76. KIIC-M Sound-on-Film Projector; right-side view with cover removed


Fig. 77. KIIC-M Sound-on-Film Projector, left-side view with cover removed

A 3 W, 4 V exciter lamp for scanning the sound track has a prefocusing socket which ensures perfect alignment of incandescent filament and sound track without adjustment when replacing lamp.

The exciter lamp is power-supplied from D. C. line.

Illuminating lamp of Projector — 8 W, 110 V.

The single-stage antimony-caesium photo-electronic multiplier of the $\Phi\Theta\Upsilon$ -1 type used in the Projector possesses great sensitivity, uniformity of characteristics, a low noise level and outstanding spectral characteristics, which ensure splendid sound reproduction of colour film sound tracks.

Power drive — single-phase, asynchronous, 50 W, 1,425 r.p.m. motor of Δ O-50 type.

Deviations in vertical steadiness of image in picture frame gate — not exceeding 0.04 mm, in horizontal steadiness — not exceeding 0.02 mm.

Magazines have film capacity of 300 metres.

Tripod, film rewinding device, tools and necessary spare parts are supplied with the KIIC-M projector.

The KIIC-M projector is supplied with either a 4-KY-12 amplifying unit with loudspeakers (Fig. 78), or KIIV-50 amplifying unit with loudspeaker (See fig. 72).


Fig. 78. KIIC-M Sound-on-Film Projector, complete set

Nominal output volume of amplifiers 10 W

Peak output volume 15 W

Reproduced range of frequencies . . . 80 to 6,000 cycles

Amplifiers 4-KY-12 and KIIV-50 are supplied with the following electronic tubes and auxiliary lamps: 6-X-7, 6H9n, 6I13, (2 pcs.), 5I14C, MH3 (neon indicator), MH-15 (6.3 V, 0.28 A incandescent lamp for KIIV-50 amplifier only).

The amplifiers have a selenium rectifier which supplies power to the Projector's 4 V, 3 W, D. C. exciter lamp.

The amplifiers enable playing of gramophone records with the aid of pick-up.

Power consumption of projector — approx. 660 W.

Amplifiers, connecting cables and accessories are provided with carrying cases for storage and transportation convenience.

**OVERALL DIMENSIONS AND WEIGHTS
OF PROJECTOR SET:**

Projector with cover (220 × 415 × 420 mm)	27 kg
Auto-transformer with cover (155 × 390 × 285 mm)	13 kg
Case with accessories (260 × 460 × 625 mm)	23 kg
Tripod (folded) (220 × 220 × 1,300 mm)	5 kg
Case with loudspeaker (4-KY-12) (250 × 340 × 540 mm)	21 kg
Case with loudspeaker and amplifier (4-KY-12) (250 × 340 × 540 mm)	27.5 kg
Case with amplifier (KIIV-50) (285 × 410 × 560 mm)	25 kg
Case with loudspeakers (KIIV-50) (245 × 460 × 550 mm)	21 kg

The Sound-on-Film Projector KIIC-M is shown packed in carrying cases (amplifier excluded) in Fig. 79.


Fig. 79. KIIC-M Sound-on-Film Projector, packed in carrying cases

**KIIT-1 35-mm STATIONARY SOUND-ON-FILM
PROJECTOR**

The modern sound-on-film projector must produce bright and sharp images and provide clear and undistorted sound.

The KIIT-1 Stationary Sound-on-Film Projector (Fig. 80) meets these requirements in full. It is intended for visual and


Fig. 80. KIIT-1 Sound-on-Film Projector

sound projection from 35-mm film with standard sound track. By special order the Projector may be supplied with KYCY-52, K3BT-3 sound reproducing devices.

Fig. 81 shows projector head, and Fig. 82 — its kinematic diagram.

The projector possesses the following advantages:

High standard visual projection. The specially-coated projection lens produces pictures uniformly sharp in screen centre as well as screen edge.

Uniform illumination of screen is ensured by a special sphero-elliptical mirror fitted in the arc lamp.

The powerful arc lamp with high-intensity carbons provides a light flood through film gate (shutter open; without film) sufficient to illuminate a screen area of 30 sq. metres at approx. 100-lux illumination.

Use of automatic carbon feed in the arc lamp eliminates any perceptible flickering in screen illumination.


Fig. 82. KIIT-1 Projector Head, kinematic diagram


Fig. 81. KIIT-1 Projector Head

Perfect steadiness of screen picture is ensured by a long film channel with lateral registration and a rigid film aperture frame.

Warping of the film in the film aperture is prevented by a heat-protection filter installed in the film channel.

High Fidelity Sound. Adequate sound optics provide sufficient brightness to the sound slit.

An oil stabilizer ensures perfect uniformity of film speed past the sound slit.

Greater Film Protection. While travelling through the Projector the film is not subjected to undue strain.

A heat protection filter placed in the path of the light beam in the film channel and a cylindrical shutter between filter and channel prevent excessive heating of the film.

Fire protection is ensured by automatic fire shutters in shutter box and film channel during breakdowns.

Ease of Operation. Precise alignment of picture frame and film aperture during film threading is made easy by an 8 W, 110 V pilot lamp. A second pilot lamp, which is automatically switched on when the lamphouse gate is open, facilitates maintenance of arc lamp.

Carbons are fed automatically. Maintenance of the arc lamp is also facilitated by projection of the positive carbon crater image on the lamphouse screen.

A wide range of feeds and feed ratios of positive and negative carbons, and independent and simultaneous adjustment of carbons by hand permit use of various types of high-intensity and plain carbons.

Fire channels can be opened to clean out the carbon deposits.

Cleaning of film channel is facilitated by the easily detachable aperture framing.

On the side of the operator the space between lamphouse and lens is closed off by the shutter box and lensmount shutter.

All projector and arc lamp controls are conveniently located and provide ease of operation.

Long Service. Automatic lubrication of the Projector Head mechanism is effected by circulating oil.

The Projector Head mechanism is hooded.

Complete Fire Safety. Fire safety is achieved by:

- fireproof magazines with fireproof channels;
- automatic safety shutter in film channel which works when film breaks in path between film channel and drum of the maltese cross mechanism;
- governor controlled automatic safety shutter in shutter box;
- heat protection filter which eliminates infra-red spectrum portion of the light beam.

SPECIFICATIONS

The KIIT-1 Projector is adapted for use in club and theatre auditoriums with seating capacity of 2,000.

Operated on high-intensity carbons with an axis brilliancy of 56,000 stillb, the Projector has an illuminating power of 3,000 lumens.

The film run is not enclosed

The film guides in the film channel are replaceable. Lateral film vibration is limited by cushioning spring insert. Dimensions of film aperture are 20.9 X 15.2 mm

Cooling of the light beam falling on the film aperture of the film channel is achieved by a heat protection filter cooled by the shutter. Light transmission factor of filter — 84 1/2%.

Intermittent movement of the film is produced by a maltese cross mechanism.

Perfect central registration of picture frame and aperture is assured

The shutter has cylindrical shape. Shutter light efficiency factor is approximately 0.5.

The sound reproducing system is equipped with a smooth-surface drum with oil speed stabilizer.

Exciter lamp — 50 W, 10 V.

Dimensions of sound slit — 2.15×0.02 mm.

Photoelectronic multiplier — type ФЭУ-1.

Type II-4 projection lens has relative aperture of 1:2. Focal lengths — 9, 10, 11, 12, 13, 14, 15, 16, 18 cm. (Projector is supplied with one of the lenses at purchaser's choice). Lenses are of specially-coated type. Their light transmission factor is 85–88%.

Power drive for Projector's mechanism is provided by three-phase asynchronous 0.25 kW, 1,400 r.p.m., 127/220 V, 50 cycle electric motor of II type.

Automatic lubrication of Projector Head mechanism is achieved by geared oil pump.

Automatic rewinding is of the dry friction type.

Projection light source — an arc lamp of the mirror type. The spheroidal mirror has a 360-mm diameter.

High-intensity carbons: the positive carbon with 8-mm diameter, negative carbon with 7-mm diameter.

Operation of lamp is achieved by direct current not exceeding 65 amps, voltage at terminals of lamp — 36–40 V.

The carbons are fed automatically. The positive and negative carbon feed ratio ranges from 1:6 to 1:1.


Fig 83. KIIT-1 Projector Arc Lamphouse

Carbon feed mechanism is driven by 36–40 V, approx. 15 W D.C. motor. Speed range of motor, controlled by rheostat, is from 3,000 to 2,000 r.p.m. Motor switches on automatically by means of relay on striking the arc.

The carbon holders are solid cast. The positive carbon holder is adjustable for carbons of various diameters. The position of the negative carbon head is also adjustable. Correct building-up of arc is ensured by energized magnetic arc control. Additional two-way shifting of the negative carbon holder allows for handling of abnormal crater formation.

The lamphouse (Fig. 83) has double lateral doors which slide up when the lamphouse opens.

Kinematic diagram of the arc lamp is given in Fig. 84. A baffle is provided to control ventilation in the hot-air exit pipe. The glowing carbon tips are projected by means of an optical system on a screen fitted in the upper part of the lamphouse. This affords visual control of the arc gap.


Fig. 84. KIIT-1 Projector Arc Lamp, kinematic diagram

The lamphouse pilot lamp operates on 100 or 220 V current and has an E-27 socket.

Upward and downward tilting range of the Projector optical axis is 6° and 17° respectively.

Capacity of reels up to 400 metres of film

Overall dimensions of Projector . . . 640 × 1,300 × 1,880 mm

Height of optical axis from floor level 1,250 mm

Weight of Projector 300 kg

K3BT-3 SOUND-REPRODUCING AMPLIFYING DEVICE

The new model of the K3BT-3 Sound-reproducing Amplifying Device reproduces sound records from 35-mm film and is adapted for work with pick-up.

The amplifying device is designed for large motion picture auditoriums with 1500 seats.

The two-band Sound-reproducing Device ensures faithful reproduction.


Fig. 85. 50Y-4 Amplifier


Fig. 86. 30A-9 Loudspeaker

The two-band amplifying system divides the sound frequency band into high-frequency and low-frequency circuits at the input and achieves subsequent amplifying of the signals in two narrow strip bands whose output is directly connected with high- and low-frequency two-way loudspeakers.

The outstanding features of the device are: practically no distortion, high output volume, low level of interference, uniform distribution of sound over entire auditorium, and highly reliable performance.

The Sound-reproducing Device consists of the following units: two 50Y-4 amplifiers (racks), (Fig. 85); one of these is a spare unit;

two 30A-9 two-way loudspeakers (Fig. 86);
three 80Y-3 photo stages fitted on projectors;
a 25A-3 monitor loudspeaker (Fig. 87);
two 6K-78 volume control resistors.

Fig. 87. 25A-3 Monitor Loudspeaker

SPECIFICATIONS

Nominal output volume of amplifier — 40 W.

Reproduced frequency range — from 40 to 10,000 cycles.

Within this range the nonlinear distortion factor does not exceed 2%.

Power supply — from 50 cycle A.C. mains with 127 or 220 rated voltage. Range of voltage regulation is from 85 to 135 V and from 170 to 220 V respectively. Connection of device to line is single-phased.

Power supply of the exciter lamp is by rectified and well flattened current. Power of exciter lamp — 50 W, 10 V.

The device operates with three stationary KIIT-1 Sound-on-Film Projectors, of which two are operated and one is a spare unit.

Correction of high and low frequencies within operating range is possible.

Supply voltage of the photoelectronic multiplier — 230 V \pm 10%.

Resistance of volume control — 10,000 ohms \pm 10%.

Types of tubes and lamps used:

in amplifier — 6XK7 (12 pcs.), 6H8C (6 pcs.), 6II3 (4 pcs.), I-807 (16 pcs.), 5I14C (6 pcs.), 5I13C (4 pcs.), BI-176 (4 pcs.), MH-7 (4 pcs.), MH-3 (4 pcs.), incandescent lamps 6.3 V, 0.28 A (4 pcs.);

in photo-cell stage — 6XK7 (3 pcs.);

in volume control panel — 5 W, 110 V, incandescent illumination lamp.

The loudspeaker units are adapted to separately reproduce high and low portions of the operating frequency range.

The loudspeaker set includes two speaker units:

a) low frequency unit, consisting of electrodynamic head of low-frequency speaker (2 pcs.), low frequency horn and directional baffles (right and left);

b) high frequency unit, consisting of electrodynamic head of high frequency speaker (2 pcs.), and high-frequency multicellular horn.

Rated output power of set 20 W

Frequency of division 550 ± 50 cycles

Average response of set not lower than 35 units

Directional characteristic variations do not exceed 5 db at 8,000 cycles frequency within a horizontal 80° angle and a vertical 50° angle.

Divergence of response does not exceed ± 6 db.

KVCY-52 SOUND-REPRODUCING AMPLIFYING DEVICE

The Stationary Sound-reproducing Amplifying Device KVCY-52 is designed to reproduce sound tracks of 35-mm film and is adapted for work with pick-up.

The amplifying unit is adapted for operation in motion picture theatres and club halls with a seating capacity of 800.


Fig. 88. 70Y-5 Amplifier


Fig. 89. 10K-4 Junction Box

The Device includes the following units:
70Y-5 Amplifier (Fig. 88) with hose and 10K-4 junction box (Fig. 89);

22B-3 Rectifier (Fig. 90);

10B-1 Rectifier (Fig. 91);

6K-16 Volume Control Panel (Fig. 92);

6Y-12 Dividing Filter;

25A-3 Monitor Loudspeaker (Fig. 93);

30A-3 Two-way Loudspeaker (Fig. 94).


Fig. 90. 22B-3 Rectifier


Fig. 91. 10B-1 Rectifier


Fig. 92. 6K-16 Volume Control Panel


Fig. 93. 25A-3 Monitor Loudspeaker


Fig. 94. 30A-3 Two-way Loudspeaker

SPECIFICATIONS

Output volume of amplifier — 40 W.
 Band of reproduced frequencies ranges from 60 to 8,000 cycles.
 Nonlinear distortion factor when operating at rated volume does not exceed 2.5% in the 70–5,000 cycles frequency band.
 Power supply — from A. C. single-phase 50 cycles mains of 127 or 220 voltage.

Regulation of the voltage supply is achieved by the 22B-3 Rectifier, which maintains normal operation of the device when voltage in mains falls from 127 to 85 V, or in 220 V mains to 170 V.

The exciter lamp is supplied with rectified and well flattened current. Exciter lamp — 10 V, 50 W.

The device is designed for use with two projectors; two devices may be used with three projectors of the stationary KIIT-1 type and the ФЭУ-1 photoelectronic multipliers.

An easy switch-over is ensured to the spare amplifying unit, and a spare KIIT-1 projector by the wiring scheme.

High and low frequency correction is provided within operating frequency range.

The 30A-3 Two-way Loudspeaker can reproduce high and low frequencies separately and represents a combination of two specific loudspeakers — of high frequency and low frequency response. Division of the operating frequency band is achieved by the 6Y-12 Filter, inserted between the 70Y-5 Amplifier and 30A-3 Loudspeaker.

Loudspeaker's maximum excitation power (at 25 V excitation voltage) is 65 W.

Photoelectronic multiplier power supply voltage — 230 V ± 10%.
 Resistance of volume control — 30,000 ohms.

The following types of tubes and lamps are used:

in 70Y-5 Amplifier — 6XK7 (2 pcs), 6H7 (2 pcs), 1-807 (4 pcs), 5L14C (2 pcs), MH-3 neon tube;

in 22B-Rectifier — 6Y-176, 6.3 V, 0.25 A tube;

in 6K-16 Volume Control Panel — 8-15 W, 110 V panel illuminating lamp.

Specifications of Two-way Loudspeaker:

Rated volume of output 20 W

Division frequency 850 cycles

Deviation of frequency characteristics within operating frequency range not exceed ± 8 db.

Average irrelative response of Loudspeaker not less than 22 units

Directional characteristics of Loudspeaker lie within an angle of ± 40° with response divergencies not exceeding 6 db. in the frequency band up to 6,000 cycles.

Overall dimensions of the Device:

70Y-5 Amplifier 418 × 398 × 200 mm

22B-3 Rectifier 418 × 398 × 200 mm

10B-1 Rectifier 232 × 280 × 130 mm

6Y-12 Dividing Filter 214 × 158 × 158 mm

6K-16 Volume Control Panel 214 × 158 × 158 mm

30A-3 Two-way Loudspeaker 730 × 1,600 × 700 mm

Weight of sound-reproducing device 187 kg

PY-65 SWITCHING RACK

The PY-65 Switching Rack (Fig. 95) is used in motion picture theatres and is designed for switching, protection and wiring of power and acoustic power supply lines of units (in projection booth and auditorium).
The design provides easy access to the terminal board, fuses and contactors.


Fig. 95. PY-65 Switching Rack

A TC-5 Light Dimmer is installed in the top part of the Switching Rack structure.

The Switching Rack provides for an independent supply of power and light from two separate inputs: "power input" and "illumination input".

The Switching Rack has switching of supply circuits for projector arcs and allows switching of an operating projector arc from one supply circuit to another.

The control panel with a voltmeter is located on the front panel of the Switching Rack at a height convenient for operation.

The electric wiring diagram provides:

- power input switch;
- three throw-over switches for three-phase on-and-off switching of the auditorium illumination circuit;
- interconnecting means of two rectifiers and three projectors;

means of quick connection of any of two operating rectifiers to the third spare projector;

possibility to substitute the third (spare) BC-60A rectifier for any of the two operating rectifiers by means of a terminal board switch;

contactors to switch the rectifiers to the power supply circuit, the contactors being controlled from the projectors; fuses.

SPECIFICATIONS

The Switching Rack is designed to switch 220/380 V, 50 cycle A. C. line circuits.

Normal operation is maintained with voltage variations of + 10% and - 20%.

Overall dimensions of Switching Rack	1,450 X 740 X 440 mm
Weight of Switching Rack	150 kg

KAT-14 AUTO-TRANSFORMER

The KAT-14 Auto-Transformer (Fig. 96) is designed for portable 16- and 35-mm sound-on-film projectors equipped with 400 W, 30 V projection lamps, the power supply coming from a single-phase 127-220 V, 50 cycle, A. C. line.


Fig. 96. KAT-14 Auto-Transformer

The Auto-Transformer's voltage is adjusted by means of a rotary switch without interrupting the electric circuit. The Auto-Transformer is supplied with voltmeter.

Protection from short-circuit is achieved by a fuse.

Switch, voltmeter, fuse and receptacles for connecting cables are mounted on the Auto-Transformer panel and have a hood to protect them from mechanical injuries during transportation.

The Auto-Transformer is adapted for use in various branches of industry for maintenance of constant voltage values of 110 V, 30 V and 5 V.

SPECIFICATIONS

Rated power 750 W
Rated output voltages 110 V, 30 V, and $5 V \pm 2.5\%$
Ranges of voltage variations:
for 127 V mains from 65 to 130 V
for 220 V mains from 165 to 230 V
Rate of voltage adjustment registered
by voltmeter 12 steps of 5-7 V each
Maximum overheating of winding $65^{\circ} C$
Overall dimensions of Auto-Transformer 390 X 285 X 155 mm
Weight of Auto-Transformer 13 kg

TPД-50 REACTIVE ARC TRANSFORMER

The TPД-50 Reactive Arc Transformer (Fig. 97) supplies motion picture projector arcs with stabilized alternating current. The Reactive Transformer cuts power consumption to half compared to consumption through a ballast (active) resistance.


Fig. 97. TPД-50 Reactive Arc Transformer

The Transformer is provided with a steel housing. To ensure cooling of Transformer the housing is perforated and its bottom is open.

Connection of arc lamp to line circuit is achieved through two single-pole arc current regulation switches and terminals which

are mounted on a panel on top of the Transformer. The panel is protected by the housing cover.

The circuit diagram of the Transformer is shown in Fig. 98; external characteristics in Fig. 99.


Fig. 98. TPД-50 Transformer Circuit Diagram

The switches regulating arc current have four steps. Apertures for the leads from "line" and "to arc" are provided in the side walls of the housing.

The Transformer is easily mounted on the projector base.


Fig. 99. TPД-50 External Characteristics of Transformer

SPECIFICATIONS

Power supply line of transformer

Rated voltage 110/220 V
 Frequency 50 cycles
 Phases 1
 Efficiency 0.8
 Power factor 0.5
 Tolerated variations in supply voltage:
 110 V line from 85 to 127 V
 220 V line from 180 to 230 V

Power supply line of arc

Rated arc voltage 20-30 V
 Rated arc current 50 A
 Rated load 1,500 W
 Stability of arc current at variations
 of arc voltage from 20 to 30 V 50 ± 5 A
 Tolerated range of arc current re-
 gulation from 40 to 55 A

Other data

Operating cycle 20-minute interval after 25-minute
 operation time
 Dimensions of transformer (without
 switch handles) 430 X 320 X 230 mm
 Overall dimensions of transformer
 (with switch handles) 450 X 320 X 285 mm
 Weight of transformer 45 kg

ЭПН PORTABLE DIFFUSE REFLECTION SCREEN

The portable hang-up screen (Fig. 100) is best suited for film showings indoors (club auditoriums, special demonstrating rooms, etc.).

The screen is hung up on the wall or on other support by rope rings at any desired height from floor level.

The screen is easily and quickly mounted before showing. It is light in weight and simple in design.


Fig. 100. ЭПН Portable Diffuse Reflection Screen

A special mechanism performs unwinding and rolling up of screen.

The screen is made of seamless linen fabric coated with special light-reflecting layer.

The screen is provided with a jacket for convenient transportation and storage.

SPECIFICATIONS

The portable screen is made available in two models:
 ЭПН-1 with efficient surface 2,000 X 1,900 mm .
 ЭПН-2 with efficient surface 2,600 X 1,900 mm
 Screen reflection factor 0.7-0.75
 Overall dimensions and weight of the folded screen:
 ЭПН-1 ЭПН-2
 Length 2,100 mm 2,700 mm
 Diameter 90 mm 100 mm
 Weight 9 kg 14 kg

TC-5 LIGHT DIMMER

The TC-5 Light Dimmer (Fig. 101) produces smooth dimming and gradual lighting of auditorium.

The Dimmer is an electromagnetic device of static action. All parts of the Dimmer are stationary and not subjected to wear. This outstanding feature places it above all other types of dimmers.


Fig. 101. TC-5 Light Dimmer

Dimming and lighting is achieved automatically on turning a switch handle.

The Dimmer is made available in two models:

- a) for connection to three-phase 127 or 220 V 50 cycle line;
- b) for connection to three-phase 220 or 380 V 50 cycle line.

SPECIFICATIONS

Rated load 5 kW
Minimum load 2 kW

Normal operation of the Dimmer is maintained at uneven inter-phase load distribution ranging up to 10% and at line voltage variations within limits of -20% and +5% of its rated value.

Lamp voltage drops from 80 to 20% of its rated value during dimming process. The process takes 25-45 sec. Lamp voltage rises from 20 to 80% of its rated value to accomplish lighting process in no more than 25 sec.

The Dimmer is designed for mounting in top part of PY-65 Switch Rack or independent installation on bracket or table. Cooling is effected by surrounding air. Air temperature must not exceed 35° C and relative air moisture must be not over 70%.

Overall dimensions of Dimmer . . . 760 X 400 X 430 mm
Weight of Dimmer 150 kg

A3C-9-10 AUTOMATIC FIRE SHUTTERS

The A3C-9-10 Automatic Fire Shutters are designed for installation in projection booth to shut off projection aperture and projection booth windows automatically during fire emergency within booth. The booth is thereby isolated from auditorium and flame and smoke are prevented from reaching the auditorium. A fire alarm and emergency lighting are switched simultaneously.

The Automatic Fire Shutter set (Fig. 102) includes equipment for two projectors.

Two automatic shutters with special protective glass for projection apertures, type A3C-9;


Fig. 102. A3C-9-10 Automatic Fire Shutters

Two automatic shutters with special protective glass for projection booth windows type A3C-10, and push-button switch;

Two automatic AB3-3 switches.

The automatic switch is fitted in the projector near the film channel and, if film ignites, it achieves automatic disconnection of the power supply line of electromagnets working four automatic shutters.

The ЭПН-1 power supply device feeds circuits of the automatic shutter electromagnets with rectified current and achieves switching of emergency lighting and fire alarm.

The ЭПН-1 power supply device has a handle for manual switching of the emergency lighting and is provided with pilot lamps.

The B-16 push-button switch is fitted at exit of booth and serves for manual switching off of supply to electromagnets of the Automatic Shutters if fire breaks out.

SPECIFICATIONS

Power supply — from single-phase 110, 120 or 220 V, 50 cycles A. C. line. Normal operation of the device is maintained with variations of line voltage ranging from -20% to +10%.

The device is designed for permanent connection to line. Surrounding air temperature must not exceed +35° and relative air moisture not over than 75%.

The design of the device's separate elements allows for combined operation of one, two or three projectors, effected by automatic shutters.

The power supply device is designed to feed the emergency lighting circuit up to 500 W lamp power.

The rectifier consists of BC-45-70 selenium bridge and TP-220 transformer. It is designed to supply electromagnets of 6 automatic shutters.

Dimensions of the shutter aperture — 150 X 150 mm.

The A3C-9 Automatic Shutter glasses are of the specular type with 85% transparency.

Overall dimensions of the Automatic Shutter:

A3C-9	480 X 240 X 50 mm
A3C-10	480 X 240 X 130 mm

Weight of set of Automatic Shutters approx. 25 kg.

J13-2 AUTOMATIC SCREEN CURTAIN WINCH

The J13-2 Automatic Winch effects opening and closing of the screen curtain used for dual purpose of decorating screen and protecting it against dust and dirt.

SPECIFICATIONS

Range of regulating curtain opening — from 1 to 10 metres.
Speed of curtain's travel — 0.2 metres per sec.
An electromechanical control device effects remote operation of the winch from one or several locations.
The starting device consists of magnetic starters with push-button controls.
Power of electric motor is 0.52 kW.
Power supply — from 220 or 380 V A. C. mains.
A hand drive is provided to adjust and check operation of the winch mechanism.
Overall dimensions of the winch . . . 600 X 500 X 308 mm

ΦC-5 FILM STORAGE CONTAINER

The ΦC-5 Film Storage Container is two-storied and has twelve separate sections provided with lugs for rolling out film reels.

A special film humidifying device in the Container ensures film preservation.

"KIEV-1" MOBILE POWER PLANT

The "Kiev-1" Mobile Power Plant (Fig. 103) is designed to feed portable projectors of the K11C-M, "Ukraina" or of similar types with single-phase alternating current.

The power plant may also be used to feed the lighting system and instruments if their power consumption does not exceed 750 W.

The power plant has small dimensions and light weight. These features add to its mobility.


Fig. 103. "Kiev-1" Mobile Power Plant

The power plant consists of a single-phase alternating current generator and a two-stroke cycle internal combustion engine, attached to the generator by a special flange.

A flexible coupling transmits rotary motion from engine to generator. The generator and engine are mounted on a framework which protects the plant from damage.

A fuel tank of approx. 5-litre capacity is fitted in the upper part of the framework.

GENERATOR SPECIFICATIONS

The 9M-3 Generator is a self-exciting single-phase A.C. electric machine, and has the following characteristics:

Voltage	115 V
Current	7 A
Power	750 W
Frequency	50 cycles
Revolutions per min.	3,000

ENGINE SPECIFICATIONS

The engine is single-cylinder, two-stroke cycle, internal gasoline combustion, with two-channel scavenging.

Constancy of speed is maintained by a centrifugal governor set for 3,000 r. p.m., which corresponds to frequency of 50 cycles.

Power (at rated 3,000 r. p.m.)	2 h.p.
Cylinder bore	52 mm
Piston stroke	58 mm
Capacity of cylinder	123 cm ³
Rate of compression	6.5
Fuel	gasoline and oil compound in 1:25 ratio
Ignition system	M-24 magneto
Spark plug	M14 X 1.25 type HA 11/10
Carburettor	K-30
Cooling	air-cooled
Overall dimensions of power plant	750 X 400 X 480 mm
Weight of power plant	60 kg

A set of accessories and tools packed in a special box are supplied with the power plant.

K3C-5 MOBILE POWER PLANT

The K3C-5 Power Plant (Fig. 104) is designed to supply 35-mm or 16-mm motion picture projectors of KIIC-M, "Ukraina" or of similar types with single-phase alternating current.


Fig. 104. K3C-5 Mobile Power Plant

The power plant may also be used to feed the lighting system and instruments if their power consumption does not exceed 750 W.

The power plant contains an internal combustion engine which drives a single-phase alternating current generator by means of flexible coupling.

Engine and generator are mounted on a rigid base frame.

To facilitate transportation the power plant may be mounted on four wheels.

Regulation of voltage of power plant is achieved by a shunt rheostat mounted in the generator panel.

ENGINE SPECIFICATIONS

The JI-3/2 engine is vertical, single-cylinder, four-stroke cycle, gasoline-operated. Centrifugal governor set for 1,500 r. p.m. maintains constancy of speed.

Power (at rated 2,200 r. p.m.)	3 h.p.
Cylinder bore	65 mm
Compression rate	4.5-5.0
Capacity	298 cm ³
Cooling	thermosyphon, capacity of tubular radiator
	5 litres
Fuel supply	by gravity
Fuel tank capacity	4.3 litres
Fuel consumption	up to 335 gr per h. p. per hour

Lubrication splash type
Lubrication system capacity 1.7 litres
Magneto type M27B
Spark plugs M-12/20 type
Carburettor K-12 T type
Overall dimensions of engine 760 X 520 X 510 mm
Weight of engine 81 kg

GENERATOR SPECIFICATIONS

The 9M-1 generator is self-exciting, single-phase A. C. electrical machine.

Voltage 115 V
Current 7 A
Power 750 W
Frequency 50 cycles
Revolutions per minute 1,500
Exciter voltage 32 V
Exciter current 3.8 A
Shunt rheostat resistance not less than 4.2 ohms
Overall dimensions of generator 480 X 300 X 260 mm
Weight of generator 42 kg
Overall dimensions of power plant 1,200 X 665 X 950 mm
Weight of power plant 180 kg

A ballast resistance to maintain load when the projector is being switched off, a connecting cord, special tools, a set of spare parts and other accessories are supplied with the Power Plant.

K9C-4 MOBILE POWER PLANT

The K9C-4 Mobile Power Plant (Fig. 105) is designed to supply projector installations with two projectors of the KIIC-M, "Ukraina" and of similar types with single-phase alternating current.

The power plant is installed in a special room fitted for operation of internal combustion engine.


Fig. 105. K9C-4 Mobile Power Plant

The power plant includes an internal combustion engine which drives a single-phase alternating current generator by flexible coupling. Engine and generator are mounted on a rigid base frame.

The power plant is provided with a switchboard installed in the power plant room, and with switchbox, mounted in proximity of projectors.

The switchboard has a voltmeter for visual control of power plant voltage, fuses, switches and plug sockets. The switchbox is also provided with fuses, switches and plug sockets.

ENGINE SPECIFICATIONS

The D8/3 engine is vertical two-cylinder, four-stroke cycle, gasoline-operated type. Constancy of speed is maintained by a centrifugal governor set for 1,500 r. p. m.


- Power (at rated 2,200 r. p. m.) 6 h. p.
- Cylinder bore 65 mm
- Compression rate 4.5-5.0
- Capacity 597 cm³
- Cooling — thermosyphon, capacity of
 - tubular radiation 7 litres
- Fan with two blades
- Carburettor type K-12Д
- Fuel supply by gravity
- Fuel tank capacity 15-18 litres
- Fuel consumption up to 335 gr per h. p. per hour
- Magneto type MMJIC-2
- Spark plugs M 12/20 type
- Lubrication splash lubricator
- Overall dimensions of engine 670 X 450 X 800 mm
- Weight of engine 100 kg

GENERATOR SPECIFICATIONS

The A11H-28.5 generator is a self-exciting single-phase alternating current electric machine.

- Voltage 120 V
- Current 15 A
- Power 1,800 W
- Revolutions per minute 1,500
- Frequency 50 cycles
- Exciter voltage 54 V
- Exciter current 3.6 A
- Shunt rheostat resistance 21 ohms + 15%
- Weight of generator 75 kg
- Overall dimensions of power plant 1,100 X 550 X 900 mm
- Weight of power plant 230 kg

The power plant is supplied with ballast resistance to maintain load when shutting off projector, connecting cord, special tools, set of spare parts, fuel tank and other accessories.

MEASURING AND INSPECTION INSTRUMENTS

16-MM-3 FILM LENGTH COUNTER

The Counter (Fig. 106) registers the length of 16-mm film. It is designed to operate on service benches with vertical rewinding devices with both forward and reverse run of film.


Fig. 106. 16-MM-3 Film Length Counter

A clearly visible dial shows length of film in hundreds, tens, ones and tenths of metres.

During forward run the counter adds the length of film, during reverse run it subtracts.

Overall dimensions of counter 119 X 106 X 68 mm
Weight of counter 1.2 kg

35-MM-3 FILM LENGTH COUNTER

The Counter (Fig. 107) registers the length of 35-mm film. It is designed to operate on service benches with horizontal rewinding devices with both forward and reverse run of film.


Fig. 107. 35-MM-3 Film Length Counter

The dial shows hundreds, tens, ones and tenths of metres during the forward run. During reverse run length of passing film is subtracted from previously registered length.

Overall dimensions of counter . . . 137 × 106 × 80 mm
Weight of counter 1.5 kg

CO-301-1 FILM INSPECTION MAGNIFIER


Fig. 108. CO-301-1 Film Inspection Magnifier

The CO-301-1 Film Inspection Magnifier (Fig. 108) is adapted to define the degree of wear of 16 and 35-mm film perforations by alignment with a standard reticule put on the lens surface.

The eye-glass of the magnifier is fitted with a device for focusing a sharp image of the reticule in the ± 2 dioptric range.

The device is of the portable pocket type.

SPECIFICATIONS

Lens magnifying power 10 ×
Overall dimensions of Magnifier . . 100 × 80 × 35 mm
Weight of Magnifier 130 g

JIIIII-16 PERFORATION PITCH INSPECTION RULE

The JIIIII-16 Rule (Fig. 109) is adapted to measure the average perforation pitch and the average percentage of shrinkage in 21 picture-frame lengths of 16-mm film.

Fig. 109. JIIIII-16 Perforation Pitch Inspection Rule

The design of the Rule enables to define shrinkage ranging from 0.5 to 1.5 %, and to inspect films of perforation pitch exceeding 7.62 mm.

Reading: of fixed scale 0.004 mm; of moving scale 0.05 %.

The Rule is of the portable pocket type.

Overall dimensions of Rule 150 × 50 × 20 mm
Weight of Rule 350 g

JIIIII-35 PERFORATION PITCH INSPECTION RULE

The JIIIII-35 Rule (Fig. 110) is adapted to measure the average perforation pitch and the average percentage of shrinkage in 3—8 picture-frame lengths of 35-mm film.


Fig. 110. JIIIII-35 Perforation Pitch Inspection Rule

The Rule is a metal plate with a film bed, pins to secure film, and scale for measuring pitch dimensions and shrinkage percentage.

The instrument is of the portable pocket type.

Overall dimensions of Rule 180 × 50 × 3 mm
Weight of Rule 50 g

ПКП-2 PERFORATION INSPECTION PROJECTOR

The Perforation Inspection Projector (Fig. 111) is portable, bench-type, and adapted to define degree of wear in 35-mm film perforation tracks.


Fig. 111. ПКП-2 Perforation Inspection Projector

The instrument projects the perforation aperture, film edge and sound track on a lens-screen.
A special reticule on the projector screen allows to define the degree of wear.

SPECIFICATIONS

- Magnifying power of Projector . . . 7.6 — 8.6 X
- Power supply 127 or 220 V A. C. line
- Illumination source 4 V, 3 W lamp
- Overall dimensions of Projector . . . 280 X 240 X 100 mm
- Weight of Projector 4.5 kg

РМФ-3 SOUND TRACK MEASURING INSTRUMENT

The РМФ-3 bench-type Measuring Instrument (Fig. 112) measures the geometric dimensions of the sound track and its position in relation to film edge.


Fig. 112. РМФ-3 Sound Track Measuring Instrument

The instrument is adapted to control:
distance from sound track axis line to film edge;
width of band to be printed;
width of zero mark lines;
width of negative film sound track.

SPECIFICATIONS

- Magnifying power of measuring instrument 20 X
- Power supply 127 or 220 V, A. C. line
- Light source 4 V, 3 W lamp
- Overall dimensions of instrument . . 300 X 245 X 110 mm
- Weight of instrument 5 kg

MH-3 INSPECTION KIT

The MH-3 Inspection Kit (Fig. 113) is designed to facilitate inspection and control of motion picture projection installations, and projection apparatus in repair shops.


Fig. 113. MH-3 Inspection Kit

The kit allows for inspection of: the optical light system; radial or end run out of sprockets; axis play; stability of picture frame in the film channel; equalization of screen illumination from two projectors; film tensioning in film channel and film tensioning on take-up frictional device; degree of perforation wear; average perforation pitch; average percentage of film shrinkage; quality of sound and image; dimensions of film channel; position of intermittent sprocket; film gate position.

The complete kit comprises:

- a) for 35-mm sound-on-film projectors:
 - device for installation of optics;
 - instrument for defining sprocket run out;
 - gauges for film channel measurements;
 - masks for light control measurements;
 - dynamometer for measuring film tensioning in film channel;

dynamometer for measuring film tensioning in take-up frictional device; .
sound test-film;
image test-film and film loop of 100 % fitness;
125-mm vernier caliper;
tape or collapsible meter;
magnifier with mount;
electric torch with battery;

- b) for inspection of 16 and 35-mm film perforations:
 - magnifier for perforation inspection;
 - 16-mm film perforation pitch inspection rule;
 - 35-mm film perforation pitch inspection rule.The kit is provided with carrying case.

Overall dimensions of case 400 X 255 X 132 mm
Weight of kit 7.5 kg

ДФЭ-2 PHOTOELECTRICAL DENSITOMETER

The ДФЭ-2 Photoelectrical Densitometer (Fig. 114) measures photographic densities in transmitted light. The instrument includes the following assemblies:

- Ferroresonance voltage stabilizer;
- optic system with lighting device;
- amplifier;
- mirror dial galvanometer;
- detachable carriage.


Fig. 114 ДФЭ-2 Photoelectrical Densitometer

The voltage stabilizer provides stability within 1% of the output voltage with line voltage variations ranging from 100 to 130 V.

The scheme of densitometer includes:

- СІВ-3 photoelectric cell;
- 6Ф6 electronic tube;
- 6 V, 15 W exciter lamp.

Measurements are facilitated by the detachable carriage construction. One measurement may be made within 2—3 sec.

The Densitometer makes possible special measurements of light intensity.

Easy access to parts of Densitometer is provided by removable lid of the instrument box.

SPECIFICATIONS

Power supply — from 120 V, 50 cycle A. C. line.
Reading of optical density values is achieved by means of galvanometer pointer deflections.
Galvanometer scale enables direct reading of optical diffusion density.
Optical density measurements from zero to 3.0 are made with an accuracy of:

- ± 0.01 with 0.0 to 1.0 density values;
- ± 0.02 with 1.0 to 2.0 density values;
- ± 0.03 with 2.0 to 3.0 density values.

Overall dimensions of Densitometer 570 X 320 X 360 mm
Weight of Densitometer 28 kg

ЭПП-4 UNIVERSAL EXPOSURE METER

The ЭПП-4 Exposure Meter (Fig. 115) is designed to determine exposure time when filming black-and-white and colour film, or photographing with day or artificial light.


Fig. 115. ЭПП-4 Universal Exposure Meter

The instrument is portable, pocket-type, supplied with special handy carrying case.

The instrument consists of the following parts:


disk-shaped selenium photoelectric cell; iris diaphragm for adjusting magnitude of the light beam; two attachments: for light measurements and measurements of average and high brilliancies respectively; metering instrument with arrangement for locking of the pointer in the deflected position; calculator with three scales bearing values of exposure-determining factors.

The Exposure Meter consists of two parts — upper and lower, connected with each other by an axis which allows the upper part to move up to an angle of 300° in relation to the lower part.

SPECIFICATIONS

Range of illumination measurings	25—100,000 lux
Range of brilliancy	30—100,000 apostilb
Viewing angle during screen illumination measurings	170°
Viewing angle when measuring high and average brilliancy by use of attachment	50°
Measuring instrument	micro-ammeter of 1.5—10 ⁻⁴ response per mm of scale division and 650 ohm frame resistance
Overall dimensions of Exposure Meter	125 × 60 × 35 mm
Weight of Exposure Meter	300 g

IIP-60 APPARATUS FOR DIFFUSED LIGHT ILLUMINATION

The IIP-60 Apparatus (Fig. 116) provides diffused light illumination for normal or newsreel filming.

The apparatus is used to illuminate backgrounds, large groups in foregrounds, sets and actors in action.


Fig. 116. IIP-60 Apparatus for Diffused Light Illumination .

The apparatus is portable, safe in operation, simple in design and is used effectively in studio and outdoor sets.

The apparatus is fitted on a tripod or suspended by special attachments.

The apparatus is equipped with a protective wire grid which is quickly and easily attached to it.


A special cable is supplied with the apparatus. Its length is 10 metres.

A switch is provided for disconnecting the electric supply line.

A set of shutters (Fig. 117) may be supplied with the apparatus on special order.

SPECIFICATIONS

Light source	incandescent 110 V. 5,000 W lamp, or similar 3,000 W lamp
Diameter of reflector	600 mm
Maximum light intensity (with 5,000 W lamp)	not less than 30,000 candles
Angle of diffusion	160°
Reflection factor	not less than 0.75
Overall dimensions of apparatus (without tripod)	250 X 300 X 500 mm
Weight of apparatus	14 kg

-Fig. 117. IIP-60 Apparatus for Diffused Light Illumination, with shutters

PJ-5 APPARATUS FOR DIFFUSED LIGHT ILLUMINATION

The PJ-5 Apparatus for Diffused Light Illumination (Fig. 118 and 119) is designed for general illumination of large sets indoors and outdoors. The Apparatus is suitable for black-and-white as well as colour filming.


Fig. 118. PJ-5 Apparatus for Diffused Light Illumination Fig. 119. PJ-5 Apparatus for Diffused Light Illumination, rear view

The Apparatus includes the following assemblies:
 apparatus body with arc lamp mechanisms and arcuated support;
 reflector installed in frame;
 diffuser with frame;
 electric motors with reduction gearing;

voltage divider;
ballast rheostat;
10-metre length of connecting cord.


Fig. 120. PII-5 Apparatus for Diffused Light Illumination, rear view with cover removed

Automatic striking and constant operation of both arc lamps are maintained with two electric motors with reduction gearing (Fig. 120).

Styling of body and structural features provide maximum operating convenience.

Heat and sound insulation is ensured by asbestos sheathing of inner surface of body casing.

The apparatus is designed for use with tripod but may be operated set on floor or suspended. Tilting may be effected to an angle $\pm 60^\circ$.

Two series-connected arc lamps and reflector with light diffuser are supplied with the apparatus.

SPECIFICATIONS

- Power supply — from 110—115 V D. C. line.
- Conditions of arc lamp operation:
 - Voltage at each arc 34—36 V
 - Current 40—42 A
 - Maximum light intensity of spotlight 15,000 candles
 - Spotlight diffusion angle 130°
- The arc lamp is designed to operate with carbons type "8-40"
 - Positive carbon with 8-mm diameter and 300-mm length burns at rate of 110 mm per hour
 - Negative carbon with 8-mm diameter and 220-mm length burns at rate of 65 mm per hour
- Distance of luminous centre of apparatus (less tripod) from floor level 520 mm
- Overall dimensions of apparatus 450 × 560 × 925 mm
- Weight of apparatus 54 kg

KIIIJ-15 ARC LAMP SPOTLIGHT

The KIIIJ-15 Arc Lamp Spotlight (Fig. 121) is designed for illumination of portrait subjects, defining details of filmed objects and producing exaggerated light contrasts during synchronous colour filming, etc.


Fig. 121. KIIIJ-15 Arc Lamp Spotlight

The spotlight apparatus includes the following assemblies: spotlight with 150-mm diameter echelon lens; ballast rheostat; tripod; 10-metre length of cable for connection to mains; shutter and set of snouts, supplied for every five devices; set of spare parts.

Design and workmanship of tripod provide for smooth lifting and lowering of the extension rod, and easy, noiseless shifting of the spotlight.

Handling of spotlight, setting of lens into mount are effected quickly and easily.

The arcuated support allows endless horizontal motion over 360°; vertical tilting is possible up to 180°. Locking devices secure the apparatus in any desired position.

SPECIFICATIONS

Power supply — from 110—115 V D. C. line.
Arc Lamp is designed to operate with:
type "8-40" positive carbon of 8-mm diameter and 220—330 mm length;
type "8-60" negative carbon of 7-mm diameter and 130-mm length.
Conditions of Arc Lamp operation:
voltage 30 V
current 30 A
Maximum light intensity 300,000 ± 10 % candles
Diffusion angle 7°
Overall dimensions of apparatus:
without tripod 300 × 400 × 470 mm
with tripod:
minimum height 1,400 mm
maximum height 2,100 mm
Weight of spotlight 13.5 kg
Weight of apparatus with tripod,
rheostat and shutters 39 kg

KПД-25 ARC LAMP SPOTLIGHT

The KПД-25 Arc Lamp Spotlight (Fig. 122) is designed for illumination for synchronous colour filming.


Fig. 122. KПД-25 Arc Lamp Spotlight

- The Spotlight includes the following main assemblies (Fig. 123):
1. body with mount and 250-mm diameter echelon lens;
 - supporting device for installation and shifting of lamp along its optical axis;
 - apparatus control instruments;
 - arcuated support for mounting of apparatus;
 - semi-automatic intensity arc lamp with instantaneous arc lighting;
 - ballast rheostat with cable;
 - shutter and set of snouts (set supplied with every five devices);

10-metre length of cable with contactor for connection with ballast rheostat and switchboards;
 set of spare parts and accessories.
 - Handling of spotlight, setting of lens in mount is effected quickly and easily.


Fig. 123. KIИ-25 Arc Lamp Spotlight, complete set

The arcuated support is designed to allow endless horizontal motion of 360°; vertical tilt is possible to an angle of 180°. Locking devices secure the apparatus in any desired position.

SPECIFICATIONS

Power supply — from 105, 115 and 125 V D. C. line.
 The Arc Lamp is designed to operate with:
 "11-75" positive carbon of 11-mm diameter and 400-mm length;
 "Extra-K" negative carbon of 9-mm diameter and 165-mm length.
 Conditions of Arc Lamp operation:
 voltage 55 V
 current 75 A
 Maximum light intensity 1,000,000 ± 20% candles
 Diffusion angle 9°
 Overall dimensions of apparatus
 (without handles) 660 × 590 × 800 mm
 Weight of apparatus 41 kg
 Weight of rheostat with cable 19 kg

KIИ-M ARC LAMP SPOTLIGHT

The KIИ-M Arc Lamp Spotlight (Fig. 124) is used for synchronous black-and-white and colour filming indoors and outdoors.

The apparatus is adapted for illumination of actors and background settings, as well as distant sets.

The apparatus includes the following assemblies:
 device consisting of arcuated support and cylindrical body with 500-mm diameter lens;


Fig. 124. KIИ-M Arc Lamp Spotlight

semi-automatic intensity arc lamp with instantaneous arc lighting;
ballast rheostat with cable;
tripod;
shutter and snouts.

SPECIFICATIONS

Power supply — from 105, 115 or 125 V D. C. line;
Power consumption — 16 to 19 kW (depending on voltage);
Mechanism of Arc Lamp driven by D. C. motor.
The Arc operates with high-intensity "ПМТ" or "ПЯ" carbons. The positive carbon with 16-mm diameter and 550-mm length. The negative carbon with 11-mm diameter and 220-mm length.
Conditions of Arc Lamp operation:
voltage 68—70 V
current 150 A
Continuous operation of Arc Lamp must not exceed 30 minutes.
Maximum light intensity 5,200,000 candles
Diffusion angle 8°
Overall dimensions of apparatus 1,655 X 800 X 2,260 mm
Weight of apparatus 95 kg

KILJI-25 INCANDESCENT LAMP SPOTLIGHT

The KILJI-25 Spotlight (Fig. 125) is used during indoor and outdoor filming for illumination of actors and sets as well as theatre and club stages.


Fig. 125. KILJI-25 Incandescent Lamp Spotlight

Fig. 126. KILJI-25 Incandescent Lamp Spotlight, with stand

The apparatus consists of a cylindrical body with a 250-mm diameter lens, a focusing device for a 2,000 W, 110 V incandescent lamp, reflector, arcuated support and 25 A switch.

A stand (Fig. 126) and 15 metres of special $2 \times 4 \text{ mm}^2$ cable are supplied with the apparatus.

SPECIFICATIONS

Maximum light intensity¹ 279,000 candles
Diffusion angle¹ 15°
Overall dimensions of apparatus $610 \times 480 \times 800 \text{ mm}$
Weight of apparatus 24 kg

¹ The above characteristics are true when a 2,000 W, 110 V incandescent lamp of 23 lm/W light efficiency is used with the apparatus.

KILJI-35 INCANDESCENT LAMP SPOTLIGHT

The KILJI-35 Spotlight (Fig. 127) is used during indoor and outdoor filming for illumination of actors and sets, as well as theatre and club stages.


Fig. 127. KILJI-35 Incandescent Lamp Spotlight

Fig. 128. KILJI-35 Incandescent Lamp Spotlight with stand

The apparatus includes: cylinder with 350-mm diameter lens; focusing device for a spotlight incandescent 5,000 W, 110 V lamp; reflector; arcuated support and a 60 A switch.

A stand (Fig. 128) 15 metres of special $2 \times 16 \text{ mm}^2$ cable and a set of spare parts and accessories are supplied with the apparatus.

SPECIFICATIONS

Maximum light intensity¹ 500,000 candles
Diffusion angle¹ 16°
Overall dimension of apparatus 690 X 570 X 930 mm
Weight of spotlight 34 kg

¹ The above characteristics are true when a 5,000 W, 110 V incandescent lamp of 26.5 lm/W efficiency is used with the apparatus.

SPECIFICATIONS

Maximum light intensity¹ 1,390,000 candles
Diffusion angle¹ 17°
Overall dimensions of apparatus 845 X 720 X 1,130 mm
Weight of apparatus 65 kg

¹ The above characteristics are true when a 10,000 W, 110 V incandescent lamp of 27.8 lm/W light efficiency is used.

KIIL-50 INCANDESCENT LAMP SPOTLIGHT

The KIIL-50 Spotlight (Fig. 129) is for indoor and outdoor illumination of actors and sets as well as theatre stages.


Fig. 129. KIIL-50 Incandescent Lamp Spotlight

Fig. 130. KIIL-50 Incandescent Lamp Spotlight, with stand

The Spotlight includes the following assemblies:
cylinder with 500-mm diameter lens; focusing device for incandescent 5,000 or 10,000 W, 110 V spotlight lamp; reflector; arcuated support and special throw-over switch for starter resistance which reduces starting current when lamp is switched on.
20 metres of special 2×25 mm² cable, stand (Fig. 130), a set of spare parts and accessories are supplied.

THE VSESOJUZNOJE EXPORTNO-IMPORTNOJE OBJEDINENIJE

"STANKOIMPORT"

EXPORTS AND IMPORTS:

Machine Tools
Wood-Working Machinery
Metal-Working Machinery (Presses, Hammers, Shears, Cold Roll Forming Machines, Punching Machines)
Rolling Mills (imports)
Measuring Instruments and Tools (for Metal Industry)
Testing Machines and Apparatus (for metals)
Optical Instruments and Equipment
Portable Electric and Pneumatic Tools (for Metal and Wood-Working)
Metal and Wood Cutting Tools
Mechanic's Tools
Lathe and Drill Chucks
Sintered Carbide and Hard-Alloy Products
Abrasives Products
Ball and Roller Bearings
Microscopes of all types
Motion Picture Equipment and Accessories
Geodetic Instruments and Equipment
Photographic Cameras, Binoculars, Magnifiers, Lenses
Crude Optical Glass and Blanks.

All inquiries and correspondence to be forwarded to:
Vsesojuznoje Exportno-Importnoje Objedinenije
"Stankoimport"
32/34, Smolenskaja-Sennaja Pl., Moscow, USSR
For cables: Stankoimport Moscow
Phone: T 4-21-32

Design and specifications of the equipment illustrated herein are subject to change without notice.

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APPLICATION

Type ВВ-15 air-circuit switch for rated voltage of 13.8 kV, rated current of 500 A and rated power of disconnection of 2000 MVA are designed for indoor operation.

DESCRIPTION

All three poles of the switch are mounted on a common bogie and controlled from a common driving device. The bogie has rotary type rollers permitting to move the switch on longitudinal and cross directions. The bogie frame forms a cupboard in which the pneumatic drive, damper, valves, etc are located.

The reservoir for a pressure air (400 l amount) is also built in the bogie frame. The rated air pressure is 20 atm excessive.

The valves and movable contacts of the switch are put in motion by a common shaft (common for all the poles) controlled by the pneumatic drive.

- The following items are fastened to the shaft:
 - a) levers for the control of the movable contacts (plates) coupled with the latter by means of insulating rods made of laminated insulation;
 - b) cams for the mechanical control of the blowing valves of three poles.
- At switching in and out the shaft is turned for 80° angle opening the blowing valves supplying the blowing out chambers (separately for each pole of the switch) with the pressure air. The air supply is accomplished at switching out as well as at switching in but in a less quantity. The blowing pipes are made of porcelain. The blowing out chambers are assembled from rectangular plates of an insulating material and they have inside fibre partitions located in the chamber in the form of a fan.
- The pressure air flow is directed across the arc. At first the arc-blowing out contacts of the

switch located outside the chambers are closed at switching in. At switching out the sequence of the contact opening is inverse. The first contacts together with the chambers are mounted at the porcelain supporting insulators. A damper is located at each of three blowing out chambers. It is necessary for accomplishing the commands "switch in" or "switch out" to put in motion the pneumatic drive of the switch. The starting valves—the valve for switching in and the valve for switching out — are designed for this aim.

The valves for switching in and out have mechanical picking up an impulse preventing from interrupting the operation of switching in and out.

The command push buttons of the local pneumatic control are located together with the switch in the distributing cupboard being located near the switch. The following items are also located in the cupboard:

- a) a locking flap and non-return valve by means of which the switch is connected with a supply main of the pressure air;
 - b) a filter for the pressure air;
 - c) two electrocontact manometer (controlling the command circuits of switching in and out in dependence on the pressure in the reservoir) each of which is equipped with an intermediate relay, an electric heater (switched in at a temperature lower +5°) accomplishes a fine operation of the electrocontact manometers at a low temperature.
- A provision is made for the adjustment of the switch by means of turning its shaft by a handle (crew) inserted into the coupling mounted on the shaft out.
- The switch weighs approximately 3,000 kg. The distributing cupboard weighs approximately 200 kg.

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
Fig. 1. Appearance of the control case with the mercury-arc rectifier:
 1 — enclosure; 2 — supporting insulators; 3 — voltmeter; 4 — ammeter of main current; 5 — ammeter of field current; 6 — rotary switch; 7 — automatic switch push-button; 8 — fan rotary switch; 9 — rectifier; 10 — ignition transformer; 11 — excitation transformer; 12 — choke of co-ordination; 13 — automatic switch; 14 — relay; 15 — excitation resistance; 16 — relay; 17 — protector; 18 — air-circuit relay

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TYPE PM300 MERCURY-ARC RECTIFIER

APPLICATION

Type PM300 metal non-pump mercury-arc rectifier is designed for supplying contact circuits of underground electric locomotive running in at mounting the rectifier:

- a) outside the mine in a closed premise at an ambient temperature of -5° up to $+35^{\circ}$ C;
- b) inside the mine in an underground compartment at an ambient temperature of $+2^{\circ}$ up to $+35^{\circ}$ C.

TECHNICAL DATA

Rated current, A	300
Voltage, V	275
Field current, A	6-8
Overloads:	
375A (25%) during 15 minutes with two-hour interval.	
450A (50%) during 2 minutes with one-hour interval.	
600A (100%) during 10 seconds.	

DESCRIPTION

Type PM300 mercury-arc rectifier is a complete set with a control case and apparatus for excitation and ignition.

The following items are mounted in this control case:
 mercury-arc rectifier,
 apparatus for excitation, ignition and cooling,
 protections of the mercury-arc rectifier, measuring instruments.

Type PM300 mercury-arc rectifier consists of a metal vacuum housing with cooling ribs both on cylindrical and taper parts of the housing.

A convex cover is welded to the housing top. Three insulated vacuum tight terminals of the main anodes and a scavenging pipe with a protective cap are all located at the convex cover. The terminals of two exciting anodes are located at the housing bottom. The anode heads made of graphite are screwed on current-carrying bars of the terminals.

The main anodes are screened with glands. A deionization filter is designed for the protection of a non-operating anode from charged particles diffusing from the arc to the nearest anode.

The cathode in the form of a metal enameled cup with a terminal is welded to the housing bottom. The current-carrying bar and bar of the

molybdenum igniter are mounted in the glass of the cathode terminal. The cathode is filled with mercury. The taper reflector protecting the main anodes from the flow of mercury vapour is located over the cathode.

A vacuum housing is strengthened to the cylindrical base inside of which the motor with a fan for cooling the rectifier is located. The cylindrical part of the rectifier housing is closed outside with a metal jacket for the direction of air flow.

The rectifier housing is filled with argon under small pressure for reducing overvoltages at starting the rectifier at low temperatures.

The rectifier is mounted in the control case by means of three porcelain insulators as at operation its housing relative to the grounded case is under voltage of 300 V.

OPERATION OF INSTALLATION CIRCUIT

At switching in the rotary switch П₂ 380 V a.c. voltage simultaneously supplies the excitation transformer ТБ and the ignition transformer ТЗ. Short-timed current of 400-450 A runs through the molybdenum igniter (immersed in mercury). A spark appears at this current in the place of the connection of the igniter with the mercury.

This spark excites discharge ignition and current runs through the exciting anodes.

The choke ДР maintains current continuity through the exciting anodes. The relay Р is switched in at running the field current and due to it the transformer ТЗ gives a voltage. Bimetallic relays ТРА are designed for the protection of the igniter terminal. The rectifier is cooled by the fan the motor of which is switched in, by the rotary switch П₁.

The rectifier protection from overloads, short circuits and reverse ignitions is accomplished by an automatic switch with manual switching in. The rectifier cannot operate without cooling.

The appearance of the control case together with the fan the motor of which is switched in, is given in fig. 1.

The principle diagram of the installation is given in fig. 2.

OVERALL DIMENSIONS AND WEIGHTS

The overall dimensions of the case with the rectifier—990 × 1040 × 1720 mm.

The rectifier without case weighs 175 kg.

The weight of the rectifier with the case 400 kg.

The whole weight of the installation is approximately 450 kg.

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SERIES ЭТ-520

APPLICATION AND CLASSIFICATION

The Series ЭТ-520 instantaneous Current Overload Relays are used in protection circuits of electrical installations (electrical machines, transformers, transmission lines) against short-circuits and overloads.

As regards connection of the sensitive element and action of the response element on the circuit breaker, these Relays belong to the group of indirect secondary relays.

As regards type of contacts, the Series ЭТ-520 Relays are divided into 3 groups:

a) ЭТ-521 — with one normally open contact;

b) ЭТ-522 — with one normally closed contact;

c) ЭТ-523 — with one normally open contact and one normally closed contact.

Each group of relays is subdivided into types differing in the tripping current setting range. The upper limit of settings for each type of relays is indicated as the denominator of a fraction, the numerator designating the corresponding group of relays.

The Type ЭТ-523/1/I Relay differs from the other Series ЭТ-520 Relays in its higher ratio of the admissible continuous current to the tripping current.

PRINCIPLE OF OPERATION AND CONSTRUCTION

As regards the construction form of the magnetic system, the Series ЭТ-520 Relays belong to the group of electromagnetic relays with a rotating armature.

The poles of the fixed core carry two coils (windings). A Z-shaped steel armature is fixed on the axle. The flux produced by the current flowing through the coils links the armature, and the latter tends to rotate clockwise towards the poles of the core. This movement is opposed by a spiral spring, the inner end of which is attached to the armature axle, while the outer end is attached to the lever of pointer.

On the axle are also fixed silver bridge contacts insulated from the axle and rocking about it (self-aligning).

At a preset tripping current flowing through the coils of the relay, the electromagnetic torque exceeds the opposing torque of the spiral spring, and the armature rotates clockwise. The bridge contacts then close the normally open (or open

the normally closed) fixed contacts. When the current flowing through the coils drops to a preset value, the armature under the action of the spiral spring returns to its initial position.

Smooth variation of the relay tripping current is accomplished by a corresponding variation of the opposing force of spring by means of pointer attached to lever.

The ends of both coils of the relay (except Type ЭТ-523/1/I Relay) are brought out to terminals and can be interconnected in series or parallel by means of jumpers. This enables to vary the tripping current setting range of the Relay.

The divisions on the scale of the Relay are marked for the coils connected in series; when the coils are connected in parallel, the current setting range is doubled.

The Type ЭТ-523/1/I Relay, as regards construction of the mechanism, differs in no way from the other Relays in the Series ЭТ-520. This

VSESOJUZNOJE OBJEDINENIJE

«MACHINOEXPORT»

type of Relay differs only in that it has an incorporated intermediate saturable current transformer; this is to ensure a higher thermal stability. The secondary winding of this transformer is connected to the coils of the Relay that are themselves connected in series. The primary winding of the transformer has four tappings connected to terminals to enable to vary the tripping current setting range.

The Series ET-520 Relays are enclosed in dustproof rectangular plastic cases, consisting of a base and body; only the Type ET-523/1/D

Relay is available in a rectangular metal case. A glass window is provided in the front of the case.

The Relays are arranged for mounting on a vertical plane and are available for either front or rear wiring connection. For front connection, the Relays are furnished with a set of current-carrying strips and screws, also with strips for fixing the Relay on a panel. For rear connection, the Relays are furnished with a set of extension screws with nuts and washers.

The Relays are arranged for sealing on site.

SPECIFICATIONS

The essential specifications of the various Types of Relays are given in Tables 1 and 2.

Table 1

TYPE	Setting range, A	Coils connected in series		Coils connected in parallel		Contacts		
		Tripping current, A	Thermal stability		Tripping current, A		Thermal stability	
			Continuous, A	1 second, A			Continuous, A	1 second, A
ET-521/0.2 ET-522/0.2 ET-523/0.2	0.05-0.2	0.05-0.1	0.3	12	0.1-0.2	0.6	24	1 NO 1 NC 1 NO; 1 NC
ET-521/0.6 ET-522/0.6 ET-523/0.6	0.15-0.6	0.15-0.3	1	45	0.3-0.6	2	90	1 NO 1 NC 1 NO; 1 NC
ET-521/2 ET-522/2 ET-523/2	0.5-2	0.5-1	4	100	1-2	8	200	1 NO 1 NC 1 NO; 1 NC
ET-521/6 ET-522/6 ET-523/6	1.5-6	1.5-3	10	300	3-6	20	600	1 NO 1 NC 1 NO; 1 NC
ET-521/10 ET-522/10 ET-523/10	2.5-10	2.5-5	10	300	5-10	20	600	1 NO 1 NC 1 NO; 1 NC
ET-521/20 ET-522/20 ET-523/20	5-20	5-10	15	300	10-20	30	600	1 NO 1 NC 1 NO; 1 NC
ET-521/50 ET-522/50 ET-523/50	12.5-50	12.5-25	20	450	25-50	40	900	1 NO 1 NC 1 NO; 1 NC
ET-521/100 ET-522/100 ET-523/100	25-100	25-50	20	450	50-100	40	900	1 NO 1 NC 1 NO; 1 NC
ET-521/200	50-200	50-100	20	450	100-200	40	900	1 NO

Table 2

ESSENTIAL SPECIFICATIONS OF TYPE ET-523/1/D RELAYS

Connection of Relay (terminal numbers)	Setting range A	Impedance, Ohms, at the following currents:			Thermal stability	
		Minimum setting current	5 A	30 A	Continuous, A	1 second, A
2-8	0.15-0.3	12	3.5	0.9	6	300
2-6	0.3-0.6	3	0.8	0.25		
2-4	0.5-1.0	1	0.3	0.08		

Reset factor 0.85 (except Type ET-521/200 Relay having a reset factor of 0.7).

Note: The above values refer to operation of the Relay on A. C. When the Relay operates on D. C. the reset factor is 10-15% less.

Tripping time is 0.15 sec at a current 1.2 times the tripping current, or 0.03 sec at double the tripping current.

The burden at a current equal to the minimum tripping current and with the coils connected in series is given in Table 3. The burden of Type

ET-523/1/D Relay is determined on the basis of the data given in Table 2.

The rupturing capacity of the contacts in a D. C. circuit with an inductive load (Time constant $T=5 \cdot 10^{-3}$ sec) is 50 W, and in a A. C. circuit 250 VA at a voltage up to 250 V and a current not exceeding 2 A.

The dielectric strength of the coils and contacts to the core is 2000 V, 50 c.p.s. for 1 minute.

Weight of Relay 1.1 kg; weight of Type ET-523/1/D Relay 3.5 kg.

Table 3

Type of Relay	ET-520/0.2; ET-520/0.6; ET-520/2; ET-520/6; ET-520/10	ET-520/20	ET-520/50	ET-520/100	ET-520/200
Burden (approx.), VA	0.1	0.2	0.8	2	8

ORDERING DATA

When ordering please state:

1. Description of Relay
2. Type of Relay
3. Front or rear wiring connection
4. Quantity of Relays required.

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UNCLASSIFIED

2301

POTENTIAL TRANSFORMERS APPLICATION

Single-phase and three-phase Potential transformers are designed for AC power plants having a rated (nominal) voltage up to 220000 V and a rated (nominal) frequency of 50 cps. They are

used for voltage and power measurement, as well as for supplying relay protection, automatic control and signalling circuits.

TECHNICAL DATA

Potential Transformers, manufactured to U.S.S.R. Standards GOST 1983-43, are designed for operation under the following rated (nominal) conditions:

a) Maximum altitude above sea level—1000 m for transformers up to and including 110000 V, and not more than 500 m for 220000 V transformers.

b) Maximum ambient temperature natural variation +40°C.

Two kinds of Potential Transformers are available:

a) Dry—i. e. self-cooled.
b) Oil immersed—i. e. enclosed in tanks filled with oil.

Dry Potential Transformers are available in the following types:

a) HOC-0.5; HOCK-3; HOCK-6 (single-phase)
b) HTC-0.5 (three-phase)

Oil immersed Potential Transformers are available in the following types:

a) HOM-6; HOM-10; HOM-15; HOM-35 (single-phase),

b) HTMK-6-48; HTMK-10 (three-phase),
c) HTMI-6, HTMI-10 (three-phase, five-core, three-winding),

d) HKΦ-110; HKΦ-220 (single-phase, cascade, three-winding, in porcelain casing).

The digit following the symbol designation denotes the highest rated (nominal) voltage for the given type of transformer, expressed in kilovolts.

Potential Transformers are capable of operating continuously at voltages 10% above (nominal) rated values.

Potential Transformers, types HOM-35, HKΦ-110, and HKΦ-220, are designed for outdoor installation, while the other types are designed for indoor installation.

Potential Transformers, types HOCK-3 and HOCK-6, are for use in switchgear boxes filled with insulating compound.

Voltage combinations, rated (nominal) and maximum power, winding connections and test voltages for Potential Transformers are indicated in Table 1.

Table 1

Type of transformer	Rated (nominal) voltage, V		Rated power in accuracy class, VA			Maximum Power, VA	Winding connections	Test voltages, kV	
	HT	LT	0.5	1	3			HT	LT
HOC-0.5	380	100	25	40	100	200	1/1-12	6	2
HOC-0.5	100	100	25	40	100	200	1/1-12	6	2
HOCK-3	3000	100	30	50	120	240	1/1-12	12	2
HOCK-6	6000	100	50	80	200	400	1/1-12	12	2
HOM-6	2100	100	30	50	120	240	1/1-12	24	2
HOM-6	3000	100	30	50	120	240	1/1-12	24	2

VSESOUZNOJE OBJEDINENIJE

«MACHINEEXPORT»

Type of transformer	Rated (nominal) voltage, V		Rated power in accuracy class, VA			Maximum Power, VA	Winding connections	Test voltages, kV	
	HT	LT	0.5	1	3			HT	LT
HOM-6	6000	100	50	80	200	400	1/1-12	32	2
HOM-10	10000	100	80	150	320	640	1/1-12	42	2
HOM-15	15000	100	80	150	320	640	1/1-12	55	2
HOM-35	35000	100	150	250	600	1200	1/1-12	95	2
HTC-0.5	380	100	50	80	200	400	Y/Yo-12	6	2
HTC-0.5	500	100	50	80	200	400	Y/Yo-12	6	2
HTMK-6-48	3000	100	50	80	200	400	Y/Yo-12	24	2
HTMK-6-48	6000	100	80	150	320	640	Y/Yo-12	32	2
HTMK-10	10000	100	120	200	480	960	Y/Yo-12	42	2
HTMH-6	3000	100-100:3	50	80	200	400	Yo/Yo-12	24	2
HTMH-6	6000	100-100:3	80	150	320	640	Yo/Yo-12	32	2
HTMH-10	10000	100-100:3	120	200	480	960	Yo/Yo-12	42	2
HKΦ-110	110000:V $\bar{3}$	100:V $\bar{3}$ -100	—	500	1000	2000	1/1-12	230	2
HKΦ-220	220000:V $\bar{3}$	100:V $\bar{3}$ -100	—	500	1000	2000	1/1-12	460	2

Notes: 1. The symbols HT and LT denote high tension and low tension, respectively.
 2. Potential Transformers, types HKΦ-110 and HKΦ-220, are available in Accuracy Classes 1 or 3, only. Other Potential Transformers can be supplied (in accordance with U.S.S.R. Standards GOST 1983-43) in the Accuracy Classes 0.5, 1 and 3 as well as in the Classes 1 and 3.

DESCRIPTION

DRY (SELF-COOLED) TRANSFORMERS

The windings are of the layer type, round winding-wire wound on cartridges of electrical cardboard. The windings are asphalt varnish impregnated. The magnetic circuit is built-up of 0.5 mm thick transformer steel stampings. The stampings are insulated with paper or water glass, and for type HOC-0.5 transformers, by annealing.

The transformers are connected to the line by means of terminals, arranged on insulating terminal boards. Transformers, HOCK-3 and HOCK-6, do not have terminal boards. The ends of windings on these transformers are brought out as flexible insulated wires. The transformers are grounded (earthed) by means of a special bolt, fastened to the transformer frame.

OIL-COOLED TRANSFORMERS (EXCEPT CASCADE TYPES HKΦ-110 and HKΦ-220)

The windings of oil-immersed Potential Transformers are of the layer type. They are round winding-wire wound on cartridges of electrical cardboard. The high tension winding is divided into two sections. The high tension winding of type HOM-35 transformer is divided into four sections and wound on a bakelite cylinder.

The transformer tanks are sheet steel welded. The tanks are circular in form, except transformers type HOM-15 which have an oval tank. The winding terminals are connected to porcelain bushings located on the tank cover. In addition to the bushings the following fittings are located on the tank cover: oil filling and "breathing" plug, lifting shackles and name plate. In the lower part of the tank (except on types HOM-6 and HOM-10) a plug is provided for draining and taking samples of the oil. The transformers are grounded (earthed) by means of special bolts fastened to the tank and suitably marked by a symbol.

All high tension windings have electrostatic shields for protection against overvoltage. The magnetic circuit is built-up of 0.5 mm thick transformer steel stampings. The stampings are insulated with paper or water glass.

CASCADE TRANSFORMERS

Cascade transformers consist of two (HKΦ-110) or four (HKΦ-220) elements connected in series and grouped in units of two elements each.

Each element comprises a two-core three-winding transformer. The secondary windings, main and additional, are located on the lower core of the lower element having a minimum potential in relation to the ground (one terminal of the primary winding is solidly grounded).

The windings are of the layer type. The high tension winding is wound on a bakelite cylinder and has an electrostatic shield on top.

The magnetic circuit is built-up of 0.5 mm thick transformer steel stampings. The stampings are paper insulated.

Unit HKΦ-110 (of two elements) is located in a special porcelain casing filled with transformer oil and mounted on a carriage truck. The upper part of the unit is closed by an expansion tank.

During installation the line terminal (from the mains) is connected to the expansion tank, while the ground lead is connected to the carriage. The ends of the secondary winding are brought out through porcelain bushings located in the lower part of transformer.

When two units are placed one on the top of the other (for cascade transformer, HKΦ-220), they are connected by a housing. The expansion tank of the upper unit serves both units.

When cascade transformers are connected in a circuit they are connected in groups in the same way as ordinary single-phase three-winding Potential Transformers.

Overall dimensions of transformers are given on Fig. 18-30.

Weights of Potential Transformers are indicated in Table 2.

Table 2

Type of transformer	WEIGHTS	
	Total weight	Weight of oil
	kg	
HOC-0.5	8	—
HOCK-3	12.7	—
HOCK-6	13.5	—
HOM-6	23	4.7
HOM-10	36.2	7.3
HOM-15	81	23
HOM-35	248	75
HTC-0.5	20	—
HTMK-6-48	47.5	15
HTMK-10	110	27
HTMH-6	105	32
HTMH 10	190	70
HKΦ-110	1360	315
HKΦ-220	2650	600

FORMULATING THE ORDER

When ordering please indicate the following data.

1) Type of transformer

2) Rated voltage of high and low tension windings (ratio)

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**ELECTRIC RESISTANCE SPOT WELDING MACHINE
TYPE МТП-300**

APPLICATION

The type МТП-300 machine is used for electric resistance spot welding of low-carbon steel parts.

BRIEF CHARACTERISTICS

Welding on the type МТП-300 machine is effected by pressing together of the parts heated by electric current.

Switching the welding transformer on and off is effected by an ignitron asynchronous contactor controlled by two electronic timers to allow multi-impulse or single-impulse welding.

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- 1) smooth adjustment of pressure between the electrodes;
- 2) constant pressure independent of electrode wear;
- 3) additional travel of top electrode;
- 4) welding current adjustment by varying the welding transformer secondary voltage;
- 5) timing adjustment of the welding cycle operations.

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SPECIFICATIONS

Rating, kVA	300
Duty cycle ("On" time percentage)	20
Primary voltage, V	380
Secondary voltage, V	from 4.87 to 9.75
Secondary voltage adjustment steps	16
Effective throat depth, mm	500, adjustable to ±50 mm

Maximum thickness of low-carbon steel parts to be welded, mm	7+7
Height of lower contact part adjustable, mm	up to 160
Operating travel of upper electrode, mm	30
Additional travel of upper electrode, mm	90
Strokes of upper electrode per minute	up to 40
Maximum operating pressure, kg	3200
Compressed air supply pressure, rated, kg/cm ²	5.5
Free air consumption, m ³ /hour	25
Cooling water required, l/hour	1200
Overall dimensions of machine, mm:	
height	2565
width	800
depth	1182
Weight, kg	1100

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1. Base with upper and lower brackets.
2. Pressure mechanism.
3. Guide with upper contact part.
4. Lower contact part.
5. Flexible connections.
6. Compressed-air system.
7. Welding transformer with step changer.
8. Cooling system.
9. Asynchronous ignitron contactor.
10. Electric timers.
11. Pedal switch.

The operation of the type МТП-300 machine is controlled by a pedal switch conveniently located for the operator.

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**TYPE MTK-2
SPOT WELDING UNIT**

A type MTK-2 welding unit is designed for spot welding ferrous and non-ferrous materials and their alloys having thickness of 0.1+0.1 mm to 0.3+0.3 mm, and also lattices whose wires have diameters of 0.6 to 1 mm.

TECHNICAL DATA

Circuit voltage, V	220
Maximum charging current, A	14
Capacitor voltage, V	500
Capacitance of capacitors:	
a) minimum # F	25
b) maximum # F	500
Number of regulating steps for the transformation ratio	5
Force on electrodes, kg	2 to 15
Output:	
a) rated, welds/min	50
b) limits of regulation, welds/min	20 to 90
Stroke of electrodes:	
a) for automatic operation, mm	5
b) for non-automatic operation, mm	20
Travel of electrodes, mm	80

CONSTRUCTION

Energy stored in the capacitors is used for welding. The amount of energy stored is regulated by changing the capacitance of the capacitors within limits from 3.12 to 62.5 W/sec.

The capacitors discharge through the primary winding of a welding transformer.

The unit may operate automatically or non-automatically. A pedal connects the unit for automatic operation. While the pedal is pressed the welding cycle is repeated at the same speed. When pressing the pedal momentarily or pressing it as far as it will go, only one weld will be made.

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
**WELDING SUSPENDED MACHINE
WITH TONGS TYPE КТГ-75**

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TYPE МТПГ-75 WELDING SUSPENDED MACHINE WITH TONGS TYPE КТГ-75

APPLICATION

The type МТПГ-75 machine is used for electric resistance spot welding of low-carbon steel parts which can be fed to stationary spot machines.

BRIEF CHARACTERISTICS

Welding on the type МТПГ-75 machine is performed by pressing together the parts heated by an electric current passing through the parts.

Switching the welding transformer on and off is effected by an ignitron asynchronous contactor. Switching the contactor on is controlled by an electronic timer.

The construction of the machine provides:

1. Smooth adjustment of pressure between the electrodes.
2. Constant pressure independent of electrode wear.
3. Welding current adjustment by varying the welding transformer secondary voltage.
4. Timing adjustment of the welding cycle operations.

These features of the machine allow to select optimum duty conditions to ensure high quality welding.

SPECIFICATIONS OF THE MACHINE

Nominal rating, kVA	75
Duty factor, %	25
Primary voltage, V	380
Secondary voltage when secondary winding turns of the transformer are connected in series, V	from 10, 12 to 19
Secondary voltage when secondary winding turns of the transformer are connected in parallel, V	from 5.05 to 9.5
Secondary voltage adjustment steps	16
Thickness of low-carbon steel parts to be welded, mm	1.5+1.5
Ratio of air circuit pressure to hydraulic pressure	1:19
Weight, kg	370
Overall dimensions of the transformer with the suspensions, mm:	
height—1805/600; width—452/400; depth—1300/280.	
Overall dimensions of the timer, mm:	
height—300; width—302; depth—195.	

SPECIFICATIONS OF THE TONGS

	Type of tongs		
	КТГ-75-1	КТГ-75-2	КТГ-75-3
Effective gab, mm	42	125	140
Maximum electrode stroke, mm	25	30	25
Distance between the arms, mm	—	94	100
Maximum pressure between the electrodes, kg	275	200	250
Air pressure in mains, atm	3	3	4.5
Air consumption m ³ /hr	9	9	13.5
Quantity of cooling water used, l/hr	600	845	600
Weight, kg	14.5	12.5	9.0
Overall dimensions, mm:			
length	460	400	315
width	75	125	255
depth	380	325	212
Electrode motion	Rectilinear	Radial	
Number of strokes per minute at duty factor—25%	to 80	to 80	to 80

The machine comprises the following essential components:

- a suspension;
 - a pneumohydraulic pressure drive;
 - a multiple-core supply cable;
 - a high-voltage hose;
 - a compressed-air system;
 - a welding transformer with a step changer;
 - a cooling system;
 - an electronic timer;
 - an ignitron asynchronous contactor.
- The operation of the machine is controlled by a knob located on the tongs.

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AUTOMATIC ELECTRIC WELDING MACHINE TYPE TC-17-M-V

APPLICATION

Type TC-17-M-V automatic welding machine is designed for arc welding under a layer of welding flux of butt joints with beveled and straight edges, lap and fillet welds with and without turning the welded article as well as circular and longitudinal welds of shells.

The machine is universal. With the aid of interchangeable parts the automatic machine can be adjusted for the needed type of weld and the corresponding welding process.

The automatic machine performs welding while travelling on the article to be welded or along a light guiding track.

The welding current is supplied to the automatic machine from any welding transformer providing a no-load voltage not lower than 60 V or a d.c. welding generator may also be used.

TECHNICAL DATA

Intensity of welding current, A	from 200 to 1,000
Diameter of electrode wire, mm	from 1.6 to 5
Rated voltage of electric equipment truck, V	220 or 380
Speed of electrode wire feed, m/min	from 52 to 403
Welding speed, m/min	from 16 to 126
Capacity of flux hopper, liters	6.5
Weight of electrode wire on drum, kg	8

SHIPPING COMPLEMENT

Each automatic machine is supplied with electric wires and a set of interchangeable and spare parts. The set consists of: a drum for electrode wire, mouthpiece nozzles, caps, feed roller and a special wrench for releasing the pressure roller, a hopper for fillet welds and contact rollers.

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UNCLASSIFIED

FLEXIBLE HOSE AUTOMATIC WELDER FOR ARC WELDING UNDER A BLANKET OF FLUX

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APPLICATION

The АДШМ-500 Automatic Welder is used:

1. For rectilinear or curved fillet automatic arc welding, under a blanket of flux, on horizontal or inclined (up to 20°) surfaces.
2. For butt or fillet manual semi-automatic arc welding, under a blanket of flux, the welds being of any shape and length, on horizontal and inclined surfaces, also welding in not easily accessible places.

1. A self-propelled automatic welding head with magnetic suction.
 2. A manual semi-automatic welding head.
 3. Electrode wire feed with push-button control.
 4. Switchgear cabinet.
 5. Pneumatic flux feeder with flux sifting screen.
- The Automatic Welder is furnished with a D. C. or A. C. supply source.

The Automatic Welder consists of:

SPECIFICATIONS

Table 1

Type of Automatic Welder	АДШМ-500-1	АДШМ-500-2	АДШМ-500-3	АДШМ-500-4
Welding	Automatic; semi-automatic			
Welding current	D. C.			
Electrode wire dia., mm	1.6 to 2			
Welding current supply source, furnished with the semi-automatic welder*	ПС-300 М	ПС-500	ПС-300 М	ПС-500
Welding current adjustment range, A	180—340	180—600	180—340	180—600
Welding current adjustment	Smooth			
Welding current, A	At 100% "On" period			
	260	400	260	400
	At 65% "On" period			
	340	500	340	500
Rated operating voltage V	35	40	35	40
No-load voltage, V	50—76	60—90	50—76	60—90
Line voltage, V	280			
Electrode wire feed rate, in per minute	1.8—7			
Flux feed	Pneumatic			
Air pressure, atm	In flux hopper			
	From 1 to 2.5			
	in line			
	From 3 to 6			
Weight of electrode wire on reel, kg	8			
Capacity of flux hopper, kg	35			
Length of flux hose (max. distance between flux feeder and operator), m	15			

*The Type of the welding current supply source, to be furnished with the Welder, should be stated when ordering (For D. C. welding with a current up to 340 A, the Type ПС-300 М welding current supply source is recommended).

VSESOJUZNOJE OBJEDINENIJE

«MACHINEEXPORT»

UNCLASSIFIED

Table 2
OVERALL DIMENSIONS OF DIFFERENT PARTS OF THE AUTOMATIC WELDER

Description	Overall dimensions, mm		Weight, kg
	Length	Width	
Automatic head, Type ГСА-2-4	Length	290	15
	Width	216	
	Height	266	
Semi-automatic head, Type ГСП-3	Length	190	0.64
	Width	66	
	Height	120	
Electrode wire feeder	Length	400	13
	Width	345	
	Height	343	
Switchgear cabinet	Length	896	122
	Width	585	
	Height	705	
Welding current supply source, Type ПС-300М	Length	1200	600
	Width	755	
	Height	1180	
Welding current supply source, Type ПС-500	Length	1400	960
	Width	770	
	Height	1110	
Welding current supply source, Type ТСА-500	Length	950	500
	Width	818	
	Height	1250	

OPERATION OF THE AUTOMATIC WELDER

The Type АДШМ-500 Automatic Welder, when the self-propelled automatic welding head with magnetic suction is used, provides automatic fillet welding. The vertical wall of the weld should then be at least 160 mm high, while its distance from the edge of the horizontal sheet should not be less than 40 mm.

When the semi-automatic welding head is used, the Automatic Welder provides butt and fillet welding, the welding head being moved by hand along the weld.

The electrical connections of the Automatic Welder permit infinite (smooth) adjustment of welding speed, when the automatic welding head is used, and of the electrode wire feed rate, when the automatic or semi-automatic welding head is used. The electrode wire feed rate is not dependent on the voltage across the arc. The welding speed and the electrode wire feed rate are adjustable by control handles arranged on the switchgear cabinet.

Radio interference suppressors are fitted on the Automatic Welder.

The Type АДШМ-500 Automatic Welder is fitted with a pneumatic flux feed hopper of a capacity sufficient for operation during a period of several hours.

Compressed air forces the flux, through a rubber hose, into the automatic or semi-automatic head. Here the air, having completed its flux transporting function, separates from the flux and escapes through a special port, while the flux flows into the arc by gravity.

The flux feeder is mounted on a carriage. An air cleaner is mounted on the same carriage to remove moisture, and oil from the compressed air supplied.

The flux feeder is fitted with a vibrating screen for sifting the flux when filling the hopper.

The rate of feed of the flux is adjustable by means of an air pressure reducing valve fitted with a pressure gauge.

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**ELECTROMETALLIZING
SPRAYER**
Model ЭМ-3А


The ЭМ-3А hand-operated electrometallizing sprayer is designed for coating metals by method of atomization.

The sprayer is used for the restoring of worn parts of machinery, for the correction of defects in castings, and for various other purposes.

SPECIFICATIONS

Operating pressure of compressed air, atm. gauge	4.5—6.0
Compressed air consumption (at pressure of 6 atm gauge), m ³ /min	1.2
Maximum feed speed of wire, m/min	2.5
Regulation limits of wire speed	1—4
Diameter of wire used, mm	1—2
Permissible value of operating current, A	130
Power consumption, kW	3—5
Weight of sprayer (less flexible hoses and wires), kg	2.2
Sprayer capacity, kg/hr	
when atomizing a 1.5-mm diameter steel wire	3.5
when atomizing a 2.0-mm diameter zinc wire	5.0

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A portable automatic welding machine that during welding travels directly on the work or along a light straight-edge guide placed on the work, is usually called a Welding Tractor.

The TC-17-M Universal Welding Tractor (Fig.1) is designed for down-hand submerged under flux (shielded) automatic butt and angle (fillet) arc welding. Welding can be carried out by means of 1.6 to 5 mm dia. electrode wire at a current of 200 to 1000 A. The rate of travel of the automatic welding machine is adjustable in the range from 16 to 126 m per hour.

UNIVERSAL WELDING TRACTOR

A portable automatic welding machine that during welding travels directly on the work or along a light straight-edge guide placed on the work, is usually called a Welding Tractor.

The TC-17-M Universal Welding Tractor (Fig.1) is designed for down-hand submerged under flux (shielded) automatic butt and angle (fillet) arc welding. Welding can be carried out by means of 1.6 to 5 mm dia. electrode wire at a current of 200 to 1000 A. The rate of travel of the automatic welding machine is adjustable in the range from 16 to 126 m per hour.

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CONSTRUCTION OF THE TRACTOR

The design of the TC-17-M Welding Tractor permits butt welding of prepared or unprepared joints, angle (fillet) welding by means of a vertical or inclined electrode, as well as lap welding. The welds can be straight or circumferential. The minimum diameter of the circumferential weld, when welding inside the container, is 1200 mm. Special devices (furnished extra) allow the tractor also to be used for double arc welding of the above types of joints by means of split electrodes.

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The construction of the tractor (Fig.2) is very simple. The tractor is driven by a single induction motor 1 and comprises a feeder 2 and a travel drive 3. Both the feeder and drive together with the motor form a single block, that carries all the principal units of the Automatic Welding Tractor: the nozzle to supply the welding current to the electrode wire 4, the bracket with control panels 5, electrical measuring instruments 6, reel for the electrode wire 7, and flux hopper 8.

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Fig.3. Welding Tractor TC-17-M arranged for gap butt welding or for prepared joints (welding by means of a "copy" jig)

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Fig. 4. Welding Tractor TC-17-M arranged for angle (fillet) welding with inclined electrode (welding by means of a „copy“ jig)

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The versatility of the tractor is achieved by means of a small set of interchangeable attachments that do not complicate the construction, the machine being of light weight and small size. The tractor can be arranged for the required type of weld by setting the feeding mechanism with the nozzle to the required angle and using the relative interchangeable attachments, a set of which is furnished with every tractor (Fig. 3, 4, and 5).

An important feature of the TC-17-M tractor, as compared with other models of the same Class, is the possibility of carrying out small size angle (fillet) welds, as well as of welding light gauge sheet structures, using small diameter electrode wire and a low welding current.

For small size angle (fillet) welds, requiring very accurate direction of the arc, the tractor is furnished with a special „copy“ jig to ensure high pre-

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direction of the electrode along the weld (Fig. 4).

Welding by means of the TC-17-M tractor is generally carried out by using „copy“ jigs, but provision is also made for hand control.

In special cases, a light straight-edge can be used to guide the tractor rollers.

The electrical equipment of the tractor, mounted in accordance with the circuit diagrams (Fig. 6 and 7), ensures raising and lowering the electrode wire and moving the tractor for auxiliary operations; switching the welding current and establishing an arc at the commencement of welding, feeding the wire and moving the tractor during welding, and finally, welding the crater and disconnecting the welding current at the end of the weld. Control of all these operations is carried out by means of conveniently arranged push-buttons.

The welding current control and switching apparatus are mounted in a special portable control box 530 × 760 × 700 mm, furnished with the machine.

The TC-17-M Welding Tractor allows to obtain high quality welds not by the use of complicated auxiliaries or by high skill of the welder, but by perfect design of the Automatic Welding Machine itself.

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Fig. 5. Welding Tractor TC-17-M arranged for open angle welding

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
ESSENTIAL SPECIFICATIONS

Welding current	200 to 1000 A
Electrode wire diameter	1.6 to 5 mm
Rated voltage of control box	220 or 380 V
Rate of feed of electrode wire	52 to 403 m/hour
Rate of welding	16 to 126 m/hour
Motor of feeder and travel drive:	
Output	0.2 kW
Speed	2900 r. p. m.
Voltage	36 V
Inclination of feeder with nozzle	up to 45°
Flux hopper capacity	10 kg
Electrode wire reel capacity	8 kg
Overall dimensions:	
Height	510 mm
Width	345 mm
Length	720 mm
Weight of tractor (excluding flux and wire)	42 kg

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Machine for single-pass electric slag welding of heavy gauge metal A-372

The Model A-372 Welding Machine is used for single-pass electric slag vertical welding of 60 to 250 mm thick metal.

Welding is accomplished by double sided positive weld shaping by means of one, two or three electrodes simultaneously, with a 20-30 mm gap between the joint edges.

The use of three electrodes considerably increases the rate of welding, and also allows to vary within large limits the chemical composition of the weld by combining various electrode wires.


Welding is possible with either direct current or alternating current; maximum efficiency is obtained with alternating current by three-phase connection of the electrodes. The electric slag process of welding, as against arc welding, permits the use of a supply source having a rigid performance curve (without inductance in the welding circuit); it is thus possible to provide balanced loading and a high power factor.

The vertical electric slag welding process is a new process carried out as follows: in the space formed by the edges of the joint to be welded and the weld shaping slides (Fig. 1), a bath of molten overheated electro-conductive slag is created, into which the electrode is immersed. The electric slag process is characterized by the absence of an arc discharge. The current flowing from the electrode to the parent metal heats the slag and maintains its temperature above steel melting point. The slag bath melts the electrode wire fed into it and fuses over the edges. The molten parent and electrode metals create a metal bath. The metal and slag baths

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are held in place between the edges of the joint by two water cooled copper slides that shape the weld.

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Principal advantages of the electric slag welding process

1. Welding is carried out in a single pass irrespective of the amount of metal to be deposited. This eliminates the necessity of multi-layer welding and removal of the slag after each pass.
2. Preparation of the edges is simplified, as it is carried out by flame cutting without bevelling.
3. The cycle of fabrication is reduced owing to the reduction in the number of intermediate operations and of operator's places required.
4. The rate of deposition is increased to 25-30g/Ah.
5. Material consumption and electric power consumption are reduced:
 - a) flux consumption 15 to 20 times;
 - b) electric power consumption down to 35%;
 - c) electrode wire consumption down to 30-40%.
6. The quality of welds and their mechanical strength are higher than in the case of multi-layer welding.

Construction

The A-372 Welding Machine (Fig. 2) is of the rail type, as it moves along a guide (rail) column parallel to the edges of the joint. A rack is provided on the column, the gear of the travel carriage 1 being in mesh with this rack (Fig. 3). As the welding bath is being filled with metal, the travel carriage moves along the weld upwards (direction of welding) at a speed of 0.4 to 5 m/hour, thus maintaining in the welding bath a constant level relative to the nozzles or the edge of the slide. A constant level maintained in the welding bath is one of the main conditions required to obtain uniform fusion of the weld edges and a high quality weld. The vertical movement speed is therefore automatically adjusted depending on the level in the bath. Speed variation within each range is by varying

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the speed of motor 2. Range variation is by means of change gears.

The travel carriage is connected to a three-electrode single motor welding head 3 feeding the electrode wire into the weld at a rate of 150 to 500 m/hour. Three-phase electric slag welding is operated from a supply source having a rigid external characteristic.

With such a supply source the welding current varies only with the rate of electrode feed, without simultaneous adjustment of the supply source. The rate of feed of each of the electrodes is varied by means of change gears. The rate of feed can also be varied during operation, simultaneously for all the electrodes, by adjusting the speed of motor 4. The welding current is supplied to the electrodes by current carrying nozzles 5. The construction of the latter is such that they allow to adjust the curvature and direction of the electrodes. The electrodes are spaced 40-50 mm apart. Welding of metal 60-90 mm thick is therefore carried out by two electrodes, while welding of metal 90 mm thick is carried out by three electrodes. Welding of metal above 110 mm thick is carried out by a horizontal reciprocating movement of the electrodes along the weld at a speed of 20 to 61 m/hour. The

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head 3 is moved for this purpose along the support 6 by means of a mechanism with an individual electric drive 7. At the extreme points of the reciprocating movement at the shaping slides, where heat dissipation is most intensive, the electrodes stop. The duration of stoppage (up to 6 sec) can be set by a special device. The points, where the electrodes are to stop and commence moving in the opposite direction, are determined by a system of terminal switches 8.

The front shaping slide 10 is suspended from bracket 9 of the travel carriage. The rear shaping slide 11 is suspended from the rod of the rear slide 12, this rod being connected to the travel carriage across the weld gap by means of strip 13.

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As flux consumption with the positive shaping welding process is considerably less than with downhill welding, the flux is fed in small doses from the flux metering device 14.

Water for cooling the slides is supplied from the works water mains or from a tank by means of a pump through the water manifold 15. The rate of water supply is adjustable by means of a valve.

For setting operations, the carriage is provided with a travel speed of 36 m/hour.

The length of rail 16 is selected to suit the maximum length of the welds. The rail can be fixed either directly on the work or on a special installation, as required.

The Machine is furnished with electrode wire reels of two types: a) For short welds—three enclosed type reels 17, having a capacity of approximately 25 kg each, for installation on the travel carriage;

b) Three open type reels, having a capacity of approximately 135 kg each, for installation separately from the Machine.

A handle is provided for loading the reels with wire; this handle must be removed before welding is started.

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Fig. 4. View of the work to be welded from the rear seam forming slide

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The welding process is controlled by means of control panel 18.

The Machine is fitted with transverse, longitudinal and radial correctors to adjust, during welding, the position of the nozzles with the electrodes in the weld gap.

Two gauges 19 are furnished with the machine to enable to set the rail parallel to the weld.

Electrical circuit

A. C. welding current is supplied to the Automatic Welding Machine from a special three-phase welding transformer, Type ТПЩ-1000-3, without choke coils; the transformer is connected "delta-star with neutral brought out" when operating with three electrodes, or "delta-V with neutral brought out" when operating with two electrodes (Fig. 6). Three or two welding transformers, Type ТСО-1000-3, designed for 220 V mains can also be used; these transformers are connected "star with neutral brought out" — star with neutral brought out", the choke coils being disconnected, for 380 V mains. D. C. welding current can be supplied from three or two welding generators, Type ПСМ-1000 or ПС-500.

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The electrode wire feeder (Fig. 7) is driven by a Type ПН-2.5, 0.7 kW motor, current to the motor being supplied from step-down transformer with sectioned secondary winding and rectifier. The rate of feed is fixed by a switch mounted on a panel and is read on an instrument connect-

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led to the rectifier which supplies the electrode wire feeder.

The vertical travel mechanism is operated by a Type ПН-2.5, 0.25 kW motor, current to this motor being supplied from a Type ЭМУ-5 A generator. The generator voltage, and thus also

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the rate of welding are automatically maintained at the values set by a special device and corresponding to a definite level of the metal bath.

Transverse adjustment of the automatic welding head is by means of a Type МАГ-2-36, 0.1 kW motor. Type МП-1 terminal switches are used


Fig. 7. Mechanism of electrode wire feeding and current-leading electrode holder of the А-372 Welding Apparatus

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for limiting the end positions of the welding head and for reversing.

The welding process controls are arranged on the control panel, while all control and switching apparatus are mounted in the control box.

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Essential specifications


Electrode wire diameter 3 mm
 Welding current, intermittent up to 1000 A
 Welding current, continuous up to 800 A
 Rated mains voltage 380 V
 Rate of electrode feed 150 to 500 m/hour
 Speed of reciprocating electrode movement 21 to 60 m/hour
 Rate of welding (vertical travel) 0.4 to 5 m/hour
 Setting speed (vertical travel) 36 m/hour
 Welding current adjustment remote stepping
 Capacity of enclosed type reel 25 kg
 Capacity of open type reel 135 kg
 Voltage on electrodes 40 to 56 V
 Slides water cooled
 Quantity of cooling water required 25 to 30 l/min

Overall dimensions: Length Width Height

Machine without rail, mm	1300	830	1660
Control box, mm	1044	644	2000
Reels (each), mm	700	295	7300

Weight:
 Machine 310 kg
 Control box 520 kg
 Rail: The weight of the rail, being dependent on its length, can be calculated as follows: $G = 71 L + 50$ kg, where L is the length of the rail in metres.

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PLEASE ADDRESS ALL ENQUIRIES IN CONNECTION
WITH PURCHASING EQUIPMENT TO:

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6201

COPPER SMELTING ELECTRIC ARC FURNACES

SERIES DM

175—250 kW, 6,000/10,000 V

II. APPLICATION

The Series DM Electric Arc Furnaces are available in two sizes and are used to smelt copper and some of its alloys. The electric supply to the furnace is through a special single-phase transformer connected to the high tension line.

II. SPECIFICATIONS

	Unit of measurement	Type of furnace	
		DM-0.25	DM-0.5
Capacity of bath	t	0.25	0.50
Allowable overload	%	20	20
Diameter of working space	mm	550	690
Length of working space	mm	680	870
Dimensions of charging door	mm	270×220	300×250
Diameter of carbon electrodes	mm	75	100
Length of electrodes	mm	1,000	1,000
Capacity of furnace transformer	kVA	175	250
Primary voltage of transformer	kV	6 or 10	6 or 10
Secondary voltage of transformer	V	100	100
Number of phases of transformer	—	1	1
Working Data for continuous (24 h) Operation:			
a) Duration of smelting	min	35—45	40—50
b) Hourly output	kg	280	450
c) Energy consumption for melting 1 t of metal	kWh	325	300
d) Electrode consumption per t of metal	kg	3.5	3.0
e) Water consumption	m ³	2.5	3.0
Weight of metal framework	t	1.4	1.5
Total weight with lining	t	2.3	2.8

The furnaces are erected on special foundations. Detailed instructions for erection, lining, drying the lining, starting and operation will be forwarded by Makers' Works.

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6132

SERIES B SALT-BATH ELECTRIC FURNACES

10-30 kW * 220/380 V * 850° C

The Series B salt-bath electric furnaces are available in three types and sizes. They are used for heating metal products in molten salts and for special kinds of thermochemical treatment.

The maximum working temperature is 850°C. The products, suspended on special devices, are introduced into the bath through a slide cover in halves.

SPECIFICATIONS

	Type		
	B-10	B-20	B-30
Working dimensions of crucible:			
diameter, mm	200	300	400
depth, mm	350	535	555
Rated power, kW	10+10%	30+5%	30+5%
Supply voltage, V	220	380 or 220	380 or 220
Number of phases	1	1/3	3
Connection of heaters	Series	Series-star	Star-delta
Number of heating zones	1	1	1
Maximum working temperature in crucible, °C	850	850	850
Operation data:			
a) no-load power, kW	6	7	9
b) industrial capacity, kg/hour	30	80	130
Weight of furnace, t	1	1.4	1.8


The type B-10 salt-bath electric furnaces are designed for direct connection to a 220 V line, while those of types B-20 and B-30 may be direct connected to a 380 or 220 V line.

DATA TO BE STATED WHEN ORDERING

When ordering please state: type of electric furnace, line voltage for connection of the furnace and number of furnaces required. Thus: Electric Furnace Type B-20, line voltage 220 V. Quantity Furnaces.

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6221

STEEL SMELTING ELECTRIC ARC FURNACES

series DCB

5,000—8,000 kW; 6,000 and 10,000 V

APPLICATION

The Series DCB Three-phase Electric Arc Furnaces are available in two sizes with melting bath of 10 and 20 ton capacities. The Series DCB Furnaces are designed for smelting steel by either the basic or acid process, as well as for refining liquid steel. The electrical supply of the furnaces is from a high-voltage station through special 3-phase furnace transformers having a choke coil built into the common shell.

SPECIFICATIONS

Description of parameter	Unit of measurement	Type of electric furnace	
		DCB-10	DCB-20
Capacity of bath	tons	10	20
Allowable overload	%	20	15
Capacity of transformer	kVA	5,000	8,000
High tension	kV	6 or	6 or 10
Low tension	V	250—228—208—144—120	261—229—204—151—132
Diameter of working space	mm	2,920	3,300
Diameter of bath at level of threshold	mm	2,830	3,180
Depth of bath to threshold	mm	475	590
Diameter of charging door	mm	680 × 940	750 × 940
Diameter of carbon electrodes	mm	350	400
Diameter of disintegration of electrodes	mm	1,100	1,200
Electrode travel	mm	2,000	2,000
Working Data for Continuous (24 hour) Operation			
Duration of smelting	hours	1.5	2
Consumption for smelting 1 ton of metal:			
a. Electrical energy	kWh	510—475	510—475
b. Carbon electrodes	kg	8—10	8—10
Water consumption	cu. m per hr	25	25
Velocity of electrode travel	m per min	1	1
Time for tilting furnace through 40°	sec	50	60
Time for lifting roof	sec	20	20
Time for rolling out the bath	sec	30	30
Weight of metal framework	tons	85	115
Total weight with lining	tons	120	175

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VOLTAGE AND CURRENTS FOR FURNACE TRANSFORMER TYPE ЭТЦК = 11000/10 FOR SYMMETRICAL VOLTAGE REGULATION OF ALL PHASES

V	A	High tension			Position of controls indicator	Low tension		Reactance of choke coil, %
		Change-over switch				V	A	
		phase A	phase B	phase C				
6,000	770	2-3	2-3	2-3	1	261	17,700	100
6,000	675	3-4	3-4	3-4	2	229	17,700	100
6,000	675	4-5	4-5	4-5	3	229	17,700	67
6,000	600	5-6	5-6	5-6	4	204	17,700	67
6,000	346	6-7	6-7	6-7	5	118	17,700	—
6,000	445	7-2	7-2	7-2	6	151	17,700	—

SHIPMENT AND ORDERING DATA

The mechanical parts of the Series DCB Electric Arc Furnaces are shipped disassembled with the high-pressure hydraulic unit for feeding the electric furnace actuating cylinders, the electric motors for lifting the movable platform, bath rotating and electrode regulating mechanisms and with the special furnace transformer. An automatic electrode regulator can be furnished on special order.

When ordering, please state:
1. Application of furnace (acid or basic process).

2. Required furnace output per 24 hours.
3. Desirable furnace capacity.
4. Design of furnace (left- or right-hand).
5. High tension voltage of Works' mains.
6. Voltage of Works' power mains.
7. Hoisting capacity of travelling crane in the pouring bay.
8. Pressure in existant compressed air lines.
9. Pressure in water line.

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HIGH FREQUENCY SETS FOR HARDENING AND HEATING TYPE МГЗ

50-100 kW • 750 V • 2500-8000 c. p. s.

I. APPLICATION

High Frequency Sets with machine generator (frequency changer) for hardening are available in three sizes: МГЗ-52, МГЗ-102 and МГЗ-108. The Sets are chiefly used for surface hardening of steel parts. They may also be used for preheating for forging or stamping, annealing, normalizing and copper brazing. The induction method for heating ensures high speed heating measured in seconds, a reduc-

tion in the amount of scale, the possibility of hardening or heating only the required portions of a product and the possibility of performing heat treatment operations in a production line sequence.

Induction heating provides best results in the mass production of one type of products.

Type МГЗ High Frequency Sets are available for power consumption 50 or 100 kW at frequencies of 2500 or 8000 c. p. s.

II. SPECIFICATIONS OF HIGH FREQUENCY SETS

Parameters of the Set	Unit of measurement	МГЗ-52	МГЗ-102	МГЗ-108
Rated consumption from mains	kW	65	130	130
Voltage of driving electric motor of frequency changer	V	380/220		
Voltage of operating circuits	V	380/220		
Phases	—	3		
Frequency of supply	c. p. s.	50		
Rated generator output	kW	50	100	100
Generator voltage	V	750		
Generator frequency	c. p. s.	2500	2500	8000
Power factor	—	0.9 (leading)		
Power of condenser bank	kVAR (reactive)	750	1000	1000
Output of high frequency transformer	kVA	500	500	500
Transformer ratio	—	30; 28; 24; 19; 16		
Number of inductor posts	pieces	2		
Maximum total time setting (heating, pause, cooling)	sec	60		
Minimum diameter of products to be hardened	mm	30	30	16
Minimum depth of hardening	mm	2	2	1
Minimum diameter of products for through heating	mm	35	35	20
Cooling water consumption	litres per min	100		
Overall dimensions of heating station:				
Width	mm	1325		
Length	mm	2280		
Height	mm	2135		
Weight of heating station	t	2		

VSESOJUZNOJE OBJEDINENIJE

«MACHINOEXPORT»

III. REGULAR EQUIPMENT AND ORDERING DATA

The complete set comprises:
1. Heating station with 2 hardening transformers.

2. Autotransformer starter, type ATI.
3. High frequency changer with electric motor voltage 380/220 V.
When ordering High Frequency Sets, please state the type of set required and its application.

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TYPES AB-2 AND AB-4 A. C. GASOLINE ELECTRIC POWER PLANTS

APPLICATION

Types AB-2 and AB-4 plants are mobile electric stations and can be used as independent sources of A.C. electric energy, having industrial frequency (50 Hz).

The plants can be used at ambient temperatures of +50° to -50° C, a relative humidity up to 98% and at an elevation above sea level up to 1000 m.

Output voltage of the plant can be set in the range of 230 to 218 V at any load from no-load up to rated current with power factor from 0.8 to 1.0.

At variations of load from no-load up to rated current and with power factor from 0.8 up to 1.0 output voltage of the plant is automatically maintained constant with an accuracy of ±4%, while the difference between the maximum and minimum values of frequency does not exceed 2 Hz.

The plants can withstand a 10% overload of nominal power for one hour at an ambient temperature not higher than +35° C and provide continuous operation at nominal power for 24 hours.

CONSTRUCTION AND TECHNICAL DATA

The plant consists of the following main units:

- 1) an internal combustion engine—type VД-1 low litre capacity gasoline engine (for AB-2) and type VД-2 (for AB-4);
- 2) one-phase or three-phase A.C. generator;
- 3) starting, control and regulation voltage equipment for the generator;
- 4) measuring and control instruments of electric output performances of the plant.

An enclosure and a tray are used for protecting the plants from mechanical damages, dirt and precipitates.

The plants are also available without enclosure and tray.

The gasoline engine rotates the generator through an elastic coupling.

The engine-generator block is located on the supports which are mounted on the plant frame by means of special rubber shock absorbers.

The instrument block is also mounted on rubber shock absorbers. The plant has a fuel tank containing gasoline which permits the operation of plant during four hours. It is also provided with a heating device for preheating the engine with a soldering lamp when starting at low ambient temperatures.

Four-cycle gasoline engine, air-cooled 8 h.p. for type AB-4 plant and 4 h.p. for type AB-2 plant is provided with a centrifugal speed regulator.

A synchronous generator is one-phase or three-phase (depending on the plant design), two-pole, with self-excitation. Nominal speed—3000 r.p.m.

The electric basic diagram of the plant has the following specific features:

- a) a special exciter for the synchronous generator is absent. The field winding is fed through the bank of selenium rectifiers mounted in the generator body;
- b) a special compounding excitation system is used for automatic output voltage stabilization of the plant with the change of load from no-load to nominal load;
- c) Protection from radio interference, overload and short-circuit protection are provided by means of special condensers and relays.

The plant has a lighting circuit of safety voltage of 12 V fed through a step-down transformer from the main plant circuit. An electric bulb illuminating the control panel and a plug socket for a portable lamp are connected in this circuit.

TECHNICAL DATA

Type of the plant	AB-2-0/230	AB-2-T/230	AB-4-0/230	AB-4-T/230
GENERATOR				
Nominal power, kW	2	2	4	4
Current	one-phase	three-phase	one-phase	three-phase
Nominal voltage, V	230	230	230	230
Power factor (at induction load)	1—0.8	1—0.8	1—0.8	1—0.8
Load current:				
at power factor, (Cos φ) = 1 A	8.7	5	17.4	10
at power factor, (Cos φ) = 0.8 A	10.9	6.25	21.8	12.5
Nominal frequency, Hz	50	50	50	50
INTERNAL COMBUSTION ENGINE				
Type	VД-1	VД-1	VД-2	VД-2
Power, h.p.	4	4	8	8
Number of cylinders	1	1	2	2
Speed, r.p.m.	3000	3000	3000	3000
Fuel	Gasoline with a 66 octane number			
Fuel consumption gm. h.p. per hour	370			
Ignition	Magneto			
Overall dimensions, mm	length	890	890	1075
	width	560	560	560
Weight without fuel, spare parts and accessories, kg	height	870	870	870
	with enclosure	not more 908	not more 870	not more 270
	without enclosure	not more 180	not more 177	not more 240

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1690T

APPLICATION

The function of a samosyn system is to transmit motion by electric means between two points. Type CC195-150T power samosyns are designed for use as transmitter samosyns in remote position type (angular displacement) self-synchronous system in metallurgical and other branches of industry.

Increased strength of electrical and mechanical construction details provide for dependable operation of the samosyn in the conditions of dust, shaking and great temperature changes.

In a complex automatic control charging system for blast furnaces, for example, a type CC195-150T samosyn is used as a remote transmitter of angular displacement of the revolving distributor in a transmitter-receiver scheme without amplifier.

SPECIFIC FEATURES

Type CC195-150T samosyn has a three-phase, distributed, double-layer stator winding. The windings may be connected in wye or in delta. The magnetic circuit of the stator is packeted from thin laminations of electrical steel. The samosyn has a salient pole rotor with a coil type, single-phase winding. The ends of the rotor

winding are brought out to slip rings which are insulated from each other and other parts of the samosyn.

Type CC195-150T samosyns are enclosed. They are designed for operation in blast furnaces and rolling mills in the tropics.

DESCRIPTION OF DESIGN

Type CC195-150T samosyns are constructed for horizontal mounting on feet. The body and bearing end shields are made of cast-iron. The terminal box for the stator and rotor windings is located on the side of the body. The stator winding can be connected in wye or in delta and the shaft can be brought out at one or both ends upon request of the customer.

Antifriction bearings are used. The slip rings are of silver, the brushes of graphite. When installed outdoors the samosyn should be protected against solar radiation by an awning.

The mechanical strength of the samosyn permits rotation speed up to 1.500 r. p. m.

BASIC TECHNICAL DATA

The samosyn is rated for continuous operation. Type CC195-150T samosyn weighs 65 kg. Technical data of a type CC195—150T samosyn are given in the table.

TECHNICAL DATA

	Ротор	Статор	Comment
Current	A. C. single-phase 50	A. C. three-phase 50	
Frequency, c. p. s.	220	appr. 350	
Rated voltage, V	6.5	appr. 200	at 45°С
Rated current, A		appr. 2.0	at 45°С
Torque for difference of 1° in angular position, kg-cm	appr. 1.3	appr. 3.5	Resistance of each connecting lead not more than 1.0 Ω
Difference in angle for rated torque, degrees	up to 9.0		Including angle of internal losses
Angle of internal losses, degrees	2.0		
Rated torque, kg-cm	8.3		
Breakdown torque, kg-cm	135		

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«MACHINOEXPORT»

FREQUENCY CHANGER SETS

constant speed frequency changer sets (with synchronous or induction motor) variable speed frequency changer sets (with D.-C. motor)

This Catalogue Issue covers:

1. Constant speed Frequency Changer Sets:
 - a) High frequency Sets;
 - b) Normal frequency Sets;
 - c) Low frequency Sets.
2. Variable speed Frequency Changer Sets.

All Frequency Changer Sets included in this Catalogue Issue are designed to operate in normal humidity enclosed rooms at an ambient temperature not exceeding $\pm 35^{\circ}\text{C}$.

The windings of the A.-C. machines and exciters (except synchronous machine field windings) have class A insulation and are suitable for a 65°C temperature rise.

The synchronous generator field windings in Sets No. 18 and No. 19 have class A insulation. The field windings of all other synchronous machines have class B insulation and are suitable for a 95°C temperature rise.

The sets are available, at Maker's option, either on common or on separate bedplates, the machines being connected together by means of either rigid or flexible couplings.

Sets designed for forced cooling air discharge or intake are not furnished with the required additional enclosures, fans, etc.

The Sets are not fitted with thermocouples or any other thermosignalling devices in the windings or cores.

All A.-C. driving motors (synchronous or squirrel cage induction motors) are suitable for direct-on-line starting.

Autotransformers or reactors for stepped down voltage starting are not furnished with the Sets. Data to enable to select a reactor or an autotransformer will be submitted on request.

The synchronous motors are suitable for operation at a reduced $\cos \phi$ (when the Sets are not

operating on full output) provided the rated field current is maintained.

All driving synchronous motors have an exciter mounted with the motor on a common bedplate, or a separate exciter set.

The dimensions and specifications of the starting and other apparatus furnished with the set will be submitted by the Maker on request.

The exciters are furnished on separate bedplates with a Single Series A driving induction motor.

The exciters, forming part of the exciter sets, are self-excited. By mutual consent, the exciters can be furnished for 110 or 220 V separate excitation.

Automatic voltage regulators and electrical machine amplifiers are not furnished with the Sets.

Generators designed for one voltage have normally 4 to 6 stator winding ends brought out. Generators designed for more than one voltage with stator winding interconnection change-over have a larger number of stator winding ends brought out, a diagram of connections being supplied with the set.

Terminal boards for stator winding interconnection change-over are fitted only on low tension (up to 1000 V) generators, for a maximum of 9 stator winding ends brought out.

The machines forming the Set are furnished with spare parts to Maker's specification.

The dimensions of the Sets in the Catalogue are approximate only and are not to be used for foundation design. Final dimension drawings will be submitted on request.

Generator type designations formerly used are given in brackets.

CONSTANT SPEED FREQUENCY CHANGER SETS

Constant speed Frequency Changer Sets are available for high frequency changing, normal frequency changing and low frequency changing.

All constant speed Sets are driven by A.-C. motors (synchronous or induction motors).

The driving motors indicated in the Specifications are those of actually manufactured Sets;

by mutual consent, they can be replaced by others of different voltage and design.

Starting of the synchronous motors is by means of a magnetic controller (Control Station) and a field discharge resistance. By mutual consent, starting can be arranged with the exciter direct connected.

SETS No. 1 and No. 2

The synchronous generator is open self-ventilated with ring lubricated pedestal bearings. The driving motor is a protected squirrel cage induction motor with ring lubricated endshield bearings.

Generator excitation is from a separate exciter set consisting of:

- 1) ПН 100 exciter, 7.6 kW, 65 V, 1000 r.p.m.;
 - 2) a driving squirrel cage induction motor, Type A62-6, 10 kW, 380 V, 970 r.p.m.
- The set is furnished complete with a Type P3B-31A field regulator for the exciter field circuit.

SETS No. 3 and No. 4

The generator is open self-ventilated, with ring lubricated pedestal bearings. The driving motor is a protected slip-ring induction motor with ring lubricated endshield bearings.

The set is furnished complete with:

- 1) a rotor magnetic controller (Control Station), Type ЧНЛ 6701-59B5, and a stator controller (Control Station), Type ЧНЛ 5101-55A5, for direct-on-line starting of the ФАМКО 1410-6 motor;
- 2) a rotor magnetic controller (Control Station), Type ЧНЛ 6701-59B, for direct-on-line starting of the ФАМКО 1410-6 motor;
- 3) a starting resistance box, Type ЯПМ-6;
- 4) a field regulator, Type PB 12L, for the exciter field circuit.

SET No. 5

The machines of this Set are open self-ventilated, with ring lubricated pedestal bearings. The generator has 6 stator winding ends brought out.

The driving motor is a synchronous motor. The set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the BC 24.5/12 exciter field circuit, and Type PB 18L for the BC 24.5/18 exciter field circuit;
- 2) a magnetic controller (Control Station), Type ЧНЛ 7301-42AO, for motor starting;
- 3) a field discharge resistance, Type CD 32B, for the ДС 1610-6 motor field circuit.

SET No. 6

The machines of this Set are open self-ventilated with ring lubricated pedestal bearings. The generator allows a single-phase intermittent load with a 20% "ON" period, 1000 kVA, 400/800 V, 1000 r.p.m., $\cos \varphi = 0.75$.

Excitation at single-phase load: 57.6 V, 247 A. The generator has 9 stator winding ends brought out to enable to change-over for 400 V or 800 V.

The driving motor is a synchronous motor. The set is furnished complete with:

- 1) a field regulator, Type PB 12L, for each exciter field circuit;
- 2) a field discharge resistance, Type CD 32B, for the motor;
- 3) a magnetic controller (Control Station), Type ЧНЛ 7301-42AO for motor starting.

SETS No. 7 and No. 8

The machines of this Set are open self-ventilated with ring lubricated pedestal bearings.

The TCO 1510-6 generator has an extra strength damper winding shorted by copper rings. The generator has a full three-phase winding for three-phase operation. The three-phase output must be settled by mutual consent. The generator has 36 stator winding ends brought out to enable to change-over for 550 or 3300 V. On single-phase operation the third phase is idle and its terminals must be insulated.

The driving motor is a synchronous motor. The set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the field circuit of each exciter;
- 2) a magnetic regulator, Type MP 120L, for the synchronous generator rotor circuit to provide extra smooth field regulation;
- 3) a rotor magnetic controller (Control Station), Type ЧНЛ 7301-32AO for starting either of the synchronous generators;
- 4) a stator magnetic controller (Control Station), Type ЧНЛ 5101-53A3 for the ДС 1405-6 motor;
- 5) a field discharge resistance, Type CH-2, for the ДС 1408-6 synchronous generator, or Type CD-32B for the ДС 1405-6 motor.

SET No. 9

The machines of this Set are available with ring lubricated pedestal bearings and in two alternative enclosures:

1. open self-ventilated, with air intake from the room through the endshields and discharge into the room through openings in the frame, the covers being then removed from these openings;
2. enclosed forced ventilated, with air intake through the endshields and discharge into the foundation pit. The forced ventilated Set must not be allowed to operate without additional suction fans.

Suction fans are not furnished with the Set. The synchronous generator can operate with frequently recurrent current peaks up to 150% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is a synchronous motor. Generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type BC 24.5/12, 11.6 kW, 50 V, 1000 r.p.m.;
- 2) an induction motor, A71-6, 14 kW, 380 V, 970 r.p.m.

The set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the BC 24.5/12 exciter field circuit, and Type PB 18L for the BC 24.5/18 exciter field circuit;
- 2) a field discharge resistance, Type CD 32B, for the motor;
- 3) a rotor magnetic controller (Control Station), Type ЧНЛ 7301-42AO, for motor starting.

SETS No. 10 and No. 11

The generator is open self-ventilated with ring lubricated pedestal bearings.

The driving motor is a squirrel cage induction motor with ring lubricated endshield bearings.

The set is furnished complete with a Type PB 12L field regulator for the exciter field circuit.

SET No. 12

The generator is open self-ventilated with ring lubricated pedestal bearings.

The generator can operate with frequently recurrent current peaks up to 150% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is a squirrel cage induction motor.

Generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type BC 24.5/11, 8.5 kW, 45 V, 1000 r.p.m.;
- 2) a driving induction motor, Type A62-6, 10 kW, 380 V, 970 r.p.m.

The set is furnished complete with a Type PB 12L field regulator for the exciter field circuit.

SET No. 13

The machines are open self-ventilated with ring lubricated pedestal bearings.

The generators allow an excess current 200% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is a synchronous motor.

The synchronous generator excitation is from a separate exciter set consisting of:

- 1) two exciters, Type BC 24.5/12, 17.5 kW, 75 V, 1500 r.p.m.;
- 2) a driving squirrel cage induction motor, Type A81-4, 40 kW, 380 V, 1450 r.p.m.

The set is furnished complete with:

- 1) a field discharge resistance, Type CH4/2, for the rotor circuit of the ДС 1707-12 motor;
- 2) a field regulator, Type PB 18L, for the field circuit of each of the BC 24.5/12 exciters, and Type PB 12L the BC 34/18 exciter field circuit;
- 3) a magnetic controller (Control Station), Type ЧНЛ 7301-42AO, for direct-on-line starting of the synchronous motor.

SET No. 14

The generator has pedestal bearings, the motor has endshield bearings, ring lubricated.

The generator is available with either of the following two enclosures:

1. Open self-ventilated, with air intake from the room through the endshields and discharge

into the room through openings in the frame, the covers being then removed from these openings;

2. Enclosed forced ventilated, with air intake through the endshields and discharge into the foundation pit.

The forced ventilated Set must not be allowed to operate without additional suction fans. Suction fans are not furnished with the Set. The synchronous generator can operate with frequently recurrent current peaks up to 150% of the rated value; the RMS value of the rotor and stator currents must then not exceed the rated value.

The driving motor is a squirrel cage induction motor.

Generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type BC 24.5/12, 11.6 kW, 50 V, 1000 r.p.m.;
- 2) an induction motor, Type A71-6, 14 kW, 380 V, 970 r.p.m.

The set is furnished complete with a field regulator, Type PB 12L, for the exciter field circuit.

SET No. 15

The machines of this Set are open self-ventilated with ring lubricated pedestal bearings.

The driving motor is a synchronous motor.

The set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the field circuit of each of the exciters;
- 2) a field discharge resistance, Type CH-2, for the motor;
- 3) a rotor magnetic controller (Control Station), Type ЧНЛ 7301-32AO, and a stator magnetic controller (Control Station), Type ЧНЛ 5101-53A3, for direct-on-line starting of the motor.

SET No. 16

The machines are open self-ventilated with ring lubricated pedestal bearings.

The driving motor is a squirrel cage induction motor.

The synchronous generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type ПН 45, 3.3 kW, 115 V, 1460 r.p.m.;
- 2) a squirrel cage induction motor, Type A 51-4, 4.5 kW, 380 V, 1440 r.p.m.

The set is furnished complete with a field regulator, Type P3B 11-B for the exciter field circuit.

Specifications are given in Tables 1, 2, 3.

CONVENTIONAL DESIGNATIONS IN THE FREQUENCY CHANGER SET SPECIFICATION TABLES

Reactances (in relative units):	Time constants (in seconds):
X_1 — Stator winding leakage;	T_{α} — Closed stator circuit time constant of field winding;
X_{σ} — Direct-axis synchronous reactance;	T_{σ} — Stator short-circuit time constant of field winding;
X'_d — Direct-axis transient reactance;	T_r — Rotor short-circuit time constant of stator winding.
X''_d — Direct-axis subtransient reactance.	

SPECIFICATIONS OF CONSTANT SPEED FREQUENCY

Set No.	Type of generator Type of motor	Output		Rating				efficiency, %	cos φ	frequency, c.p.s.	current, A	voltage, V	Field rating ¹	Stator phase resistance at 15°C, ohm	Rotor winding resistance at 15°C, ohm	Short-circuit ratio
		kVA	kW	voltage, V	stator current, A	speed, r.p.m.	stator current, A									
Sets, 150 c.p.s.																
1	Generator ГС4 1508-12	750	600	3300	131	1500	94.4	0.83	150	54	29.5	0.103	0.384	0.52		
	Motor ДАМСО 1410-4	—	680	3000	156	1480	94.0	0.89	50	—	—	0.344	—	—	—	
2	Generator ГС4 1508-12	750	600	3300	131	1500	94.4	0.83	150	54	29.5	0.103	0.384	0.52		
	Motor ДАМСО 1508-4	—	680	6000	79	1480	93.0	0.88	50	—	—	0.51	—	—	—	
Sets, 133 c.p.s., for supply																
3	Generator ГС4 1506-16 (ВГ440-1000-133)	440	352	550	462	1000	90.8	0.8	133	136	41.7	0.01	0.216	0.76		
	Motor ФАМСО 147-6	—	390	500	572	985	92.5	0.85	50	525	470	0.0205	0.007	—	—	
4	Generator ГС4 1506-16 (ВГ440-1000-133)	440	352	550	462	1000	90.8	0.8	133	136	41.7	0.01	0.216	0.76		
	Motor ФАМСО 1410-6	—	380	600	46	985	92.0	0.87	50	385	615	1.36	0.0082	—	—	

¹ In the case of slipping motors, rotor data are given in this column
² cos φ leading.
³ Excitation from a separate exciter set.

SPECIFICATIONS OF CONSTANT SPEED FREQUENCY

Set No.	Type of generator Type of motor	output		Rating				efficiency, %	cos φ	frequency, c.p.s.	current, A	voltage, V	Field rating	Stator phase resistance at 15°C, ohm	Rotor winding resistance at 15°C, ohm	Short-circuit ratio
		kVA	kW	voltage, V	stator current, A	speed, r.p.m.	stator current, A									
Low-voltage Set																
5	Generator ГС 1510-6(С940-1000)	940	750	400	1360	1000	94.7	0.8	50	192	45	0.0015	0.164	0.77		
	Motor ДС 1610-6	1250	945	6000	120.5	1009	94.0	0.8	50	285	40	0.305	0.099	0.77		
Set with change-over for																
6	Generator ГС 1510-6(С600-1000)	600	450	400	855	1000	94.1	0.75	50	148	34.5	0.00165	0.164	1.0		
	Motor ДС 1508-6	650	483	6000	62.5	1000	93.0	0.8	50	158.5	37.6	0.817	0.161	0.76		
Single phase Sets with stator winding																
7	Generator ГСО 1510-6(СО550-1000)	550	165	3300	167	1000	84.8	0.3 ¹	50	45	11	0.234	0.187	1.4		
	Motor ДС 1408-6	360	263	6000	34.7	1000	91.4	0.8	50	138	36	1.43	0.181	1.3		
8	Generator ГСО 1510-6(СО550-1000)	550	165	3300	167	1000	84.8	0.3 ¹	50	45	11	0.234	0.187	1.4		
	Motor ДС 1405-6	330	242	380	502	1000	91.6	0.8	50	146	35.5	0.0079	0.176	0.65		

¹ cos φ leading

CHANGER SETS (HIGH FREQUENCY)

Table 1

X _l	X _d	X _d ^m	X _d ^m	Time constants			Starting data			Flywheel effect of rotor, kg·m ²	Cooling air required, m ³ /min	Weight of set complete, kg	Exciter		
				T _{do} , sec	T _d , sec	T _o , sec	I _{start} /I _{rated}	M _{start} /M _{rated}	M _{max} /M _{rated}				Type	output, kW	voltage, V
for test transformers															
0.213	2.06	0.44	0.44	0.97	0.207	0.09	—	—	—	390	1.2	8150	s		
0.213	2.06	0.44	0.44	0.97	0.207	0.09	—	—	—	390	1.2	9100	s		
—	—	—	—	—	—	—	—	—	—	4.8	1.0	1.9	170	1.7	
to spinning machine spindle motors															
0.174	1.61	0.41	0.41	0.65	0.165	0.043	—	—	2.05	420	1.1	6650	BC 24. 5/11	8.5	45 1000
—	—	—	—	—	—	—	—	—	1.8	125	1.2				
0.174	1.61	0.41	0.41	0.65	0.165	0.043	—	—	2.05	420	1.1	7200	BC 24. 5/11	8.5	45 1000
—	—	—	—	—	—	—	—	—	2.4	170	1.2				

CHANGER SETS (NORMAL FREQUENCY)

Table 2

X _l	X _d	X _d ^m	X _d ^m	Time constants			Starting data			Flywheel effect of rotor, kg·m ²	Cooling air required, m ³ /min	Weight of set complete, kg	Exciter		
				T _{do} , sec	T _d , sec	T _o , sec	I _{start} /I _{rated}	M _{start} /M _{rated}	M _{max} /M _{rated}				Type	output, kW	voltage, V
for test Station															
0.082	1.54	0.25	0.21	2.34	0.32	0.125	—	—	2.1	370	1.4	14000	BC24. 5/12	11.6	50 1000
0.103	1.53	0.21	0.154	2.19	0.36	0.038	5.2	2.5	2.1	600	2.1		BC24. 5/18	16.0	50 1000
400 and 800 V, for test Station															
0.0575	1.23	0.14	0.14	2.7	0.31	0.135	—	—	2.6	370	0.9	8940	BC29. 5/14	20	65 1000
0.093	1.65	0.218	0.165	2.3	0.28	0.035	5.2	2.2	2.0	300	1.3		BC24. 5/11	8.5	45 1000
change-over for 550 and 3300 V, for test Station															
0.05	0.75	0.125	0.083	1.58	0.26	0.019	—	—	—	390	1.0	9100	BC24. 5/11	8.5	45 1000
0.062	1.07	0.125	0.1	1.74	0.2	0.02	7.3	3.5	2.8	160	0.8		BC24. 5/11	8.5	45 1000
0.05	0.75	0.125	0.083	1.58	0.26	0.019	—	—	—	390	1.0	8500	BC24. 5/11	8.5	45 1000
0.123	1.88	0.232	0.189	1.34	0.165	0.032	4.0	1.9	2.0	100	0.9		BC24. 5/11	8.5	45 1000

SPECIFICATIONS OF CONSTANT SPEED FREQUENCY

Set No.	Type of generator Type of motor	Rating										Field rating		Short-circuit ratio
		Output		voltage, V	rator current, A	speed, r.p.m.	efficiency, %	cos φ	frequency, c.p.s.	current, A	voltage, V	Stator phase at 15° C. min	Rotor winding resistance at 15° C. min	
kVA	kW	kw	amp											rpm
Set, 30 c.p.s., for supply														
9	Generator ГСЧ 1510-6 (НГ560-600)	560	224	240	1345	600	88.7	0.4	30	169	41	0.00148	0.17	0.57
	Motor ДС 1508-10	450	328	6000	43.3	600	91.1	0.8	50	162.5	41	1.36	0.177	0.85
Set, 29 c.p.s., for supply														
10	Generator ГСЧ 1510-6	410	123	307	770	585	84.2	0.3	29	140	34.5	0.00293	0.171	0.73
	Motor ДАМСО 147-10	—	200	3000	52	585	90.5	0.82	50	—	—	0.61	—	—
11	Generator ГСЧ 1510-6	410	123	307	770	585	84.2	0.3	29	140	34.5	0.00293	0.171	0.73
	Motor ДАМСО 1410-10	—	200	6000	27	590	89.0	0.8	50	—	—	3.0	—	—
12	Generator ГСЧ 1407-6 (НГ220-600)	220	88	232	548	585	83.3	0.4	29	143	35	0.00565	0.172	0.7
	Motor ГАМ-6-128-10	—	180	380	270	585	91.5	0.82	50	—	—	0.0164	—	—
Set, 25 c.p.s., for supply														
13	Generator ГСЧ 1510-6*	470	188	200	1355	500	87.0	0.4	25	169	41	0.00148	0.17	0.57
	Motor ДС 1707-12	675	500	6000	65	500	93.0	0.8	50	186	51	0.64	0.192	1.09
Set,														
14	Generator ГСЧ 1510-6 (НГ400-428)	400	160	172	1345	428	86.0	0.4	21.5	169	41	0.00148	0.17	0.57
	Motor ДА 1510-14	—	260	6000	35	418	89.7	0.8	50	—	—	1.9	—	—
Set;														
15	Generator ГСЧ 1507-6 (НГ200-375)	200	150	150	770	375	90.5	0.75	18.7	161	28	0.0026	0.122	0.87
	Motor ДС 1506-16	230	163	380	380	375	88.5	0.80	50	124	41.2	0.0171	0.234	0.84
Set,														
16	Generator ГСЧ 1407-6 (НГ67-333)	67	50	138	290	333	88.0	0.75	16.6	15.2	93	0.0079	3.82	1.05
	Motor ДАМСО 1504-18	—	60	380	144	324	87.0	0.78	50	—	—	0.0467	—	—

* Excitation from a separate exciter set.
 † The Set consists of two identical generators and a common driving motor.

CHANGER SETS (LOW FREQUENCY)

Table 3

X ₁	X _d	X' _d	X'' _d	Time constants			Starting data			Flywheel effect of motor, kg·m ²	Cooling air required m ³ /sec	Weight of set complete, kg	Exciter			
				T' _d , sec	T'' _d , sec	T _s , sec	M _{start} /M _{rated}	M _{start} /M _{rated}	M _{max} /M _{rated}				Type	output kW	voltage V	speed, r.p.m.
to roll table motors																
0.089	2.27	0.109	0.109	3.68	0.176	0.109	—	—	4.0	370	1.1	9300	BC24.5/18	7.9	45	600
0.115	1.62	0.24	0.165	1.74	0.26	0.025	4.5	2.5	2.3	390	1.2	—	—	—	—	—
to roll table motors																
0.068	1.74	0.153	0.153	3.33	0.293	0.126	—	—	6.1	370	0.9	8950	BC29.5/14	12.5	50	600
—	—	—	—	—	—	—	4.2	1.0	1.8	125	1.1	—	—	—	—	—
0.068	1.74	0.153	0.153	3.33	0.293	0.126	—	—	6.1	370	0.9	9450	BC29.5/14	12.5	50	600
—	—	—	—	—	—	—	5.6	1.1	2.5	170	1.2	—	—	—	—	—
0.1	1.72	0.2	0.2	1.76	0.205	0.079	—	—	4.4	140	0.7	4750	1	—	—	—
—	—	—	—	—	—	—	5.0	1.2	2.2	61	0.65	—	—	—	—	—
to roll table motors																
0.09	2.27	0.109	0.109	3.68	0.176	0.109	—	—	4.3	370	1.1	17700	1	—	—	—
0.121	1.16	0.236	0.171	1.7	0.344	0.04	5.0	1.8	2.6	2200	1.4	—	BC34/18	21	65	500
21.5 c.p.s.																
0.089	2.27	0.109	0.109	3.68	0.176	0.109	—	—	4.3	370	1.1	9850	1	—	—	—
—	—	—	—	—	—	—	4.0	0.9	1.8	475	1.4	—	—	—	—	—
18.7 c.p.s.																
0.078	1.44	0.17	0.17	2.7	0.32	0.112	—	—	2.5	300	0.7	8500	BC29.5/14	7.5	32	375
0.145	1.48	0.33	0.215	1.05	0.23	0.021	2.5	1.7	2.2	420	0.8	—	BC29.5/14	7.5	50	375
16.6 c.p.s.																
0.0575	1.21	0.1	0.1	2.17	0.18	0.06	—	—	2.65	140	0.3	4520	1	—	—	—
—	—	—	—	—	—	—	3.5	0.9	1.8	285	0.35	—	—	—	—	—

VARIABLE SPEED FREQUENCY CHANGER SETS

The variable speed and variable frequency Frequency Changer Sets consist of synchronous generators and D.-C. driving motors.

As regards method of D.-C. driving motor speed regulation, the Frequency Changer Sets can be grouped as follows:

1. Speed regulation by motor field weakening (Sets Nos. 17, 18, 20, 21, and 24);
2. Speed regulation by supply voltage variation and by motor field weakening (Sets Nos. 19 and 22);
3. Speed regulation by motor supply voltage variation (Set No. 23).

On request, the variable speed Sets can be furnished with tachogenerators and centrifugal switches arranged on the shaft of the set. The centrifugal switches are adjusted to trip at 10% overspeed. The tachogenerators are either of the permanent magnet Type M9T-7/10, or for 220 V D.-C. separate excitation, Type PH-5.

In the variable speed and variable frequency sets, for driving the synchronous generators are used D.-C. motors of individual design, series MI and PI00.

The series MI motors are with a compensating winding and are used in heavy duty sets (large excess currents, wide range speed regulation by field weakening). The maximum operating cut-out excess current of the series MI 2.25 motors is 2.5 times the rated value.

The series PI00 motors have no compensating winding, but have commutating poles and are used in small output sets of lighter duty. The maximum operating cut-out excess current of the series PI00-1.5 motors is 2.0 times the rated value. The RMS value of the motor operating current must not exceed the rated value.

Motor starting is either by gradual increase of the supply voltage or by means of starting resistances, provided the starting current is not more than 2.5 times the rated value. Heavier starting currents are admissible by mutual consent with the Maker. The motor is started with the full field current flowing.

Starting apparatus for the motors are not furnished by the motor Maker. All data required for ordering starting apparatus will be submitted on request.

The motors are furnished complete with direct hand operated field regulators, or remote hand operated, or motor operated. The driving motor is either 220 or 110 V D.-C.

The motor excitation, on Frequency Changer Sets with speed regulation by means of main field weakening only, can be independent of the external supply source or of the mains supply.

The motors on Frequency Changer Sets, with speed regulation by means of supply voltage variation are available with separate excitation. The field voltage of the motors is given in Table 4.

The series PI00 motors are designed with light series windings to provide a stable (slightly drooping) speed curve.

SET No. 17

The motor and generator are forced-ventilated and have an outside casing. Air intake is from the machine room and air discharge into the bedplate.

The Set is mounted on a bedplate, with ring lubricated pedestal bearings.

The synchronous generator can operate with current peaks up to 200% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is D.-C. Type MI250-333. Motor speed regulation is by main field weakening.

Generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type BC 24.5/12, 11.6 kW, 50 V, 1000 r.p.m.;
- 2) a squirrel cage induction motor, Type A71-6, 14 kW, 380 V, 970 r. p. m.

The Set is furnished complete with a field regulator, Type P3B-31B, for the motor field circuit, and Type PB 12L for the exciter field circuit. The Set is fitted with a Type PH-5 tachogenerator and a Type PLIC-1000 centrifugal switch.

SET No. 18

The machines of this Set are open self-ventilated, with pedestal bearings, mounted on a common bedplate.

The synchronous generator can operate with frequently recurrent current peaks up to 200% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is D.-C., Type MI190-333. Speed regulation of the motor is by main field weakening.

The generator has separate excitation, 110 V.

The Set is furnished complete with:

- 1) a tachogenerator, Type PH-5;
- 2) a field regulator, Type P3M-1A, to regulate the motor speed between 400 and 1000 r.p.m.;
- 3) a centrifugal switch, Type PLIC-1000.

SET No. 19

The machines of this Set are open self-ventilated, mounted on a common bedplate.

The FC4 1407-6 generator is suitable for the following duties:

- 1) 260 kVA, 78 kW, 480 V, 312 A, 1200 r.p.m., 60 c.p.s., $\cos \phi = 0.3$. Excitation: 95 V, 20 A. Intermittent;
- 2) 16.6 kVA, 5 kW, 96 V, 100 A, 240 r.p.m., 12 c.p.s., $\cos \phi = 0.3$. Excitation: 60 V, 13 A.

The generator has separate excitation, 110 V. The driving motor is D.-C., Type PI04-1. Motor speed regulation is between 750 and 1200 r.p.m. by main field weakening and from 240 to 750 r.p.m. by supply voltage variation.

The motor is furnished with a Type PB 24L field regulator.

The Set is fitted with a Type PH-5 tachogenerator.

SET No. 20

The motor and generator are forced ventilated and have an outside casing. Air intake is from the machine room and air discharge into the bedplate.

The set is mounted on a bedplate, with ring lubricated pedestal bearings.

The synchronous generator can operate with current peaks up to 200% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is D.-C., Type MI700-500. Motor speed regulation between 500 and 1200 r.p.m. is by main field weakening.

Generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type BC 24.5/12, 17.5 kW, 75 V, 1500 r.p.m.;
- 2) a squirrel cage induction motor, Type A71-4, 20 kW, 380 V, 1450 r.p.m.

The Set is furnished complete with:

- 1) a field regulator, Type PB 18L, for the exciter field circuit and Type P3B-41B for the motor field circuit;
- 2) a tachogenerator, Type PH-5;
- 3) a centrifugal switch, Type PLIC-1000.

SET No. 21

The machines of this Set are open self-ventilated, with antifriction bearings, mounted on a common bedplate.

The driving motor is D.-C., Type PI04-1. Motor speed regulation is by main field weakening.

Generator excitation is from a separate exciter set consisting of:

- 1) an exciter, Type BC 24.5/11, 8.5 kW, 45 V, 1000 r.p.m.;
- 2) a squirrel cage induction motor, Type A62-6, 10 kW, 380 V, 970 r.p.m.

The Set is furnished complete with:

- 1) a field regulator, Type PBM-1A for speed regulation of the PI04-1 motor;
- 2) a field regulator, Type PBM-1, for voltage regulation of the BC 24.5/11 exciter.

SET No. 22

The machines of this Set are open self-ventilated, with ring lubricated pedestal bearings, mounted on a common bedplate.

The synchronous generator can operate with frequently recurrent current peaks up to 150% of the rated value; the RMS value of the stator and rotor currents must then not exceed the rated value.

The driving motor is D.-C., Type MI270-750. Motor speed regulation between 750 and

1350 r.p.m. is by main field weakening, between 750 and 150 r.p.m. by supply voltage variation.

Generator excitation is from an exciter set consisting of:

- 1) an exciter, Type BC 24.5/11, 8.5 kW, 45 V, 1000 r.p.m.;
- 2) a squirrel cage induction motor, Type A62-6, 10 kW, 380 V, 970 r.p.m.

The Set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the exciter field circuit, and Type PB 24L for the motor field circuit;
- 2) a centrifugal switch, Type PLIC-1000;
- 3) a tachogenerator, Type PH-5.

SET No. 23

The machines of this Set are open self-ventilated, with pedestal bearings, mounted on a common bedplate.

The driving motor is D.-C. Type MI700-1000. Motor speed regulation is by supply voltage variation.

Generator excitation is from an exciter set consisting of:

- 1) an exciter, Type PH-100, 7.6 kW, 65 V, 1000 r.p.m.;
- 2) a squirrel cage motor, Type A62-6, 10 kW, 380 V, 970 r.p.m.

The Set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the exciter field circuit, and Type PBM-1 for the motor field circuit;
- 2) a tachogenerator, Type PH-5;
- 3) a centrifugal switch, Type PLIC-1000.

SET No. 24

The FC4 1508-12 generator is open self-ventilated with ring lubricated pedestal bearings, mounted on a bedplate.

The motor is protected self-ventilated, with antifriction bearings.

The generator has a sectioned stator winding with 36 ends brought out for interconnection change-over to allow operation at 6300/600 V at both frequencies.

The driving motor is D.-C., Type PI03. Motor speed regulation is by main field weakening.

Generator excitation is from an exciter set consisting of:

- 1) an exciter, Type PH-100, 7.6 kW, 65 V, 1000 r.p.m.;
- 2) a squirrel cage induction motor, Type A62-6, 10 kW, 380 V, 970 r.p.m.

The Set is furnished complete with:

- 1) a field regulator, Type PB 12L, for the exciter field circuit;
- 2) a field regulator, Type P3B-31B, for the motor field circuit.

SPECIFICATIONS OF VARIABLE SPEED AND VARIABLE

Set No.	Type of generator Type of motor	Output		Rating					Field rating				Short-circuit ratio
		kVA	kW	voltage, V	stator (rated) current, A	speed, r.p.m.	efficiency, %	cos φ	frequency, c.p.s.	current, A*	voltage, V	Short-circuit ratio	
Frequency Changer													
17	Generator ГЧ4 1507-6	400	160	400	578	1000	88.4	0.4	50	152	26	0.0026	
	Motor МП 250-333	183	53.2	133	975	333	82.3	—	16.6	—	—	—	
Frequency Changer													
18	Generator ГЧ4 1407-6 (C200-1000/400)	200	120	400	289	1000	89.6	0.6	50	19.1	98	0.0079	
	Motor МП 190-333	80	48	160	727	400	86.6	—	20	—	—	—	
Frequency Changer													
19	Generator ГЧ4 1407-6*	200	100	400	289	1000	87.5	0.5	50	19.1	98	0.0079	
	Motor П 404-1	—	110	500	240	1200	90.5	—	—	3.2	220	—	
Frequency Changer													
20	Generator ГЧ4 1510-6	1130	452	480	1355	1200	91.2	0.4	60	169	41	0.00148	
	Motor МП 700-500	470	188	200	2600	500	87.0	—	25	—	—	—	
Frequency Changer													
21	Generator ГЧ4 1404-6	202	60.6	97.5	1200	1500	78.2	0.3	75	116	24.5	0.000326	
	Motor П 104-1	67.5	20.2	32.5	390	500	72.7	—	25	—	—	—	
Frequency Changer													
22	Generator ГЧ4 1508-8	540	178	414	755	1350	85.5	0.33	90	113	28	0.00134	
	Motor МП 270-750	60	27	46	290	150	79.0	0.45	10	—	—	—	
Frequency Changer													
23	Generator ГЧ4 1508-12	600	450	800	434	1000	93.0	0.75	100	102	56	0.00645	
	Motor МП 700-1000	60	45	164	217	200	86.6	—	20	74	40	—	
Frequency Changer													
24	Generator ГЧ4 1508-12*	250	50	6300	23	1000	68.0	0.2	100	112	61	2.91	
	Motor П 103	—	85	440	216	500	89.5	—	50	—	—	—	

* Output at other duties—see description of Set.
 * The generator has a sectioned stator winding with 36 ends brought out to change over for 6300 or 600 V and

FREQUENCY CHANGER SETS

Table 4

Set No.	Short-circuit ratio	Reactances				Time constants			Flywheel effect of rotor, kg·m²	Cooling air required, m³/sec	Weight of Set complete, kg
		X ₁	X _d	X' _d	X'' _d	T' _d , sec	T'' _d , sec	T _a , sec			
Set 50 to 16.6 c.p.s.											
0.12	1.1	0.059	1.1	0.132	0.132	2.65	0.322	0.11	260	0.6	8800
Set 50 to 20 c.p.s.											
3.83	1.05	0.059	1.2	0.105	0.105	2.36	0.206	0.063	140	0.45	7300
Set 60 to 12 c.p.s.											
3.83	1.05	0.059	1.2	0.105	0.105	2.36	0.206	0.063	140	0.45	4950
Set 60 to 25 c.p.s.											
0.17	0.57	0.09	2.27	0.109	0.109	3.68	0.176	0.109	370	1.1	12600
Set 75 to 25 c.p.s.											
0.148	1.1	0.083	1.045	0.135	0.135	1.17	0.15	0.063	80	0.45	4600
Set 90 to 10 c.p.s.											
0.172	1.0	0.068	1.22	0.13	0.13	2.1	0.225	0.1	280	1.0	7650
Set 100 to 20 c.p.s.											
0.384	0.9	0.123	1.45	0.25	0.25	1.22	0.21	0.137	390	1.0	9800
Set 100 to 50 c.p.s.											
0.384	0.54	0.166	2.0	0.32	0.32	2.44	0.392	0.045	390	0.9	5950

allows operation at 6300/600 V at both frequencies.

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UNCLASSIFIED

GAMMA-RAY APPARATUS FOR INDUSTRIAL RADIOGRAPHIC SCANNING

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I. APPLICATION

Type ГП-Со-0,5-1 Gamma-ray apparatus is intended for laboratory or shop use as an industrial gamma-ray defectoscope, and may also be used for investigating highly inaccessible elements of structural assemblies.

II. CHARACTERISTICS OF THE APPARATUS

The source of gamma-ray radiation in the type ГП-Со-0,5-1 apparatus is a radioactive isotope of cobalt, Co^{60} , with a radiation hardness of the order of 125 mega R, and an radioactivity of 0.5 grams of eqvt. Radium. The radioactivity of the isotope decreases with time and reaches a value of one half of the initial activity in 5.3 years.

III. BRIEF DESCRIPTION OF THE DESIGN

The cobalt isotope is placed in a protective casing fastened on a standard. The casing proper consists of a pear-shaped cast iron body with a cast lead lining and is provided with a removable part in which the isotope may be carried to places or objects of limited accessibility.

When not in use, the isotope is stored in the central, spherical part of the casing but, when necessary, may be brought into the necessary working position with the aid of a control mechanism which has three fixed positions.

When working with the protective casing secured to the standard, two working positions are available, one of which corresponds to a position in which the isotope is located in the lighter section of the casing opposite a window which allows a conical beam of gamma-rays to be emitted (apex angle 53 degrees). The second working position corresponds to a position of the isotope in which it is brought out of the pro-

TECTIVE casing, at the end of the lighter section, to a distance of 50 mm from the casing along its longitudinal axis (where it emits a ring beam of gamma-rays).

When working with the removable, portable part of the protective casing, only the first working position is allowable for use and, in which case, the radiation is emitted as a ring beam of gamma-rays. The removable part of the protective casing is carried to places of investigation by the aid of two special and removable handles.

The standard has a base mounted on three, swivelling castors. The travelling carriage assembled on the column of the standard is controlled by a chain drive. It carries a horizontal tubular frame with a fork in which the protective casing is secured such that it may be displaced with respect to two mutually perpendicular axes.

The available angular displacements of the protective casing are:

- 135° about the axis of the casing,
- 0 to 90-degree displacement of the axis of the casing in a vertical plane.

The apparatus occupies 1.4x1.4 m of the surface and is 2.1 m height.

It weighs 165 kg.

IV. SHIPPING COMPLEMENT OF THE APPARATUS

The shipping complement of the ГП-Со-0,5-1 Apparatus includes:

- 1 The castor-mounted standard,
- 2 The protective casing, complete with a charge of the cobalt isotope and the fork
- 3 Isotope position-control mechanism
- 4 Two handles for carrying the removable part of the protective casing.

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