
50X1-HUM

Page Denied

STAT

STAT

MODERN LABORATORY SUPPLY LABORATORIUMSAPPARATE

STAT

FRIEDRICH GEYER

**LABORATORIUMSAPPARATE
KOMMANDITGESELLSCHAFT
ILMENAU-THÜRINGEN**

Fernsprecher / Telephone: 2683/2684

Drahtanschrift / Cable Address: GEYERLABOR ILMENAU

Kto. 4030, K.-Nr. 11 1030 Deutsche Notenbank Ilmenau

ÜBERSICHT:

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SUMMARY:

	SEITE/PAGE
INDEX	1 ... 14
LABORATORY FURNITURE	1 ... 14
MODERN LABORATORY SUPPLY	15
HEATING UNITS	16 ... 27
ELECTRIC SWITCHING DEVICES	28 ... 31
MOTOR DRIVEN APPARATUS	32 ... 39
ELECTROCHEMICAL APPARATUS	40 ... 47
MEDICAL LABORATORY APPARATUS	48 ... 63
SPECIAL TESTING APPARATUS	64 ... 74
METALLOGRAPHIC OUTFIT	75 ... 80
PORTABLE LABORATORY UNITS	81 ... 86
VACUUM ENGINEERING	87 ... 102
FACILITIES FOR RADIOACTIVE WORKS	103 ... 130

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Illustrations and specifications in this catalogue are as accurate as possible at the time of publication, but we are constantly modifying and improving designs, and hence they are subject to alteration in detail.

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Table listing various laboratory equipment and their page numbers, including items like Nadelventil, Netzschlußgeräte, and various pumps and filters.

INDEX

Table listing various laboratory equipment and their page numbers, including items like Aceton determination, Ammonia determination, and various containers and pumps.

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
For the equipment of laboratories of all types, we recommend our special construction, described on the following pages. This construction allows to combine laboratory benches and fume hoods by use of

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

Our system of construction has been approved of in hundreds of laboratories furnished by us of home and abroad, for having the following advantages:
Quick and simplified mounting and dismounting of benches and fume-cupboards, being combined of standardized parts

Low weight of only 175 kg per square meter of bench base, considering all tops, pipes, and structural parts.
Each carrier frame is provided with 4 adjustable levelling devices, by means of which any inclination of the floor can be compensated, even considering the

longest laboratory benches.

Pipes and connections are easy installable, for only a minimum of manipulations is necessary to take out the cupboard and drawer units standing under the bench top.

Adaptability to the conditions given by the building and by the purpose, which the laboratory is intended for


Abb. links zeigt den Zusammenbau eines Tisches aus Normbauelementen
Left hand illustration shows the mounting of a bench combined of standardized structural parts

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
movable, standardized and serial manufactured structural parts.

Our system of construction has been approved of in hundreds of laboratories furnished by us of home and abroad, for having the following advantages:
Quick and simplified mounting and dismounting of benches and fume-cupboards, being combined of standardized parts

Low weight of only 175 kg per square meter of bench base, considering all tops, pipes, and structural parts.
Each carrier frame is provided with 4 adjustable levelling devices, by means of which any inclination of the floor can be compensated, even considering the

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Laboratory in the National Research Institute Cairo
benches combined of standardized structural parts

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Through a considerable number of projects and equipments for large institutes, we have collected experiences, which are for the benefit of our customers, furnishing their laboratories. So we have worked out standards for room and building arrangements, the consideration of which guarantees a suitable and well-projected laboratory. (drawings are given on page 4/5)

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We have formed new preceptions for the arrangement of plumbing units and electric installations, which are to be considered before the details of the building are completed. Moreover we are planning the ventilation of laboratory rooms and fume extraction plants.

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We place at our customer's disposal our experiences regarding flooring and wall finishes, arrangement of heating and lighting, windows and doors, as well as precautionary measures against fire and accident.

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We shall be able to produce a better and more representative layout of the laboratory facilities the sooner your projecting orders reach us, considering the latest trends of the technique.

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
On the following pages of this list only Standard Structural Elements are mentioned. Apart from these constructions, we are continually delivering special equipments of all types, suitable for laboratories of factories, research and educational institutions.

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Ideal Measurements for Laboratory Rooms:

Distance between the middle of the windows 3,00 to 3,125 m
 Depth of the room 6,75 and 8,15 m
 Height of the room under consideration of German Industrial Standards DIN 5034 regarding the lighting of working rooms 3,50 to 4,00 m
 Height of window-sill 1,00 m

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

A: Laboratory well-equipped, due to its purpose.
 Sufficient accommodation on the walls, writing and titration possibilities in front of the windows and with convenient liberty to move.

B: Laboratory equipped inappropriately with insufficient space on the walls, inconvenient moving limit and lack of writing places.

C: Chemical College Room
 Left side: Fume hoods in a straight line with the benches, windows easy to access.
 Right side: Fume hoods built in the frame of the windows, in this way insufficient natural lighting and mutual hindering of the persons working in this part of the room.

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


LAB BENCHES WITH WOODEN CARRIER FRAMES

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

Reagenzgefäße
made of light metal, acid proof lacquered, with
fluorescent tube lighting, glass plates, plugs,
sockets and switches.
Number of sections according to choice.
Size of each section: a) 0,82 m, b) 0,94 m

Plumbing System
for gas, water, compressed air, steam, vacuum,
waste, with taps made of brass, acid proof var-
nished, sinks made of acid proof stoneware in
various dimensions.

Desktop Tops, 45 mm thick
in segments suitable for the carrier frames.
a) made of plywood with a cover of oak,
b) made of plywood with a cover of plastic
material.

Größe:
0,65 m wide, 0,55 m deep and 0,75 m high
*FNK with 5 drawers
*FNK with folding door
*FNK with 2 drawers and one folding door

Dimensions:
1,29 m wide, 0,55 m deep and 0,75 m high
*AK with 2 sliding doors
*BK with 2 drawers and 2 sliding doors
*EK with 5 drawers, writing plate, sliding doors
*EK with 5 drawers and a folding door
*EK with 5 large and 11 small drawers

*Refrigerator Units 998
with motor 3-phased 220/380 Volt, capacity of
the freezing compartment 130 liter. For tem-
peratures up to - 4° C for making ice cubes.

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

K1

Carrier frame
made of wood for taking up the cap-
board units, with 4 levelling screws, bolts
for the bench top, the socket, the plum-
bing system as well as to join units
throughout the line.

Dimensions:
K2 0,71 m wide, 0,59 m deep, 0,75 m high
K3 1,00 m wide, 0,59 m deep, 0,75 m high
K1 1,35 m wide, 0,59 m deep, 0,75 m high

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Bench Socket
made of wood to be fixed to the floor.
Bench sockets can also be used, corre-
sponding with the carrier frames.


Mikroskopieraal in einem Hygiene-Institut - Microscopical Laboratory in an Institute of Public Health

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LAB BENCHES WITH STEEL CONSTRUCTED CARRIER FRAME

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

Standard Stoneware Parts for Ventilations Arrangements for Laboratory
 Flame Cupboards
 acid proof, all round brown glazed, in various sizes
Exhaust Fans
 free from vibration, direct-coupled with three phase motor 220/380 Volt
Extract Sliders
 for regulating passage of air, made of porcelain
Flame Cupboard, Cabin Type
 suitable for mounting of large apparatus such as stillizers, desulfating columns etc.
 Dimensions
 1,40 m wide, 0,80 m deep, 2,50 m high

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

Flame Cupboard, Cabin Type
 suitable for mounting of large apparatus such as stillizers, desulfating columns etc.
 Dimensions
 1,40 m wide, 0,80 m deep, 2,50 m high
Cupboard Units MK for hoods
 made of 20 mm plywood, treated with acid resisting lacquer for setting in corner frames of wood or steel tubing.
MKK as above, but with inside coating of plastic material and junctions to the fume extraction system.

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Fume Hoods with sash window
 made of wood, glazed and with fluorescent tube lighting.
 made of steel, for working with perchloric acid and for use in tropical climates. Rear and side walls made of wire-glass.
 made of plastic material, for working with hydrofluoric acid or radio-active substances. Glass panes are made of transparent plastic material "Plexiglas".

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Dimensions
 MK3 1,40 m wide, 0,78 m deep, 1,60 m high
 MK4 1,00 m wide, 0,78 m deep, 1,60 m high
Laboratory Bench Tops
 of concrete in iron frames, covered with tiles, acid-resistant finished. Gaps for sinks, dripwaste and taps are taken into consideration.
Switch Cases
 installed according to VDE (German Association of Electro-Technicians)
 2 plugs and sockets with earth-contact
 1 switch for lighting

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
Carrier Frame K4 and K5
 as shown on page 7 but with switch-board for the electrical installation as well as for the front controlled service taps.
Dimensions:
 MK 1,20 m wide, 0,59 m deep, 0,75 m high
 K5 0,89 m wide, 0,59 m deep, 0,75 m high
Frame Sockle
 made of wood, to be fixed to the floor

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


Sink Unit II
with waste box, closed by a flap, 1 cupboard on the
left side with folding door. The unit is completed by
a draining board and a peg draining rod, 50x50
cm. A sink of glazed stoneware
is provided for, to be mounted on the wall

Desk Room Bench
made of steel tube with leveling feet, bench top
made of concrete coated with plastic material, having
hole cuts for 2 sinks, also made of plastic material.
F17 1,40 m lang, 0,78 m wide, 0,90 m high
Laboratory Writing Desks and Benches for Micro-
scopic Work on page 12

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Titrertisch Bench
for working in large scale
(described on page 33)

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




Standard Steinzeug Parts for Ventilation
Arrangements for laboratory
fume cupboards
acid proof, all round brown glazed, in
various sizes.

Exhaust Fans
free from vibration, direct-coupled with
three phase motor 220/380 Volt

Lüftungsschieber
for regulating passage of air, made of
porcelain

Tiefabzug
suitable for mounting of large apparatus
such as stillizers, distilling columns etc.
Dimensions
1,40 m wide, 0,80 m deep, 2,50 m high

Cupboard Units MK for hoods
made of 20 mm plywood, treated with acid
resisting lacquer for setting in carrier frames
of wood or steel tubing.






MNK as above, but with inside coating of
plastic material and junctions to the fume
extraction system.

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Fume Hoods with sash window
made of wood, glazed and with fluorescent tube
lighting,
made of steel, for working with perchloric acid
and for use in tropical climates. Rear and side-
walls made of wire-glass,
made of plastic material, for working with hydro-
fluoric acid or radio-active substances. Glass
panes are made of transparent plastic material
"Pirocyl"

Dimensions
MK3 1,40 m wide, 0,78 m deep, 1,60 m high
MK4 1,00 m wide, 0,78 m deep, 1,60 m high

Support frames
of concrete in iron frames, covered with tiles,
laboratory bench tops
of concrete in iron frames, covered with tiles,
acid-proof joints. Cops for sinks, drip-paste
and taps are taken into consideration


Switch Case
installed according to VDE (German Association
of Electro-technicians)
2 poles and sockets with earth-contact
1 switch for lighting
1 safety-switch with pilot-lamp for the exhaust fan

Carrier Frame K4 and K5
as shown on page 7 but with switch-board for
the electrical installation as well as for the front
controlled service taps.
Dimensions
K4 1,29 m wide, 0,59 m deep, 0,35 m high
K5 0,89 m wide, 0,59 m deep, 0,35 m high
Frame Sockel
made of wood, to be fixed to the floor

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Belagte Table
made of steel tube with levelling flanges, vibra-
tion-damping.


Dimensions:
PS1 1,40 m long, 0,60 m wide, 0,80 m high
PS4 1,00 m long, 0,60 m wide, 0,80 m high

Sink Unit II
with waste box, closed by a flap, 1 cupboard on the
left side with folding door. The unit is completed by
a draining board and a peg draining rack 50x50
cm. A sink of glazed stoneware
is provided for, to be mounted on the wall


Dark Room Bench
made of steel tube with levelling feet, bench top
of concrete coated with plastic material, having
hole cuts for 2 sinks, also made of plastic material
F17 1,40 m long, 0,78 m wide, 0,90 m high

laboratory Writing Desks and Benches for Micro-
scopic Work on page 12

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Titration Bench
for working in large scale
(described on page 34)

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




Histological Bench HT
Height of bench: 78 cm, bench top: 140x63 cm, black stained and varnished oak. On the left side 2 rows of 7 drawers, on the right 5 drawers and a removable plate to put up a microtome. The middle drawer has a safety central lock closing all drawers. The drawers contain numbered partitions for storage of 3910 microscope slides.

Bench for Microscopic Work XT
Height of bench: 78 cm, bench top: 140x63 cm, black stained and varnished oak. On the left side 2 drawers and a cupboard with shelves, on the right a writing plate can be drawn out.


Bench for Small Laboratories AT
2350 mm length, 670 mm wide, 900 mm high (1890 mm including back) varnished white. Equipment: bench top covered with tiles and acid-proof joints, reagent bottle shelf combined with stands for burettes and pipettes. On the left side a sink of glazed stoneware with back skirting and fittings, as well as a peg draining rack. On the right a shelf for filtering, weighing and microscopic work. Nickel-plated water tap, two-way bench stand, two-way gas tap as well as 2 plug-boxes.

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


ALL PROJECTS for fitting out laboratories, institutes and lecture halls demand experience in many technical disciplines. Your laboratories and the equipment of your institute will be improved by consultation with us and utilization of our designs.

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Laboratory Infrared Heater

Technical Data
 Voltage: 220 V
 Models of Infra-Red Lamp:
 small model BGW 80/250
 large model BGW 120/250
 Wattage: 250 W
 Largest thermal effect:
 inside of the rays at 8-10 cm diam.
 Base: 35x20 cm
 Over-all height support rod incl.: 700 mm
 Weight: small model 7.5 kg
 large model 9.0 kg

Application.

The Infrared Heater is suitable as a heating source for fractional distillations, fine-evaporations, extractions, melting, and for many other purposes. The resistor added to this device makes possible the adjustment of reproducible quantities of energy and also the fixing of an exact regulation of temperature without any regulating apparatus. This possibility is made easier when extracting with ether. Vapours of oil and water usually arising in water baths, oil baths, and vapour baths cease by using the Infrared Heater "ILS". If one intends to regulate temperatures by means of contact thermometers or other regulating devices, there will be obtained such a high constancy of temperature, as has been obtained by no heater ever before. In tests, this device was found to be useful for regulating temperatures in the range of $\frac{1}{10}^{\circ}\text{C}$ without any difficulty.

Description of this Device

A robust cast base carries the radiator. Upon the base are mounted:
 the circuit close, the circuit breaker, and the slide resistor the support for the radiator case and other implements.

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
Laboratory Infrared Heater

The radiator (heating lamp) itself is mounted into a case of light metal. It can be adjusted in height and is surrounded by a projector, reflecting the heating rays in one direction. The radiator case is covered by a set of porcelain concentric rings and can be accommodated to the largeness of the employed vessel (flask, tray etc.). In this case the radiator heats from below. After swinging the heater into the contrary direction the infrared lamp may be used for irradiating substances, lying in trays or tumblers, which are to be heated, dried or evaporated. The heating of liquids or substances in flat or round bottom trays, beakers, and glasses in contrast to other heating sources, will take place simultaneously in the whole of the vessel, because the absorption of the incident infra-red rays provokes thermal effects in the radiated medium on account of their great wave length. The thermal efficiency of the Infrared Heater is so high, that in comparison to a heating plate of the same wattage, there will be a much faster heating result.

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
This illustration shows a switching schema. The relay units NGN 55/NG 55 and the electronic relay ER 55 being described on p. 28-31 are suitable as switch apparatus together with adjustable contact thermometers.

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This graph explains the efficiency of both the heating sources on the same consumption of energy

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Laboratory Heating Units

A varnished case of cast aluminium serves as a covering for the face-ground heating plate. The external radiation of the plate is reduced to a minimum by a radiation plate and insulating air cushions.

with sand-bath tray 135 mm diam., single unit with a support rod Code: EH 1

Devices in series on a base plate
with 2 heating positions Code: EH 2
with 3 heating positions Code: EH 3
with 4 heating positions Code: EH 4

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Technical Data:

Stepped electric heating
switch step III 400 W
switch step II 250 W
switch step I 150 W
switch step 0 switched off



Wattage:
Number of heating positions 1 2 3 4
Wattage: 400 800 1200 1600

Voltage: 220 V
Connection to the network: socket with a protecting cover

Accessories for each heating position: support rod and sand-bath tray 50 mm of height 135 mm diam.

Base: with 1 heating position 20x20 cm
with 2 heating positions 45x25 cm
with 3 heating positions 65x25 cm
with 4 heating positions 85x25 cm

Height of the rod: 80 cm
Weight: 3 kg each heating position

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Laboratory Heating Plates

Heating by electricity. Side covering by oblong perforated plates for ventilation and protection from disadvantage of radiation. With 3-step switch.

Technical Data:
Voltage: V ≈ 220 220
Wattage: W 1600 2200

Three-step heating: W 600/1000/1600 1000/1200/2200
Temperature: max. + 400° C
Plate: ground gray cast iron
Base mm 350x260 480x320
Weight: approx kg 7.0 14.0
Size of the Plate mm 300x225 450x300
Codes: LP 30 LP 45

Sand Baths





Heating by electricity. With 3 step-switch
Technical Data: the same as LP 30 and LP 45
Weight: approx. kg 7.5 14.0
Plates mm 300x225 450x300
Codes: SB 30 SB 45

Sand Baths

Heating by electricity. With automatic regulation of temperature by means of contact thermometer and relay
Technical Data: the same as LP 30 and LP 45
Base mm 420x260 600x320
Weight: approx. kg 11.0 18.0
Size of Plate: mm 300x225 450x300
Codes: SB 30 R SB 45 R

Laboratory Heating Plates

Heating by electricity. Intended for installation into laboratory benches. With connecting terminals and not fixed three-step cylindrical switch, for the installation
Technical Data: the same as LP 30 LP 45
Weight: approx. kg 3.0 9.0
dimensions mm 300x225 450x300
Codes: LP 30 T LP 45 T

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Water Bath

W 20 Standard model, heating by electricity, with nickleed container of copper, handles, porcelain concentric rings, constant level and support rod. Three-stage switching.

Technical Data:
 Bath container inner width 200 mm
 depth 105 mm
 Wattage: 1100 W
 Voltage: 220 V \approx
 Base: 30x30 cm
 Overall height: 75 cm
 Weight: 5 kg
 Code: W 20

Water Bath
 W 20 R with automatic regulation of temperature by adjustable contact thermometer and relay.

Technical data the same as W 20
 Weight: 6 kg
 Code: W 20 R

Oil Bath
 O 20 R for temperatures up to +300° C, bath container of copper, with automatic regulation of heat


Weight: 1100 W
 Code: O 20 R

Test Tube Stand Water Bath Insert for 24 glasses 160x16 mm code: RG 24 for 16 glasses 100x30 mm code: RG 16

Water Baths in Series
 Heating by electricity. Especially suited to serial working, intended for several simultaneous operations in the common bath container or in separate containers, singly connectable and regulable in three steps. With porcelain concentric rings, support rod for each heating position and common constant level. Voltage: 220 V \approx

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Block Water Baths of copper With 6 Openings each of 150 mm of inner width, 6 sets of concentric Porcelain Rings

Nickeled Cover Plate
 Water Container 450x300x150 mm, 10 l capacity
 With Constant Level on a Four-footed Stand with Water Top


Block Water Bath
 Heating by electricity with 2 interchangeable heaters each with a heating power of 600 W, Voltage: 220 V \approx , altogether 1200 W, With a three-step switch, Base 45x30 cm Height 38 cm Weight 10 kg Code: W 6 E

Block Water Bath
 Heating by gas
 Dimensions the same as W 6 E
 Base: 45x30 cm
 Total height: 45 cm
 Weight: 6 kg
 Code: W 6 G

Block Water Bath
 Heating by steam, by means of a perforated steam tube in the water bath (so called sniffing steam)
 Base: 45x30 cm
 Overall height: 38 cm
 Weight: 8 kg
 Code: W 6 D

Block Water Bath
 Heating by electricity, intended for the installation into laboratory benches. Opening required for the installation 45x30 cm
 Voltage: 220 V \approx , Wattage: 1200 W, Three-stage switch
 Depth: 15 cm
 Weight: 6 kg
 Code: W 6 E

Block Water Bath
 Heating by steam in the water bath. Intended for the installation into laboratory benches. Opening required for the installation 45x30 cm
 Depth: 15 cm
 Weight: 7 kg
 Code: W 6 D-E

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

Apparatus for Nitrogen Determination acc. to Kjeldahl
 The Complete Apparatus Consist of:
A Heating Digestion Shelf
B Distilling Apparatus with vertical metal-condenser
C All the Necessary Glass and Rubber Parts
 The apparatus are delivered as Mikro-Models to hold 100 ml capacity flasks or Macro-Models to hold 500 ml capacity flasks. The units are gas or electrically heated on a voltage of 220 V_~.
 Base of the complete unit . . . 100x50 cm

Apparatus, Kjeldahl, Micro-Model
 with 100 ml flasks, Parts A, B, C. Electrically heated, voltage 220 V_~. Three step switch for each heating position
 for 4 determinations:
 Wattage 800+800=1600 W
 weight 16 kg Code: KE 41
 for 6 determinations:
 Wattage 1500+1500=3000 W
 Weight 20 kg Code: KE 61
 Heating by gas:
 for 4 determinations:
 weight 10 kg Codes: KG 41
 for 6 determinations:
 weight 14 kg Code: KG 61

Apparatus, Kjeldahl, Macro Model
 with 500 ml flasks, Parts A, B, C. Electrically heated, voltage 220 V_~. Three step switch for each heating position
 for 4 determinations:
 Wattage: 1400+1400=2800 W
 weight 21 kg Code: KE 45
 for 6 determinations:
 Wattage: 2100+2100=4200 W
 weight 30 kg Code: KE 65

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
Apparatus, Kjeldahl, Macro Model
 Heating by gas
 for 4 determinations
 weight 15 kg Code: KG 45
 for 6 determinations.
 weight 23 kg Code: KG 65

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Apparatus for determination of residual nitrogen
 acc. to FARNAS-WAGNER are described on p. 57

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Digestion Shelves for Kjeldahl Apparatus (A)
Micro-Model
 electrically heated, 220 V_~, for flasks 100 ml
 4 determinations, wattage 800 W, weight 5,6 kg Codes AE 41
 6 determinations, wattage 1500 W, weight 7,0 kg Codes AE 61

Gas heated, for flasks 100 ml
 4 determinations, weight 3 kg Codes AG 41
 6 determinations, weight 5 kg Codes AG 61

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Macro-Model
 electrically heated, 220 V_~, for flasks 500 ml
 4 determinations wattage 1400 W, weight 10 kg Codes AE 45
 6 determinations wattage 2100 W, weight 18 kg Codes AE 65

Gas heated, for flasks 500 ml
 4 determinations, weight 3 kg Codes AG 45
 6 determinations, weight 6 kg Codes AG 65

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Fume Ducts of Lead for the Digestion Shelves
 for 4 determinations KAB 4
 for 6 determinations KAB 6

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Electric Water Still WD 12, made of Copper. Intended for continuous working, fitted with automatic feeding of the boiling chamber by preheated cooling water. Boiling chamber and condenser tube are internally tinned. The apparatus is electrically heated by interchangeable heating elements. The still is equipped with a three-step switch, pilot lamp and a mains plug with protective cover.

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
PIPE FURNACE for Carius Tubes

Consisting of an aluminum cast block with bores to take the interchangeable iron protective sheaths for glass combustion tubes.

The heating elements are directly installed in the block, thus obtaining a fast and uniform distribution of heat all over the oven.

Technical Data:	
External dimensions:	length 950 mm, width 330 mm
Height:	560 mm
Wattage:	2,1 Kw
Switch steps:	220 V ~
	III = 2,1 Kw, II = 1,4 Kw, I = 0,7 Kw
Number of pipes:	3 1 2
Inner width of bores:	29 46 30
Protective sheaths outside/inside diam.:	28/24 44,5/39,5 38/33
Glass tubes, max. outs. diam.:	23 38 32
Useful length mm:	620 620 620
Codes:	RO 3 RO 1 RO 2

The regulation of temperature is effected by adjustable dial contact thermometer and relay T 30.

Working temperature: max. + 450° C

Equipped with splitter protective hood. Base: 960 x 350 mm Weight: about 34 kg

The metal block is insulated by glass wool, the case constructed of sheet metal with oblong slots. On account of the good insulation and heat radiation, the external temperature of the oven remains very low, even when working at maximum temperature.

The oven itself is mounted on an iron stand. The front can be closed by a swivelled steel plate. On the opposite side there are guide rails into which the splitter protective hood can be placed. This attachment is necessary, because it may happen that the combustion tubes crack on account of the increase of pressure within the pipes while working.


The dial contact thermometer is mounted on the top of the furnace housing. Another bore is closed by a screw and serves for setting in the control thermometer. The dial thermometer has a measuring pointer for indicating the temperature and a driving contact. By means of a separate key a further red pointer can be adjusted to the temperature desired.

The driving contact and the red pointer cause the switching impulse for the relay T 30. The latter is fixed on the right side of the furnace and underneath the plug with protective cover for the mains connection.

The electric equipment is installed acc. to the prescriptions of VED and is tested by the DAMV.

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Drying Oven and Hot Air Sterilizer ET 48
Electrically heated with automatic regulation of temperature by a rod expansion regulator and relay, fitted with a three-step switch, double-walled and insulated by glass wool. Air intake from the bottom. Adjustable air vent and glass control thermometer on the top of the apparatus. A special latch used on the door of the oven prevents explosions by releasing the door, should an undue pressure develop within the oven.

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Drying Oven, small size, ET 248 Electrically heated, double walled, in the operation chamber there are two perforated and removable shelves. Automatic regulation of temperature by an adjustable contact thermometer and a relay. Three step switch and pilot lamp for the heating are installed into the stand. Mains connection made possible by the plug of the apparatus.

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
All switch- and regulating devices are built in the socket, and the furnace is ready to work after connecting it to the mains. The socket is separated from the oven by an air cushion, so that switch, regulator, pilot lamp will not undergo a notable rise in temperature.

When opening the door, the hot inside of it will flap downwards thus avoiding molestation by heat. Moreover the door in this position forms a suitable platform for inserting the crucibles.

By adjusting the lever of the regulator the desired temperature is obtained and kept constant automatically. Independent from this a temperature indicator shows the respective temperature in the glowing room. The distribution of heat has been arranged to ensure a uniform temperature all over the muffle room. The maximum temperature of +1050° C is reached in nearly an hour.

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Adjustable Contact Thermometers

Adjustable contact thermometers with rotary magnet have been successfully installed in many of our laboratory devices as contactors for regulation of temperatures. They have been also used in the laboratory practice and for technical purposes. Such contact thermometers may be loaded up to three Watt only. By using our relay units and our electronic relay circuit breaking capacities up to 1500 W are obtainable. To obtain still higher capacities additional relays should be used.

Contact Thermometers with rotary magnet, adjustable

Measuring range	Immers. tube up to 300 mm in length	Immers. tube up to 600 mm in length
0...+ 50° C; 1/10°	KT 50.300	KT 50.600
0...+ 100° C; 1/10°	KT 100.300	KT 100.600
0...+ 250° C; 1/10°	KT 250.300	KT 250.600
0...+ 400° C; 1/10°	KT 400.300	KT 400.600

Protective Frames for contact thermometers with screw thread R^{3/4} / Immers. tube: 300 mm KTS.300 600 mm KTS.600 Range for screwing on, 60 mm diameter KTSF

Relay Units and Electronic Relay

They have the task of transferring small control impulses into large electric switching powers. These switching units are well approved in the laboratory practice. The apparatus are ready for service immediately and they are able to be composed quite easily into units. By taking into use these apparatus, the improvised and unclear conducting wires in the laboratory range can be done away with. Furthermore the protective earth, which answers to the prescriptions is guaranteed. Hereby the prescriptions, which comply to the working protection for laboratories are fulfilled.

Situated on the contact terminals of the apparatus are low voltages by which, absolute safety against accidents is obtained, and the connection of sensible contactors as for example instruments with contact switches, adjustable contact thermometers and also of photoelectric cells etc. is also made possible. The apparatus are delivered for working and closed circuit systems and with reversing installations for both of these methods of operation.

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Relay Unit NGN 55


Normally this device is forwarded for working current with a circuit breaking capacity up to 1500 Watt, in accordance to the load table.

When the adjusted temperature is obtained the NGN switches off the consumer. When specially ordered the device can be forwarded according to the closed circuit system, at which when the contact point is obtained, the device is switched on, for example in connection with minimum contact thermometers, signalling devices, the ventilator of an air conditioning plant etc.

Description of the Device

The device NGN is installed in a varnished housing and can be placed on benches. Located on the right front side is the overflowing socket for the possibility of the mains connection. A plug and socket with earth connection is located on the left side.

In the upper cover plate, the following are installed: The switch for switching off and on of the device itself and the switch for the consumer, the pilot lamps "Voltage" and "Heating". A mercury switch tube is located in the device for the operating current. By the application of a contact thermometer, NGN is switched off when attaining the adjusted temperature of the consumer, for example, heating devices, infrared heaters, motors, magnets or the like.

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

Technical Data:
 Voltage: 220 V ac
 Circuit breaking capacity according to the load capacity table
 Self consumption: ca. 3 W
 Contact load: (switching capacity) 220 V, 14 mA = 3 W
 Automatic connection through the already mentioned protective earth.
 Base: 25 x 12 cm
 Weight: 1.6 kg
 Code: NGN 55

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Relay Unit NG 55

This device is most suitable for working and closed circuit arrangement for alternating current. It obtains a circuit breaking capacity up to 1500 W according to the load index. This device contains a reversing unit. The relay unit requires alternating current.

A special point of interest of the NG 55 is the contactor circuit, which operates on a low voltage. A low voltage is situated on the terminals of the contactor which is derived from a transformer with protective network by the way of a rectifier.

Description of the Device:

NG 55 is enclosed in a varnished case which is able to be placed on a bench. On the right-hand side the socket is located for the mains connection. To the left there is a plug and socket with ground connection for the connection to the consumer. The following are installed into the top cover plate: The switch for switching on or off the device itself and the switch of the consumer being controlled, reversing switch for adjusting the device onto the operating current or closed circuit, pilot lamps for "Voltage" and "Heating". The connections for the contactor are also found on the front cover, which are accordingly marked. Installed in the inside of the casing are the following: The low tension transformer with protective network for driving the contactor circuit, the rectifier with smoothing device, the reversing relay, and the relay T 15 as a controlling device for the consumer.

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Load Index NG 55
the same as NGN 55 p. 29
but only for Alternating Current

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Technical Data:
Operating Voltage: 220 V ~
Circuit breaking capacity according to the load table: 1500 W
Self-consumption: 5 W
Contact load (reaction capacity): 8 V 20 mA = 0,16 W
Base: 25,12 kg
Weight: 25,12 kg
Code: NG 55

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Electronic Relay ER 55

The electronic relay contains a reversing unit for working and closed circuit system. The circuit breaking capacity obtains 1500 W, according to the load index and it is driven by alternating current.

By the application of the electron tube in the relay arrangement, the device is suitable for the most sensitive contactors.

The contactor operates practically without current, while the current flowing through the contactor circuit amounts to 10⁻⁴ A. The control of the electronic relay can be realized by a contact galvanometer or by a photo-electric cell. On account of the extreme min. load of the contact, the contactor is to be treated very cautiously.


Description of the Device:

ER 55 is enclosed in a varnished case. On the right side the socket of the device is located. On the left a plug and socket with ground connection for the connection of the consumer. Enclosed in the cover plate are the following: The switch to switch off and on the device and the consumer to be controlled, the reversing switch for adjusting the device on working or closed circuit, pilot lamps for "Voltage" and "Heating", and the connections for the contactor.

The sensitivity regulator is fixed up on the back of the device. By the use of this sensitivity regulator it is possible for instance to regulate the switching moment of a photo-electric cell and to adjust onto various intensities of light.

The following are built in the device: the electron tube with the corresponding resistor, the sensitivity regulator, the reversing relay and the relay T 15 with a 1500 W circuit breaking capacity proposed as a device to control the consumer.

Load Index ER 55
the same as NGN 55 on p. 29
but only for Alternating Current

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
Technical Data:
Operating Voltage: 220 V ~
Circuit breaking capacity according to the load index: 1500 W
Self-consumption: 20 W
Contact load (reaction capacity): 8 V 20 mA = 0,16 W
Tube complement: 1 electron tube UEL 51
Base: 25,12 kg
Weight: 25,12 kg
Code: ER 55

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Magnetic Stirrer with Electric Heating Plate

Advantages of the Magnetic Stirrer Technique: The modern method of the stirring technique makes it possible to stir in containers upwards from the bottom and at the same time heating is also carried out. The stirring bars of which are sealed into the glass or bedded into the rubber are freely movable in the stirring container. There exists the possibility of stirring in entirely closed containers and arrangements, in liquid gases, in pressure containers or in vacuum. In this device therefore, no impeller shafts and consequently no stuffing boxes, packings are included. Furthermore no corrosion of the motor on the chuck, or delay in boiling in the stirring container will take place.

Description of the Device: A case which is made of porcelain serves as a support for the stirrer, and at the same time supports the motor, as well as the on and off switch and the stop switch for the regulation of r. p. m. Adapted on the case is the cooling jacket and very close to this the hot plate. The support rod is fixed on the case, on which the heating plate and other parts, which are necessary for the examinations, are clamped.

The stirrer is equipped with a 220 V A.C./D.C. motor, the shaft of which bears a permanent bar magnet. The magnetic field is closed by the stirring bar which lays in the stirring vessel. When the chief magnetic bar is driven by the motor, then the magnetic field will rotate and with this the stirring bar which is situated in the beaker. By carefully regulating the r. p. m. it can be attained quite easily, that the stirring bar rotates without trouble and the desired stirring efficiency will be obtained.

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Technical Data:
 Stirring speed regulable: 150--1000 r.p.m.
 Hot plate: 150 mm diam., 500 W
 Operating voltage: 220 V ~
 Total connected load: 550 W
 Overall height: 350 mm
 Stator height: 300-370 mm
 Weight: 6,0 kg
 Code: R 50

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R 50

When the plate is heated, then the water cooling under the hot plate must be put into operation. A small amount of water is quite sufficient to cool the device.

The hot plate is fixed on the cooling cap with 1 mm space between them and then screwed securely on the support rod. The stirrer can be operated without having to be heated. The parts which are included in the case need no particular attention. The are however to be protected against dampness.

Accessories:
 Water Bath with constant level for R50
 Weight: 0,8 kg
 Code: WA
 Stirring Bar, of cylindrical shape, magnelized, glass sealed
 30x7 mm
 15x7 mm
 Code: RA 1
 Code: RA 2
 bedded into rubber, of rectangular stone shape
 27x17x10 mm
 Code: RA 3
 Permanent magnet bar for R50 (spare part) Code: RHM

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
Magnetic Stirrer R 57 simplified, small type in a case of varnished sheet with support rod and working plate made of plastic, without heating plate

Technical Data:
 Stirring speed, regulable 150--1000 r.p.m.
 Working plate, diameter 16 cm
 Operating voltage 220 V ~
 Wattage 40 W
 Height, working plate 20 cm
 Over-all height with rod 55 cm
 Base 16x16 cm
 Weight 4 kg
 Code: R 57

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Titration Bench T 50

with four titrating positions, with magnetic stirring from below.

The illumination by luminiscence lamps, which are installed into the working bench and in the top of the bench behind frosted glass, makes the carrying out of titrations easy.

Description of the Device:

The bench frame consisting of steel tube construction is fitted out with levelling feet. In the edge of the bench, the electric motor is located as well as 4 permanent bar magnets with the appliances to drive them.

On the front side of the edge the following are located: 1 switch for the electric motor, 1 rotating knob for the potentiometer for the continual regulation of the stirring devices as well as two switches for the illumination. All enclosed parts are very easy to access only by opening the front flap.

A white flagged bench board, with an opening for the illumination and covered with frosted glass, serves as a working place.

The top-piece of the table is also glazed by panes of frosted glass and bears the supports for the titration burettes, as well as 4 supply bottles each having a capacity of 5 l for the titration liquids.

The Standard Equipment:

4 Schellbach burettes with zero point adjustment and overflow, 50 ml $\frac{1}{10}$ ml with clamps, 4 white bottles, each 5 l capacity with connecting tubes, 8 cylindrical stirring bars, glass sealed, 30x7 mm, 4 stirring bars rubber embedded 27x17x10 mm large.

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Titration Support T 51 with magnetic stirrer

This device, in connection with pH meters and glass electrodes, serves specially for all titration work. An electrode support allows a rapid interchanging of the titration vessel.

Fixed on two support rods, are two titration burettes, a bench for titrating, provided with a spout of plastic, and also removable illumination. The magnetic stirrer, which is driven by the 220 V A.C./D.C. motor, is enclosed in the case, just as the speed regulation and the pilot lamp for the magnetic stirrer.

Standard equipment: 2 titration burettes 50 ml graduated $\frac{1}{10}$ ml, with stopcock with holder, 3 stirring bars glass-sealed 30x7 mm, 1 support for the electrodes.

Technical Data:
Voltage: 220 V ~
Wattage: 50 W
Base: 20x25 cm
Over-all Height: 45 cm
Weight: 13 kg
Code: T 51



Accessories:
removable illumination
Code: T 51 B

Titration Support TIT titration bulb of JENA glass

The double walled tempered titration bulb according to Gerstenberger and Sack makes it possible to carry out titrations in various temperatures. Titrations in an inert gas atmosphere can also be carried out. The support and the ground plate are equipped with motor and speed control.

Standard equipment:
1 titration beaker, JENA glass, with glass cover, colored reference electrode, of low resistance, bulb electrode 9221, burette of 50 ml capacity $\frac{1}{10}$ ml, condenser, 100 wire, thermometer -10+130°C, gas inlet pipe G 3, screwed joint for the cover plate, 2 spiral springs.
When requested the equipment can also be made of other electrode outfit.


Technical Data:
Voltage: 220 V ~
Base: 40x40 cm
Weight: 9 kg
Wattage: 50 W
Over-all Height: 70 cm
Code: TIT

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

Universal Shaking Machine
with motor, resistance, gearing, interchangeable shaking board and ground plate for 60 - 285 shaking movements per minute.
With: 1 board with a fixing device for 4 bottles each having a capacity from 250 to 1000 ml and 1 shaking board, for 1 bottle having a capacity of 1 to 10 litre

Technical Data:
Voltage 220 V ≈
Power 120 W
Base 78 cm
Height 95 cm
Weight 32 kg
Code SMF

Accessories:
Shaking board for 180 test tubes, 16 to 18 mm in diam., with 6 racks
Weight: approx. 4,0 kg Code: SMF/R

Shaking Machine according to Kahn
The technical data and the design correspond to that of the Universal Shaking Machine.
Equipm:
6 test tube racks each for 30 Kahn test tubes 75x12 mm. With 275/285 shaking movements/min.
Weight: approx. 27 kg Code: SMK

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Shaker according to Wagner
For 10 Stohmann flasks each 1 Lt. capacity. With electro motor and sliding resistance.

Technical Data:
Operating Voltage 220 V ≈
Power 120 W
Base and space required for the rotating movement Length ca. 78 cm, width 95 cm, height 110 cm
Weight 32 kg
Code SW 10

Accessories:
10 reducing caps for Stohmann flasks each 1/2 l capacity.
Weight: approx. 0,2 kg Code: SW 10R

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
Siebrüttel-Maschine
Sieve Shaking Machine

Description of the Device:

The structural elements of the device are mounted on a wooden plate, which is carried by a steel tube frame. The sieve cage is suspended on three pillars by using rubber connections. The shaking movement is developed by an adjustable eccentric load. According to the structure of the material to be tested, the shaking movement can be increased or weakened by adjusting the eccentric load.

The sieve cage is joined with the eccentric wheel by means of rubber/metal blocks and ball-bearings, through which the machine runs almost noiseless. It allows the setting-in up to 8 sieves by a cross-formed clamp with the receiver and cover having a diameter of 150 or 200 mm.

The driving results by an A.C./D.C. motor and V-belt. All rotating parts run in ball-bearings.

The shaking speed depends on the number of revolutions per minute of the electric motor, the running of which can be regulated by a sliding resistance which is fixed up under the working plate.

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Technical Data:
Diameter of sieve rings 150/200 mm
Motor 220 V ≈
Wattage 120 W
Height of the bench 80 cm
Base ca. 95x60 cm
Weight ca. 45 kg
Code SR

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MOTOR DRIVEN APPARATUS

Zubehör zur Siebrüttelmaschine

inserted by means of a counter ring. On each ring the meshes respect, the width of the hole is specified. The gauze can be changed quite easily.

Spare sieve linings and plates can be had at any time. The sets of sieves of 200 mm in diam., with rings of black plastic material (code "P") are delivered as a standard design. Rings which are made from cast light metal bear the code "LM".

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Accessories for the Sieve Shaking Machine Sieves acc. to DIN 1171 (German Industrial Standards): consisting of wire gauze made of brass or phosphor-bronze. Sieves acc. to DIN 1170: with perforated metal sheets.

They consist of individual sieve rings of black plastic, bakelite or cast aluminium, having a height of 55 mm, 150 or 200 mm diameter. Inside are conical gaskets, the sieve plates are neatly

Siebe DIN 1171 Sieves acc. to DIN 1171		Siebe DIN 1170 Sieves acc. to DIN 1170	
Preßstoffring 150 oder 200 mm Durchmesser Rings of black plastic "Bakelite" 150 or 200 mm diameter	Leichtmetall- güßring 150 oder 200 mm Durchmesser Rings of cast light metal 150 or 200 mm diameter	Preßstoffring 150 oder 200 mm Durchmesser Rings of black plastic "Bakelite" 150 or 200 mm diameter	Leichtmetall- güßring 150 oder 200 mm Durchmesser Rings of cast light metal 150 or 200 mm diameter
Idöhe Maschenweite width of meshes in mm	Idöhe Lochweite width of openings in mm	Idöhe Lochweite width of openings in mm	Idöhe Lochweite width of openings in mm
Kurzbez. - Code 0,04 P 0,045 P 0,05 P 0,055 P 0,06 P 0,07 P 0,1 P 0,125 P 0,16 P 0,2 P 0,25 P 0,315 P 0,4 P 0,5 P 0,63 P 0,8 P	Kurzbez. - Code 0,04 LM 0,045 LM 0,05 LM 0,055 LM 0,06 LM 0,07 LM 0,1 LM 0,125 LM 0,16 LM 0,2 LM 0,25 LM 0,315 LM 0,4 LM 0,5 LM 0,63 LM 0,8 LM	Kurzbez. - Code 2 P 3 P 4 P 5 P 6 P 7 P 8 P 9 P 10 P	Kurzbez. - Code 2 LM 3 LM 4 LM 5 LM 6 LM 7 LM 8 LM 9 LM 10 LM
Deckel „D“ und Aufhängesätze „A“ für die Siebsätze Covers "D" and Receivers "A" for the sets of sieves			
		Preßstoff black plastic "P"	Leichtmetall Light Metal "LM"
		Kurzbez. - Codes D 150 P A 150 P D 200 P A 200 P	Kurzbez. - Codes D 150 LM A 150 LM D 200 LM A 200 LM

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
4,5 l
300 V
250 W
130 x 16 cm
126 cm
ca. 105 kg
KM 2
KM2/P20
KM2/S18

Ball Mill

This machine consists of a case made of cast aluminium and two hardened porcelain drums which are sealed dust proof each having a capacity of 4,5 l. The mill is mounted on a special stand, which has rubber feet and two boards which can be drawn out to place things on and also a drawer. The electric driving is released by geared motor with a speedometer to perform reproducible crushing for two various speeds. Two perforated spouts hold back the porcelain balls in the crushing drums, when the substance, which has been ground, is emptied. 3 kg of porcelain balls, which are 20 mm in diameter belong to the equipment.

Technical Data:
2 Crushing Drums, each having a capacity
Operating Voltage: three-phase current
Wattage
Base
Overall Height
Weight
Code
Accessories:
Porcelain balls of 20 mm diam.
Steel balls of 8 mm diam.

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
Portable Blower KLB


The apparatus, which is able to be transported quite easily, can be applied in the laboratory for many purposes. The blower is provided with a safety valve and a stopcock on the pressure connections. An independent oiling device which is regulable, with a glass oil container and an oil inspection glass for the conveyance of the oil, is provided for the lubrication.

Technical Data
Blower Output
Pressure Capacity
Operating Voltage
Wattage
Base
Overall Height
Weight
Code

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Technical Data:
Operating Voltage: 220 V ~
Wattage: 50 Watt
Height: 16 cm
Base: approx. 30,30 cm
Weight: approx. 9,5 kg
Code: T 55

Electronic Titrimeter

Description of the Device:
The titrimeter combines a burette and an electrode support including a magnetic stirrer and an electrical measuring part. In connection with the common electrodes for pH measurements, the variation of the pH value of the solution which is to be titrated can be measured by means of electricity.

An absolute pH measurement is not able to be carried out with this apparatus, because the titration is based on an usual adjusted zero point and the measuring instrument is used for determination of the potential threshold. Other electrodes of a small resistance made of JENÄ glass can be applied, also our electrodes of silver iodide. This device is especially suitable for oxidation-reduction-titrations, because the potential differences are especially large here, while on the otherhand, halogen titrations require a little experience with these operations.

The following are enclosed in the front panel: The measuring instrument of 0-100 µA, with zero point adjustment, and sensitivity regulation, switch for the magnetic stirrer with a rotating knob for the regulation of the r.p.m., also 2 pilot lamps for the stirrer and the measuring component.

The following are located on the back side of the case:
The valve fuse of 500 mA, the ground bushing as well as the mains plug with protective cover, valves, rectifier and stabilizers etc. all are easy to reach behind the back plate of the device, which can be removed.

Standard Equipment:
Electronic Titrimeter T 55
Special Burettes, with selfabob strips and stopcock
One-armed clamps for the burettes with clamp holders
Magnetic stirring bars, glass sealed the same, bedded in rubber.
Holder for the electrode.
As a reference electrode: 1 silveriodide electrode
As a measuring electrode: 1 JENÄ glass electrode 9130

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Determinator of Ash
Contents "ASZ"
for Sugar Solutions

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This apparatus operates according to the principle of the Wheatstone bridge arrangement as a conductivity meter, especially to determine the ash contents from sugar solutions, particularly in the case of serial analyses. Moreover the apparatus can be applied for the measuring of conductivity in liquids.

This apparatus has the shape of a desk which is quite easy to survey, the enclosed measuring instrument is fitted out with a regulator for the sensitivity. Enclosed in the front panel are the following: Switch and pilot lamp for the mains component as well as a circular shaped scale which is subdivided in one guide dial and five measuring ranges, furthermore the corresponding measuring range switch for the measuring ranges.

The five measuring ranges contain divisions in "Percentage of Ash". The specific conductivity and the resistance value of liquids can be ascertained by the tables which are included.

Measuring Ranges:
Ash Contents: 3,0 to 0,0015 %
Resistance W: 205 to 380.000 Ohm
Conductance G ($G = \frac{1}{W}$): 4856 to $2,6 \cdot 10^{-6} S$
Spec. Cond. κ ($\kappa = \frac{1}{W \cdot C}$): 1553 to $0,84 \cdot 10^{-6} S$
Constancy of the Conductivity cell C = 0,32
Measuring according to the 100-ml-method, samples 5 g

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Determinator of Ash Contents in Sugar Solutions

To the equipment belong:

- 1 conductivity cell with platinum electrodes (constancy of the cell C = 0,32).
- 1 precision thermometer +15...+25° C: 1/10.
- 2 sealing stoppers, rubber tube and pinch cock for the outlet.
- 1 support for the conductivity cell and a set of connecting cables with pole shoes.
- 2 tables to correct the temperature.
- 5 conversion tables to ascertain specific values of resistance and conductivity.

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Technical Data:

Voltage by transformation efficacious at the bridge 18 V, 50 c/s
 Wattage Base including cond. cell Weight Code

Meßbereiche Measuring Ranges	Aschengehalt Ash Contents	spez. Leitfähigkeit Specific Conductivity	Widerstand Resistance
1	ca. 3-0,6 ‰	1533-336,0 (10-4S)	203,9-9520 Ohm
2	.. 0,8-0,12 ‰	448-67,2 (10-4S)	714-4761 ..
3	.. 0,16-0,025 ‰	89,6-14,00 (10-4S)	3571-22857 ..
4	.. 0,033-0,005 ‰	18,48-2,80 (10-4S)	17316-114285 ..
5	.. 0,0095-0,0015 ‰	5,32-0,84 (10-4S)	60150-380952 ..

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With the tables which are included, the conductivities of alkaline, acid and saline solutions can be ascertained. More detailed directions for operation are added to each apparatus (8 pages).

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Reference to Literature:

Dr. E. A. Lange: The electric conductivity of impure sugar solutions and the application of the same for the determination of the ash contents in products of sugar factories. Journal German Soc. of Sugar Industry 1931, P. 267/412.

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
Electrolyser Outfit ELY 2

This is a modern laboratory outfit, having two working positions for the rapid quantitative determination of metals in solutions. It contains electric hot plates and magnetic stirrers for the electrolytic vessels. Moreover the outfit can be used as a battery charging device and as a direct current source. The measuring instruments can also be applied for direct current measurements also outside of the electrolysis.

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The case is covered with a strong plate of light metal the enlarged rims or borders and outlets prevent spilled liquids to penetrate into the outfit.

A plate of transparent plastic protects the board for instruments and the switching grips. This plate can be pushed up, as long as the apparatus is operating.

The working plate supports two magnetic stirrers, the speed of which can be adjusted and two hard chromium-plated hot plates, which can be regulated in three stages.

The supports of the electrolyser outfit consist of hard rubber, the tops of which can be turned. The electrode supports can be adjusted in a forward direction, the angle between the two supports can be varied quite easily.

The electrolyser outfit is arranged for a mains connection of 220 Volt alternating current and for a battery current supply of 12 Volt.

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
ELY 2
The electric equipment exhibits devices which offer the following operating advantages:
 The unit can be used for charging batteries (of six cells) only by a simple manipulation of the switches. The instruments are able to be switched in the charging current to control the charging process.
 This equipment is also able to be used as a low-volt direct current source for all kinds of laboratory work.
 The electrolyzer current is coarsely tapped on a step transformer. The fine intervals between the individual steps of the transformer are accurately regulated by a fine resistance, to 1/10 Volt.
 In the range from 0-7 Volts, voltage steps from 0,5 to 0,5 Volt can be tapped and in the range from 7-10 Volt voltage steps from 1 Volt to 1 Volt.
 All the parts of the electric equipment and the stirring devices are enclosed in the case of the electrolyzer outfit. After removing the bottom plate and the side and rear sheet coverings, then all the installed parts are quite easy to access.

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Spezial-Bench ELY 2 T
 for the electrolyzer outfit. The bench top, which is covered with P.V.C., is set up on the right hand side as to be able to support the electrolyzer. At the left the raised top of the table is covered with tiles. There are five drawers - two of which are fitted with safety locks for the storage of the electrodes - and a cupboard with two flap doors. To right hand: a cupboard for the battery with two flap doors and a special switch for the single cells of the battery.
 Base: 1,4 m wide, 0,62 m in depth
 Height of the benches: 47 respect. 85 cm
 Weight: approx. 116 kg
 Code: ELY 2 T

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ELECTRO-ANALYSER

Technical-Data:
Electrolysis Voltage: variable 0-10 Volt
Electrolysis Current: variable 0-10 Amp.
Steps of Voltage: dependent on the voltage and the resistances in the electrolyzer vessels, in the range from 0-7 Volt in stages from 0,5 Volt in the range from 7-10 Volt in stages from 1 Volt 1/10 Volt by crank resistance
Fine Regulation: 220 Volt alternating Voltage, or 12 Volt battery
General Charging Voltage: 0-14 Volt
Source of Energy: 750 Watt
Measuring Instruments: 1 voltmeter 0-10 Volt, 1/10 graduation, 1 ammeter 0-10 A, 1/10 grad.
Control Lamps: for heating, stirring, electrolysis
1 Fuse: 1 Ampere
2 Heat Plates: each having three heating steps 250, 150, 100 Watt with regulation of speed
2 Magnetic Stirrers: with regulation of speed
Accessories: 2 stirring bars, bedded in rubber
Dimensions of the Case: 4 stirring bars, bedded in glass (without electrolyzer-support) 700x400x50 mm
Base with Bench: 1400x600 mm
Total Weight with bench: Net approx. 165 kg, gross weight 300 kg (2 cases)
Code: ELY 2

We recommend:

ELECTRODES Produced by: VEB FREIBERGER BLEIHOTTEN HALSBROCKE/SA.

Electrode made of platinum according to Winkler
 Size: The cylinder is 35 mm diam. and 50 mm high
 Over-all Height: 150 mm, gauze consisting of wire 0,12 mm in diam.
 Weight of the wire gauze electrode approx. 16 g
 at 0,25 mm wire gauze approx. 24 g
 Coil: 0,9 mm wire diam., 140 mm high.
 Wire gauze electrodes are also delivered in gold platinum
 Weight: approx. 6 g

Double Wire Gauze Electrode according to Fischer
 For rapid electrolysis, dulled
 External Cylinder: 38 mm diam., 50 mm high
 Internal Cylinder: 34 mm diam., 40 mm high, with glass armature
 Wire gauze of a 0,12 mm wire
 Weight: approx. 38 g
 Unperforated or perforated cylindrical electrodes, as also spiral electrodes and disc electrodes can be supplied.

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
Electrolytic Polisher
"Electropol" for Metallographic Purposes

The basis of the electrolytic polishing method is as follows: into an electric circuit, consisting of a current source, anode and cathode, a sample of metal is inserted and a circulating electrolyte produces the leading connection between cathode and anode. During the polishing, the sample itself is abraded and the metal ions from the sample,

smoothed through the emission of metal ions from the sample, which serves as an anode. A tubular V2A electrode is applied as the cathode through which flows the electrolyte fluid. The pumping system provides for the rapid rinsing of the sample by the electrolyte fluid.

The outfit consists of: a varnished case of sheet metal with carrying grips and detachable rear plate, the retractable polishing chamber made of plastic material with built in pumping plant driven by motor, the electrical equipment and the switch clock.

The current density can be adjusted according to the sample and to the electrolyte fluid. The voltmeter and the ammeter can be reversed from 10 to 100 Volt and 200 Volt, respect, from 1 A to 10 A, so that the current density and voltage, which are necessary for the etching can be read exactly with 1/10 of the polishing value.

The time of polishing is controlled by a switch-clock, having a dial-plate 0-60 sec. By means of the clock-switch the polishing begins. The electrodes are only under electric current as long as the switch-clock runs and as the adjusted time has been obtained. In special circumstances, when longer polishing times are required, the switch-clock can be bypassed by means of the "Selection Time" and "Stop Time" switches.

The current of the electrolyte fluid in the chamber can be regulated by a rotary resistor.

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
CHEMICAL APPARATUS

It has been concluded, the sample polishing current. This is carried out by "Polishing" to "Etching", without voltage again. Immer in connection with a metal able polishing current, after the end to a 220 V AC mains. ver, which separates both electrodes that no tension to earth can take

face which is being polished and d. The insertion of samples of sible, by removing the electrode

Data sheet Density Slopes Immer ammeter variable Motor Saltpetre nitrate Weight Code

Application of suitable electrolytes, the recommended are electrolytes on the acid, density 1.6

yes, ground samples of steel were polished

g Chamber made of plastic material, with in-sucking system and electric motor easily removed from the complete hour having to loosen a wire connector. A multiple-plug serves for connection to the polishing unit. The is provided with an overflow and an outlet which can be

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
Steel tempered superheated
 Polishing current 0.3 A
 Voltage 10
 Polishing time 10 s
 Cathode-gap 6 mm
 The surface of the sample was aftertreated with 1% alcoholic nitric acid

Steel ordinary tempered
 Polishing current 0.5 A
 Voltage 10
 Polishing time 10 s
 Cathode-gap 6 mm
 The surface of the sample was mechanically rough-polished and aftertreated with 1% alcoholic nitric acid

Austenitic steel
 Polishing current 0.22 A
 Voltage ca. 22 V
 Polishing time 25 s
 Cathode-gap 3 mm
 The surface of the sample was after-treated with V2A-etchant

High-speed steel
 Polishing current 0.55 A
 Voltage ca. 22 V
 Polishing time 15 s
 Cathode-gap 4 mm
 Sufficient corrosion of the sample was already effected by the electrolytic polishing

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„Electropol“

After the polishing procedure has been concluded, the sample can be etched by 1/10 of the polishing current. This is carried out by changing the switch from „Polishing“ to „Etching“, without having to adjust the voltage value again.

An inserted variable transformer in connection with a metal rectifier, delivers the adjustable polishing current, after the complete machine is connected to a 220 V AC mains. A further protective transformer, which separates both electrodes from the mains, ensures that no tension to earth can take place on the sample.

The space between the surface which is being polished and the cathode can be varied. The insertion of samples of various thickness is also possible, by removing the electrode support.

Technical Data:
 Surface capable of being polished 4 Amp/cm²
 Adjustable Voltages 0-200 V
 Measuring Range of the Voltmeter 0-100 V
 Measuring Range of the Ammeter 0-10 A
 Input of the Pump Motor 20 W
 Operating Voltage 220 V
 Weight 3.3 kg
 Dimensions of Case 750 x 550 x 550 mm
 Code E POL

As regards the possibilities of application of suitable electrolytes, the referring literature gives this information. Recommended are electrolytes on the basis of perchloric acid or of alcohol in various constitutions, e.g. 35% perchloric acid, density 1.6 24% ethyl alcohol 12% glycerine
 By the application of these electrolytes, ground samples of steel were polished within 20 seconds.

Polishing Chamber


which is made of plastic material, with installed pumping system and electric motor can be easily removed from the complete outfit, without having to loosen a wire connection. A multiple-plug serves for coupling the chamber to the polishing unit. The chamber is provided with an overflow attachment and an outlet which can be closed.
 Code: E POL

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
Paraffin Embedding Unit B 25

Standard model, Made of brass, double walled for waterfilling. The working chamber is provided with corrugated inner walls to enlarge the surface of the heat transfer, and also with two removable shelves. The outside covering consists of plated white plastic material. The door is heat insulated. The unit can be fixed to wall by the aid of two loops. The ventilation of the working chamber is caused by an adjustable tube and by the thermometer tube which is supplied with openings.

The electric heating has a three step switch. The automatic temperature regulation is carried out by an adjustable contact thermometer and a relay. Furthermore it is provided with a pilot lamp, a switch and a mains plug with cord.

A water level indicator with float, a drain tap and a control thermometer complete the equipment.

Paraffin Embedding Unit B 25 T with Drying Compartment

The construction of this unit is in accordance to that of the Standard Model B 25, this unit however is additionally equipped with an incorporated drying compartment 14x16x24 cm. which is closed by a glass door.

Technical Data:	B 25	B 25 T
Working Temperature	+ 30...+ 70° C	
Working Chamber	25 x 25 x 25 cm	
Capacity of Water	14 l	
Operating Voltage	220 V ~	
Wattage	350 W	
Switch Steps	250/200/150/0 W	
Wall Surface	40 x 70 cm	
Weight	22 kg	28 kg
Code	B 25	B 25 T

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Bacteriological Incubators

The Standard Models B 30 and B 40 are fitted out with one-leaf doors, whereas the larger models B 60 and B 80 with folding doors. However the working chamber is closed extra on the inside by glass doors.

The incubators have electric heating with coarse regulation for the heat by a three step switch, and automatic temperature regulation by adjustable contact thermometer and relay.

The working temperature attains values in the range from + 30...+ 60° C, the temperature constancy amounts to ± 1/2° C. The electric heating elements can be changed quite easily.


Besides the electric heating, the Special Model B 40K contains additionally a water cooling with automatic controlled water valve. The working temperature can be adjusted in the range from + 15...+ 60° C. The temperature of the cooling water must be 5° C less than lowest required working temperature. The regulation of the temperature proceeds automatically. For operation above 30° C the water valve can be switched off.

The incubators B 30P and B 40P are to be heated with petroleum. The automatic regulation of temperature between + 30 and + 37° C proceeds by a special thermostat.


General Equipment:

The bacteriological incubators are made of copper. The walls of the working chamber are corrugated to increase the reflection of heat. The water jacket is externally insulated. The outer walls are covered with white plates of plastic material. All incubators have adjustable tubes for ventilation.

The following belong to the equipment. Temperature regulation with contact thermometer, a control thermometer, a water level indicator and a water draining cock. In the working chamber are located two perforated removable shelves.

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
Coagulator B 40 S for Inspissating Sera

The construction and the external equipment of the coagulator correspond to that of the Standard Model B 40. The working temperature can be adjusted in the range from +30...+85°C. The apparatus has an enlarged insulation, corresponding to the raised working temperature, electric heating and automatic regulation of temperature.

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Slanting Racks for the coagulator, for the inclined arrangement of the culture tubes 180x18 mm, made of nickel-plated brass 21x23 cm large.
Code B 40 S/EG

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Bacteriological Incubator B 41 E. This is a portable model especially suited for expeditions. It is heated by petroleum and is provided with automatic temperature regulation, by means of a special thermostat. The four feet can be folded up. The inside of the apparatus consists of copper and is covered by brown plates of plastic material, framed by copper. A control thermometer in a wooden case also belongs to the equipment. When being transported, all individual parts, including the petroleum heating device, and as also the frame can be put into the working chamber.

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
Technical Data:

Working Temperature +30...+37°C
Working Chamber 42 x 32 x 22 cm
Outside Dimensions 55 x 42 x 37 cm
Weight 24 kg
Code B 41 E

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Cooling Plate

To cool fresh cast preparations of paraffin. The plate is made of nickleled brass and is fitted out with hose nipples, so that the tap water can flow through. The outer covering consists of white plates made of plastic material. Water level and levelling screws make it possible to mount the plate in a horizontal position.

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Warming Plate


Made of nickleled brass. Is used to dry micro slides, cover glasses and microscopic preparations. The water filling is electrically heated and the temperature can be regulated automatically. The device has a working plate made of high-polished nickleled brass. To this belong: a water level and levelling screws to mount the device in a horizontal position, water filling nipples and a water level tube. The side walls are covered with white plates of plastic material.

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Technical Data:

Working Plate Capacity of Water Voltage Wattage Weight Code

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Water Bath according to Tröster


Serological water baths for inactivations, electrically heated. These large thermostatic water baths are suitable for tests in series and in bulk. The baths consist of double insulated walls, and outside they have a lasting varnished case of sheet metal with metal mountings of nickleled brass. On top the baths are closed by glazed coverings, which are put up in slanting positions, so that the condensed water can flow off. By the aid of holding supports, the glazed covers can be held in an open position. For the capacity of water, a level regulator and an outlet cock are provided for. The proper test tube racks are nickle plated and stand on a sieve plate in the bath. The working chamber is tinned. The electric heating proceeds through the heating elements which are easily interchangeable, and is fitted out with a three-step switch for coarse regulation and a pilot lamp. The automatic regulation of temperature is effected by an adjustable contact thermometer and a relay. The baths can be delivered in two different sizes.

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
WT 144

WT 144
External Dimensions mm 680 x 235 x 380
Working Chamber mm 638 x 200 x 180
Capacity of Water 20
Test Tubes 144
Test Tube Racks 6
Working Temperature +30...+60°C
Voltage 220 V ≈
Wattage Watt 1200
Base approx. mm 880 x 250
Weight approx. kg 21

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
Water Bath for Serological Tests

heated by electricity for inactivations in series and in bulk. The bath has double insulated walls and nickeled frames. Test tube racks in one or two rows can be hung up or removed very quickly. A level regulator and an outlet tap are provided for the water charge. The electric heating is carried out by interchangeable heating elements which have a three-step switch for coarse regulation. An adjustable contact thermometer with relay is provided for the automatic regulation of the temperature. Two pilot lamps make it possible to control the current supply and the heating.

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Technical Data:
 Dimensions 630x430x290 mm
 Working Chamber 570x355x100 mm
 Capacity of Water 20 l
 Test Tube Sticks 15
 Working Temperature +30...+57°C
 Voltage 220 V ~
 Power 1200 W
 Base 65x45 cm
 Weight 25 kg
 Code **WR 725**

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Water Bath according to Griffin for gas heating.

Copper bath basin with two openings of 130 mm in diameter and 1 set of concentric porcelain rings. Useful bath depth is 100 mm, and is suitable for the insertion of glass or porcelain beakers. On the cover is a tube with the proper control thermometer in the range from 0-100°C.

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Dimensions: 28 cm Height
 33 cm Length
 18 cm Width
Weight approx. 4 kg
 Code **WGR**
Accessories: Double flame
 burner
 Code **WGR/B**

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Special Baths for Melting Paraffin for Histological Purposes



Water Bath WOB, Model München, electrically heated, with a copper jacket and a nickelated cover plate. It has 10 hard soldered-in copper beakers to take up the proper cups, and also a water jacket with levelling and drain off cock.

The following parts belong to the equipment:
 6 cups 100 ml and 6 cups 200 ml with handles.

The electric heating is carried out by two easy interchangeable heating elements, the automatic regulation of the temperature is effected by means of an adjustable contact thermometer and a relay. The following parts are on the front: three-step switch for the coarse regulation and two pilot lamps for tension and heating.

Water Bath WLM, Model Leipzig Heated by electricity, with a copper jacket and a nickelated cover plate. It has six soldered-in hollows to take up six semi-globular containers of copper, each having a capacity of 240 ml. The following parts belong to the equipment:
 6 containers 240 ml with handles, 1 level and 1 drain cock.

The electric heating is carried out by the way of two easy changeable heating elements, with a three-step switch for coarse regulation. Furthermore, a switch as well as two pilot lamps for the current supply and the heating are provided for. An adjustable contact thermometer with relay serves for the automatic regulation of the temperature.

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Hot Water Funnel
On tripod, electrically heated, with a copper jacket and easy changeable electric heating element, furthermore with regulation for the water level, with inserted glass funnel, 130 mm in diameter.

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Steam Sterilizer, according to Koch
Proposed for flowing steam. The heating protective jacket is varnished and has a nickel-plated double walled cover with control thermometer.
A reflux for the condensed water is warranted by the shape of the cover. A water level regulator, water level indicator and drain-off cock complete the equipment.

The working chamber is tinned. The sterilizer is heated by electricity, it has a three-step switch and a pilot lamp.

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
Device for the Determination of Ammonia (Micro Method)
Especially for rest-N-determinations according to Parnas-Wagner, improved nitrogen determination by Kjeldahl. It can be applied for serial examinations, one determination is carried out in 5-6 minutes. The device operates with steam distillation. The glass parts have standard ground joints.

Parts of the Device

Micro distilling bulb, made of JENA glass of 100 ml capacity, condenser with bent piece and glass-spiral, reflux condenser with cock, T-shaped tube with funnel and cock, round bottom flask 2 l for the development of steam, beaker 250 ml, canted flask 100 ml. Metal parts: 2 pairs of springs, 2 supports with three clamps and clamp holders, 1 radial gas burner.

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
Heating Digestion Shelves, micro size for 100 ml flasks
electrically heated 220 V ≈
for 4 determinations (800 Watt) Code AE 41 Weight 5,6 kg
for 6 determinations (1500 Watt) Code AE 61 Weight 7 kg
with gas heating Code AG 41 Weight 2,5 kg
for 6 determinations Code AG 61 Weight 4,5 kg

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
Apparatus for the Determination of Acetone, according to Krainick, Micro size for 6 determinations.
 The paraffin bath is electrically heated by interchangeable heating elements and has a three-step switch for coarse regulation, pilot lamp and mains plug.
 The automatic regulation of the temperature is carried out by an adjustable contact thermometer and relay.
 The separate ice bath has a cover and a water drain cock.

The outfit: 6 sets of glass parts, each consisting of Micro distilling bulb made of JENA glass with funnel and cock, Condenser-bridge, with JENA immersion filter-stick and glass frit, cylindrical graduated tube. The glass parts have fastening hooks, fastening springs and rubber connections.

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
Barcroft-Warburg-Apparatus for blood-gas-analysis and for the measurement of cell respiration.
 The thermostatic water bath has a stirrer and is heated by electricity by means of interchangeable heating elements, a three-step switch and pilot lamp. The regulation of temperature is effected by an adjustable contact thermometer with relay.

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Barcroft-Warburg-Apparatus

The shaking mechanism for the reaction flasks consists of an electric drive with motor and regulation for the speed of the shaking process, as well as a coupling to manipulate the manometers while the motor is running.

The 8 manometers are connected by standard ground joints with the reaction flasks, a holding attachment for 8 manometers belongs to the equipment.

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W 3 III

Accessories:

Manometer acc. to BARCROFT, on a wooden frame, with mirror scale, and with two reaction flasks, for smaller charges (0,1)

Spare reaction flasks, with standard ground joint

Manometer acc. to BARCROFT-ROBERTS, mounted on a wooden board, with two reaction flasks for larger charges.

Spare Manometer with standard ground joint

Manometer acc. to WARBURG, with stop-cock, graduated test tube and standard ground joints, without reaction flask.

Reaction Flask, acc. to WARBURG with ground in glass stoppers and standard ground joints,

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


van Slyke Apparatus for Blood Gas Analysis,
volumetric method.

The apparatus consists of a base with a frame made of wood, with lighting, glass parts ready for work, and the shaking mechanism.

The following parts belong to the equipment:

Graduated manometer tube, tubing system, regulating and ventilating cocks, an elastic rubber connection, reaction pipette with water jacket, thermometer and mercury levelling vessel, and also the shaking device with electric motor and sliding resistance.

Technical Data:
Operating Voltage 220 V
Wattage 50 W
Base 65 x 45 cm
Over-all Height 125 cm
Weight, approx. 25 kg
Code SLVM

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

van Slyke Apparatus for Blood Gas Analysis,

volumetric method. The glass parts and the shaking device are attached to a stand. The reaction vessel and the absorption pipette are suspended and able to be moved.

Gas burette of 150 ml and graduated 0-40 ml: 1/10, reaction vessel 0-10 ml: 1/10, absorption pipette, levelling bulb, shaking mechanism, with adjustable connecting gear, as well as electric motor and sliding resistance for the regulation of the r. p. m.

Technical Data:
Operating Voltage 220 V
Wattage 30 W
Base 40 x 30 cm
Over-all Height 80 cm
Weight, approx. 6 kg
Code SLV

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Macerating Plant,


- 1. the proper Macerating Plant MZ
- 2. the Degreasing Apparatus ME
- 3. the Bleaching Apparatus MB and MS
- 4. the Shaking Board

In order to obtain clean anatomical preparations, JENA glass is preferred for the material of the glass containers, which has a capacity of 150 litres. The glass containers are suspended in steel frames which are especially constructed for this purpose and are free from stress. They are heated by the aid of glycerine baths. The copper glycerine baths are in close contact with the glass containers with special armatures and have heating elements.

The regulation of the temperature is effected by an adjustable contact thermometer and relay.

The apparatus has proved its efficiency in numerous pathological institutes. The individual layout parts consist of the following structural components:

- Bleaching Apparatus „MB“** (Illustr. pag 63) is proposed to bleach macerated and degreased bones. Steel frame with three legs, 1500 mm high, Base 850x850 mm, with a two part flanged ring and hinged cover, brake cylinder, rubber padding and holding screws, as well as the necessary holding devices for all parts of the apparatus. **MBa** Container made of JENA glass, capacity 150 litres, height 108 cm, diameter 51 cm, width of the neck 41 cm **MBb** Cap made of JENA glass with 1 tube, height 35 cm, diameter 41 cm **MBc**
- Degreasing Apparatus „ME“** to degrease macerating preparations Steel frame with 3 legs, 1500 mm in height, Base 850x850 mm, with a one piece flanged ring, rubber padding, fixing screw and holding device. **MEa** Container made of JENA glass **MEb**

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
MEDICAL LABORATORY APPARATUS

10-“ME“

Cap made of lined copper, with a doubled walled flange cover, adapted to this is an open worm condenser.
Fixed to the cover the following are located
Copper condenser with a danga-over cock "Dissolvent Distillation and Discharge Return", each containing 1 water feeding valve for the condenser and the worm condenser, safety tube, hand grips and detachable distillation tube.
Drain cock made of glass, with 32 mm boring
1 flange connection
Glycerine bath
An additional safety contact thermometer, intended for a constant cut-off temperature of +120° C.
2 two-part perforated bottoms made of V2A plate.


Macerating Apparatus "MZ" (Illustr. pag. 61)

for the preparation of bones by extracting the remains of the meat.
Steel frame, 1500 mm in height
Container made of JENA glass, having a capacity of 150 litres
Cap made of JENA glass with tube, having a diam. of 41 cm and a height of 35 cm.
The following are suited to the cap. 1 control thermometer 0 to 100° C, length of the shaft 90 cm with rubber stoppers.
1 bent glass tube with rubber stoppers for the connection of a suction line.
1 glass tube with ground fit, 100 cm in length for venting.
1 drain cock, JENA glass with 32 mm boring and glass flanges
1 connector for the flange 32 b 405.
1 glycerine bath made of copper, diam 56 cm, height 48 cm, full capacity approx 40 litres, with special armature for the carrying out of the flowing-out connection through the bottom of the bath, with two heating elements of 400 and 1200 Watt and stuffing box screwing to convey the current, inserted contact thermometer 0-100° C, with socket, length of the immersion tube 45 cm, armature for the water level, filling container of 4 litres capacity with cock, drain cock.
2 two-part perforated bottoms made of V2A plate
1 water jet injector made of artificial material with built-in return valve, connecting sleeve and 1 m of vacuum tubing.
When requested, a pit which is made from JENA glass with the necessary flange connections to hold the sludge, can be delivered, as illustrated on p. 61.
However when the frame is put 35 cm higher, one must be sure to put an corresponding pit under the apparatus.

Switch Board "MS"

65 cm in width 90 cm high 15 cm in depth.
Consisting of a varnished wooden frame with front panel made of Melakort and serving lid for connecting terminals.
The following are on the front
4 switching relays T 15, switching capacity 1.500 Watt, 220 V~ with 4 pilot lamps.
3 automatic safety fuses 15 A.
3 three-stage switch for the heating, which is marked, "Maceration" "Degreasing" "Bleaching"
as well as three pilot lamps.

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
Baader Test Apparatus

The Baader Test Apparatus serves for the determination of sludge forming tendency of insulating oils and for testing the saponification number.

The oil samples are kept at a constant temperature of +95° C and during the period of 48 hours, glass or metal-coils, which are suspended on glass stirrers, are dipped 72000 times (25 ± 1 stroke/min) in the test tubes, containing the oil to be examined.

Thereby thorough mixing of the oil samples with air, also a contact of all oil particles with the metal of the immersed coils (iron, copper or lead) is attained.

Coils of glass can be used, when the influence of the temperature and air only is to be tested.

The oil samples are tested with regard to: turbidity, change of colour, settled water, sediment, insoluble sludge and acid reaction. The ageing inclination is taken into account by the saponification value (Vz Cu) of the copper. Please refer to the German Standards DIN 51554 and 53659.

The Baader Test Apparatus may also be used for testing the attack of oil on insulating material, for example.

The temperature constancy of 1/10° C which is required for the Baader Test Apparatus is obtained with the double thermostat which consists of a water bath and an oil bath.

An immersion-stirrer, which runs in the oil bath, is driven by an electric motor with the necessary regulating resistance. Noiseless running is obtained by an intermediate pinion gear and toothed wheel made from plastic material (Novotex). A counter records the number of strokes.

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Description of the Apparatus:


- A Base to take up B.
- A2 A second base plate to take up the support stand when assembling the apparatus.
- B Support stand.
- C Immersion-stirrer with drive.
- D A double thermostat consisting of 1 water bath, 1 oil bath to transform the heat, with electric heating and control device.
- E Test tubes (2 items) with standard ground joint, Liebig condenser and glass stirrer.

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Technical Data:

- Test Tubes: for eight simultaneous determinations, (Standard Ground NS 40), each sample requires 60 ml.
- Voltage: 220 Volt ± 5% Wattage: 700 Watt
- Thermostat: electrically heated water bath, which has a capacity of 9.5 litre, 600 watt, with water level and drainage.
- Regulation of the temperature: by an adjustable contact thermometer and relay 115.
- Temperature Range: +20...+100° C.
- Oil Bath: suspended in the water bath as a transformer of heat, oil required approx. 3.8 litre.

Means of Cooling: by tap-water. 2 metal reflux condensers on the water bath, 8 Liebig condensers on the test tubes.

Immersion-stirrer: with 5 glass stirrers.

Drive: electric motor 70 Watt, speed reduced by worm gear to 131 r.p.m.

Regulation: sliding resistance, switch and pilot lamps.

Control of stroke: counter with zero adjustment.

Accessories: plate with stand to put away the immersion stirrer, 1 ceiling socket.

Over-all Height: approx. 75 cm. Base: 52x53 cm. Weight: approx. 45 kg. Code: BT 8

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Apparatus for the Determination of Evaporation Residue of Fuels similar to the ASTM and GOST Standards

The apparatus is mounted on a wooden base and consists of the following parts:

- 1 Small blower with electromotor delivering compressed air of approx. 600 mm pressure
- 2 Gas washing bottles, to purify the compressed air
- 1 Flowmeter to control the passage of air, measuring range 10-60 l/min
- 1 Oil Bath with electrical heating and automatic temperature regulation by an adjustable contact thermometer and relay, a device to take the oil samples as well as a copper serpentine pipe, for the constant heating of the air current, with a two-piece air outlet.

Two test tubes, 100 ml capacity, filled with the sample are inserted in the container which is fixed to the cover of the apparatus. This container immerses into the electrically heated oil bath which is brought to a constant temperature of +160°C. Located in the oil bath is a copper serpentine pipe with a total length of 4.5 m and a diameter of 6 mm.


10 to 60 litres of air per minute are conducted through this serpentine pipe. The air passing through attains the temperature of the oil bath, and is supplied to the test tubes by the way of a distributor head located on the cover.

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Technical Data:

Working temperature Number of tests Oil required Voltage Wattage Base Height Weight Code

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Elutriating Apparatus

The elutriating apparatus according to SCHULZE-HARKORT is suitable for rapid and routine analysis of elutriation for process control in factories producing ceramics. The apparatus is supplied for one determination, or also consisting of 3 parts on a stand, with a common water container. The nickel plated elutriating pipe with funnel and overflow to regulate the water level, is delivered with 11 various jets: Size of bore 1-6 mm diameter. The glass elutriating funnel, which has a mark for filling at the top and overflow, contains a stop cock, which is situated below and is of wide bore for the rapid removal of the sludge residue after each analysis. The elutriating pipe may be adjusted in height, and due to the changeability of the elutriating boring it can be used for various elutriating speeds without having to alter the water pressure.

Technical Data:

Apparatus for one determination: brass elutriating pipe, nickel plated, adjustable, 11 jets with bores of 1, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5 and 6 mm. One elutriating funnel with filling mark made of glass with overflow. 1 heavy stand with clamps and holders.

Base approx. 20x15 cm overall height 80 cm weight approx. 5 kg Code: SH 1

Apparatus consisting of three parts: 3 complete elutriating pipes, 3 sets of jets, 3 elutriating funnels made of glass, water container of brass (nickel plated), with approx. 15 l capacity with inlet and outlet as well as clamps and holders to fix the elutriating pipe and funnel on the stand. Base approx. 90x30 cm overall height approx. 120 cm weight approx. 30 kg Code: SH 3

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
SPECIAL TESTING APPARATUS

Apparatus for Determining the Melting Point of Ash, Micro-method

Description of the Standard Apparatus:

The following parts are incorporated in a varnished case of sheet metal: 1 galvanometer with scale from +20...+1600°C, electric stopping mechanism and key switch, 1 blast burner with pipes for the supply of gas and oxygen.

Assembled on the housing is a support for fastening a platinum-platinum-rhodium thermo-couple with ceramic end, shovel shaped for the reception of the samples of ash. The thermo-couple can be adjusted in vertical and horizontal position.

On a second support the following parts are adapted in a swivel manner:

- 1 protective cylinder, chimney shaped, and made of ceramic.
- 1 feeding funnel, 1 blowing-out device for the cleansing of the thermo-couple, when the sintering point is to be determined.
- The gas burner head, with a protective cylinder which juts out of the housing.

Apart from the scale of the pyrometer which is adapted on the front panel of the housing there are fixed:

- 2 connecting terminals for the key switch, 1 yellow handle for the regulation of the gas, 1 blue handle for the fine regulation of the oxygen supply.

The optical outfit which belongs to the standard apparatus is attached to a swivel support and is adjustable for the various distances of the objects. It consists of:

- 1 Binocular stereomicroscope SM XVI manufactured by VEB Carl Zeiss, Jena, with the following equipment:
 - 1 Double objective 2.5
 - 1 Orthoscopic eye-piece 12.5 x
 - 1 ditto - adjustable
 - 2 Intermediate tubes, 23.2/30 mm
 - 2 Eye shells
 - 1 Eyepiece crossline micrometer 10 mm²
 - 2 Damping filters NG 3

Operating voltage for the electrical stopping mechanism: 220 V~
Base: 30x30 cm. Over-all Height: 35 cm. Weight: 10 kg

Code Apparatus for determining the melting point of ash, micro-system, in the previously mentioned description Standard outfit for visual examinations MAB

ACCESSORIES:

Illustration see page 71

To the Apparatus for Determining the Melting Point of Ash

Photomicrographic Device for Camera:

- Milieu-Adapter with rapid changer and adjustable telescope 3.2 x with adjustable eye-piece as well as built-in lined plate. Code MAB 30 60 11
- MF adaption for SM XVI. Code MAB 30 60 30

Furthermore for Choice according to the type of camera which is at hand:

- Adapter for Exakta-Vorex, Kine-Exakta and Exa. Code MAB 30 60 28
- or for Contax S+D and Praktica. Code MAB 30 60 27
- or for Praktina. Code MAB 30 60 29

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Apparatus for Determining the Melting Point of Ash, - Micro-method -

The apparatus serves for the rapid micro-determination of the melting characteristic of combustible ashes according to Dold's Pöchlmler modified by Reinhardt, ascertaining the characteristic points in the temperature melting range of combustible ashes.

The softening point and the flowing point can be determined and thereby rapid decisions in firing of steam boilers and orientating tests in reference to the qualification of solid combustion materials, can be taken with regard to the residue in the combustion chamber of boiler plants.

With the apparatus it is possible to carry out individual, visual observations, and with the appropriate accessories photographic recordings with a camera (Contax, Exakta, Exo, Praeflix or Praktina) can also be performed.

In principle the examinations proceed in the following way: By means of a charging device, the sample of ash is brought to a small scoop-formed platinum-platinum-rhodium thermo-couple, which is heated from below by a blast burner operating on oxygen and gas.

The thermo-couple is connected to the galvanometer with temperature scale (+20...+1600° C), which is incorporated in the housing of the apparatus.

A binocular stereomicroscope serves to observe the samples of the ash and the changes while temperature is increasing. During the observation, and by actuating a key, the pointer of the pyrometer can be stopped at that temperature, which characterizes a certain change of state (softening or melting point), of the sample located on the thermo-couple.

A very clear and plastic observation of the ash sample is made possible with the standard outfit of this apparatus up to an enlargement of 30 fold.


Accessories, which are deliverable make possible the use of the stereomicroscope separately for general laboratory purposes up to a 125 fold enlargement.

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Apparatus for Determining the Melting Point of Ash, with Photomicrographic Outfit acc. to page 69

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Further Possibilities of Application of the Optical Outfit.

The standard optical equipment may also be applied for other microstereoscopic examinations in the laboratory. For this purpose a convenient stand is obtainable as well as another double objective and eyepieces to obtain various scales of enlargement.

- Additional accessories:**
- 1 Stand for SM XVI with round base and interchangeable plate with black and a white surface, as well as two springs. **MAB 30 11 26**
 - 1 Double objective 5,0 x **MAB 30 24 04**
 - 2 Eyepieces 10 x **MAB 30 31 51**
 - 2 Eyepieces 16 x **MAB 30 31 52**
 - 2 Eyepieces 25 x **MAB 30 31 53**

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
Stability Test Apparatus according to Bergmann-Junk, nitrogen oxide fission method.
 The apparatus consists of a copper heating bath to take a glycerine-water mixture, with a massive cover made of brass; screwed into this are 10 protective sheaths also made of brass for the stability tubes. For protection against vaporization from the heating bath a Sahlstedt reflow condenser is fixed into the cover. The exterior jacket consists of varnished steel.
 The heating is electrical by 2 interchangeable heating elements, the electrical temperature regulation by an adjustable contact thermometer and relay. A control thermometer serves for further supervision, measuring range $\pm 100 \dots \pm 160^\circ \text{C}$.
 An electrically driven stirrer with a flexible shaft is employed to guarantee a constant temperature in the whole heating bath.

Electrical outfit:
 A three stage switch with a control lamp for "Voltage" and "Heating";
 Switch and pilot lamp, as well as a sliding resistance for the stirrer.

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Accessories:
Test-Tubes (Stability-Tubes), according to Bergmann-Junk with ring mark, standard ground joint and safety tube. Code: St BI

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
Safety Cabinet for the Stability Apparatus acc. to Bergmann-Junk
 Width 650 mm, depth 510 mm, over-all height 1550 mm, consisting of a wooden stand.

In the lower part: One rear partition made of plywood, two side sashes with wire-glass, two front sashes from which one is made of clear glass and the other of wire glass.

In the top part: Two folding doors with wire-glass, when in an open position, one of the two lower sashes can close the opening of the top part.

On the top of the safety cabinet 10 lifting devices with cords made from Perlon are located and also a suspension frame to lift out the stability tubes from the stability apparatus, when the protective cabinet is closed. Weight 65 kg Code: St BI

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Electrophoresis Apparatus EPH 1

The apparatus is designed to the principle according to which particles of different nature show different velocity of migration in electric fields. It serves for electrophoretic examination in chemical and clinical laboratories. Two different apparatus are at hand: **IEPH 1** generating 1000 V D.C. or **IEPH 2** generating 300 V D.C. The units consists of:

- 1 mains supply unit EPH 1 or EPH 2 and
- one or more electrophoresis vessels EPH K for 4 determinations.

The mains supply units are contained in varnished plate case with:

- 1 measuring instrument from 0-6 mA, reversible to 60 mA
- 16 connecting terminals for the electrophoresis vessel, type EPH 1 or 8 connecting terminals for the electrophoresis vessel, type EPH 2
- 1 enclosed voltage selector with 10 step-switch
- 1 fine regulation between the coarse stages
- 1 voltage stabilisation $\pm 2\%$ pilot lamp, cord and plug
- 1 set of high tension proof cables with connection screws and protected stainless steel carbon electrodes.

The electrophoresis vessels are made of glass and plastic. They are only delivered on special order.

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Technical Data:

	EPH 1	EPH 2
Voltage, mains supply	220 V	220 V
secondary	1000 V	300 V
Range	6/60 mA	6/60 mA
Dimensions	60/32/28 cm	40/23/20 mA
Weight	20 kg	14 kg
To connect directly to electrophoresis vessels	two EPH K	one EPH K
EPH K:	each 4 = 8	determinations
Size of one vessels EPH K:	43x43 cm	
Weight:		9 kg

Spare Parts:
Stability tube **GR 150** Rectifier tube **AZ 12**


4-5 electrophoresis vessels for 15 to 20 determinations can be connected to the mains supply unit by means of terminal blocks.

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Grinder for Metallographic Specimens with two or more shafts acc. to specification of types page 76

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The machine consists of a steel frame, which support all moving and rotating parts. It is provided with a wooden cover and built-in partition for emery papers and tools. The bench top is coated with plastic material, well polished and easily cleaned.

The speed may be varied from 300 to 900 r.p.m. The grinding shafts are mounted on plain bearings and carry the polishing heads of 250 mm diam, which can quickly be removed. Each shaft can be released and regulated separately by a hand lever.

The grinding heads rotate in bowls, closed by a cover to protect from dust when not in use.

A locking ring fits into a groove around the rim of the head and holds the emery papers firmly in place.

Each grinding head can be fitted out with mechanical sample-holder, thus avoiding holding of the samples by hand. Large specimens are provided with a boring for the point of a holding needle. For samples of an unsuitable shape, special holders can be manufactured.

The grinder is driven by a three-phase motor 220/380 Volt and supplied with "off" and "on" switch. The bench top is provided with a drip-waste of acid-proof stoneware and a water tap.

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Accessories and additional Parts for the Metallographic Grinders: Sample Holders

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for simultaneous attachment of several specimens, quickly to be removed by conical locking ring. The holder itself is not worn out by the grinding process.
Type SH 1 for 4 specimens upto 30 mm diam. Code SH 14
Type SH 3 for 2 specimens upto 20 mm thickness or 30 mm diam. with removable jaws Code SH 22

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Spare Friction plates for the grinding shafts
Spare Friction plates for the automatic drive

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
Polishing or Grinding Machine with one shaft,

mounted on ball bearings to support either polishing or grinding head, connected to the shaft by a cone and running in a bowl, which is closed by a cover.
The emery paper is fixed by a spring and the polishing cloth by a locking ring, both easily interchangeable.

A three-stage drive makes it possible to operate on 420, 675 or 825 r.p.m. by changing the V-belt.
The implement is installed in a robust wooden frame.
Design with polishing head of 300 mm diam. Code: MP 1P
Design with grinding head of 250 mm diam. Code: MP 1S

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
Metallographic Polisher with 2 shafts

with polishing heads of 300 mm diam. and friction drive. The speed may be varied stepwise from 300 to 900 r.p.m.
Each shaft can be released and regulated separately.
The machine is equipped with a three-phase motor 220/380 V, the current consumption is 1000 Watt.
The remaining features correspond to the grinder MS with 2 shafts. Code: MP 2

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
Polishing or Grinding Machine with 1 Shaft, to be mounted directly to the wall. Universal speed adjustment from 350 to 1000 r.p.m. by means of a hand lever. The machine is driven by a three-phase motor 220/380 Volt, the current consumption is 300 Watt.

Alternatively this type is fitted out with grinding or polishing heads. Code MWP
Design with polishing head 300 mm diam. Code MWS
Design with grinding head 250 mm diam. Code MWS

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Accessories and Spare Parts for Metallographic Polishers

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Standard Accessories for Metallographic Works

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Kurzbez. AET

Etching Bench

1,40 m long, 0,75 m wide, 0,90 m high acid-proof varnished, cupboard unit with 2 drawers and sliding doors, reagent bottle shelf made of wood, etching sink of acid-proof stoneware, two-way water tap, two-way gas tap and plug-box on standard

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
Kurzbez. B1

Fume Hood

1,40 m long, 0,78 m wide, 0,90/2,5 m high, acid-proof varnished, cupboard unit with sliding doors, bench top covered with tiles, combustion hood with sash window, fume duct, lighting and fume extraction system. The unit is provided with plumbing for gas, water, waste and electric current

Code M3

Storage-Dispensing Unit

for metallographic samples 1,30 m long, 0,60/0,40 m wide, 0,90/2,10 m high, varnished, lower section with 2 drawers and sliding doors, upper section with 2 glazed sliding doors and shelves

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Kurzbez. ET 248

Drying Oven ET 248 illustration on page 26

range of temperature from 40 to 180° C. Working compartment: 23 x 25 x 30 cm, automatic regulation of temperature, 1500 Watt, Supplied for 110 or 220 Volt consumption of current: 1500 Watt.

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Kurzbez. SPP

Metallographic Specimen Mounting Press

designed for mounting metallographic specimens in order that they can be ground, polished and examined under magnification. With interchangeable moulds 22 and 40 mm diam.

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- Corrosive reagent, Oberhoff Nitric acid, Hydrochloric acid, Copper Ammonium chlor. sol.

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- Picric Acid, Iron (II) chloride solution, Sodium Picrate, Solution A acc. to Fry, Solution B acc. to Fry

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- Diluted Sulfuric Acid

Reagent Bottles for Etching Liquids

brown glass, stoppered, with sand-blast labels

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- Copper Ammonium chloride Heyn, Cadmium Acetate solution, Hydrochloric Acid, dil., Diluted Sulfuric Acid, Potassium Ferricyanid sol. A, Potassium Ferricyanid sol. B, Corrosive reagent for stainless steel

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- Corrosive reagent, Benedicks, Corrosive reagent, Villela, plastic flask for Hydrofluoric acid

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
Camera Microscope for direct light "Mikrophot A 21" for metallographic tests and photomicrography

designed for the examination of opaque objects such as metals, both for visual observation and photomicrography. The camera is incorporated in the stand of the microscope, and serves for plates 4.5 x 6 cm or films 24 x 36 mm.

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Quick-Test Units for steel and iron as well as for non-ferrous metals

These portable quick-test units are proved for rapid orientating tests of raw materials such as steel and non-ferrous components.

They are very convenient and most simple in construction, so that tests may be carried out everywhere without any special preparations. This type of apparatus cannot replace exact analysis of steel and other metals in a well-equipped laboratory. They are however the means for

Rapid Testing

and particularly suitable for the prevention of mistakes when assaying metals in the store-room.

Two portable implements are at hand

Model EST 8000 for iron and steel
Model NE 8001 for non-ferrous metals
based on the Spark and Spot Analysis.

The Spark Test serves for ascertaining the essential components of alloys as well as the carbon contents of steel and iron. When grinding the sample, there appear specific phenomena of combustion in various colours and forms, such as bulbs, spikes, drops, bulbs and the like. An assortment of test samples with exact analysis is included for comparing the kind of sparks with each other. This method, of course, requires a certain measure of experience.


A full set of reagents for carrying out the Spot Test is comprised. By adding a drop of the spot reagent, there appear various colour reactions on the tested sample, by which the components are easily determined. For the same purpose test-papers are provided.

A so-called "Poldi-hammer" serves for carrying out the Hardness-Test. Two leaflets EST 8000 and NE 8001 with well outlined test plans and tables give all information with regard to the carrying out of these tests.

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Quick-Test Unit for Iron and Steel
Portable, cased apparatus with handle and lock
Dimensions: 47x32x20 cm.
Weight approx. 8 kg
Code EST 8000

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Table of Contents EST 8000

- 1 Hand grinding machine with 3 different disks
- 1 Chisel
- 1 Hammer
- 1 Blunt scraper
- 1 File, triangular
- 1 Magnet
- 1 Set of spark-test samples with certificate
- 1 Wooden handle to hold the spark-test samples
- 1 Book of nickel-test-paper
- 1 Book of molybdenum-test-paper
- 1 Poldi-hammer for hardness test

1 Set of Spot-Test-Reagents

- 250 ml Washing solution (Nitric acid 1:2)
- 250 ml Hydrogen peroxide 3%
- 250 ml Hydrochloric acid 1.16-1.19
- 250 ml Ammonium persulphate saturated

- 7 Pipette bottles of 50 ml with etched names
- Silver nitrate
- Hydrogen peroxide
- Hydrochloric acid
- Ammonium persulphate
- Nitric acid
- Ammonia
- Washing solution
- 2 Pipette bottles without inscriptions
- 6 Test tubes
- 10 Glass rods
- 2 Bottles with wide mouth, 50 ml
- 1 Graduated cylinder 100 ml
- 11 Stopper-holders


- 250 ml Nitric acid 1.40
- 250 ml Ammonia conc. 0.910
- 250 ml Sulphuric acid 1.84
- 250 ml Solution of Silver nitrate 0.2%
- 250 ml Solution of copper nitrate 0% copper nitrate on 50 ml dist. water

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Quick-Test Unit for non-ferrous Metals
Portable, cased apparatus with handle and lock
Dimensions: 37x16x37 cm
Weight approx. 7 kg
Code NE 8001

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Quick-Test Unit for non-ferrous Metals
Nitric acid, hydrochloric acid
Ammonia, Silver nitrate
Ammonium persulphate
Sodium hydroxide solution 20%
Potassium iodide solution
Sulphurous acid 3%
Ammonium-mercury-thiocyanate
Potassium hydroxide solution
Ammonium thiocyanate sol. sat.
Cadmium sulphate 5%
Bottles, narrow mouth, brown, with stopper, 50 ml and etched names
Hydrogen peroxide 3%
Acetic acid
Tetraoxanthronchin in sodium hydroxide solution, Dist. water

Bottles, wide mouth, brown, with stopper, 50 ml and etched names
Sodium bismuth
Sodium chloride
6 Test tubes 160x16 mm
6 Glass rods
1 File
1 Hammer 120 mm long
1 Hammer 250 g
1 Graduated cylinder 100 ml
1 Funnel 55 mm diam.
100 Filter papers 7 cm diam.
4 Sheet of emery paper different grades
1 Hardness-tester "Poldi"

1 Set of Spot-Test-Reagents
250 ml Nitric acid conc.
250 ml Hydrochloric acid conc.
100 ml Potassium persulphate 3%
250 ml Ammonia conc.
100 ml Potassium iodide solution, saturated
100 ml Sulphurous acid 3%
100 ml Ammonium-mercury-thiocyanate
250 ml Potassium hydroxide sol. conc.
50 ml Acetic acid 98/100% pure
100 ml Silver nitrate solution, saturated

1 Set of Spot-Test-Reagents
100 ml Ammonium persulphate solution, saturated
100 ml Ammonium thiocyanate solution, saturated
250 ml Sodium hydroxide sol. 25%
25 ml Sodium bismuth, pure, p. 25
100 ml Cadmium sulphate sol. 5%
150 g Sodium chloride pure p. 25
1 book of nickel-test paper
1 book of titanium-test paper
1 book of potassium-ferrocyanide test-paper
Code NE 8001 CH

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
Apparatus for Testing of Boiler-Feed-Water acc. to
Prof. Dr. Spittigerber and K. H. Ulrich

The control of feed-water in steam boilers, hot- and warm-water-plants etc. is outermost important in view of the safety of working for the durability and of a rational fuel consumption. In these plants. In many countries this control is even stipulated by law.

This is an up-to-date portable laboratory unit installed in a wooden cabinet 58 cm wide, 57 cm high and 15 cm deep. The most important tests of feed-water can be carried out with regard to:

- degree of hardness in crude water, soft water, condensed water as well as in feed-water
- degree of lime-saturation
- total salt-content (density), p- and m-value, (alkalinity)
- total chloride-content
- colorimetric test of phosphate
- pH-value

The implement comprises 4 special burettes with storage tanks, which are ready for use only after swiveling out the appropriate holders. Moreover the practical construction permits leaving the burettes and reservoirs in filled state even if the apparatus is transported. All parts of the unit are standardized and can be easily substituted.

Outfit:
4 special burettes with reservoir and deeply cut sand-blasted labels "Soap Solution B+B" .. burette 0-5 ml/0-260 dH
"Soap Solution B+B 1+9" .. 0-2 ml/0-10 dH
"Hydrochloric Acid n/10" .. 0-10 ml/Schellbach
"Silver Nitrate" .. 0-10 ml/Bronnglass
and all sundries necessary for the mentioned tests.
Weight: 14 kg
1 Set of Reagents

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Kurzbes. **NWUZ**
Kurzbes. **NWUZ CH**

Supplementary Outfit for the control of feed-water in steam boilers, determining **oxygen (O₂)** and **free carbon dioxide (CO₂)**

The apparatus comprises two oxygen determination flasks according to Winkler, two swiveling mounted special burettes with reservoirs for "Sodium Thiosulfate Sol. n/100", "Sodium Hydroxyde n/50" and the necessary sundries for the tests.

Dimensions of case: 51x37x16 cm, weight 7 kg
1 set of Reagents weight 1 kg

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Kurzbes. **NWU 2**
Kurzbes. **NWU 2 CH**

Simplified Apparatus (small model) for the control of feed-water in steam boilers, determining **degree of hardness, carbonate hardness, total salt content (density) p- and m-value (alkalinity)**

The implement comprises two swivel mounted special burettes with reservoirs for "Soap Solution Boutron-Boudier", "Hydrochloric Acid n/10" and sundries for the tests.

Dimensions of case: 51x37x16 cm, weight 7 kg
1 set of Reagents weight 1 kg

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
Apparatus for Testing of Drinking Water portable in a wooden casing for the determination of **degrees of hardness, carbonate hardness, iron, potassium permanganate number and pH-value.**

Outfit: 6 swiveling mounted special burettes and reservoirs with deeply cut sand-blasted labels and the necessary sundries: automatic quick acting pipette 5+10 ml, pipette- and dropping flasks, washing bottle, spot plate, test tubes, spirit alcohol Bunsen burner with hose and reservoir etc.

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

We manufacture suitable vacuum pumps for all working ranges of vacuum technology. These pumps are described on pages 90 to 91.

They are classified as follows:

Rotary Vacuum Pumps driven by Electric Motors and Diffusion Pumps

by which mercury or silicone oil is used as pumping medium and which are electrically heated.

For the selection of the pump it is decisive in which range of the vacuum technology one wishes to operate

In the range of 760 to about 5 Torr and above all for pressures below 20 Torr, the single stage rotary vacuum pump ED 3 has proved to be economical and reliable in service, especially since this pump is fitted out with gas ballast device and is therefore also suitable for the pumping out of condensable vapours.

The gas ballast device was installed in this model since in the range of coarse vacuum there may be large quantities of condensable vapours e.g. in case of vacuum drying and filtration work.

The double stage rotary vacuum pumps ZD 4, ZD 8 and Z 30 should on no account be employed in the coarse vacuum range. Apart from the high costs, these pumps require a more attentive servicing and do not offer any advantages as against single stage pumps for coarse vacuum

In the fine vacuum range, up to a necessary working pressure from 1 Torr, single stage pumps are also used with advantage

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Only in the case of lower pressures it is justified and economical to use the double stage rotary vacuum pumps ZD 4, ZD 8 or Z 30.

In the case of fine vacuum work also the gas ballast device is of decisive importance as very often condensable vapours must be pumped off e.g. in the case of fine vacuum freezing, drying or fine vacuum impregnation.


When working in the high vacuum range (pressures below 10⁻³ Torr) mercury or oil diffusion pumps should be used, as, for example, our models Q 10 and Q 12 or O 21 and O 22 which are made from Jena Rosotherm glass. These diffusion pumps require a pre-vacuum for operation. The following are quite suitable as pre-vacuum pumps:

- the rotary vacuum pump ED 3 for the mercury diffusion pumps Q 10 and Q 12
- the rotary vacuum pumps ZD 4 and ZD 8 for the oil diffusion pumps O 21 and O 22.

If mercury vapours which arise during use of mercury diffusion pumps (vapour pressure of mercury at 20° C ~ 10⁻³ Torr) should adversely affect the working chamber these vapours must be frozen out in a cooling trap (which is to be fed with liquid air).

It should be remembered that a cooling trap strongly throttles the sucking speed of a diffusion pump. For this reason, and due to the poisonous mercury, for many applications plants with oil diffusion pumps are preferable, although the costs of these are a little higher.

The silicone oil which is used as a pumping medium for these pumps has a vapour pressure of 10⁻³ Torr at room temperature, so that there is no need, to apply a cooling trap. It has been proved that the application of oil diffusion pumps, results in shorter evacuation times, and above all the supply of the pump plant with liquid air is also done away with.

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

Single Stage Rotating Vacuum Pump Model ED 3
 Oil pump for the suction of gases and vapours. Suitable as backing pump for mercury diffusion pumps. This pump has a gas ballast device for reliably pumping off condensable vapours. A reflux of the pump oil into the vacuum apparatus will not occur.

In the top of the housing, baffles are incorporated on the suction and pressure sides, which prevent impurities from entering the chamber of the pump.

The oil charge amounts from 30 to 70 ml. It is recommended to work with the lowest amount, when for example, aggressive mediums are to be sucked off (as for instance, acid vapours) and a frequent change of oil is necessary. This pump is arranged for a rapid changing of the oil.

Displacement-Diagram of the Rotary Vacuum pumps see page 89

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

Double Stage Rotary Vacuum Pumps, Model ZD 4 and ZD 8
 These oil-air pumps are intended for the pumping out of gases and vapours They serve as backing pumps for diffusion and vapour jet pumps of all kinds. Special distinguishing marks of the apparatus are the separate arrangement of the pre- and high vacuum stages and the automatic oil regeneration, according to DAP 7566 and DBP 958690. (The principle of the oil regeneration is described on page 93.) Due to this the working oil of the high vacuum stage remains fit for use without limits when properly treated. The pre-stage of the pumps is arranged for a quick change of the oil. On the pre-stage the pumps are fitted with a gas ballast device. The oil consumption of the pumps is limited. A return flow of the pump oil into the vacuum apparatus does not take place. Furthermore these pumps are recommendable due to their high working reliability when in continuous operation, their small dimensions, light weight and their high end vacuum in the case of gas ballast operations. Bases: ZD 4 wood, ZD 8 cast lightmetal.

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
VACUUM ENGINEERING

Double Stage Rotary Vacuum-Pump Z30

The construction of the pump warrants absolute plant safety during simple operating. In consideration to the rational capacity, this pump is pre-eminently intended as backing pump for large oil diffusion pumps and as high vacuum roots-type blowers as well as for all fine vacuum work e. g.: freeze drying plants, impregnation, metal sintering, vacuum soldering etc. The two pump-stages are kept separately. The shafts of the pump are arranged vertical so that the pump driving can be advantageous located in the allround closed base-plate of the pump. The driving of the pump results by Vee-belt which is stretched by idler pulleys with ball bearings. The pre-stage of the pump has a gas ballast device and the high vacuum stage is equipped with automatic oil regeneration acc. to DAP 7566 and DBP 958 690. The high ultimate vacuum of the pump during gasballast operations, makes possible to drive the pump with gasballast device permanently and so to utilize the advantages of the gasballast operations. Noteworthy are the small dimensions and the low weight of the pump in consideration to the displacement of 30 m³/h.

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
rotary vacuum pump with automatic oil regeneration according to DAP 7566 and DBP 958 690

1. Pump out of operation
2. High vacuum stage
3. Junction channel
4. Boring
5. Oil container
17. Connecting tube

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
VAKUUM-TECHNIK

Double Stage Rotary Vacuum-Pu
The construction of the pump was in consideration to the rational and as high vacuum roots-type of silencing, vacuum soldering etc. The two pump-stages are kept separate and are located in the oil reservoir stretched by idler pulleys with ball bearings. The pre-stage of the pump has access to DAP 7566 and DBP 958.6. The high ultimate vacuum of the device permanently and so to utilize the small dimensions.

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The Gas Ballast Device

For the reliable pumping out of condensable vapours, the rotary vacuum pumps ED 3, ZD 4, ZD 8 and Z 30, which operate on principle of Goeda should be used. In this case a condensation in the pump is avoided from the very beginning as there is a certain amount of fresh air (gas ballast) continuously running into the scoop chamber through an opening fitted with a regulating screw.

As the drawing in the margin illustrates, the gas ballast valve is adapted in such a way, that the gas ballast only flows into that part of the scoop chamber, which is separated from the suction side of the pump.

The operation of the gas ballast valve is based on the fact, that the supply of the gas ballast into the scoop chamber of the pump effects an increase of pressure. Hereby the valve closing the exhaust side of the pump, can be opened in time, so that the vapours which are in the scoop chamber need not be compressed to condensation.

Any amount of condensable vapours can be pumped with a gas ballast pump. It is recommended to take care, that the vapour should not enter into the pump with a temperature which is higher than 20° C

Eventually a water cooled condenser should also be mounted into the suction line of the pump.

The gas ballast device does not offer sufficient protection when sucking off oil-soluble vapours, e. g. petroleum. The same refers to vapours and also gases, which attack the material of the pump and the working oil mixture.

For the sucking of aggressive material, oil pumps are suitable only within certain limits. In such cases devices (e. g. absorption towers) must be mounted in the suction line, which bind the aggressive mediums. However they only offer a limited protection.

In such operating conditions, only single stage oil-air pumps with oilsuperposed pressure valve should be used, for example our model ED 3.

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Automatic Oil Regeneration, Method of Operation
 The new oil supply for the double stage rotary vacuum pumps ZD 4, ZD 8 and Z 30 according to DAP 7566 and DBP 958.690 warrants the regeneration of the working oil of the high vacuum stage II. The pre-stage I is kept separate from the high vacuum stage II. The channel 13 serves to join the two stages.

During the pumping process, a low pressure is generated in housing 12, whereby gases and vapours which are in the working oil are freed and sucked through 13 into the pre-stage I.

The boring 15, affords a plentiful lubrication of the scoop chamber and the moving pump parts.


If the pump is no longer in motion, then atmospheric air, enters into the housing 12. The working oil of the high vacuum stage is pressed into the reservoir 16 which is under vacuum. Hereby the oil-charge is protected against the influence of the atmospheric air.

Mechanical impurities of the oil are deposited on the bottom of the reservoir 16, and do not arrive anymore in the pump chamber of the high vacuum stage, into which the clean regenerated oil flows back as soon as the pump is once more set running and a certain vacuum is obtained in the pump housing 12.

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Double stage rotary vacuum pump with automatic oil regeneration according to DAP 7566 and DBP 958.690
 a Pump in operation b Pump out of operation
 1 Pre-vacuum stage II High vacuum stage
 2 Suction connection 13 Junction channel
 11 High vacuum pump 15 Boring
 12 Housing 16 Oil container
 17 Connecting tube

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
High Vacuum Pumping Plant P Hg 1
This pump stand is primarily intended for high vacuum evaporation of metals, for the evacuation of vacuum tubes on pump gages and similar physical high vacuum operations.

The complete plant easily viewable it is incorporated in a steel tube frame and is fitted out with a switch board for the servicing and supervision of the pumps. A vacuum plate of Jena Duron glass is inserted in the rubber coated bench plate for high vacuum work under a receiver bell jar. Due to the standard ground joints in the receiver plate, vacuum tight inlets from the working devices are thus made possible. Separate experimental systems, such as pump gages, test valves etc. can be directly adopted to the connections of the pumps.

The freezing out of the mercury and grease moisture takes place in the cooling trap, which is to be fed with liquid air. After attaining the high vacuum in the receiver, the backing pump may be switched off. A large acceleration of the pumping process can be noted whilst there is no overcharge of the pumps.

Outfit: Steel tube frame, switch board, bench plate, vacuum plate, bell jar, single stage rotary pump ED 3 with electric motor, Mercury diffusion pump Q 12 with electric heating, valve block, cooling trap, Dewar vessel with protective jacket, rotary vacuummeter ML 11,3, U-tube manometer MU, connecting joints for the electric vacuummeter, 500 g rotary pump oil, 100 g Ramsay grease.

The Model P Hg 1 is simple in service and insensitive against the intrusion of air during operation.

Technical Data

Ultimate vacuum	10 ⁻² Torr
Operating time for evacuation	4-5 Min.
Vacuum plate with 2 sockets and 2 sockets	275 mm Ø NS 29/20 NS 26/22
Bell jar diam.	230 mm
Height	300 mm
Voltage	1000 W
Cooling water is required	
Base/Over-all height	65/100/140 mm
Weight, Ede	75 kg P Hg 1

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
High Vacuum Pumping Plant PO 1
PO 1 is rational in operation, and represents a suitable plant for all high vacuum operations. The device is incorporated into a steel tube frame which is quite easy to survey. For the supply of the high vacuum serves a double-stage rotary vacuum pump and an oil diffusion pump which has a high speed of suction, and which is driven with silicone oil. This driving medium is neither poisonous nor injurious to health. The vapour pressure of the silicone oil is less than 10⁻⁴ Torr. The pumping plant operates without liquid air. A water cooled oil-moisture trap holds tight the last traces of oil-moisture from the high vacuum receiver. Particular attention has been paid to construction of a short vacuum line, securing a high displacement. A vacuum-gauge serves for the control of the pressure.

The following parts belong to the Apparatus:
Steel tube frame, switch board, bench plate, vacuum plate, bell jar, a double-stage rotary vacuum pump ZD 8 with electric block with stopcock for the diffusion pump, change-over cock for the backing pump as well as a cock for the intake of dry pre-vacuum container, water cooled oil-moisture trap, 1 Pirani-gauge with 2 measuring heads, 1 ionisation gauge according to Penning, with one measuring head, 500 g rotary pump oil, 200 g silicone oil, 100 g Ramsay grease.

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
Technical Data

Ultimate vacuum	10 ⁻⁴ Torr
Operating time for evacuation	4-5 Min.
Vacuum plate with 2 sockets and 2 sockets	275 mm Ø NS 29/20 NS 26/22
Vacuum connection - ground socket	230 mm
Bell jar, diam.	300 mm
Voltage	500 W
Cooling water is required	
Base/Over-all height	130/65 cm
Weight	ca. 140 kg
Code	PO 1

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
High Vacuum Freeze Drying Plant HVG 1

The device is suitable for refrigerated dryings in laboratories in the temperature range from -30 to -70° C on histological preparations, surgical implantation materials, living cells, sensitive sera, bacterio, virus cultures, antibiotics and other pharmaceutical preparations. Water vapours are sucked from the substances direct by the pump plant, without the use of liquid air traps or absorption adapters.

No liquid air is necessary. The plant operates with a high suction speed and short drying times. The pre-pump can temporarily be switched off without interruption of the drying process. An electric gauge with two measuring instruments and measuring heads serves for the control of pressure.

Outline:

Steel tube frame, switch board, bench plate, rotary vacuum pump ZD 8 with electro-motor, oil diffusion pump O 22 with electric heating, water cooled oil-vapour trap, pre-vacuum vessel, 3 drying rooms for cooling with dry-ice acetone mixture and 1 Pirani gauge with two measuring heads.

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Special Accessories to the High Vacuum Freeze-Drying Plant HVG 1

So as to be able to dry larger quantities of material at relatively high temperatures (0 to -30° C) in an economical way in vials, (the direct sucking of the water vapour through the pump plant is only advantageous in the case of low drying temperatures), the following can be additionally delivered on request:

1 Pump rake for the connection of 16 vials, with liquid



air condenser for the precipitation of the water vapour.

The liquid air condenser is cooled with liquid air or a dry-ice acetone mixture. The Dewar vessel, which belongs to the standard equipment, serves for the reception of the cooling medium. The liquid air condenser with pump rake, is connected to the plant by a standard ground joint NS 45/40 where otherwise the drying chamber is fixed.

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
The illustration in the margin shows the drying chamber 250 X 40 mm belonging to the standard equipment HVG 1 and the drying chamber 150 X 60 mm.

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MERCURY AND OIL DIFFUSION PUMPS
Made of Heat Resistant JENA "Rosotherm"-Glass

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Electric heating device for O 21 and 22, 220 V, 450 W
Code: O 21/22/H

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Gauge according to Mc-Leod ML 113 is a precision vacuummeter. A water jet or a rotary pump is necessary for the operation. To increase the measuring range, a swivelable vacuummeter can be additionally connected.

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
Measuring range: 0,05 ... 10⁻³ Torr, with ML 3,5 = 10 ... 10⁻³ Torr, capacity of the measuring chamber 113 ml, measuring capillary: 1,0 mm diam., vacuum connection standard ground socket 19/26, necessary Hg (chemically pure) 250 ml, over-oil height: approx. 60 cm, on a steel tube frame, weight approx. 4 kg

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Gauge, swivelable type smallest dimensions. The measuring process requires only turning up to 90°. A high measuring accuracy is obtained by a wide testing capillary tube.

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Modell ML 11,3 range 0,1 ... 10⁻¹ Torr Connection by standard ground socket 14,5/35, Hg charge 16 ml, measuring tube 1 mm diam.

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U-Tube Manometer Model MU can be filled very comfortable under high vacuum. Access of air bubbles into the measuring tube is impossible, measuring range 1 ... 60 Torr, connection standard ground cone 14,5/35, Hg 3 ml, measuring tube 3 mm diam.

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Electric vacuum measuring instruments for the continuous measurement of the total pressure of gases and vapours.

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Pirani Gauge. The pressure dependence of the heat conductivity in rarefied gases and vapours serves as measuring principle.

The Pirani gauge is fitted with 2 measuring heads, and 2 dial instruments, so that the pressure at two measuring positions, can be read simultaneously without any change-over. The dial of the measuring instruments is furnished with a pressure scale calibrated in Torr. The measuring heads are interchangeable and all their physical data are normalized.

The Pirani gauge is specially suitable for the vacuum measuring at fine vacuum plants and for the pre-vacuum measuring of oil diffusion pumps. It is used in high vacuum plants in conjunction with an ionisation gauge.

This apparatus is also supplied with 1 measuring head.

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Ionization Gauge acc. to Penning. The ionization capacity of vapours is used for the measuring. The apparatus is extraordinarily reliable in operation and entirely insensitive against air intrusion. The servicing of the apparatus is quite simple.

The measuring head, with a conical ground standard socket NS 29/32 and contains the electrode system, on which a voltage of 1500 V is applied. Two cermet permanent magnets for the amplification of the ion current, are adopted externally to the measuring tube which is made from glass. Measuring tube and magnets are contained in a protective housing made of metal. It is connected to the apparatus by a high tension cable.

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Accessories for High Vacuum-Plants

Needle Valve made of steel, for the inlet of minute amounts of gas and air with:

2 standard ground sockets 14,5/23, possibilities of regulation 0 ... 50 Torr X l/sec. Length 150 mm, Weight approx. 0,4 kg. Code **NV**

Current Lead-In made from steel with standard ground joints, NS 29/32 suitable to the receiver plates and cover plates of the pump stands and vacuum apparatus, loads up to 50 Amp. Weight approx. 1,0 kg. Code **STD**

Turning Handle made of steel, high vacuum tight, for the receiver plates and cover plates with standard ground socket 29/32, Weight approx. 0,7 kg. Code **VDD**

Dewar Vessel made from Jena glass 20, of cylindrical shape, double-walled and silver-plated.

Internal Height	200	200	250	mm
Internal Diam.	40	60	80	mm
Code	VDG 1	VDG 2	VDG 3	


Protective Jacket for the Dewar Vessels, varnished sheet metal, with a holder for fixing to stand rods.

For Dewar Vessels of:	200/40	200/60	250/80	mm
Code	VDG 1 S	VDG 2 S	VDG 3 S	

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



Pipetter tongs with syringe and interchangeable pipettes
 Code: PIZ

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Tongs, flexible, with digitate and ungulate heads
 Code: Z b

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

Pipetter "Mikro" with adjustable syringe and interchangeable pipettes
 0,1, 1,0, 2,0 ml
 Code: PIM

Pipetter "Pantograph" with interchangeable syringes and bulb pipettes
 1, 2, 5, 10 and 20 ml
 Code: PIZ


Transparent Plastic Shield made of 10 mm "Plexiglyl", 700 mm high, 600 mm long, sidewall 900 mm long, with which can be closed, hinges made
 Code: SP 700

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




Pipetter „Mikro“
with adjustable syringe and
interchangeable pipettes
0,1, 1,0, 2,0 ml
Code: PIM

Pipetter „Pantograph“
with interchangeable syringes
and bulb pipettes
1, 2, 5, 10 and 20 ml
Code: PIZ

Transparent Plastic
Shield
made of 10 mm „Piacryl“,
700 mm high, 600 mm long,
sidewall 200 mm long, with
two openings for manipulation which can be closed, hinges made
of „Perlon“
Code: SP 700

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The establishing of a labora-
tory for radioactive work
requires an accurate project
observing all prescriptions
regarding the safety of the
isotope workers and the ge-
neral public.

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Lead Brick A
160 x 80 x 50 mm, 4-sided with prismatical mortises and tenons (material: lead-antimony alloy, lacquered with a yellow acid-proofed chlorinated rubber lac)
Code: Pb A

Lead Brick B
bottom closed, basis unit, only 3-sided with prismatical mortises and tenons
Code: Pb B

Corner Brick C
all-sided with prismatical mortises and tenons
Code: Pb C

Corner Brick D
bottom closed, other sides with prismatical mortises and tenons
Code: Pb D

Ball Brick E
160 x 80 x 50 mm, as above "A"
Code: Pb E

Ball Brick F
160 x 160 x 50 mm.
Code: Pb F

Window brick G with lead glass of 50 thickness, 160 x 160 x 50 mm.
Code: Pb G

Cast iron Bricks
lacquered with a yellow acid proofed chlorinated rubber lac

- Form A - like lead bricks Code: Fe A
- Form B - Code: Fe B
- Form C - Code: Fe C
- Form D - Code: Fe D
- Form E - Code: Fe E
- Form F - Code: Fe F
- Form G - Code: Fe G

Glass-Mirror
300 x 400 mm, mounted on an adjustable support, with handles for manipulation on the right hand
Code: SP

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

Sink unit made of plastic
material
size of the bench top 1000 x 600 mm
height of the bench 900 mm
with rear wall 1600 mm
with arm operated taps for cold and
warm water, knee operated waste
outlet and a sink made of black plastic
material (600 x 400 x 200 mm) for
direct connection to the canalisation.

Set up on the bench plate is the
rear wall of plastic material with two
shelves, rods and drain
Code: INA

alternatively:
Sink unit suitable for connection to a
waste container, with switch for the
water level indicator, pilot-lamp and
movable support for the stainless steel
container
Code: INO

Liquid waste container
made of stainless steel
contents 10 Liter
height 400 mm
Diameter 250 mm
cover flanged on the container, with
tubule for water-level-indicator
Code: V2A

Automatic water
level-indicator
with swimmer and switching valve
Code: WSA

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Solid waste container

made of sheet iron lacquered with chlorinated rubber lac, with foot
operated segment closing and separate pail made of plastic material
or paraffined paper
Code: AEK

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Plastic pail

315 mm high, 300 mm in diam. (diam. with handle turned down: 380 mm)
Code: EKK
Pail made of paraffined paper
Code: EK

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Abchlußstein L, 120 x 140 x 200 mm
Kurzbz.: Baryt L


Baryte-concrete bricks

dovetail-form, angle 120°, for pro-
tection against scattered radiation
Basic brick H, 120 x 140 x 200 mm
Code: Baryt H

Structural brick I
120 x 140 x 200 mm
Code: Baryt I

Structural brick, half size, K
120 x 140 x 190 mm
Code: Baryt K

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
Lead Container
for transports and storage only within the laboratory with cover and aluminum tube

interior	height 130 mm
Weight 21,5 kg	diameter 50 mm
	Code: BB


basket of hoop-iron with handle
Code: BBT

stand with devices for holding the lead container in several positions
Code: BBS

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ACTIVE WORK

with 2 fixed rubber
d and plug
interior and exterior

Code: GB

Code: GB 5t

glove box with the

Code: GB/AV

Code: BGP


GERÄTE FÜR


Equipments for air-suction and auxiliary air plants of radiochemical fume hoods

1. Air filter for filtering of the suction air stream, interchangeable, in PVC (plastic) casing, dimensions 600 x 500 mm, filled with folded asbestos paper of 3.8 m² surface for an airflow of 500 m³/h, cased in a housing of 400 mm overall height, to be installed on the top of the fume hood, with connection flange.
2. Air ducts and fittings made of PVC (plastic) or acid proof stoneware, in different dimensions.
3. Exhaust fan made of PVC (plastic) of different capacity, accommodated to requirements, with direct coupled three-phase electromotor for 380/380 V with ground plate, vibration dampers and flexible parts.
4. Auxiliary air fan made of acid proof lacquered iron sheet, of different capacity, accommodated to requirements with direct coupled three-phase electromotor 220/380 V with ground plate, vibration dampers and flexible connection parts.
5. Air filters interchangeable for filtering the auxiliary air, cased type, with connecting flanges.
6. Heating unit for 220 or 380 V, switched in two steps, with connecting flanges, wattage for one hood compartment: 6 kw
7. Flanged pipe with airflow-switch, designed to switch-off the heating unit when the auxiliary air-supply is interrupted or choked.
8. Flanged pipe with two installed adjustable contact thermometers for regulation of auxiliary air temperature i. acc. with the outside air temperature.
9. Switchboard „A“ is to be installed near by the fume hood, with 4 double push-button switches to switch on and off the two steps of the heating unit as well as the suction and auxiliary air fans, with pilot lamps and fuses.
10. Switchboard „B“ to be installed on the suction-fan, with protective switch and relay, terminals and fuses.
11. Switchboard „C“ to be arranged on the auxiliary air fan, with motor-protective-switch, 2 relays for the heating unit, 2 relays type T 50, 2 relays type T 1, fuses and terminals.
12. Fixing appliances and bases for pipings and other parts of the plant.

Wiring diagram of ventilation plant

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
Bleiglasplatten, 300 x 400 x 10 mm

Glove-box

- 1 removable tray of stainless-steel, 860 x 510 x 15 mm
- 1 throttle valve for air intake
- 1 cable opening, which can be closed
- 1 gas-tap
- 1 water-tap
- 1 waste funnel made of plastic
- 1 sluice on the exterior side closed by flap door, in the interior by sliding door
- 1 suction tube with throttle flap, removable filtering box and one set of filters made of glass-wool tissue

- 2 manipulating openings 200 mm diameter with 2 fixed rubber gloves
 - 1 low tension lamp with throttle, switch, cord and plug thereunder
 - 1 Placryl window tightly welded with the interior and exterior covering of plastic material
 - 1 flange with cover 160 mm diameter
 - 3 support rods on the rear-wall
 - 2 slides of wood
- Stand of steel tube 920 x 550 x 800 mm
- Set of flexible junction parts (P. V. C.), 90° cornered for connecting the glove box with the exhaust duct
- Lead glass plates, 300 x 400 x 10 mm

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

FACILITIES FOR RADIOACTIVE WORK

Radiochemical fume hood

- Size: 1410 mm length, 850 mm width, 3100 mm total height
 Work surface: 1360 mm length, 760 mm width, height of the bench 900 mm
 Outoff: × 2 gas-taps and
 × 2 water taps mounted on the switch-board, outlets in the interior
 of the hood
 × 2 waste-funnels of plastic material for nonactive, liquid waste,
 mounted appr. 200 mm over the work surface
 × 1 active waste funnel, which can be closed, made of stainless
 steel, in the center of the work surface
 2 throttle-valves for the auxiliary air
 1 handle for the movable section of the baffle
 1 lighting device
 1 filtering box mounted over the roof in the suction duct
 × 1 light metal box supporting on its switchboardlike cover:
 2 electricity sockets 230 Volts with protecting flap and earth-wire
 1 electricity socket 380 Volts with flap and earth-wire
 1 switch for lighting
 1 switch for the water-level-indicator, pilot lamp, buzzing signal
 device
 Stand: steel-tubes with levelling feet
 Hood with baffle: plywood of different thickness, completely covered with plastic
 foil 2 to 6 mm thick and massive black and white plastic
 material, corners rounded
 Sash: transparent plastic material "Plexiglas" 10 mm thick, counter-
 balanced, with two manipulating openings
 Joints: between hood and work surface must be filled by the special putty:
 Phialopal G, after having mounted the hood
 Bendtop: concrete, covered with 1 mm stainless steel, bordered
 Thickness: 50 100 150 mm
 Maximal loading capacity: 400 700 900 kg
 Codes: IM 1/50 IM 1/100 IM 1/150

Sink hood

- same measurements and materials as above
 Bendtop: wood covered with black plastic material, bottom covered with lead
 Outoff: 1 black plastic sink 500 x 400 x 250 mm, for connection to the
 drainage
 1 ditto, with treaddoorkite cover, for connection to liquid waste
 container. Both sinks inside covered with lead.
 2 knee-operated water taps and sink outlets
 1 switch for lighting
 2 plastic shelves on the sidewalls
 Without electricity sockets, gas and water taps, marked by × Code: IM 2

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


FRIEDRICH GEYER / LABORATORIUMSAPPARATE KG / ILMENAU

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
FACILITIES FOR RADIOACTIVE WORK

The illustrations are demonstrating two different manners of mounting the suction fan and the auxiliary air device. In the auxiliary air duct a filtering and an electrical heating device, a regulation of temperature and a safety-switch are built-in. The specialhoods for radioisotope works have an air conducting, which it makes possible to let pass over the full front of an uniform air-flow. Its velocity has been fixed at 150 linear feet per minute with the cabinet front open 1 foot. The withdraw of both light and heavy fumes may be achieved by the proper adjustment of a movable section of the back baffle. The hood is fitted with auxiliary air devices for working without hindering the air flow behind the closed sash-window. The regulation of the auxiliary air flow is effected by throttle valves to right and to lefthand. The inlet behind the sash contains an air by-pass opening, which becomes exposed, when the sash is lowered. This "by-pass" prevents excessive velocity of air at the work surface when the sash is nearly closed. The exhaust blowers and auxiliary fans must be so dimensioned, that it becomes possible to work always with the same air quantities, whose flow paths are regulated by the throttle valves. In such a manner an over or underpressure can be prevented. The electrical outfit must be constructed in such a manner, that the heating unit could be only switched on, if the auxiliary air fan would be running and this device should be only switchable, if the suction-fan would be already switched on. After having given to us your order to project your radiochemical laboratory pay attention please, specially to our prepositions concerning the air-circulating plant.


Abzug-Schema / Radiochemical Hood, Scheme

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
Storage-cabinet for beta-emitters
 constructed of sheet iron, inside and outside covered with plastic material (PVC). Each compartment is closed by a downward movable flap and consists of two cells, 190 mm wide, 920 mm high, 960 mm deep, each with "Placryl" door and stainless steel tray. If desired, the combination of two small cells to a greater one of 190 mm width, 530 mm height, 60 mm depth, can be arranged for storing long vessels, pipettes etc.

storage cabinet
 size 1 with 3 compartments and 6 small cells respectively with 4 small cells and 1 big cell
 875 mm width
 990 mm height
 300 mm depth
 75 kg weight
 STOR 1
 Codes: STOR 1 G

size 2 with 6 compartments and 10 small cells respectively with 8 small cells and 2 big cells
 1710 mm width
 990 mm height
 300 mm depth
 150 kg weight
 STOR 2
 Codes: STOR 2 G

Vault for Gamma-emitters
 consisting of separate cells made of sheet iron, inside covered with PVC, with massive cast iron door and one tray of stainless steel for each cell. The mounting of the vault is to be carried out by local masons on charge to the orderer


Vault:
 size 1 with 2 small cells each 920 mm wide, 940 mm high, 300 mm deep with 2 big cells each 920 mm wide, 380 mm high, 300 mm deep
 total measuring:
 size 1 875 mm wide 990 mm high 300 mm deep total weight 800 kg
 size 2 with 4 small cells each 920 mm wide, 940 mm high, 300 mm deep with 4 big cells each 920 mm wide, 380 mm high, 300 mm deep
 total measuring:
 size 2 1710 mm wide 990 mm high 300 mm deep total weight 1600 kg
 Code: BET 2

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FACILITIES FOR RADIOACTIVE WORK

Island plumbing unit

steel constructed, covered with plastic material, cross section: 60 x 60 cm, with the necessary gas, water, vacuum and waste pipes, fuse-board and electrical installation, all parts acid-proof lacquered

size 1: height 1.0 m

Codes IFS 1

Quilfit:


- 2 columns 200 mm high each with two cold water taps and one warm water tap
 - 2 sinks of acid-proof stoneware, 540 x 350 x 180 mm, for direct connection to the waste system
 - 2 columns, 90 mm high, each with two gas taps
 - 2 switch boards each with two electric sockets 220 Volt with flap, 1 socket 380 Volt with flap and ground connection
- cover plates: concrete with red tiles
plastic sheathing: black Melakart

size 2: height 2,2 m

Codes IAS 2

Quilfit:

- Sink-side:
- 1 sink of acid-proof stoneware 540 x 350 x 180 mm for direct connection to the waste system
- over that:
- 1 cold water tap in a height of 30 cm
 - 1 warm water tap in a height of 30 cm
 - 1 cold water tap in a height of 55 cm
 - 1 drainboard of plastic with pins
 - 1 board for 5 litre flasks (aqua destillata)
 - 2 shelves made of light metal, glass covered
 - 1 fluorescent tube with switch and throttle
- other sides:
- 3 cold water taps
 - 3 two-way gas taps
 - 2 vacuum taps

IAS 2

- 6 light metal shelves, glass covered
 - 3 fluorescent tubes with throttles and switches
 - 1 voltmeter 0-400 Volt
 - 2 sockets with flap and ground connection 220 volt
 - 2 sockets with flap and ground connection 380 volt
 - 2 adjustable ventilation openings for glove-box connection
 - 6 support rods
- in the interior:
- 1 central suction duct
 - 1 vacuum pump ED 3 with safety vessel
- top part: covered with white plastic Melakart
- lower part: black plastic Melakart with door and flap

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
Kurzbez.: IAG 1


Plumbing unit for wall mounting

made of wood, acid-proof lacquered, tapboard covered allround with plastic, with the necessary gas, water and waste-pipes, fuse-board and electric installation

Outfit:

- size 1: 2,8 m
- 2 luminiscent lighting tubes with choke unit
 - 4 support rods, 13 mm diam.
 - 4 water taps on columns of 300 mm height
 - 2 plastic waste funnels
 - 3 gas-columns, 90 mm high
 - 2 switches for lighting
 - 3 sockets, 220 volt, with flap
 - 1 socket, 380 volt, with flap

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Kurzbez.: IAG 2

Größe 2: 4,2 m

- 3 luminiscent tubes with choke unit
- 4 support rods, 13 mm diam.
- 6 water-taps on columns of 300 mm height
- 3 plastic waste funnels
- 4 gas-columns, 90 mm high
- 3 switches for lighting
- 3 sockets, 220 volt, with flap
- 1 socket, 380 volt, with flap

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Kurzbez.: G1 F

Standard lab benches

for radio-isotopes work consisting of a steel constructed carrier frame, 900 mm high, acid-proof lacquered, with levelled feet and connection parts for drawer and cupboard units. Max. permissible load 1800 kg with bench-top made of concrete, 1400 x 630 mm, 60 mm thick, with 1 mm stainless steel covering, with 15 mm border and edge-protection of 60 mm height. Max. permissible load 380 kg

Code: G1 V 60

bench-top like above, but 100 mm thick. Max. permissible load 700 kg

Code: G1 V 100

bench-top like above, but 150 mm thick. Max. permissible load 900 kg

Code: G1 V 150

bench-top, 1400 x 630 x 60 mm thick made of wood, covered with 4 mm black plastic material (P.V.C.), with edge-protection of plastic 60 mm high. Max. permissible load 500 kg



Code: G1 K

bench-top like above, but with 10 mm border. Max. permissible load 500 kg

Code: G1 KR

bench-top made of concrete, covered with red acid-proof tiles, joints filled by acid-proof putty "Combindal", with special border-stones. Max. permissible load 380 kg

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



H2A

Removable cupboard units

made of wood, allround acid-proof lacquered, for steel constructed carrier frames

Drawer unit

with two drawers and writing plate
Code: ZG 1

Cupboard unit

with sliding doors and adjustable wooden shelf
1,20 m length, 0,54 m high, 0,50 m deep
Code: H1A

Cupboard unit

with sliding doors, adjustable wooden shelf and two drawers
1,20 m length, 0,54 m high, 0,50 m deep
Code: H1B

Cupboard unit

with sliding doors, adjustable wooden shelf, two drawers and writing plate
1,20 m length, 0,54 m high, 0,50 m deep
Code: H1C

Cupboard unit

with trap door, adjustable wooden shelf and 5 drawers
1,20 m length, 0,54 m high, 0,50 m deep
Code: H1D

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
H2A

Removable cupboard units

made of wood, allround acid-proof lacquered, for steel constructed carrier frames

Cupboard unit

with 8 small and 2 big drawers
1,20 m length, 0,54 m high, 0,50 m deep
Code: H2E

Cupboard unit

with 4 drawers
0,50 m length, 0,54 m high, 0,50 m deep
Code: H2Fa

Cupboard unit

with trap door and movable intermediate shelf
0,50 m length, 0,54 m high, 0,50 m deep
Code: H2Fb

Cupboard unit

with 3 drawers
0,50 m length, 0,44 m high, 0,50 m deep
Code: H2Fbb

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
The illustrations on page 124/125 may explain the inclusion of a radio isotope laboratory in an existing building. The sketches are founded on the supposed existence of a corridor, running centrally or laterally through the building, at which the labs may be connected.


The suitability of the proposed designs is only guaranteed, if the proposed level is not exceeded

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

- Six-room lab
- I. hot lab IV. low level lab
 - II. active storage room V. nonactive annex
 - III. sluice with WC, shower bath, monitors etc. VI. counting room

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

Sheet iron hood
with swivelled mounted holders for plastic pails suitable for waste
ste age tanks
700 mm long, 600 mm wide, 600 mm high
applicable for openings of 590x690 mm

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

Plant for preparation of liquid waste in radio-chemical laboratories:

1. Water meter
2. Control counter
3. Acoustic signal level
4. Light signal level
5. Waste container for cover water level cap. with filter of stoneware, 52 ltr.
6. Preparation and level drain system 500 ltr. tank made of stoneware.
7. GM-tube for control of sewer water activity in the slime pit
8. Concrete base
9. Concrete protective wall
10. Distilling vessel of acid proof steel
11. Water gage
12. Filler neck for chemicals
13. Light signal unit for water level
14. Acoustic signal unit
15. Distilling apparatus with gas heating
16. Neutralization and dilution tank

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Water Still of V2A used as „Concentrator“
for the preparation of radioactive liquid waste. Our
“liquid waste container” (Page 109) is recommended
as distilling bulb and must be ordered separately.

Technical Data:
Gas consumption approx. 6 m³/h
Gas connection R 1/2"
Cooling water connection R 1/2"
Base approx. 130x60 cm
Over-all height approx. 195 cm
Weight approx. 40 kg
Code KONZ

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FACILITIES FOR RADIOACTIVE WORK

Commendable Materials for Radio Chemical Laboratories

Special synthetic resin cement „Phtalopal G“

for the production of flawless cement connections between stainless steel and plastic material, tiles, wood, glass etc. which are necessary for example in the hoods so as to warrant connection from bench top and side wall (made from plastic) without joints. It represents a thermo-plastic resin from the phthalic acid poly-alcoholic type, which has a low ability for absorbing and retaining active materials and which is gluing nearly indissolubly steel and plastic material.

Producer:
VEB Chemische Werke Buna, Schkopau

„PCD 13“

PC Glue Solution „PCD 13“ to produce flawless adhesive connections from parts made of Ekadur (polyvinylchloride hard) or Ekallit (polyvinylchloride soft).

Producer:
VEB Elektrochemisches Kombinat, Bitterfeld

Silicone Oil Emulsion „OE E 6008“

for conferring hydrophobicity on test pipettes and other measuring devices. By the application of „OE E 6008“ the liquids which were used in the measuring apparatus, loose their surface tension and do not adhere to the glass walls. The application from „OE E 6008“ is extraordinary practical in the case of radio-chemical working, because it restricts the contamination from glass apparatus to a minimum.

Producer:
VEB Chemiewerk, Nünchritz/Sa.

Stripable Lacquer „Ilmofol“ on the basis of plastic material

It is applied for the preparation of trays, special endangered places of bench tops, walls, apparatus parts etc. The material becomes hard in a very short time and can be removed as a film, from a hard ground.

Producer:
VEB Lackfabrik, Ilmenau

Dekontamination paste „SEP“

serves for elimination of contamination on rubber parts, measuring instruments, fittings, bench tops etc., caused by radioactive substances. The paste can be laid on easily by spatula or duster and prevents from smear over and spread out radioactive substances. It is not poisonous, effects no skin disease, can easily be washed off and is soluble in water.

Producer:
Pharmaz.-chem. Laboratorium „DORMED“
Geve & Co. KG., Leipzig C 1, Paul-Cruner-Str. 24

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When small amounts of these materials are demanded, then they can be forwarded by us. For larger amounts, please apply directly to the producer.

FRIEDRICH GEIER / LABORATORIUMSAPPARATE KG / ILMENAU